



Micro™ Power Options



Item#
CAPSOM505001

Transcend Micro PowerAway® Battery

Extended Run Time

- **Long-Lasting Power:** Offers 2-3 nights or 17.7 hours of run time on a single charge at a pressure setting of 9cm, ensuring restful sleep without frequent recharging.
- **Quick Charge Time:** Recharges in less than 4 hours, so it's ready to go when you are.

Reliability and Safety

- **FAA-Approved:** Cleared for in-flight use, making it an essential for travelers
- **Backup Power:** Provides a permanent backup power source during outages, ensuring uninterrupted sleep therapy.
- **Long Life Span:** Capable of 1,000+ full charge/discharge cycles for long-term reliability



Item#
CAPSOM505049

Transcend Micro Sun Watt Solar Panel Charger

Utilize solar energy to power your CPAP machine, eliminating the need for traditional power sources.


The foldable design of the solar panel and compact size of the battery make transportation hassle-free.

Key Features:

- 40W Portable Solar Panel
- 3 Charging ports: DC output, USB-A, and USB-C
- Charges fully in under 4 hours from < 5%
- If charged every day, you could replenish the battery in less than 2 hours.
- 1 year warranty

Item#
CAPSOM505035

Transcend Micro Plug Adapters (AU/EU/UK)



Travel with ease and confidence with the Transcend Plug Adapters. Designed for use with the Transcend Micro CPAP machine and its PSA5 power supply, these adapters ensure you can continue your sleep therapy wherever your journey takes you.

- **Universal Compatibility:** The Transcend Plug Adapters are designed to work with the Transcend PSA5 power supply that comes with the Transcend Micro.
- **International Use:** These adapters are designed for international use, fitting outlets in the EU, AU and UK. They allow you to stay powered up wherever your travels take you.
- **All-in-One Solution:** This package combines multiple plugs into one compact unit, eliminating the need to worry about compatibility.
- **Package Includes:** Each package contains three Transcend Plug Adapters, offering great value.

Item#
CAPSOM505010

Transcend Micro Power Supply



Replacement power supply set for Transcend Micro. This unit is a direct replacement for the original accessory, and it can quickly and easily charge the powerful P8 battery. This replacement set comes with the AC power brick and US wall plug.