

# Product Catalog



Contract PP-FA-2138  
**Medical Gas  
Equipment**



Contract PP-NS-1753  
**Medication Delivery**



Contract PP-NS-1741  
**Respiratory  
Disposables**



Contract PP-NS-1762  
**Resuscitation  
Devices**



Contract PP-OR-2217  
**Laryngeal Masks**



Contract BT-OR-019  
**Disposable  
Anesthesia  
Products**



## Sales Representatives



**Egle Sventnickaite**  
872-258-1654  
esventnickaite@sunsethcs.com



**Mason Moss**  
262-788-1862  
mmoss@sunsethcs.com



**Owen Ruth**  
312-533-2451  
owenr@sunsethcs.com



**Cindy Geighes**  
312-533-2478  
cgeighes@sunsethcs.com




**Hunter Baudino**  
262-788-1855  
hbaudino@sunsethcs.com




**Wyatt Kanemoto**  
262-788-8941  
wkanemoto@sunsethcs.com

 *Director, Strategic Accounts*  
**Lisa Iannazzone** RN, BSN  
636-222-3134  
lisa@sunsethcs.com

 *Manager, Product & Business Development*  
**Tom Munar**  
312-533-2476  
tmunar@sunsethcs.com



 *Contract Manager, Strategic Accounts*  
**Gina Rosenberg**  
914-837-5215  
grosenberg@sunsethcs.com

 *RRT Clinical Specialist*  
**Krystle Bodine** BS, RRT  
kbodine@sunsethcs.com

## Corporate Office

Phone 872-275-3210

Fax 312-997-9985

141 W Jackson Blvd Ste 1950  
Chicago IL 60604

## Central Returns/Distribution Center

2725 Alft Ln  
Elgin IL 60124

## West Returns/Distribution Center

6750 W Van Buren St  
Phoenix AZ 85043

# Content

## Contract PP-FA-2138

### Medical Gas Equipment Page

Microflow Microdial Flow Meter	4
--------------------------------	---

## Contract PP-NS-1741

### Respiratory Disposables Page

Lynx Venturi Oxygen Delivery	5
Spirometers	6
Drainage Bags	6
Bacteria Filters	7
Breathing Circuits	8-10
Pediatric NIV Masks	10
CPAP and Oxygen Filters	11-12
Oxygen Cannula and Accessories	13-14
Oxygen Masks	15-16
Air Entrainment Masks	16
Oxygen Supply Tubing	17
Firesafe	17
Oxygen Connectors	18
Oxygen Accessories	19
Berman + Guedel Airways	20
HMEs	21
Trach Masks and Adapters	21
Suction Catheter Kits	22
Suction Catheters	23
Trach Ties	23
Connectors	23

## Contract PP-NS-1753

### Medication Delivery Page

Inhaler Chamber	24
Nebulizer Compressors	24
Peak Flow Meter	24
Adapters	24
Nebulizer Kits	25
Aerosol Masks	25

## Contract PP-NS-1762

### Resuscitation Devices Page

Resuscitators	26
Manometers	27
PEEP Valves	27
Reservoir Bags	28
Pocket Resuscitator Set	28
Breathing Masks	28

## Contract PP-OR-2217

### Laryngeal Masks Page

Silicone Disposable	29
Silicone Reusable	29
Flexible Silicone Disposable	29
Flexible Silicone Reusable	29
PVC Disposable	30
Flexible PVC Disposable	30

## Contract BT-OR-019

### Disposable Anesthesia Products Page

PeDIA	31
-------	----

## Microflow Microdial Flowmeter

2024 Innovation Celebration Winner



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure Range	Flow Range*
<b>RES8160048</b>	345 – 500 kPa	<b>A</b>
<b>RES8160050</b>	350 – 550 kPa	<b>E</b>



DISS Inlet • DISS Outlet

Item#	Inlet Pressure Range	Flow Range*
<b>RES8160049</b>	345 – 500 kPa	<b>A</b>
<b>RES8160051</b>	250 – 450 kPa	<b>E</b>

\*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
<b>Range A</b>	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
<b>Range E</b>	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

## Assured, repeatable accuracy to enable safe, effective weaning.

### What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.

### Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned



## Lynx Venturi Oxygen Delivery


### Low Concentration

Item# **RES2400**

Each

Adapter Only **RES240LC**

- Design replaces multiple venturi adapters
- Connects to oxygen tubing and regulates flow
- Provides precision FiO<sub>2</sub> delivery from 24-60% with less than 2% variance during mild, moderate, and severe respiratory distress
- Excellent CO<sub>2</sub> washout



### Unprecedented simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

### Do no harm

Quicker de-escalation of O<sub>2</sub> levels and prevention of excessive CO<sub>2</sub> in the patient.

### Reduction in waste

Remains with patient throughout EMS and hospital journey, replacing multiple adapters which may get lost or contaminated during transition.

### High Concentration

Item# **RES2401**


Each

Adapter Only **RES241HC**

- Allows rapid de-escalation of FiO<sub>2</sub> in the EMS ambulance or ED to prevent CO<sub>2</sub> buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia (PaO<sub>2</sub>/FiO<sub>2</sub>) allows early intervention and better outcomes
- Excellent CO<sub>2</sub> washout

Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Capabilities could include replacing non-rebreathers and/or partial rebreather products



### Unmatched precision

Provides specific levels of oxygen with less than 2% variance. Comparable products vary by as much as 30%.

## Volumetric Incentive Spirometer

12/case

Item#  
**RES055**  
2500ml

**RES056**  
5000ml

**RES055P**  
Pediatric



SCAN FOR VIDEO



- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise
- Dual-sided, clear and legible identification of volumes
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration
- Compact and ergonomic design

## Three Ball Spirometer

Item#  
**RES057-10**  
10/pack

**RES057**  
12/case

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient



- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

## Drainage Bag with Y-Adapter

Item# **RES071**  
750cc  
50/case



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

## Drainage Bag with T-Adapter

Item# **RES070**  
750cc  
50/case



### Final Bacteria Filter

Item# **BF516**

Final Bacteria Filter  
1/4in Barbed

**10/pack**

OEM# P0750X030



### Final Bacteria Filter — 3/16in

Item# **BF316**

3/16 Stepped Barb and 1/4in  
Straight Barb  
on Both Sides

**10/pack**

**BF316-1** 1/pack



### Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item#

**BF2200**

22mm Bacteria  
Disk Filter

**5/pack**



Attach between the  
machine and tube for  
additional protection  
against bacteria and  
allergens

#### Filter Specifications

Bacterial Filtration Efficiency: (BFE) 99.99+%

Viral Filtration Efficiency: (VFE) 99.9+%

NaCl Particulate Filtration Efficiency: 96+%

Flow Resistance: Approximately 8.8 mmH2O at 30 LPM

Internal Volume: 66.3 mL without fittings / 78.3 mL with fittings

Connections:

• Patient End: 22 mm OD / 15 mm ID

• Machine End: 22 mm ID

### Bacterial/Viral HME Filter – Adult

Item# **BF104**

**Each**



- Bacterial efficiency 99.9999%
- Viral efficiency 99.999%
- Humidification efficiency 37.5 mg/l H2O (500 ml tidal volume)
- Coupled with a hydrophobic bacterial viral filter membrane



## Passive Non-Heated Breathing Circuits

For use with Respiration Trilogy  
**Each**

Exhalation Port	Item#	OEM#
Swivel	<b>VC1000</b>	1073228 & 1073236
Fixed	<b>VC1100</b>	1069467 & 1069468



## Adult Passive Heated Circuit with Exhalation Port

**30/case**

For use with Philips  
Respiration Trilogy

Item#  
**VC1003**


**VC1003P**   
Pediatric Size

Comparable to OEM# 1076046,  
1073227, 1073235



### Includes:

- 1 VC1114 Passive Heated Breathing Circuit
- 1 VEN002 Swivel Exhalation Port



## Passive Circuit w/Exhalation Port and Bacterial/Viral Filter

With 22mm cuffs

Item#  
**VC1005**  
**Each**

OEM# 1147378, 1132340, 7615



## Passive Heated Breathing Circuit

Item#  
**VC1114**  
For use with  
MR850 Humidifier  
**10/case**

**VC1114P**  
Pediatric Size



## Swivel Connector Silencer

Item#  
**VEN002**  
**Each**





## Non-Heated Breathing Circuits



For use with Respiration Trilogy  
22mm cuffs

**30/case**

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
<b>VC1001</b>	1073223 & 1073232	Active	No	Adult	22mm
<b>VC1004</b>	1069443 & 1069444	Passive	Yes	Adult	22mm
<b>VC1004P</b>	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#  
**VC1001**  
Active Adult

OEM#  
1073223  
1073232




Item#  
**VC1004**  
Passive Adult

OEM#  
1069443  
1069444




To use active circuits with Respiration Trilogy, simply cut ends off pressure and exhalation lines.



Item#  
**VC1004P**  
Passive Pediatric

OEM#  
1073224  
1073233



## Adult NIV Patient Circuit


Item# **VC1160**

Compatible with V60, No Filters

**10/case**

Comparable to

OEM# 1069210



## Adult NIV Single Limb Patient Circuit with Main Flow Filter

Item# **VC1161**

Includes 1 BF2200 22mm Bacteria Disk Filter

**10/case**




## Adult NIV Single Limb Patient Circuit with Main Flow Filter and Exhalation Port Filter

Item# **VC1162**

Includes 2 BF2200 22mm Bacteria Disk Filters

**10/case**



## MiniMe Pediatric Nasal Mask



Item# **CM004M**

HCPCS# A7034

With Medium Headgear

**Each**

OEM# 60214



## MiniMe 2 Pediatric Nasal Mask

Item#

**CM027M**

Mask with Headgear

• Non-vented  
version available

OEM# 60218

HCPCS# A7034



## ResMed S9 and AirSense 10 Series Filters

### Item# **CF2107**

Disposable Filter

**3/pack** • 2<sup>1</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>8</sub>"

HCPCS# A7038

**CF2107-1** 1/pack OEM# 36850

**CF2107-2** 2/pack OEM# 36851

**CF2107-50** 50/pack OEM# 36853

### Item# **CF2107HA**

Hypoallergenic Filter

**3/pack** • 2<sup>1</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>8</sub>"

HCPCS# A7038

**CF2107HA-1** 1/pack OEM# 36855

**CF2107HA-2** 2/pack OEM# 36856

**CF2107HA-50** 50/pack  
OEM# 36858



## ResMed AirSense 11 Series Filters

### Item# **CF2111**

Disposable Filter

**3/pack**

HCPCS# A7038

**CF2111-1** 1/pack OEM# 39300

**CF2111-2** 2/pack OEM# 39301

**CF2111-50** 50/pack OEM# 39303

### Item# **CF2111HA-3**

Hypoallergenic Filter

**3/pack**

1<sup>1</sup>/<sub>4</sub>" x 1<sup>3</sup>/<sub>16</sub>"

HCPCS# A7038

**CF2111HA-1** 1/pack

**CF2111HA-2** 2/pack



## 3B Luna G3 Style Filters

### Item# **CF8006-6**

Disposable Filter

**6/pack**

HCPCS# A7038

**CF8006-1** 1/pack

**CF8006-2** 2/pack

### Item# **CF8006F-1**

Reusable Filter

**1/pack**

HCPCS# A7039



## 3B Luna and Luna II Style Filters

### Item# **CF8005**

Disposable Filter

**12/pack** • 2" x 1<sup>3</sup>/<sub>16</sub>"

HCPCS# A7038

### Item# **CF8005F-1**

Reusable Filter

**1/pack** • 2" x 1<sup>3</sup>/<sub>16</sub>"

HCPCS# A7039

**CF8005F-2** 2/pack




## 3B Stratus 5 Style Filter

### Item# **OF6030**

Foam Cabinet Filter

**10/pack**

OEM# STR1030



## Respironics DreamStation 2 Series Filters

Item# **CF1012**  
Disposable Filter  
**6/pack**

OEM# 1142831

HCPCS# A7038

**CF1012-1** 1/pack OEM# 1142828

**CF1012-2** 2/pack OEM# 1142829

**CF1012-3** 3/pack

**CF1012-50** 50/pack

Item# **CF1012F**  
Reusable Filter  
**2/pack**

OEM# 1142687

HCPCS# A7039

**CF1012F-1** 1/pack



## Respironics DreamStation Series Filters

Item# **CF1010**  
Disposable Filter  
**6/pack**

OEM# 1122519

HCPCS# A7038

**CF1010-1** 1/pack

**CF1010-2** 2/pack

Item# **CF1010F**  
Reusable Filter  
**2/pack**

OEM# 112446 (1-pack)

HCPCS# A7039

**CF1010F-1** 1/pack

**CF1010F-6** 6/pack



## EverFlo Kit

Item# **OF1003SP**  
1 Intake Filter  
1 Final Bacteria Filter  
**Each**



## SolAiris V Kit

Item# **OF5005SP**  
2 Foam Cabinet Filters  
1 Extended Life Bacteria Filter  
w/Muffler  
1 Final Bacteria Filter  
**Each**



## Mobilaire 3 and 5 Kit

Item# **OF9001SP**  
2 Foam Cabinet Filters  
1 HEPA Filter 6in Capsule  
1 Final Bacteria Filter  
**Each**



## Perfecto2 Kit

Item# **OF9004SP**  
2 Foam Cabinet Filters  
1 HEPA Filter 6in Capsule  
1 Final Bacteria Filter  
**Each**



## Filter for ResMed Astral Ventilator

Item#  
**VF003-4**  
**4/pack**



Item# **VF003-H**  
Filter Housing  
**Each**



Item# **VF003-KIT**  
4-Filter Kit  
w/Filter Housing  
OEM# 27939



## Humidoflo High Flow Nasal Cannula

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respironics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)



### 10/case

Size	Item#
Small	<b>HFN001AS</b>
Medium	<b>HFN001AM</b>
Large	<b>HFN001AL</b>

### With F&P Adapter

Size	Item#
Small	<b>HFN001ASA</b>
Medium	<b>HFN001AMA</b>
Large	<b>HFN001ALA</b>

## Adult Oxygen Cannula

Item#	Tubing Length	Connector	Qty
<b>RES1100</b>	None	Standard	50
<b>RES1100N</b>	None	Standard	50
<b>RES1104</b>	4 ft	Standard	50
<b>RES1107</b>	7 ft	Standard	50
<b>RES1107N</b>	7 ft	Standard	50
<b>RES1107NU</b>	7 ft	Universal	50
<b>RES1107U</b>	7 ft	Universal	50
<b>RES1110</b>	10 ft	Standard	50
<b>RES1114N</b>	14 ft	Standard	50
<b>RES1114U</b>	14 ft	Universal	50
<b>RES1115</b>	15 ft	Standard	50
<b>RES1115U</b>	15 ft	Universal	50
<b>RES1125</b>	25 ft	Standard	25
<b>RES1135</b>	35 ft	Standard	20
<b>RES1150U</b>	50 ft	Universal	20



- 1-6 LPM
- 6-channel supply tubing for safety

## Soft Cannula

Item#	Tubing Length	Connector	Qty
<b>RES1100S</b>	None	Standard	50
<b>RES1107S</b>	7 ft	Standard	50
<b>RES1107SU</b>	7 ft	Universal	50
<b>RES1114S</b>	14 ft	Standard	50
<b>RES1114SU</b>	14 ft	Universal	50
<b>RES1115S</b>	15 ft	Standard	50
<b>RES1125S</b>	25 ft	Standard	25
<b>RES1125SU</b>	25 ft	Universal	25



- 1-6 LPM
- 6-channel supply tubing for safety



## Pediatric, Infant, and Neonate Cannula

50/case

Item#	Size	Tubing Length	Connector	Touch	Ear Cushions
<b>RES1100P</b>	Pediatric	None	Standard	Standard	No
<b>RES1107I</b>	Infant	7 ft	Standard	Standard	No
<b>RES1107NE</b>	Neonate	7 ft	Standard	Standard	No
<b>RES1107P</b>	Pediatric	7 ft	Standard	Standard	No
<b>RES1107PEC</b>	Pediatric	7 ft	Standard	Standard	Yes
<b>RES1107PECU</b>	Pediatric	7 ft	Universal	Standard	Yes
<b>RES1107SI</b>	Infant	7 ft	Standard	Soft	No
<b>RES1107SIU</b>	Infant	7 ft	Universal	Soft	No
<b>RES1107SNE</b>	Neonate	7 ft	Standard	Soft	No
<b>RES1107SP</b>	Pediatric	7 ft	Standard	Soft	No
<b>RES1107SPU</b>	Pediatric	7 ft	Universal	Soft	No
<b>RES1114P</b>	Pediatric	14 ft	Standard	Standard	No
<b>RES1114PU</b>	Pediatric	14 ft	Universal	Standard	No
<b>RES1114SP</b>	Pediatric	14 ft	Standard	Soft	No
<b>RES1114SPU</b>	Pediatric	14 ft	Universal	Soft	No

- 1-6 LPM
- 6-channel supply tubing for safety

## Adult Cannula with Ear Cushions

Item#	Tubing Length	Connector	Qty
<b>RES1100EC</b>	None	Standard	50
<b>RES1104EC</b>	4 ft	Standard	50
<b>RES1107EC</b>	7 ft	Standard	50
<b>RES1107ECU</b>	7 ft	Universal	50
<b>RES1114EC</b>	14 ft	Standard	50
<b>RES1115EC</b>	15 ft	Standard	50
<b>RES1125EC</b>	25 ft	Standard	25

## Ear Cushions for Cannula

Item# **RES020**  
Cushions Only  
50 Pairs/case





- Add cushions to existing cannula

## Secure Cannula Stickers

Item# **RES080**  
25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs





## Simple O<sub>2</sub> Mask

With 7ft Safety Tubing  
Elongated

Item# **RES2100**  
Adult Size  
50/case



- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

## Medium Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector
<b>RES2100NT</b>	Adult	None	Standard
<b>RES2100P</b>	Pediatric	7 ft	Standard
<b>RES2100U</b>	Adult	7 ft	Universal
<b>RES2203I</b>	Infant	None	Standard
<b>RES2203P</b>	Pediatric	None	Standard
<b>RES2203A</b>	Adult	None	Standard
<b>RES2203A7</b>	Adult	7 ft	Standard
<b>RES2203I7</b>	Infant	7 ft	Standard
<b>RES2203P7</b>	Pediatric	7 ft	Standard

- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



## High Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length
<b>RES2202A</b>	Adult	None
<b>RES2202A7</b>	Adult	7 ft
<b>RES2202I</b>	Infant	None
<b>RES2202I7</b>	Infant	7 ft
<b>RES2202P</b>	Pediatric	None
<b>RES2202P7</b>	Pediatric	7 ft



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

## Non-Rebreather Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	Safety Vent
<b>RES2102</b>	Adult	7 ft	Standard	No
<b>RES2102I</b>	Infant	7 ft	Standard	No
<b>RES2102IM</b>	Infant	7 ft	Standard	No
<b>RES2102P</b>	Pediatric	7 ft	Standard	No
<b>RES2102PS</b>	Pediatric	7 ft	Standard	Yes
<b>RES2102PSM</b>	Pediatric	None	Standard	Yes
<b>RES2102PSU</b>	Pediatric	7ft	Universal	Yes
<b>RES2102PU</b>	Pediatric	7 ft	Universal	No
<b>RES2102S</b>	Adult	7 ft	Standard	Yes
<b>RES2102SM</b>	Adult	None	Standard	Yes
<b>RES2102U</b>	Adult	7 ft	Universal	No

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



## Adult Air Entrainment Mask

With Flex Tube

50/case

Item#	Size	Kit	Tubing Length	Connector
<b>RES2400A</b>	Adult	Mask Only	7 ft	Standard
<b>RES2400AK</b>	Adult	Yes	7 ft	Standard
<b>RES2400AKU</b>	Adult	Yes	7 ft	Universal
<b>RES2400AU</b>	Adult	Mask Only	7 ft	Universal
<b>RES2400P</b>	Pediatric	Mask Only	7 ft	Standard
<b>RES2400PK</b>	Pediatric	Yes	7 ft	Standard
<b>RES2400PKU</b>	Pediatric	Yes	7 ft	Universal
<b>RES2400PU</b>	Pediatric	Mask Only	7 ft	Universal



## Oxygen Supply Tubing

Item#	Length	Connector	Color	Qty
<b>RES3007</b>	7 ft	Standard	None	50
<b>RES3007G</b>	7 ft	Standard	Green	50
<b>RES3007U</b>	7 ft	Universal	None	50
<b>RES3015</b>	15 ft	Standard	None	50
<b>RES3025</b>	25 ft	Standard	None	25
<b>RES3025G</b>	25 ft	Standard	Green	25
<b>RES3025U</b>	25 ft	Universal	None	25
<b>RES3040G</b>	40 ft	Standard	Green	20
<b>RES3050G</b>	50 ft	Standard	Green	20

- 6-channel
- PHT/DEHP-free



**Universal Connector**

## FireSafe Nozzle

Item#


**RES010FND**

**10/pack**



- A thermal fuse designed to stop the flow of oxygen in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire

**Firesafe™  
Nozzle**



## Christmas Tree Connector

### Item# RES001

Green without Swivel  
10/pack

**RES001-30**  
30/pack



### Item# RES002

Green with Swivel  
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

## Oxygen Tubing Connector — Non-Swivel

### Item# RES012

Fits 3/16in and 1/4in  
ID Tubing  
10/pack



- Barbed ends ensure secure, leak-free connections

## Tubing Reducer Connector

### Item# RES011

Fits 3/16in, 1/4in,  
5/16in, 3/8in, and 7/16in  
ID Tubing  
10/pack



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

## Humidifier Tubing Adapter

### Item# RES016

22mm ID x 22mm OD  
with Lip  
10/pack



## Circuit Connectors

### Item# RES013

22mm OD x  
15mm ID/  
22mm OD  
10/pack



### Item# RES014

22mm OD x  
22mm OD  
With Lip  
10/pack



## Oxygen Tubing Connector — Swivel

### Item# RES018

Male/Male  
10/pack

**RES018-1** 1/pack  
**RES018-25** 25/pack



## Oxygen Bleed-In Adapter

### Item# RES019

22mm with 6mm Rotating Stem  
10/pack

**RES019-1** 1/pack



## Oxygen Tubing Y-Connector

### Item# RES017

Fits 3/16in and 1/4in ID Tubing  
(All 3 sides are the same)  
10/pack

- Connects two oxygen sources for higher flows



## Connector with 6mm Oxygen Stem

### Item# RES058

22mm with 6mm Stem  
10/pack



## Finger Pulse Oximeter

Item# **RES5102**  
Each



- Oversized 1.3" OLED display face
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- Includes lanyard and manual

## Wrenches

Item# **RES031**  
With Chain  
10/pack  
**RES031-1** 1/pack



## Water Traps

Item# **RES037**  
Push On End Caps  
Each



- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

## Oxygen Liter Meter

Item# **RES034**  
0-10 Liter/Minute  
Each



- Oxygen pen
- Works with tanks and concentrators

## Humidifier Connectors

Item#  
**HUM008**  
L-Connector 12in  
10/pk  
**HUM008-1** 1/pack

**HUM008L**  
L-Connector 21in  
10/pk

**HUM008S**  
Straight Connector 9in  
10/pk



## Covers for Concentrators



Item# **EC251530C**  
25×15×30"  
1.5mil • Clear  
250 bags/roll

## Color Coded Dual Channel Berman Oral Airway

50/pack

Item#	Size	Color
<b>A0A40BC-50</b>	40mm	Pink
<b>A0A50BC-50</b>	50mm	Blue
<b>A0A60BC-50</b>	60mm	Black
<b>A0A70BC-50</b>	70mm	White
<b>A0A80BC-50</b>	80mm	Green
<b>A0A90BC-50</b>	90mm	Yellow
<b>A0A100BC-50</b>	100mm	Red
<b>A0A110BC-50</b>	110mm	Orange



- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

## Oropharyngeal Guedel Airway

10/pack

Item#	Size	Color
<b>A0P500GC</b>	50mm	Blue
<b>A0P601GC</b>	60mm	Black
<b>A0P702GC</b>	70mm	White
<b>A0P803GC</b>	80mm	Green
<b>A0P904GC</b>	90mm	Yellow
<b>A0P1005GC</b>	100mm	Red



- Designed to overcome associated risks of loose or detached bite blocks
- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free



## Flex HME






Item# **RES027FLEX**  
40/case



<b>TV Range (mL)</b>	300-1500
<b>Moisture Output (mg H<sub>2</sub>O/L)</b>	33 mg/L
<b>Weight (g)</b>	17
<b>Resistance (cm H<sub>2</sub>O)</b>	0.5 @ 30 LPM 1.34 @ 60 LPM
<b>Dead Space (mL)</b>	75
<b>Connection (mm)</b>	22M,15F

## Heat Moisture Exchangers

30/case

					
<b>Item#</b>	<b>RES027</b> T-HME	<b>RES027OP</b> T-HME with Oxygen Port	<b>RES027IL</b> Inline HME with Gas Sampling Port	<b>RES027NE</b> Inline HME — Neonatal	<b>RES027HV2</b> Foam HME — Adult/Child
<b>TV Range (mL)</b>	250-1000	150-1250	70-600	15-50	250-1000
<b>Moisture Output (mg H<sub>2</sub>O/L)</b>	30, Vt = 0.1L	36.3 mg/L	34.5 at VT = 250 ml 32.5 at VT = 500 ml	30	29
<b>Weight (g)</b>	4.5	7.4	10	4.5	21
<b>Resistance (cm H<sub>2</sub>O)</b>	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.2 @ 30 LPM	0.8 @ 30 LPM 2.1 @ 60 LPM	0.5 @ 10 LPM	0.52 @ 30 LPM 1.32 @ 60 LPM 2.54 @ 90 LPM
<b>Dead Space (mL)</b>	12.5	7	21	2.5	29
<b>Connection (mm)</b>	15F	15F 5.5mm O <sub>2</sub> port	22M/15F-15M	15M,15F	22F/15M, 22M/15F

## Adapter for T-HME

Item#  
**RES027A**  
Each

**RES027A-30** 30/pack



## Flexible Connector

50/case

Item# **TRAA004P**  
Pediatric Size  
15mm F x 15mm M  
Length 4.9-7.5cm  
OEM# 3215



## Tracheostomy Mask

50/case

Item#  
**RES2130**  
Adult Regular Size

**RES2130L**  
Adult Large Size

**RES2230**  
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O<sub>2</sub> supply tubing

## Suction Catheter Kits

### 50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	<b>TRATK305</b>
6fr	Y	DeLee	Pediatric	<b>TRATK306</b>
8fr	Y	DeLee	Pediatric	<b>TRATK308</b>
8fr	Y	Whistle	Adult	<b>TRATK108</b>
8fr	Chimney	Whistle	Adult	<b>TRATK108CV</b>
10fr	Y	Whistle	Adult	<b>TRATK110</b>
10fr	Chimney	Whistle	Adult	<b>TRATK110CV</b>
12fr	Y	Whistle	Adult	<b>TRATK112</b>
12fr	Chimney	Whistle	Adult	<b>TRATK112CV</b>
14fr	Y	Whistle	Adult	<b>TRATK114</b>
14fr	Chimney	Whistle	Adult	<b>TRATK114CV</b>
16fr	Y	Whistle	Adult	<b>TRATK116</b>
16fr	Chimney	Whistle	Adult	<b>TRATK116CV</b>

- Sterile
- Includes Catheter and Gloves

## Suction Catheter Kits with Cup

### 50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	<b>TRATK405</b>
6fr	Y	DeLee	Pediatric	<b>TRATK406</b>
8fr	Y	DeLee	Pediatric	<b>TRATK408</b>
8fr	Y	Whistle	Adult	<b>TRATK208</b>
8fr	Chimney	Whistle	Adult	<b>TRATK208CV</b>
10fr	Y	Whistle	Adult	<b>TRATK210</b>
10fr	Chimney	Whistle	Adult	<b>TRATK210CV</b>
12fr	Y	Whistle	Adult	<b>TRATK212</b>
12fr	Chimney	Whistle	Adult	<b>TRATK212CV</b>
14fr	Y	Whistle	Adult	<b>TRATK214</b>
14fr	Chimney	Whistle	Adult	<b>TRATK214CV</b>
16fr	Y	Whistle	Adult	<b>TRATK216</b>
16fr	Chimney	Whistle	Adult	<b>TRATK216CV</b>


- Sterile
- Includes Catheter and Gloves



## Suction Catheters

Size	Valve Type	Tip	Intended Use	Packaged	Item#	Qty
5fr	Y	DeLee	Pediatric	Coiled	<b>TRASC305C</b>	50
5fr	Y	DeLee	Pediatric	Straight	<b>TRASC305S</b>	100
6fr	Y	DeLee	Pediatric	Coiled	<b>TRASC306C</b>	50
6fr	Y	DeLee	Pediatric	Straight	<b>TRASC306S</b>	100
8fr	Y	DeLee	Pediatric	Coiled	<b>TRASC308C</b>	50
8fr	Y	DeLee	Pediatric	Straight	<b>TRASC308S</b>	100
10fr	Y	Whistle	Adult	Coiled	<b>TRASC110C</b>	50
10fr	Y	Whistle	Adult	Straight	<b>TRASC110S</b>	100
12fr	Y	Whistle	Adult	Coiled	<b>TRASC112C</b>	50
12fr	Y	Whistle	Adult	Straight	<b>TRASC112S</b>	100
14fr	Y	Whistle	Adult	Coiled	<b>TRASC114C</b>	50
14fr	Y	Whistle	Adult	Straight	<b>TRASC114S</b>	100
16fr	Y	Whistle	Adult	Coiled	<b>TRASC116C</b>	50
16fr	Y	Whistle	Adult	Straight	<b>TRASC116S</b>	100

- Sterile
- Single patient use
- Latex-free



## Tracheostomy Ties with Fastener Tabs

### 20/case


#### One-piece

- Easy on

#### Item#

**RES029** Adult

**RES029P** Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton


#### Two-piece

- Easy off

#### Item#

**RES030** Adult

**RES030P** Pediatric



## Tracheostomy Swivel Elbow Connectors

#### Item# **RES021**

With suction port  
15mm ID x 15mm OD  
**10/pack**



## SperaChamber

### Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler

Item#  
**RES195-10**  
10/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free

**RES195-50**  
50/case



### SperaChamber with Mask 10/case

Item#	Size
<b>RES195MS</b>	Small (Infant)
<b>RES195MM</b>	Medium (Pediatric)
<b>RES195ML</b>	Large (Adult)

### Mask Only

Item#	Size	Qty
<b>RES195MASK-S</b>	Small (Infant)	Each
<b>RES195MASK-S-10</b>	Small (Infant)	10/box
<b>RES195MASK-M</b>	Medium (Pediatric)	Each
<b>RES195MASK-M-10</b>	Medium (Pediatric)	10/box
<b>RES195MASK-L</b>	Large (Adult)	Each
<b>RES195MASK-L-10</b>	Large (Adult)	10/box

Adult Mask

## Nebulizer Compressor

Item#  
**NEB100**

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Foam Filters (BF018 x5)

**NEB100-R**

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



### Features

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

### Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

SCAN FOR VIDEO



## Peak Flow Meter

Item# **RES052**  
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**  
Pediatric Mouthpiece for Peak Flow Meter  
10/case

## Valved T-Adapter, Aerosol

Item# **NEB9001T**  
22mm ID x 22mm OD  
10/case



## Valved T-Adapter — Adult

Item#  
**VEN001**  
22F-18F-22M/15F  
Each



## Plastic Intake Filter with Head For Respironics Inspiration 626 & 929

Item# **BF016**  
10/pack



## Disposable Nebulizer Kit

Item#

**RES090**

Adult Size

**50/case**

- 6cc jet nebulizer bottle
- Mask
- 7ft supply tube

**RES090A**

No Accessories

**50/case**

**RES090NM**

- 6cc jet nebulizer bottle
- 7ft supply tube

**50/case**

**RES092**

Pediatric Size

**50/case**

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube



## Disposable Nebulizer Kit with T-Piece

**50/case**

Item#

**RES091**

Disposable  
Single Patient

**RES091U**

Disposable  
Single Patient,  
Universal Connector

**RES091NA**

Disposable  
Single Patient,  
No Aerosol

HCPCS# A7003



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

### Includes:

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube

## Reusable Nebulizer Kit with T-Piece

Item# **RES093**

Reusable —  
Single Patient

**10/case**

HCPCS# A7005

- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°



### Includes:

- 10cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube

## Disposable Nebulizer Kit with Angled Mouthpiece

**50/case**

Item#

**RES096**

Disposable  
Single Patient

HCPCS# A7003

**RES096U**

Disposable  
Single Patient,  
Universal Connector



### Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube

## Aerosol Mask

**50/case**

Item#

**RES2110**

Adult Size

**RES2110F**

Adult Size  
w/ Flex Tube

**RES2210**

Pediatric Size

- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



## Adult Resuscitator

12/case

Item#	Pop-Off	Mask
<b>RES8100</b>	No	None
<b>RES8100H</b>	Yes	None
<b>RES8101</b>	Yes	#4
<b>RES8102</b>	Yes	#5



## Child Resuscitator

12/case

Item#	Pop-Off	Mask
<b>RES8110H</b>	Yes	None
<b>RES8111</b>	Yes	#2
<b>RES8112</b>	Yes	#3



## Infant Resuscitator

12/case

Item#	Pop-Off	Mask
<b>RES8120H</b>	Yes	None
<b>RES8121</b>	Yes	#0
<b>RES8122</b>	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve



### Disposable Pressure Manometer

Item# **RES8201**  
25/case



- Pressure range 0-60 cm H2O
- Indicator strip
- All models have the pressure range emboss printed on the housing
- Some users like to also have a color strip / indicator

### Disposable Pressure Manometer, L-Shaped

Item# **RES8201L**  
25/case



- Pressure range 0-60 cm H2O
- Indicator strip
- All models have the pressure range emboss printed on the housing
- Some users like to also have a color strip / indicator

### Reusable PEEP Valve, Orange

2-10cm H2O  
10/case

Item#  
**RES8203R**  
15mm ID



**RES8203RL**  
30mm ID

- Available in orange color in 2-10cm H2O and blue color in 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Materials: Polycarbonate, silicone, stainless steel for spring

### Reusable PEEP Valve, Blue

5-20cm H2O  
10/case

Item#  
**RES8204B**  
15mm ID



**RES8204BL**  
30mm ID

- Available in orange color in 2-10cm H2O and blue color in 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Materials: Polycarbonate, silicone, stainless steel for spring

### Disposable PEEP Valve, Orange

2.5-10cm H2O  
10/case

Item#  
**RES8205R**  
15mm ID



**RES8205RL**  
30mm ID

- Available in 2.5-10cm H2O and 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Single patient use
- Materials: Polycarbonate, silicone, stainless steel

### Disposable PEEP Valve, Blue

5-20cm H2O  
10/case

Item#  
**RES8205B**  
15mm ID



**RES8205BL**  
30mm ID

- Available in 2.5-10cm H2O and 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Single patient use
- Materials: Polycarbonate, silicone, stainless steel

## Reservoir Bag, PVC

25/case

Item#	Size
<b>RES8300P</b>	0.5L
<b>RES8301P</b>	1L
<b>RES8302P</b>	2L



## Reservoir Bag, Silicone

25/case

Item#	Size
<b>RES8300S</b>	0.5L
<b>RES8301S</b>	1L
<b>RES8302S</b>	2L



Suitable for reusable resuscitator.

## CPR Pocket Resuscitator Set

Item# **RES8200**

12/case

- High transparency for better visualization
- Non-reusable, single use only
- 100% Medical level PVC material
- Latex-free components
- Economical choice
- Blow molding



## Breathing Mask

20/case

Item#	Size	Hooks
<b>SBC01N</b>	Neonatal (#0)	No
<b>SBC01I</b>	Infant (#1)	No
<b>SBC01P</b>	Pediatric (#2)	No
<b>SBC01S</b>	Adult Small (#3)	No
<b>SBC01M</b>	Adult Medium (#4)	No
<b>SBC01L</b>	Adult Large (#5)	No
<b>SBC01XL</b>	Adult Extra Large (#6)	No
<b>SBC02N</b>	Neonatal (#0)	Yes
<b>SBC02I</b>	Infant (#1)	Yes
<b>SBC02P</b>	Pediatric (#2)	Yes
<b>SBC02S</b>	Adult Small (#3)	Yes
<b>SBC02M</b>	Adult Medium (#4)	Yes
<b>SBC02L</b>	Adult Large (#5)	Yes
<b>SBC02XL</b>	Adult Extra Large (#6)	Yes



## Silicone Disposable Laryngeal Masks

Item#	Size
<b>ALM1001</b>	#1
<b>ALM1015</b>	#1.5
<b>ALM1002</b>	#2
<b>ALM1025</b>	#2.5
<b>ALM1003</b>	#3
<b>ALM1004</b>	#4
<b>ALM1005</b>	#5

- Medical grade silicone for super compatibility and patient comfort
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters
- Available with seven sizes for all patient use



## Silicone Reusable Laryngeal Masks

Item#	Size
<b>ALM1001R</b>	#1
<b>ALM1015R</b>	#1.5
<b>ALM1002R</b>	#2
<b>ALM1025R</b>	#2.5
<b>ALM1003R</b>	#3
<b>ALM1004R</b>	#4
<b>ALM1005R</b>	#5

- Reusable Silicone Laryngeal Masks are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Serial no. is printed on each connector for easy tracking
- Additional user card for record of device usage times and patient information is supplied in each pouch



## Flexible Silicone Disposable Laryngeal Masks

Item#	Size
<b>ALM3001</b>	#1
<b>ALM3015</b>	#1.5
<b>ALM3002</b>	#2
<b>ALM3025</b>	#2.5
<b>ALM3003</b>	#3
<b>ALM3004</b>	#4
<b>ALM3005</b>	#5

- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness



## Flexible Silicone Reusable Laryngeal Masks

Item#	Size
<b>ALM3001R</b>	#1
<b>ALM3015R</b>	#1.5
<b>ALM3002R</b>	#2
<b>ALM3025R</b>	#2.5
<b>ALM3003R</b>	#3
<b>ALM3004R</b>	#4
<b>ALM3005R</b>	#5

- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness



### PVC Disposable Laryngeal Masks

Item#	Size
<b>ALM2001</b>	#1
<b>ALM2015</b>	#1.5
<b>ALM2002</b>	#2
<b>ALM2025</b>	#2.5
<b>ALM2003</b>	#3
<b>ALM2004</b>	#4
<b>ALM2005</b>	#5

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



### Flexible PVC Disposable Laryngeal Masks

Item#	Size
<b>ALM4001</b>	#1
<b>ALM4015</b>	#1.5
<b>ALM4002</b>	#2
<b>ALM4025</b>	#2.5
<b>ALM4003</b>	#3
<b>ALM4004</b>	#4
<b>ALM4005</b>	#5

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



### Size Patient Weight and Age

Size	Patient Weight and Age	
#1	Up to 5 kg	Pediatric
#1.5	5-10 kg	Pediatric
#2	10-20 kg	Pediatric
#2.5	20-30 kg	Pediatric
#3	30-50 kg	Pediatric
#4	50-70 kg	Adult
#5	70-100 kg	Adult

## PeDIA - Pediatric Device for Induction of Anesthesia

Item#

**ANEP021Y**

With Mask Adapter Elbow


**12/pack**



*2025 Innovation Celebration Winner*

Use with  
mask

Or use  
without  
mask



### Turn panic into play

For children, undergoing surgery is already unsettling, and the process of induction is often the most traumatizing part. Often kids will kick and scream, believing they're being trapped by the gas or suffocated by the mask. The experiences create negative emotional scarring that can last weeks, months or even years.

Specially designed for children undergoing surgery, the Pediatric Device for Induction of Anesthesia device, or "PeDIA," for short, includes a colorful whistling balloon, providing kids with something fun to focus on and turning a stressful moment into a lighter one.

Whether they're squeezing the balloon and focusing on its whistle or directly putting themselves to sleep by breathing the gas through balloon delivery, the PeDIA balloon captivates children's attention, allowing their active participation in the induction process.

It's a simple way to add a brighter touch to the medical experience, ease anxiety and bring smiles to patients, parents and staff.

This unique approach—presenting an interactive, familiar toy-like device—not only diverts children's focus from the impending procedure but also significantly reduces the anxiety and resistance commonly associated with mask induction. The aim is to decrease the anxiety felt by over 50% of pediatrics as well as parents and staff.

The PeDIA can be used alone as a mask-free device or, for younger children, with an anesthesia mask.

**SCAN FOR VIDEOS**

