



# Product Catalog

Contract MS1087

## Respiratory Disposables



Contract MS1049

## Non-Invasive Ventilation



Contract MS1097

## Medication Delivery



Contract MS1486

## Open Suction



Contract MS7709

## PeDIA Anesthesia Induction



**vizient**<sup>®</sup>  
Awarded Supplier



## Sales Representatives



**Egle Sventnickaite**  
872-258-1654  
esventnickaite@sunsethcs.com



**Mason Moss**  
262-788-1862  
mmoss@sunsethcs.com



**Owen Ruth**  
312-533-2451  
owenr@sunsethcs.com



**Cindy Geighes**  
312-533-2478  
cgeighes@sunsethcs.com




**Hunter Baudino**  
262-788-1855  
hbaudino@sunsethcs.com



**Wyatt Kanemoto**  
262-788-8941  
wkanemoto@sunsethcs.com

*Director, Strategic Accounts*  
**Lisa Iannazzone** RN, BSN  
636-222-3134  
lisa@sunsethcs.com

*Manager, Product & Business Development*  
**Tom Munar**  
312-533-2476  
tmunar@sunsethcs.com



*Contract Manager, Strategic Accounts*  
**Gina Rosenberg**  
914-837-5215  
grosenberg@sunsethcs.com

*RRT Clinical Specialist*  
**Krystle Bodine** BS, RRT  
kbodine@sunsethcs.com

### Corporate Office

Phone 872-275-3210

Fax 312-997-9985

141 W Jackson Blvd Ste 1950  
Chicago IL 60604

### Central Returns/Distribution Center

2725 Alft Ln  
Elgin IL 60124

### West Returns/Distribution Center

6750 W Van Buren St  
Phoenix AZ 85043

# Content

## Contract MS1049

### Non-Invasive Ventilation Page

---

|                             |       |
|-----------------------------|-------|
| Adult / Pediatric NIV Masks | 4-8   |
| CPAP Tubing                 | 9     |
| Connectors                  | 9     |
| Headgear                    | 10    |
| Chinstraps                  | 10-11 |

## Contract MS1087

### Respiratory Disposables Page

---

|                                   |       |
|-----------------------------------|-------|
| Lynx Venturi Oxygen Delivery      | 12    |
| Oxygen Masks                      | 13    |
| Oxygen Cannula and Accessories    | 14-17 |
| Oxygen Supply Tubing              | 18    |
| Pulse Oximeters                   | 18    |
| Firesafe                          | 19    |
| Humidifier                        | 20    |
| Oxygen Connectors and Accessories | 20-21 |
| Suction Machine Filters           | 22    |
| Oxygen Filters                    | 22    |
| Bacteria Filters                  | 22    |
| HEPA Filters                      | 22    |
| Ventilator Filter                 | 23    |
| Air Entrainment Masks             | 23    |
| Breathing Circuits                | 24-27 |
| CPAP Tubing                       | 27    |
| CPAP Filters                      | 28-29 |
| HMEs                              | 30    |
| HME Accessories                   | 30    |
| Tracheostomy Masks                | 30    |
| Trach Ties                        | 31    |
| Berman and Guedel Airways         | 31    |
| Spirometers                       | 32    |
| Aerosol                           | 32    |

## Contract MS1097

### Medication Delivery Page

---

|                      |    |
|----------------------|----|
| Inhaler Chamber      | 33 |
| Nebulizer Compressor | 33 |
| Peak Flow Meter      | 33 |
| Adapters             | 33 |
| Nebulizer Kits       | 34 |

## Contract MS1486

### Open Suction Page

---

|                       |       |
|-----------------------|-------|
| Suction Catheters     | 35    |
| Trach Care Kits       | 35    |
| Suction Catheter Kits | 35-36 |

## Contract MS1486

### Closed Suction Page

---

|                 |       |
|-----------------|-------|
| Standard Switch | 37-38 |
| Deluxe Switch   | 39-40 |

## Contract MS7709

### PeDIA Anesthesia Induction Page

---

|  |    |
|--|----|
| PeDIA Pediatric Device for Induction of Anesthesia | 41 |
|--|----|

## Ray Pediatric Masks Comparable to Nonny masks

### Mask with Headgear, Tube, and Replaceable Cushion

| Size        | Item#          |
|-------------|----------------|
| Small       | <b>CM061S</b>  |
| Medium      | <b>CM061M</b>  |
| Large       | <b>CM061L</b>  |
| Extra Large | <b>CM061XL</b> |


### Replacement Cushions

| Size        | Item#            |
|-------------|------------------|
| Small       | <b>CUCM061S</b>  |
| Medium      | <b>CUCM061M</b>  |
| Large       | <b>CUCM061L</b>  |
| Extra Large | <b>CUCM061XL</b> |

### Headgear Only

| Size  | Item#         |
|-------|---------------|
| Small | <b>HG061S</b> |
| Large | <b>HG061L</b> |

SCAN CODE to see  
Ray Full Face Mask  
on YouTube



Available Size Guide

### Nasal Mask with Headgear, Tube, and Replaceable Cushions

| Comes with    | Item#             |
|---------------|-------------------|
| S&M Cushions  | <b>CM060-KITS</b> |
| L&XL Cushions | <b>CM060-KITL</b> |

### Replacement Cushions

| Size        | Item#            |
|-------------|------------------|
| Small       | <b>CUCM060S</b>  |
| Medium      | <b>CUCM060M</b>  |
| Large       | <b>CUCM060L</b>  |
| Extra Large | <b>CUCM060XL</b> |

### Headgear Only

| Size  | Item#         |
|-------|---------------|
| Small | <b>HG060S</b> |
| Large | <b>HG060L</b> |



SCAN CODE to see  
Ray Nasal Mask  
on YouTube




## Sol Full Face CPAP Mask

### Mask with Headgear and Replaceable Cushion

10/case

| Size   | Item#         |
|--------|---------------|
| Small  | <b>CM205S</b> |
| Medium | <b>CM205M</b> |
| Large  | <b>CM205L</b> |



### A premium mask without the premium price

#### Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

#### Quiet vent

Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

#### Quick release connector

Rapidly and easily separates or connects mask and the tube

#### Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment



## Full Face Mask

### With Standard Headgear

#### Standard Elbow, Non-Vented

| Size   | Item#         |
|--------|---------------|
| Small  | <b>CM041S</b> |
| Medium | <b>CM041M</b> |
| Large  | <b>CM041L</b> |

#### Single Limb Elbow, Non-Vented, AAV

| Size   | Item#           |
|--------|-----------------|
| Small  | <b>CM041SEE</b> |
| Medium | <b>CM041MEE</b> |
| Large  | <b>CM041LEE</b> |

#### Single Limb Elbow, Vented, AAV

| Size   | Item#            |
|--------|------------------|
| Small  | <b>CM041SEEV</b> |
| Medium | <b>CM041MEEV</b> |
| Large  | <b>CM041LEEV</b> |

**novaplus<sup>+</sup>**

#### Standard Elbow, Non-Vented

| Size   | Item#          |
|--------|----------------|
| Small  | <b>VCM041S</b> |
| Medium | <b>VCM041M</b> |
| Large  | <b>VCM041L</b> |

#### Single Limb Elbow, Non-Vented, AAV

| Size   | Item#            |
|--------|------------------|
| Small  | <b>VCM041SEE</b> |
| Medium | <b>VCM041MEE</b> |
| Large  | <b>VCM041LEE</b> |

#### Single Limb Elbow, Vented, AAV

| Size   | Item#             |
|--------|-------------------|
| Small  | <b>VCM041SEEV</b> |
| Medium | <b>VCM041MEEV</b> |
| Large  | <b>VCM041LEEV</b> |

#### Replacement Cushions

| Size   | Item#              |
|--------|--------------------|
| Small  | <b>CUCM005S-10</b> |
| Medium | <b>CUCM005M-10</b> |
| Large  | <b>CUCM005L-10</b> |



### With Cap Style Headgear

#### Standard Elbow, Non-Vented

| Size   | Item#          |
|--------|----------------|
| Small  | <b>CM041SC</b> |
| Medium | <b>CM041MC</b> |
| Large  | <b>CM041LC</b> |

#### Single Limb Elbow, Non-Vented, AAV

| Size   | Item#            |
|--------|------------------|
| Small  | <b>CM041SEEC</b> |
| Medium | <b>CM041MEEC</b> |
| Large  | <b>CM041LEEC</b> |

#### Single Limb Elbow, Vented, AAV

| Size   | Item#            |
|--------|------------------|
| Small  | <b>CM041SEEV</b> |
| Medium | <b>CM041MEEV</b> |
| Large  | <b>CM041LEEV</b> |

- Unobtrusive design offers a clear line of sight
- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone cushion increases comfort while maintaining optimal seal

## Sunset Deluxe Nasal CPAP Masks

### Mask with Headgear and Replaceable Cushion

| Size   | Item#         |
|--------|---------------|
| Small  | <b>CM006S</b> |
| Medium | <b>CM006M</b> |
| Large  | <b>CM006L</b> |

HCPCS# A7034



Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.



## Nasal Mask

### With Standard Headgear

| Size   | Item#         |
|--------|---------------|
| Small  | <b>CM008S</b> |
| Medium | <b>CM008M</b> |
| Large  | <b>CM008L</b> |

**novaplus<sup>+</sup>**

| Size   | Item#          |
|--------|----------------|
| Small  | <b>VCM008S</b> |
| Medium | <b>VCM008M</b> |
| Large  | <b>VCM008L</b> |



- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

## Total Face Mask

### Standard Elbow

| Size   | Item#         |
|--------|---------------|
| Small  | <b>CM040S</b> |
| Medium | <b>CM040M</b> |
| Large  | <b>CM040L</b> |

### EE Elbow

| Size   | Item#           |
|--------|-----------------|
| Small  | <b>CM040SEE</b> |
| Medium | <b>CM040MEE</b> |
| Large  | <b>CM040LEE</b> |



## Veraseal 2 Disposable Gel Mask with Headgear



### Standard Elbow, Non-Vented

**Sleepnet®**

| Size        | Item#          |
|-------------|----------------|
| Small       | <b>CM50933</b> |
| Medium      | <b>CM50934</b> |
| Large       | <b>CM50935</b> |
| Extra Large | <b>CM50936</b> |

**novaplus<sup>+</sup>**

| Size        | Item#           |
|-------------|-----------------|
| Small       | <b>VCM053S</b>  |
| Medium      | <b>VCM053M</b>  |
| Large       | <b>VCM053L</b>  |
| Extra Large | <b>VCM053XL</b> |



### Single Limb Elbow, Non-Vented, AAV

**Sleepnet®**

| Size        | Item#          |
|-------------|----------------|
| Small       | <b>CM50937</b> |
| Medium      | <b>CM50938</b> |
| Large       | <b>CM50939</b> |
| Extra Large | <b>CM50940</b> |

**novaplus<sup>+</sup>**

| Size        | Item#           |
|-------------|-----------------|
| Small       | <b>VCM054S</b>  |
| Medium      | <b>VCM054M</b>  |
| Large       | <b>VCM054L</b>  |
| Extra Large | <b>VCM054XL</b> |



### Single Limb Elbow, Vented, AAV

**Sleepnet®**

| Size        | Item#          |
|-------------|----------------|
| Small       | <b>CM50827</b> |
| Medium      | <b>CM50828</b> |
| Large       | <b>CM50829</b> |
| Extra Large | <b>CM50830</b> |

**novaplus<sup>+</sup>**

| Size        | Item#           |
|-------------|-----------------|
| Small       | <b>VCM052S</b>  |
| Medium      | <b>VCM052M</b>  |
| Large       | <b>VCM052L</b>  |
| Extra Large | <b>VCM052XL</b> |

- The Veraseal® 2 disposable full face mask brings AIR°gel® Technology to the acute care environment for a giant step forward in patient comfort. Three color coded elbow designs make choosing the right mask for CPAP, single, or dual limb circuit NIV systems simple. The Veraseal® 2 fits quickly and easily while enhancing patient comfort. Both vented (with leak ports) and non-vented (without leak ports) versions are available for CPAP, Bilevel, and ventilator therapy.
- The AIR°gel® cushion provides the utmost in gentle cushioning to protect from pressure points and irritation. AIR°gel® promotes advanced skin protection, superior performance, and unbeatable comfort, all in a lightweight mask.
- The active headgear connector moves as the patient moves for a reliable seal.
- Breathable headgear lets patients enjoy soft, lightweight stability they can count on for a comfortable fit.

## Sleepnet Innova

### Nasal Mask with Headgear

#### Vented

| Size         | Item#          |
|--------------|----------------|
| Small/Medium | <b>CM50180</b> |
| Medium/Large | <b>CM50181</b> |



| Size         | Item#           |
|--------------|-----------------|
| Small/Medium | <b>VCM034SM</b> |
| Medium/Large | <b>VCM034ML</b> |

#### Non-Vented

| Size         | Item#          |
|--------------|----------------|
| Small/Medium | <b>CM50594</b> |
| Medium/Large | <b>CM50595</b> |



| Size         | Item#             |
|--------------|-------------------|
| Small/Medium | <b>VCM034NVSM</b> |
| Medium/Large | <b>VCM034NVML</b> |



Innova® CPAP masks feature patented AIR°gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA

## Sleepnet Ascend Mask Kit

### Full Face System

Item# **CM037F**

#### Each

Vented

#### Mask Only

| Size   | Item#         |
|--------|---------------|
| Small  | <b>CM037S</b> |
| Medium | <b>CM037M</b> |
| Large  | <b>CM037L</b> |

#### Replacement Cushions

| Size   | Item#          |
|--------|----------------|
| Small  | <b>CM037CS</b> |
| Medium | <b>CM037CM</b> |
| Large  | <b>CM037CL</b> |



### Nasal System

Item# **CM037N**

#### Each

Vented

#### Mask Only

| Size   | Item#          |
|--------|----------------|
| Small  | <b>CM037NS</b> |
| Medium | <b>CM037NM</b> |
| Large  | <b>CM037NL</b> |

#### Replacement Cushions

| Size   | Item#           |
|--------|-----------------|
| Small  | <b>CM037NCS</b> |
| Medium | <b>CM037NCM</b> |
| Large  | <b>CM037NCL</b> |



## Sleepnet Phantom®2 and iQ®2

### Phantom®2

Item#  
**CM050**

Replacement  
Cushion  
**CUCM050**



Replacement  
Frame/Tubing  
Assembly  
**CM051**

Sleepnet's patented AIR°gel® technology provides maximum comfort while minimizing skin irritation and pressure points on the nasal bridge. Unlike many other gel masks, the unique qualities of AIR°gel® allow the mask to seal with minimal tension. By design AIR°gel® is twice as soft as competing gel products.

### iQ®2

Item#  
**CM049**

Replacement  
Cushion  
**CUCM049**



Ask us about  
our 30-day  
mask guarantee

- Sleepnet has the only mask that you can shape to each individual face. Flex, Fit, Forget!
- Cushions may be removed and replaced without purchasing a new mask.

- Magnetic Clips — Easy on and off!
- 360° Swivel — Less tube drag for active sleepers.
- One size cushion — Larger fitting range.

## MiniMe 2 Pediatric Nasal Mask



### Vented

| Size        | Item#          |
|-------------|----------------|
| Extra Small | <b>CM60216</b> |
| Small       | <b>CM60217</b> |
| Medium      | <b>CM60218</b> |
| Large       | <b>CM60221</b> |

### Non-Vented

| Size        | Item#          |
|-------------|----------------|
| Extra Small | <b>CM60256</b> |
| Small       | <b>CM60257</b> |
| Medium      | <b>CM60258</b> |
| Large       | <b>CM60261</b> |



### Vented

| Size        | Item#           |
|-------------|-----------------|
| Extra Small | <b>VCM055XS</b> |
| Small       | <b>VCM055S</b>  |
| Medium      | <b>VCM055M</b>  |
| Large       | <b>VCM055L</b>  |

### Non-Vented

| Size        | Item#           |
|-------------|-----------------|
| Extra Small | <b>VCM056XS</b> |
| Small       | <b>VCM056S</b>  |
| Medium      | <b>VCM056M</b>  |
| Large       | <b>VCM056L</b>  |



- Maximum comfort, minimum weight. Patented AIR°gel makes MiniMe masks some of the most gentle fitting masks on the market.
- The shapeable shell mask can be easily molded to fit individual needs.
- 100% designed and manufactured in the USA

## Pediatric SleepWeaver® CPAP Mask

Item#  
**CM101325**



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.


- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

## MiniMe Pediatric Nasal Mask



| Headgear Size | Item#         | OEM#  |
|---------------|---------------|-------|
| Small         | <b>CM004S</b> | 60213 |
| Medium        | <b>CM004M</b> | 60214 |
| Large         | <b>CM004L</b> | 60215 |

HCPCS# A7034



| MASK SIZE | HEAD CIRCUMFERENCE |
|-----------|--------------------|
| SMALL     | 14-16in            |
| MEDIUM    | 17-18in            |
| LARGE     | 19-21in            |

Measure around head above eyebrows for head circumference

## SleepWeaver® Advance Nasal Mask

Item#  
**CM017B**

Mask with Adjustable Headgear

- Easy to store and pack for travel

OEM# 100274  
HCPCS# A7034





## NIV Tubing

Item#

**TUB06** Each  
**SUH0610** 10/pack

19mm with 22mm cuffs



- Sunset's most popular item
- ISO standard 22mm cuffs
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

## Standard Elbow, Non-Vented

Item# **CAP009**  
10/case



## Single Limb Elbow, Non-Vented, AAV

Item# **CAP008**  
10/case



## CPAP Exhalation Port

Item# **VEN002-10**  
10/case



## Single Use Exhalation Port

Item# **VEN003**  
10/case



## Standard 4-Point Headgear

Item#  
**HG003S** Small  
**HG003M** Medium  
**HG003L** Large



## Cap Style Headgear

Item# **HG002**



## Universal Neoprene Chinstrap

Item#  
**CS001**  
 Universal Chinstrap  
 HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows



## Deluxe Chinstrap

| Size | Item#         |
|------|---------------|
| 26in | <b>CS004S</b> |
| 28in | <b>CS004L</b> |

Compare to  
 Respironics Deluxe  
 OEM# 302425  
 HCPCS# A7036

- Latex-free
- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort



## Comfort Chinstrap

Item# **CS005**  
 Compare to  
 ResMed Sullivan  
 OEM# 16015  
 HCPCS# A7036

- Latex-free
- One size fits most
- Easy to use —attaches to itself at the top of the head
- Soft cloth chin cup adds comfort



## Premium Chinstrap

Item# **CS006**  
 Compare to  
 Respironics Premium  
 OEM# 1012911  
 HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.



## Ruby Style Adjustable Chinstrap




| Size        | Item#              |
|-------------|--------------------|
| Regular     | <b>CS007-ADJ</b>   |
| Extra Large | <b>CS007-ADJXL</b> |

HCPCS# A7036

- Latex-free
- Fully adjustable
- Super soft materials for extra comfort



## Ruby Style Chinstrap



For fitting, measure completely around the head, as indicated by the dashed line above.

| Size        | Item#         |
|-------------|---------------|
| 21.5-22.5in | <b>CS007S</b> |
| 22-23in     | <b>CS007M</b> |
| 23-24in     | <b>CS007L</b> |

Compare to Puresom Ruby

OEM# TMS-09 • HCPCS# A7036

- Latex-free
- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

## Halo Style Chinstrap

Item# **CS025**

Compare to BreatheWear Halo Chinstrap

OEM# H3407

HCPCS# A7036



## Lynx Venturi Oxygen Delivery


### Low Concentration

Item# **RES2400**

Each

Adapter Only **RES2401C**

- Design replaces multiple venturi adapters
- Connects to oxygen tubing and regulates flow
- Provides precision FiO<sub>2</sub> delivery from 24–60% with less than 2% variance during mild, moderate, and severe respiratory distress
- Excellent CO<sub>2</sub> washout



### Unprecedented simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

### Do no harm

Quicker de-escalation of O<sub>2</sub> levels and prevention of excessive CO<sub>2</sub> in the patient.

### Reduction in waste

Remains with patient throughout EMS and hospital journey, replacing multiple adapters which may get lost or contaminated during transition.

### High Concentration

Item# **RES2401**


Each

Adapter Only **RES241HC**

- Allows rapid de-escalation of FiO<sub>2</sub> in the EMS ambulance or ED to prevent CO<sub>2</sub> buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia (PaO<sub>2</sub>/FiO<sub>2</sub>) allows early intervention and better outcomes
- Excellent CO<sub>2</sub> washout

Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Capabilities could include replacing non-rebreathers and/or partial rebreather products



### Unmatched precision

Provides specific levels of oxygen with less than 2% variance. Comparable products vary by as much as 30%.

## Medium Concentration Oxygen Mask

50/case

| Item#            | Size      | Tubing Length | Connector |
|------------------|-----------|---------------|-----------|
| <b>RES2100NT</b> | Adult     | None          | Standard  |
| <b>RES2100U</b>  | Adult     | 7 ft          | Universal |
| <b>RES2100</b>   | Adult     | 7 ft          | Standard  |
| <b>RES2203P</b>  | Pediatric | None          | Standard  |
| <b>RES2100P</b>  | Pediatric | 7 ft          | Standard  |
| <b>RES2203A</b>  | Adult     | None          | Standard  |
| <b>RES2203A7</b> | Adult     | 7 ft          | Standard  |
| <b>RES2203I</b>  | Infant    | None          | Standard  |
| <b>RES2203I7</b> | Infant    | 7 ft          | Standard  |
| <b>RES2203P7</b> | Pediatric | 7 ft          | Standard  |

- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



## High Concentration Oxygen Mask

50/case

| Item#            | Size      | Tubing Length |
|------------------|-----------|---------------|
| <b>RES2202A</b>  | Adult     | None          |
| <b>RES2202A7</b> | Adult     | 7 ft          |
| <b>RES2202I</b>  | Infant    | None          |
| <b>RES2202I7</b> | Infant    | 7 ft          |
| <b>RES2202P</b>  | Pediatric | None          |
| <b>RES2202P7</b> | Pediatric | 7 ft          |

## Non-Rebreather Oxygen Mask

50/case

| Item#             | Size      | Tubing Length | Connector | Safety Vent |
|-------------------|-----------|---------------|-----------|-------------|
| <b>RES2102I</b>   | Infant    | 7 ft          | Standard  | No          |
| <b>RES2102IM</b>  | Infant    | None          | Standard  | No          |
| <b>RES2102P</b>   | Pediatric | 7 ft          | Standard  | No          |
| <b>RES2102PS</b>  | Pediatric | 7 ft          | Standard  | Yes         |
| <b>RES2102PSM</b> | Pediatric | None          | Standard  | Yes         |
| <b>RES2102PSU</b> | Pediatric | 7ft           | Universal | Yes         |
| <b>RES2102PU</b>  | Pediatric | 7 ft          | Universal | No          |
| <b>RES2102S</b>   | Adult     | 7 ft          | Standard  | Yes         |
| <b>RES2102SM</b>  | Adult     | None          | Standard  | Yes         |
| <b>RES2102U</b>   | Adult     | 7 ft          | Universal | No          |
| <b>RES2102</b>    | Adult     | 7 ft          | Standard  | No          |



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



## Humidoflo High Flow Nasal Cannula

10/case

| Item#            | Size   | F&P Adapter |
|------------------|--------|-------------|
| <b>HFN001AS</b>  | Small  | No          |
| <b>HFN001AM</b>  | Medium | No          |
| <b>HFN001AL</b>  | Large  | No          |
| <b>HFN001ASA</b> | Small  | Yes         |
| <b>HFN001AMA</b> | Medium | Yes         |
| <b>HFN001ALA</b> | Large  | Yes         |

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respironics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)



With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerable due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O<sub>2</sub> concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.



## Adult Oxygen Cannula

| Item#            | Tubing Length | Connector | Qty |
|------------------|---------------|-----------|-----|
| <b>RES1100</b>   | None          | Standard  | 50  |
| <b>RES1104</b>   | 4 ft          | Standard  | 50  |
| <b>RES1107N</b>  | 7 ft          | Standard  | 50  |
| <b>RES1107U</b>  | 7 ft          | Universal | 50  |
| <b>RES1107NU</b> | 7 ft          | Universal | 50  |
| <b>RES1107</b>   | 7 ft          | Standard  | 50  |
| <b>RES1115U</b>  | 15 ft         | Universal | 50  |
| <b>RES1115</b>   | 15 ft         | Standard  | 50  |
| <b>RES1125</b>   | 25 ft         | Standard  | 25  |
| <b>RES1135</b>   | 35 ft         | Standard  | 20  |
| <b>RES1150U</b>  | 50 ft         | Universal | 20  |
| <b>RES1150Z</b>  | 50 ft         | Standard  | 20  |



## Soft Cannula

| Item#            | Tubing Length | Connector | Qty |
|------------------|---------------|-----------|-----|
| <b>RES1100S</b>  | None          | Standard  | 50  |
| <b>RES1107SU</b> | 7 ft          | Universal | 50  |
| <b>RES1107S</b>  | 7 ft          | Standard  | 50  |
| <b>RES1114S</b>  | 14 ft         | Standard  | 50  |
| <b>RES1114SU</b> | 14 ft         | Universal | 50  |
| <b>RES1115S</b>  | 15 ft         | Standard  | 50  |
| <b>RES1125S</b>  | 25 ft         | Standard  | 25  |
| <b>RES1125SU</b> | 25 ft         | Universal | 25  |

## Flared Cannula

| Item#            | Tubing Length | Connector | Qty |
|------------------|---------------|-----------|-----|
| <b>RES1100F</b>  | None          | Standard  | 50  |
| <b>RES1107F</b>  | 7 ft          | Standard  | 50  |
| <b>RES1107FU</b> | 7 ft          | Universal | 50  |
| <b>RES1114F</b>  | 14 ft         | Standard  | 50  |
| <b>RES1125F</b>  | 25 ft         | Standard  | 25  |
| <b>RES1150F</b>  | 50 ft         | Standard  | 20  |

## Adult Cannula with Ear Cushions

Standard Touch, Standard Tip

| Item#             | Tubing Length | Connector | Qty |
|-------------------|---------------|-----------|-----|
| <b>RES1100EC</b>  | None          | Standard  | 50  |
| <b>RES1104EC</b>  | None          | Standard  | 50  |
| <b>RES1107ECU</b> | 7 ft          | Universal | 50  |
| <b>RES1107EC</b>  | 7 ft          | Standard  | 50  |
| <b>RES1114EC</b>  | 14 ft         | Standard  | 50  |
| <b>RES1115EC</b>  | 15 ft         | Standard  | 50  |
| <b>RES1125EC</b>  | 25 ft         | Standard  | 25  |



## Cannula Ear Cushions

Item# **RES020**  
50 Pairs/case

- Add cushions to existing cannula



## High Flow Soft Cannula

25/case

| Length | Item#             |
|--------|-------------------|
| 4ft    | <b>RES1104SHF</b> |
| 7ft    | <b>RES1107SHF</b> |
| 15ft   | <b>RES1115SHF</b> |
| 25ft   | <b>RES1125SHF</b> |




- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

## 7 ft Cannula With CO<sub>2</sub> Sampling Line

50/case

| Item#            | Size      | Connector | Luer   | Filter |
|------------------|-----------|-----------|--------|--------|
| <b>RES8007</b>   | Adult     | Standard  | Male   | No     |
| <b>RES8007F</b>  | Adult     | Standard  | Male   | Yes    |
| <b>RES8007I</b>  | Infant    | Standard  | Male   | No     |
| <b>RES8007IF</b> | Infant    | Standard  | Male   | Yes    |
| <b>RES8007IU</b> | Infant    | Universal | Male   | No     |
| <b>RES8007P</b>  | Pediatric | Standard  | Male   | No     |
| <b>RES8007PF</b> | Pediatric | Standard  | Male   | Yes    |
| <b>RES8007PU</b> | Pediatric | Universal | Male   | No     |
| <b>RES8007U</b>  | Adult     | Universal | Male   | No     |
| <b>RES8107</b>   | Adult     | Standard  | Female | No     |
| <b>RES8107F</b>  | Adult     | Standard  | Female | Yes    |
| <b>RES8107I</b>  | Infant    | Standard  | Female | No     |
| <b>RES8107IF</b> | Infant    | Standard  | Female | Yes    |
| <b>RES8107P</b>  | Pediatric | Standard  | Female | No     |
| <b>RES8107PF</b> | Pediatric | Standard  | Female | Yes    |



## Pediatric, Infant, and Neonate Cannula

50/case


| Item#              | Size      | Tubing Length | Connector | Touch    | Ear Cushions |
|--------------------|-----------|---------------|-----------|----------|--------------|
| <b>RES1100P</b>    | Pediatric | None          | Standard  | Standard | No           |
| <b>RES1107I</b>    | Infant    | 7 ft          | Standard  | Standard | No           |
| <b>RES1107NE</b>   | Neonate   | 7 ft          | Standard  | Standard | No           |
| <b>RES1107P</b>    | Pediatric | 7 ft          | Standard  | Standard | No           |
| <b>RES1107PEC</b>  | Pediatric | 7 ft          | Standard  | Standard | Yes          |
| <b>RES1107PECU</b> | Pediatric | 7 ft          | Universal | Standard | Yes          |
| <b>RES1107SI</b>   | Infant    | 7 ft          | Standard  | Soft     | No           |
| <b>RES1107SIU</b>  | Infant    | 7 ft          | Universal | Soft     | No           |
| <b>RES1107SNE</b>  | Neonate   | 7 ft          | Standard  | Soft     | No           |
| <b>RES1107SP</b>   | Pediatric | 7 ft          | Standard  | Soft     | No           |
| <b>RES1107SPU</b>  | Pediatric | 7 ft          | Universal | Soft     | No           |
| <b>RES1114SP</b>   | Pediatric | 14 ft         | Standard  | Soft     | No           |
| <b>RES1114SPU</b>  | Pediatric | 14 ft         | Universal | Soft     | No           |

## Secure Cannula Stickers

Item# **RES080**

25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



## Oxygen Supply Tubing

| Item#            | Length | Connector | Color | Bore      | Qty |
|------------------|--------|-----------|-------|-----------|-----|
| <b>RES3007SB</b> | 7 ft   | Standard  | None  | Smooth    | 50  |
| <b>RES3007U</b>  | 7 ft   | Universal | None  | Channeled | 50  |
| <b>RES3007</b>   | 7 ft   | Standard  | None  | Channeled | 50  |
| <b>RES3007G</b>  | 7 ft   | Standard  | Green | Channeled | 50  |
| <b>RES3015G</b>  | 15 ft  | Standard  | Green | Channeled | 50  |
| <b>RES3015</b>   | 15 ft  | Standard  | None  | Channeled | 50  |
| <b>RES3025G</b>  | 25 ft  | Standard  | Green | Channeled | 25  |
| <b>RES3025U</b>  | 25 ft  | Universal | None  | Channeled | 25  |
| <b>RES3025</b>   | 25 ft  | Standard  | None  | Channeled | 25  |
| <b>RES3030</b>   | 30 ft  | Standard  | None  | Channeled | 20  |
| <b>RES3030G</b>  | 30 ft  | Standard  | Green | Channeled | 20  |
| <b>RES3040G</b>  | 40 ft  | Standard  | Green | Channeled | 20  |
| <b>RES3050G</b>  | 50 ft  | Standard  | Green | Channeled | 20  |
| <b>RES3050</b>   | 50 ft  | Standard  | None  | Channeled | 20  |

Universal Connector



## Finger Pulse Oximeter

Item# **RES5102**  
Each



- Oversized 1.3" OLED display face
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- Includes lanyard and manual





## Firesafe Thermal Shut-Off Valve

Item#  
**RES010-10**  
Bidirectional  
10/case




- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire
- Smaller
- Laser etched
- Expiration date
- Human readable UDI
- Designed to last 5 years



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

### Fatalities per 100,000 patients



## FireSafe Nozzle

Item#  
**RES010FND**  
20/pack



## FireSafe Tubing Connector

Item#  
**RES010TC** 100/case



## Disposable Bubble Humidifier 500 ml

Item#

**RES050**

Up to 15 LPM • 6 psi • 500 ml

**50/case**

HCPCS# E0560

OEM# HUM001



- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient

## Bacterial/Viral HME Filter

Item#

**BF2200****5/pack**

## Oxygen Tubing Connector

Item# **RES012**

Fits 3/16in and 1/4in ID Tubing

**10/pack**

- Barbed ends ensure secure, leak-free connections

## Christmas Tree Connector

Item# **RES001**

Green without Swivel

**10/pack****RES001-30 30/pack**

- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Item# **RES002**

Green with Swivel

**10/pack**

## Circuit Connectors

Item# **RES014**

22mm OD x 22mm OD

With Lip

**10/pack**

## Humidifier Tubing Adapter

Item# **RES016**

22mm ID x 22mm OD

with Lip

**10/pack**

## Bacterial/Viral HME Filter - Adult

Item#

**BF104****Each**

- Bacterial efficiency 99.9999%
- Viral efficiency 99.999%
- Humidification efficiency 37.5 mg/l H2O (500 ml tidal volume)
- Coupled with a hydrophobic bacterial viral filter membrane



## Wrenches

Item# **RES031-1**

With Chain

**1/pack**

## Tubing Reducer Connector

Item# **RES011**

Fits 3/16in, 1/4in,  
5/16in, 3/8in, and 7/16in  
ID Tubing

**10/pack**



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

## Circuit Connectors

Item# **RES013**

22mm OD x 15mm ID/22mm OD

**10/pack**



## Oxygen Bleed-In Adapter

Item# **RES019**

22mm with 6mm  
Rotating Stem

**10/pack**

**RES019-1 1/Pack**



## Oxygen Tubing Y-Connector

Item# **RES017**

Fits 3/16in and 1/4in ID Tubing  
(All 3 sides are the same)

**10/pack**

- Connects 2 oxygen sources for higher flows



## Connector with 6mm Oxygen Stem

Item# **RES058**

22mm with 6mm Stem

**10/pack**



## Humidifier Connectors

Item# **HUM008**

L-Connector 12in

**10/pk**

**HUM008-1 1/Pack**

Item# **HUM008L**

L-Connector 21in

**10/pk**

Item# **HUM008S**

Straight Connector 9in

**10/pk**



## Oxygen Tubing Connector — Swivel

Item# **RES018**

Male/Male

**10/pack**

**RES018-1 1/Pack**

- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing




## Swivel Connector Silencer for Ventilator Circuits and CPAP

Item# **VEN002**

**10/pack**



### Suction Machine Inline Filters



Item# **BF1412**  
Final Bacteria Filter  
1/4in to 1/2in Stepped Barb  
on Both Sides  
**10/pack**



Item# **BF3812**  
3/8in to 1/2in  
**10/pack**



Item# **BF1312**  
1/3in to 1/2in  
**1/pack**



Item# **BF516TS**  
5/16in Barbed with  
Threaded NPT  
**10/pack**

### HEPA Filter for 3B Stratus

Item# **BF6020**  
HEPA Filter  
**Each**  
OEM# STR1020



### 6in HEPA Capsule with Removable Sticker

Item# **BFH6S**  
Threaded  
**Each**



### 3B RESmart Style Filter

Item# **CF8004**  
Disposable Filter  
**3/pack** • 3<sup>5</sup>/<sub>16</sub>" x 1<sup>5</sup>/<sub>8</sub>"  
HCPCS# A7038



### Oxygen Filters



Item# **BF316**  
Final Bacteria Filter • 3/16 Stepped  
Barb and 1/4in Straight Barb on  
Both Sides  
**10/pack**  
**BF316-1** 1/pack



Item# **BF516**  
1/4in Barbed on Both Sides  
**10/pack**



Item# **BF1438**  
1/4in to 3/8in  
**10/pack**

### Invacare Platinum 5 and 10 Style Filter

Item# **OF9002**  
Foam Cabinet Filter  
9" x 27/8" x 1/2"  
**10/pack**  
OEM#1107412



### Invacare Perfecto2 Style Filter

Item# **OF9004**  
Foam Cabinet Filter  
67/8" x 33/4" x 1/2"  
**10/pack**  
OEM# 1143492



### Extended Life Intake Bacteria Filter

Item# **OF1025D**  
**Each**



### DeVilbiss/Drive SolAiris III, 303 Series Style Filter

Item# **OF5004**  
Foam Cabinet Filter  
33/4" x 33/4" x 1/2"  
**10/pack**  
OEM# 30302605



## Filter for ResMed Astral Ventilator

Item# **VF003-4**  
4/pack



Item# **VF003-H**  
Filter Housing  
Each



Item# **VF003-KIT**  
4-Filter Kit  
w/Filter Housing  
OEM# 27939



## Adult Air Entrainment Mask

50/case

| Item#             | Size      | Kit       | Tubing Length | Connector | Flex Tube |
|-------------------|-----------|-----------|---------------|-----------|-----------|
| <b>RES2400A</b>   | Adult     | Mask Only | 7 ft          | Standard  | Yes       |
| <b>RES2400AK</b>  | Adult     | Yes       | 7 ft          | Standard  | Yes       |
| <b>RES2400AKU</b> | Adult     | Yes       | 7 ft          | Universal | Yes       |
| <b>RES2400AU</b>  | Adult     | Mask Only | 7 ft          | Universal | Yes       |
| <b>RES2400P</b>   | Pediatric | Mask Only | 7 ft          | Standard  | Yes       |
| <b>RES2400PK</b>  | Pediatric | Yes       | 7 ft          | Standard  | Yes       |
| <b>RES2400PKU</b> | Pediatric | Yes       | 7 ft          | Universal | Yes       |
| <b>RES2400PU</b>  | Pediatric | Mask Only | 7 ft          | Universal | Yes       |



## 90-Degree CPAP Tubing Elbow

Item# **CAP1027**  
Each

The Sunset 90-Degree CPAP Tubing Elbow Adapter connects standard 22mm or slim 15mm CPAP tubes to most CPAP machines at a 90-degree angle, relieving stress on the tubing connection and creating a more convenient pathway to the bed.

Comparable to OEM# 37394






## Passive Non-Heated Breathing Circuits

### Each

For use with Respiration Trilogy

| Exhalation Port | Item#         | OEM#              |
|-----------------|---------------|-------------------|
| Swivel          | <b>VC1000</b> | 1073228 & 1073236 |
| Fixed           | <b>VC1100</b> | 1069467 & 1069468 |



## Non-Heated Breathing Circuit


For use with Respiration Trilogy

Item# **VC1001**

Active Adult, No Water Trap  
22mm cuffs

**30/case**

OEM#  
1073223 & 1073232



## Adult Passive Heated Circuit with Exhalation Port

For use with Philips  
Respiration Trilogy

Item#  
**VC1003**

**VC1003P**  
Pediatric Size


**30/case**

Comparable to OEM# 1076046,  
1073227, 1073235



### Includes:

- 1 VC1114 Passive Heated Breathing Circuit
- 1 VEN002 Swivel Exhalation Port



## Passive Circuit w/Exhalation Port and Bacterial/Viral Filter

With 22mm cuffs

Item#  
**VC1005**  
**Each**

OEM# 1147378, 1132340, 7615





## Non-Heated Breathing Circuits

For use with Respiration Trilogy  
22mm cuffs


**30/case**

| Item#          | OEM#              | Mode    | Water Trap | Size      | Tube Diameter |
|----------------|-------------------|---------|------------|-----------|---------------|
| <b>VC1001</b>  | 1073223 & 1073232 | Active  | No         | Adult     | 22mm          |
| <b>VC1004</b>  | 1069443 & 1069444 | Passive | Yes        | Adult     | 22mm          |
| <b>VC1004P</b> | 1073224 & 1073233 | Passive | Yes        | Pediatric | 15mm          |


Item#  
**VC1001**  
Active Adult  
OEM#  
1073223  
1073232




To use active circuits with Respiration Trilogy, simply cut ends off pressure and exhalation lines.



Item#  
**VC1004**  
Passive Adult  
OEM#  
1069443  
1069444



Item#  
**VC1004P**  
Passive Pediatric  
OEM#  
1073224  
1073233



## Passive Heated Breathing Circuit

**10/case**


For use with MR850 Humidifier

Item#

**VC1114**

**VC1114P**

Pediatric Size



## Adult NIV Patient Circuit


Item# **VC1160**

Compatible with V60, No Filters

**10/case**

Comparable to

OEM# 1069210



## Adult NIV Single Limb Patient Circuit with Main Flow Filter

Item# **VC1161**

Includes 1 BF2200 22mm Bacteria Disk Filter

**10/case**




## Adult NIV Single Limb Patient Circuit with Main Flow Filter and Exhalation Port Filter

Item# **VC1162**

Includes 2 BF2200 22mm Bacteria Disk Filters

**10/case**



## Non-Heated Dual Limb Breathing Circuits

30/case

| Item#          | Size      | Humidifier Limb | Water Trap |
|----------------|-----------|-----------------|------------|
| <b>VC2001P</b> | Pediatric | No              | No         |
| <b>VC2001</b>  | Adult     | No              | No         |
| <b>VC2002P</b> | Pediatric | Yes             | No         |
| <b>VC2002</b>  | Adult     | Yes             | No         |
| <b>VC2003</b>  | Adult     | Yes             | Yes        |



## Slim CPAP Tube

Item#

**TUB006SS**

6ft x 15mm with  
22mm cuffs

**Each**

HCPSC# A7037



### Highly durable

Strain relief cuff decreases separation and tearing

### Virtually no memory

Designed to be the most flexible tube available

### Minimal weight

One of the lightest tubes on the market at 2.6 ounces

### Great value

Affordable and sustainable

## CPAP Tubing

Item#

**TUB008** 8ft

**TUB010** 10ft

19mm with 22mm cuffs

**Each**

HCPSC# A7037



- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

## Heated Tube for Respirationics

Item#


**TUB06HPR**

Heated Tube for Philips Respirationics  
Dreamstation Series and System One

**Each**

OEM# HT15

HCPSC# A4604



- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off

### ResMed S9 and AirSense 10 Series Filters

**Item# CF2107**

Disposable Filter

**3/pack** • 2 1/8" x 1 3/8"

HCPCS# A7038

**CF2107-1** 1/pack OEM# 36850

**CF2107-2** 2/pack OEM# 36851

**CF2107-50** 50/pack OEM# 36853

**Item# CF2107HA**

Hypoallergenic Filter

**3/pack** • 2 1/8" x 1 3/8"

HCPCS# A7038

**CF2107HA-1** 1/pack OEM# 36855

**CF2107HA-2** 2/pack OEM# 36856

**CF2107HA-50** 50/pack

OEM# 36858



### ResMed AirSense 11 Series Filters

**Item# CF2111**

Disposable Filter

**3/pack**

HCPCS# A7038

**CF2111-1** 1/pack OEM# 39300

**CF2111-2** 2/pack OEM# 39301

**CF2111-50** 50/pack OEM# 39303

**Item# CF2111HA-3**

Hypoallergenic Filter

**3/pack**

1 1/4" x 1 3/16"

HCPCS# A7038

**CF2111HA-1** 1/pack

**CF2111HA-2** 2/pack



### 3B Luna G3 Style Filters

**Item# CF8006-6**

Disposable Filter

**6/pack**

HCPCS# A7038

**CF8006-1** 1/pack

**CF8006-2** 2/pack

**Item# CF8006F-1**

Reusable Filter

**1/pack**

HCPCS# A7039



### 3B Luna and Luna II Style Filters

**Item# CF8005**

Disposable Filter

**12/pack** • 2" x 1 3/16"

HCPCS# A7038

**Item# CF8005F-1**

Reusable Filter

**1/pack** • 2" x 1 3/16"

HCPCS# A7039

**CF8005F-2** 2/pack




### 3B Stratus 5 Style Filter

**Item# OF6030**

Foam Cabinet Filter

**10/pack**

OEM# STR1030





## Respironics DreamStation 2 Series Filters

### Item# **CF1012**

Disposable Filter

**6/pack**

OEM# 1142831

HCPSC# A7038

**CF1012-1** 1/pack OEM# 1142828

**CF1012-2** 2/pack OEM# 1142829

**CF1012-3** 3/pack

**CF1012-50** 50/pack



### Item# **CF1012F**

Reusable Filter

**2/pack**

OEM# 1142687

HCPSC# A7039

**CF1012F-1** 1/pack



## Respironics DreamStation Series Filters

### Item# **CF1010**

Disposable Filter

**6/pack**

OEM# 1122519

HCPSC# A7038

**CF1010-1** 1/pack

**CF1010-2** 2/pack



### Item# **CF1010F**

Reusable Filter

**2/pack**

OEM# 112446 (1-pack)

HCPSC# A7039

**CF1010F-1** 1/pack

**CF1010F-6** 6/pack



## Flex HME






Item# **RES027FLEX**  
40/case



|                                   |                               |
|-----------------------------------|-------------------------------|
| <b>TV Range (mL)</b>              | 300-1500                      |
| <b>Moisture Output (mg H2O/L)</b> | 33 mg/L                       |
| <b>Weight (g)</b>                 | 17                            |
| <b>Resistance (cm H2O)</b>        | 0.5 @ 30 LPM<br>1.34 @ 60 LPM |
| <b>Dead Space (mL)</b>            | 75                            |
| <b>Connection (mm)</b>            | 22M,15F                       |

## Heat Moisture Exchangers

30/case

| Item#                             |  |  |  |  |  |
|-----------------------------------|---|---|--|---|---|
| <b>Item#</b>                      | <b>RES027</b><br>T-HME  | <b>RES0270P</b><br>T-HME with<br>Oxygen Port                                      | <b>RES0271L</b><br>Inline HME —<br>Large Volume                                    | <b>RES027NE</b><br>Inline HME —<br>Neonatal   | <b>RES027HV2</b><br>Foam HME —<br>Adult/Child                                       |
| <b>TV Range (mL)</b>              | 250-1000  | 150-1250  | 50-250   | 15-50   | 250-1000  |
| <b>Moisture Output (mg H2O/L)</b> | 30, Vt = 0.1L   | 36.3 mg/L   | 31, Vt = 0.1L  | 30  | 29  |
| <b>Weight (g)</b>                 | 4.5   | 7.4   | 10   | 4.5   | 21  |
| <b>Resistance (cm H2O)</b>        | 0.32 @ 30 LPM<br>0.81 @ 60 LPM<br>1.62 @ 90 LPM                                   | 0.2 @ 30 LPM  | 0.4 @ 10 LPM   | 0.5 @ 10 LPM  | 0.52 @ 30 LPM<br>1.32 @ 60 LPM<br>2.54 @ 90 LPM                                     |
| <b>Dead Space (mL)</b>            | 12.5  | 7   | 13   | 2.5   | 29  |
| <b>Connection (mm)</b>            | 15F   | 15F<br>5.5mm O2 port  | 15M,22M/15F  | 15M,15F   | 22F/15M, 22M/15F  |

## Adapter for T-HME

Item# **RES027A**  
Each



## Flexible Connector

Item# **TRAA004**  
15mm F x 15mm M  
Length 4.9-7.5cm  
50/case  
OEM# 3215



## Tracheostomy Mask

50/case

Item# **RES2130**  
Adult Regular Size

**RES2130L**  
Adult Large Size

**RES2230**  
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O2 supply tubing

## Trach Ties

### 20/case


#### One-piece

- Easy on

##### Item#

**RES029-17** 2.5x43cm (17in)

**RES029P-11** 2.5x30cm (11.5in)



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft 60% Polyester, 20% Cotton, 20% Nylon

#### Two-piece

- Easy off


##### Item#

**RES030P-10** 2x25.5cm (10in)

**RES030P-14** 2x35.5cm (14in)

**RES030P-18** 2x45.5cm (18in)

**RES030-20** 2.5x51cm (20in)



## Color Coded Dual Channel Berman Oral Airway

### 50/pack

| Item#              | Size  | Color  |
|--------------------|-------|--------|
| <b>AOA40BC-50</b>  | 40mm  | Pink   |
| <b>AOA50BC-50</b>  | 50mm  | Blue   |
| <b>AOA60BC-50</b>  | 60mm  | Black  |
| <b>AOA70BC-50</b>  | 70mm  | White  |
| <b>AOA80BC-50</b>  | 80mm  | Green  |
| <b>AOA90BC-50</b>  | 90mm  | Yellow |
| <b>AOA100BC-50</b> | 100mm | Red    |
| <b>AOA110BC-50</b> | 110mm | Orange |



- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

## Oropharyngeal Guedel Airway

### 10/pack

| Item#            | Size  | Color  |
|------------------|-------|--------|
| <b>AOP500GC</b>  | 50mm  | Blue   |
| <b>AOP601GC</b>  | 60mm  | Black  |
| <b>AOP702GC</b>  | 70mm  | White  |
| <b>AOP803GC</b>  | 80mm  | Green  |
| <b>AOP904GC</b>  | 90mm  | Yellow |
| <b>AOP1005GC</b> | 100mm | Red    |



- Designed to overcome associated risks of loose or detached bite blocks
- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

## Volumetric Incentive Spirometer

12/case

Item# **RES055**  
2500ml

**RES056**  
5000ml

**RES055P**  
Pediatric

- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise.
- Dual-sided, clear and legible identification of volumes.
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration.
- Compact and ergonomic design.



SCAN FOR VIDEO



## Three Ball Spirometer

Item# **RES057**  
12/case

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient



- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

## Valved T-Adapter, Aerosol 22mm ID x 22mm OD

Item# **NEB9001T**  
30/case



## Aerosol Mask

50/case

Item# **RES2110**  
Adult Size

Item# **RES2110F**  
Adult Size  
w/ Flex Tube

Item# **RES2210**  
Pediatric Size



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

## SperaChamber

### Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler

Item#  
**RES195-10**  
10/case

**RES195-50**  
50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



Adult Mask



### SperaChamber with Mask 10/case

| Item#           | Size               |
|-----------------|--------------------|
| <b>RES195MS</b> | Small (Infant)     |
| <b>RES195MM</b> | Medium (Pediatric) |
| <b>RES195ML</b> | Large (Adult)      |

### Mask Only

| Item#                  | Size               | Qty    |
|------------------------|--------------------|--------|
| <b>RES195MASK-S</b>    | Small (Infant)     | Each   |
| <b>RES195MASK-S-10</b> | Small (Infant)     | 10/box |
| <b>RES195MASK-M</b>    | Medium (Pediatric) | Each   |
| <b>RES195MASK-M-10</b> | Medium (Pediatric) | 10/box |
| <b>RES195MASK-L</b>    | Large (Adult)      | Each   |
| <b>RES195MASK-L-10</b> | Large (Adult)      | 10/box |

## Nebulizer Compressor

Each

Item#  
**NEB100**

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Foam Filters (BF018 x5)

**NEB100-R**

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



### Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

### Features

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

SCAN FOR VIDEO




## Peak Flow Meter

Item# **RES052**  
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**  
Pediatric Mouthpiece for Peak Flow Meter  
10/case

## Plastic Intake Filter with Head For Respiration Inspiration 626 & 929

Item# **BF016**  
10/pack



## Valved T-Adapter

Item# **VEN001**  
22F-18F-22M/15F  
10/pack





## Disposable Nebulizer Kits with Mouthpiece

**50/case**

**All include 6cc jet nebulizer medicine cup**

**Accessory options**

- T-piece
- Mouthpiece
- Corrugated tube
- 7ft supply tube
- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal



| Item#           | Size  | Accessories                   |
|-----------------|-------|-------------------------------|
| <b>RES091</b>   | Adult | All                           |
| <b>RES091NA</b> | Adult | No corrugated tube            |
| <b>RES091U</b>  | Adult | All, Plus Universal Connector |

## Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**

Disposable  
Single Patient  
**50/case**

**RES096U** With Universal Connector

Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube



## Disposable Nebulizer Kits with Mask

**50/case**

**All include 6cc jet nebulizer medicine cup**

**Accessory options**

- Mask
- 7ft supply tube

| Item#         | Size      | Accessories |
|---------------|-----------|-------------|
| <b>RES090</b> | Adult     | With Mask   |
| <b>RES092</b> | Pediatric | All         |



- Disposable single patient use for long term delivery procedures
- 6cc bottle is easy to assemble and disassemble


Medicine Cup  
Item# **RES090A**  
**50/case**



## Suction Catheters - Calibrated

50/case

| Item#              | Size | Tip     | Valve   | Type     | Sterile |
|--------------------|------|---------|---------|----------|---------|
| <b>TRASC106C</b>   | 6fr  | DeLee   | Chimney | Coiled   | No      |
| <b>TRASC108C</b>   | 8fr  | DeLee   | Chimney | Coiled   | No      |
| <b>TRASC108S</b>   | 8fr  | DeLee   | Chimney | Straight | No      |
| <b>TRASC110C</b>   | 10fr | Whistle | Y       | Coiled   | No      |
| <b>TRASC110S</b>   | 10fr | Whistle | Y       | Straight | Yes     |
| <b>TRASC112C</b>   | 12fr | Whistle | Y       | Coiled   | No      |
| <b>TRASC112S</b>   | 12fr | Whistle | Y       | Straight | Yes     |
| <b>TRASC114C</b>   | 14fr | Whistle | Y       | Coiled   | No      |
| <b>TRASC114CVS</b> | 14fr | Cloude  | Chimney | Straight | No      |
| <b>TRASC114S</b>   | 14fr | Whistle | Y       | Straight | Yes     |
| <b>TRASC116C</b>   | 16fr | Whistle | Y       | Coiled   | No      |
| <b>TRASC118C</b>   | 18fr | Whistle | Y       | Coiled   | No      |
| <b>TRASC118S</b>   | 18fr | DeLee   | Chimney | Straight | No      |



## Trach Care Kits

50/case

### Sunset Trach Care Kit

Item# **TRATK003**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- Gauze
- 2 Cotton Tips
- Brush,
- Twill Tape
- Trach Dressing

### Sunset Mini Trach Kit

Item# **TRATK004**

Includes:

- 3 Pipe Cleaners
- 4 Gauze
- 2 Cotton Tips
- Brush
- Twill Tape
- Dressing

### Sunset Trach Care Kit with Forceps

Item# **TRATK005**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- 3 Gauze
- 3 Cotton Tips
- Brush
- Trach Dressing

## Suction Catheter Kits with Cup and Gloves

50/case

| Item#           | Catheter Size | Tip     | Valve   | Sterile |
|-----------------|---------------|---------|---------|---------|
| <b>TRATK206</b> | 6fr           | DeLee   | Chimney | Yes     |
| <b>TRATK208</b> | 8fr           | Whistle | Y       | No      |
| <b>TRATK210</b> | 10fr          | Whistle | Y       | No      |
| <b>TRATK212</b> | 12fr          | Whistle | Y       | No      |
| <b>TRATK214</b> | 14fr          | Whistle | Y       | No      |
| <b>TRATK218</b> | 18fr          | Whistle | Y       | Yes     |



## Pediatric Suction Catheter Kits

50/case

| Item#           | Catheter Size | Cup |
|-----------------|---------------|-----|
| <b>TRATK306</b> | 6fr           | No  |
| <b>TRATK308</b> | 8fr           | No  |
| <b>TRATK408</b> | 8fr           | Yes |

- Sterile
- Y-Valve DeLee Catheter
- Includes Gloves

## Mini Catheter Tray with Gloves

50/case


| Item#               | Catheter Size      |
|---------------------|--------------------|
| <b>TRATK310</b>     | 10fr               |
| <b>TRATK312</b>     | 12fr               |
| <b>TRATK312VCVN</b> | 18fr, Vinyl Gloves |
| <b>TRATK314</b>     | 14fr               |
| <b>TRATK316</b>     | 16fr               |
| <b>TRATK318</b>     | 18fr               |

- Coiled Catheter
- Calibrated

## Sleeved Suction Coiled Catheter Kits with Cup

50/case

| Item#           | Catheter Size |
|-----------------|---------------|
| <b>TRATK410</b> | 10fr          |
| <b>TRATK412</b> | 12fr          |
| <b>TRATK414</b> | 14fr          |
| <b>TRATK418</b> | 18fr          |



## Closed Ventilation Suction System

# Standard Switch

Sterile and latex-free

ISO color coded for instant recognition of catheter sizes.




Once the airway is clear, withdraw the catheter completely into the plastic sheath. A black mark on the catheter indicates full withdrawal.



With the unique isolation valve, 72hr closed suction system features safe and simple access for mechanical ventilated patients for 72 hours of continuous use.

Smooth tubing surface and atraumatic rounded catheter tip for minimal tracheal wall injury.



# Closed Ventilation Suction System

## Standard Switch

Sterile and latex-free

### 24 Hour Adult

#### Endotracheal 530mm With Swivel Connector

| Item#        | Description    |
|--------------|----------------|
| VEN72453010M | 10Fr, MDI Port |
| VEN72453012M | 12Fr, MDI Port |
| VEN72453016M | 16Fr, MDI Port |

#### Endotracheal 530mm With T-Piece

| Item#        | Description    |
|--------------|----------------|
| VEN82453010M | 10Fr, MDI Port |
| VEN82453012M | 12Fr, MDI Port |
| VEN82453014M | 14Fr, MDI Port |
| VEN82453016M | 16Fr, MDI Port |

#### Tracheostomy With Swivel Connector

| Item#        | Description           |
|--------------|-----------------------|
| VEN72430510M | 305mm, 10Fr, MDI Port |
| VEN72430512M | 305mm, 12Fr, MDI Port |
| VEN72430514M | 305mm, 14Fr, MDI Port |
| VEN72430516M | 305mm, 16Fr, MDI Port |

#### Tracheostomy With T-Piece

| Item#        | Description           |
|--------------|-----------------------|
| VEN82430510M | 305mm, 10Fr, MDI Port |
| VEN82430512  | 305mm, 12Fr           |
| VEN82430512M | 305mm, 12Fr, MDI Port |
| VEN82430514  | 305mm, 14Fr           |
| VEN82430514M | 305mm, 14Fr, MDI Port |
| VEN82430516  | 305mm, 16Fr           |
| VEN82430516M | 305mm, 16Fr, MDI Port |
| VEN82440512M | 405mm, 12Fr, MDI Port |
| VEN82440514M | 405mm, 14Fr, MDI Port |
| VEN82440516M | 405mm, 16Fr, MDI Port |

### 72 Hour Adult

#### Endotracheal 530mm With Swivel Connector

| Item#        | Description    |
|--------------|----------------|
| VEN77253010M | 10Fr, MDI Port |
| VEN77253012M | 12Fr, MDI Port |
| VEN77253014M | 14Fr, MDI Port |
| VEN77253016M | 16Fr, MDI Port |

#### Tracheostomy 305mm With Swivel Connector

| Item#        | Description    |
|--------------|----------------|
| VEN77230510M | 10Fr, MDI Port |
| VEN77230512M | 12Fr, MDI Port |
| VEN77230514M | 14Fr, MDI Port |
| VEN77230516M | 16Fr, MDI Port |

### 72 Hour Neonatal/Pediatric

#### With Swivel Connector

| Item#       | Description |
|-------------|-------------|
| VEN77230505 | 305mm, 5Fr  |
| VEN77230506 | 305mm, 6Fr  |
| VEN77230507 | 305mm, 7Fr  |
| VEN77230508 | 305mm, 8Fr  |



### 24 Hour Neonatal/Pediatric

#### With Y-Adapter

| Item#        | Description           |
|--------------|-----------------------|
| VEN66000500  | 305mm, 5Fr            |
| VEN66000500M | 305mm, 5Fr, MDI Port  |
| VEN66000600  | 305mm, 6Fr            |
| VEN66000600M | 305mm, 6Fr, MDI Port  |
| VEN66000700  | 305mm, 7Fr            |
| VEN66000700M | 305mm, 7Fr, MDI Port  |
| VEN66000800  | 305mm, 8Fr            |
| VEN66000800M | 305mm, 8Fr, MDI Port  |
| VEN66001000  | 405mm, 10Fr           |
| VEN66001000M | 405mm, 10Fr, MDI Port |

#### With Swivel Connector

| Item#       | Description |
|-------------|-------------|
| VEN72430505 | 305mm, 5Fr  |
| VEN72430506 | 305mm, 6Fr  |
| VEN72430507 | 305mm, 7Fr  |
| VEN72430508 | 305mm, 8Fr  |

## Closed Ventilation Suction System

# Deluxe Switch (one-handed operation)

Sterile and latex-free

### Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

### Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

### MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

### Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.


### Isolation Valve\*

Seals the irrigation chamber to provide a better cleaning of the catheter. Allows for continuous ventilation and prevents backflow of saline during cleaning to ensure patient comfort and safety.

\*Only available with 72 Hour catheter

### Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.



# Closed Ventilation Suction System

## Deluxe Switch (One-handed operation)

Sterile and latex-free

### 24 Hour Adult

#### Endotracheal 540mm

| Item#       | Description                     |
|-------------|---------------------------------|
| VEN160100A  | 10Fr                            |
| VEN160100AM | 10Fr, MDI Port                  |
| VEN162100A  | 10Fr, Wedge + Adapter           |
| VEN162100AM | 10Fr, MDI Port, Wedge + Adapter |
| VEN160120   | 12Fr                            |
| VEN160120M  | 12Fr, MDI Port                  |
| VEN162120   | 12Fr, Wedge + Adapter           |
| VEN162120M  | 12Fr, MDI Port, Wedge + Adapter |
| VEN160140   | 14Fr                            |
| VEN160140M  | 14Fr, MDI Port                  |
| VEN162140   | 14Fr, Wedge + Adapter           |
| VEN162140M  | 14Fr, MDI Port, Wedge + Adapter |
| VEN160160   | 16Fr                            |
| VEN160160M  | 16Fr, MDI Port                  |
| VEN162160   | 16Fr, Wedge + Adapter           |
| VEN162160M  | 16Fr, MDI Port, Wedge + Adapter |

#### Tracheostomy 305mm

| Item#       | Description                     |
|-------------|---------------------------------|
| VEN161100A  | 10Fr                            |
| VEN161100AM | 10Fr, MDI Port                  |
| VEN163100A  | 10Fr, Wedge + Adapter           |
| VEN163100AM | 10Fr, MDI Port, Wedge + Adapter |
| VEN161120   | 12Fr                            |
| VEN163120   | 12Fr, Wedge + Adapter           |
| VEN161140   | 14Fr                            |
| VEN163140   | 14Fr, Wedge + Adapter           |
| VEN161160   | 16Fr                            |
| VEN161160M  | 16Fr, MDI Port                  |
| VEN163160   | 16Fr, Wedge + Adapter           |

### 72 Hour Adult

#### Endotracheal 540mm

| Item#       | Description                     |
|-------------|---------------------------------|
| VEN365100A  | 10Fr                            |
| VEN365100AM | 10Fr, MDI Port                  |
| VEN367100A  | 10Fr, Wedge + Adapter           |
| VEN367100AM | 10Fr, MDI Port, Wedge + Adapter |
| VEN365120   | 12Fr                            |
| VEN365120M  | 12Fr, MDI Port                  |
| VEN367120   | 12Fr, Wedge + Adapter           |
| VEN367120M  | 12Fr, MDI Port, Wedge + Adapter |
| VEN365140   | 14Fr                            |
| VEN365140M  | 14Fr, MDI Port                  |
| VEN367140M  | 14Fr, Wedge + Adapter           |
| VEN365160   | 16Fr                            |
| VEN365160M  | 16Fr, MDI Port                  |
| VEN367160   | 16Fr, Wedge + Adapter           |
| VEN367160M  | 16Fr, MDI Port, Wedge + Adapter |

#### Tracheostomy 305mm

| Item#       | Description                     |
|-------------|---------------------------------|
| VEN366100A  | 10Fr                            |
| VEN366100AM | 10Fr, MDI Port                  |
| VEN368100A  | 10Fr, Wedge + Adapter           |
| VEN368100AM | 10Fr, MDI Port, Wedge + Adapter |
| VEN366120   | 12Fr                            |
| VEN366120M  | 12Fr, MDI Port                  |
| VEN368120   | 12Fr, Wedge + Adapter           |
| VEN368120M  | 12Fr, MDI Port, Wedge + Adapter |
| VEN366140   | 14Fr                            |
| VEN366140M  | 14Fr, MDI Port                  |
| VEN368140   | 14Fr, Wedge + Adapter           |
| VEN368140M  | 14Fr, MDI Port, Wedge + Adapter |
| VEN366160   | 16Fr                            |
| VEN366160M  | 16Fr, MDI Port                  |
| VEN368160   | 16Fr, Wedge + Adapter           |
| VEN368160M  | 16Fr, MDI Port, Wedge + Adapter |

### Neonatal/Pediatric

| Item#        | Description           |
|--------------|-----------------------|
| VEN26000500  | 305mm, 5Fr            |
| VEN26000500M | 305mm, 5Fr, MDI Port  |
| VEN26000600  | 305mm, 6Fr            |
| VEN26000600M | 305mm, 6Fr, MDI Port  |
| VEN26000700  | 305mm, 7Fr            |
| VEN26000700M | 305mm, 7Fr, MDI Port  |
| VEN26000800  | 305mm, 8Fr            |
| VEN26000800M | 305mm, 8Fr, MDI Port  |
| VEN26001000  | 405mm, 10Fr           |
| VEN26001000M | 405mm, 10Fr, MDI Port |



### Swivel Straight Connector

Item#  
**VEN005-20**

Disposable  
22mm ID, 15mm OD  
**20/pack**





## PeDIA - Pediatric Device for Induction of Anesthesia

Item#

**ANEP021Y**

With Mask Adapter Elbow

**12/pack**



### Turn panic into play

For children, undergoing surgery is already unsettling, and the process of induction is often the most traumatizing part. Often kids will kick and scream, believing they're being trapped by the gas or suffocated by the mask. The experiences create negative emotional scarring that can last weeks, months or even years.

Specially designed for children undergoing surgery, the Pediatric Device for Induction of Anesthesia device, or "PeDIA," for short, includes a colorful whistling balloon, providing kids with something fun to focus on and turning a stressful moment into a lighter one.

Whether they're squeezing the balloon and focusing on its whistle or directly putting themselves to sleep by breathing the gas through balloon delivery, the PeDIA balloon captivates children's attention, allowing their active participation in the induction process.


It's a simple way to add a brighter touch to the medical experience, ease anxiety and bring smiles to patients, parents and staff.

This unique approach—presenting an interactive, familiar toy-like device—not only diverts children's focus from the impending procedure but also significantly reduces the anxiety and resistance commonly associated with mask induction. The aim is to decrease the anxiety felt by over 50% of pediatrics as well as parents and staff.

The PeDIA can be used alone as a mask-free device or, for younger children, with an anesthesia mask.

Use with  
mask

Or use  
without  
mask



**SCAN FOR VIDEOS**

