



Product Catalog

Contract MS1087

Respiratory Disposables



Contract MS1049

Non-Invasive Ventilation



Contract MS1097

Medication Delivery



Contract MS1486

Open Suction



Contract MS7709

PeDIA Anesthesia Induction



vizient[®]
Awarded Supplier



Sales Representatives



Egle Sventnickaite
872-258-1654
esventnickaite@sunsethcs.com



Mason Moss
262-788-1862
mmoss@sunsethcs.com



Owen Ruth
312-533-2451
owenr@sunsethcs.com



Cindy Geighes
312-533-2478
cgeighes@sunsethcs.com




Hunter Baudino
262-788-1855
hbaudino@sunsethcs.com




Wyatt Kanemoto
262-788-8941
wkanemoto@sunsethcs.com



Director, Strategic Accounts
Lisa Iannazzone RN, BSN
636-222-3134
lisa@sunsethcs.com



Manager, Product & Business Development
Tom Munar
312-533-2476
tmunar@sunsethcs.com



Contract Manager, Strategic Accounts
Gina Rosenberg
914-837-5215
grosenberg@sunsethcs.com



RRT Clinical Specialist
Krystle Bodine BS, RRT
kbodine@sunsethcs.com

Corporate Office

Phone 872-275-3210

Fax 312-997-9985

141 W Jackson Blvd Ste 1950
Chicago IL 60604

Central Returns/Distribution Center

2725 Alft Ln
Elgin IL 60124

West Returns/Distribution Center

6750 W Van Buren St
Phoenix AZ 85043

Content

Contract MS1049

Non-Invasive Ventilation Page

Adult / Pediatric NIV Masks	4-8
CPAP Tubing	9
Connectors	9
Headgear	10
Chinstraps	10-11

Contract MS1087

Respiratory Disposables Page

Lynx Venturi Oxygen Delivery	12
Oxygen Masks	13
Oxygen Cannula and Accessories	14-17
Oxygen Supply Tubing	18
Pulse Oximeters	18
Firesafe	19
Humidifier	20
Oxygen Connectors and Accessories	20-21
Suction Machine Filters	22
Oxygen Filters	22
Bacteria Filters	22
HEPA Filters	22
Ventilator Filter	23
Air Entrainment Masks	23
Breathing Circuits	24-27
CPAP Tubing	27
CPAP Filters	28-29
HMEs	30
HME Accessories	30
Tracheostomy Masks	30
Trach Ties	31
Berman and Guedel Airways	31
Spirometers	32
Aerosol	32

Contract MS1097

Medication Delivery Page

Inhaler Chamber	33
Nebulizer Compressor	33
Peak Flow Meter	33
Adapters	33
Nebulizer Kits	34

Contract MS1486

Open Suction Page

Suction Catheters	35
Trach Care Kits	35
Suction Catheter Kits	35-36

Contract MS1486

Closed Suction Page

Standard Switch	37-38
Deluxe Switch	39-40

Contract MS7709

PeDIA Anesthesia Induction Page

PeDIA Pediatric Device for Induction of Anesthesia	41
--	----

Ray Pediatric Masks Comparable to Nonny masks

Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	CM061S
Medium	CM061M
Large	CM061L
Extra Large	CM061XL


Replacement Cushions

Size	Item#
Small	CUCM061S
Medium	CUCM061M
Large	CUCM061L
Extra Large	CUCM061XL

Headgear Only

Size	Item#
Small	HG061S
Large	HG061L

SCAN CODE to see
Ray Full Face Mask
on YouTube



Available Size Guide

Nasal Mask with Headgear, Tube, and Replaceable Cushions

Comes with	Item#
S&M Cushions	CM060-KITS
L&XL Cushions	CM060-KITL

Replacement Cushions

Size	Item#
Small	CUCM060S
Medium	CUCM060M
Large	CUCM060L
Extra Large	CUCM060XL

Headgear Only

Size	Item#
Small	HG060S
Large	HG060L



SCAN CODE to see
Ray Nasal Mask
on YouTube




Sol Full Face CPAP Mask

Mask with Headgear and Replaceable Cushion

10/case

Size	Item#
Small	CM205S
Medium	CM205M
Large	CM205L



A premium mask without the premium price

Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

Quiet vent

Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

Quick release connector

Rapidly and easily separates or connects mask and the tube

Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment



Full Face Mask

With Standard Headgear

Standard Elbow, Non-Vented

Size	Item#
Small	CM041S
Medium	CM041M
Large	CM041L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM041SEE
Medium	CM041MEE
Large	CM041LEE

Single Limb Elbow, Vented, AAV

Size	Item#
Small	CM041SEEV
Medium	CM041MEEV
Large	CM041LEEV

novaplus⁺

Standard Elbow, Non-Vented

Size	Item#
Small	VCM041S
Medium	VCM041M
Large	VCM041L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	VCM041SEE
Medium	VCM041MEE
Large	VCM041LEE

Single Limb Elbow, Vented, AAV

Size	Item#
Small	VCM041SEEV
Medium	VCM041MEEV
Large	VCM041LEEV

Replacement Cushions

Size	Item#
Small	CUCM005S-10
Medium	CUCM005M-10
Large	CUCM005L-10



- Unobtrusive design offers a clear line of sight
- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone cushion increases comfort while maintaining optimal seal

With Cap Style Headgear

Standard Elbow, Non-Vented

Size	Item#
Small	CM041SC
Medium	CM041MC
Large	CM041LC

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM041SEEC
Medium	CM041MEEC
Large	CM041LEEC

Single Limb Elbow, Vented, AAV

Size	Item#
Small	CM041SEEV
Medium	CM041MEEV
Large	CM041LEEV

Sunset Deluxe Nasal CPAP Masks

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM006S
Medium	CM006M
Large	CM006L

HCPCS# A7034



Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.

Nasal Mask

With Standard Headgear

Size	Item#
Small	CM008S
Medium	CM008M
Large	CM008L

novaplus⁺

Size	Item#
Small	VCM008S
Medium	VCM008M
Large	VCM008L



- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

Total Face Mask

Standard Elbow

Size	Item#
Small	CM040S
Medium	CM040M
Large	CM040L

EE Elbow

Size	Item#
Small	CM040SEE
Medium	CM040MEE
Large	CM040LEE



Veraseal 2 Disposable Gel Mask with Headgear



Standard Elbow, Non-Vented

Sleepnet®

Size	Item#
Small	CM50933
Medium	CM50934
Large	CM50935
Extra Large	CM50936

novaplus⁺

Size	Item#
Small	VCM053S
Medium	VCM053M
Large	VCM053L
Extra Large	VCM053XL



Single Limb Elbow, Non-Vented, AAV

Sleepnet®

Size	Item#
Small	CM50937
Medium	CM50938
Large	CM50939
Extra Large	CM50940

novaplus⁺

Size	Item#
Small	VCM054S
Medium	VCM054M
Large	VCM054L
Extra Large	VCM054XL



Single Limb Elbow, Vented, AAV

Sleepnet®

Size	Item#
Small	CM50827
Medium	CM50828
Large	CM50829
Extra Large	CM50830

novaplus⁺

Size	Item#
Small	VCM052S
Medium	VCM052M
Large	VCM052L
Extra Large	VCM052XL

- The Veraseal® 2 disposable full face mask brings AIR°gel® Technology to the acute care environment for a giant step forward in patient comfort. Three color coded elbow designs make choosing the right mask for CPAP, single, or dual limb circuit NIV systems simple. The Veraseal® 2 fits quickly and easily while enhancing patient comfort. Both vented (with leak ports) and non-vented (without leak ports) versions are available for CPAP, Bilevel, and ventilator therapy.

- The AIR°gel® cushion provides the utmost in gentle cushioning to protect from pressure points and irritation. AIR°gel® promotes advanced skin protection, superior performance, and unbeatable comfort, all in a lightweight mask.

- The active headgear connector moves as the patient moves for a reliable seal.

- Breathable headgear lets patients enjoy soft, lightweight stability they can count on for a comfortable fit.

Sleepnet Innova

Nasal Mask with Headgear

Vented

Size	Item#
Small/Medium	CM50180
Medium/Large	CM50181



Size	Item#
Small/Medium	VCM034SM
Medium/Large	VCM034ML

Non-Vented

Size	Item#
Small/Medium	CM50594
Medium/Large	CM50595



Size	Item#
Small/Medium	VCM034NVSM
Medium/Large	VCM034NVML



Innova® CPAP masks feature patented AIR°gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA

Sleepnet Ascend Mask Kit

Full Face System

Item# **CM037F**

Each

Vented

Mask Only

Size	Item#
Small	CM037S
Medium	CM037M
Large	CM037L

Replacement Cushions

Size	Item#
Small	CM037CS
Medium	CM037CM
Large	CM037CL



Nasal System

Item# **CM037N**

Each

Vented

Mask Only

Size	Item#
Small	CM037NS
Medium	CM037NM
Large	CM037NL

Replacement Cushions

Size	Item#
Small	CM037NCS
Medium	CM037NCM
Large	CM037NCL



Sleepnet Phantom®2 and iQ®2

Phantom®2

Item#
CM050

Replacement
Cushion
CUCM050



Replacement
Frame/Tubing
Assembly
CM051

Sleepnet's patented AIR°gel® technology provides maximum comfort while minimizing skin irritation and pressure points on the nasal bridge. Unlike many other gel masks, the unique qualities of AIR°gel® allow the mask to seal with minimal tension. By design AIR°gel® is twice as soft as competing gel products.

iQ®2

Item#
CM049

Replacement
Cushion
CUCM049



Ask us about
our 30-day
mask guarantee

- Sleepnet has the only mask that you can shape to each individual face. Flex, Fit, Forget!
- Cushions may be removed and replaced without purchasing a new mask.

- Magnetic Clips — Easy on and off!
- 360° Swivel — Less tube drag for active sleepers.
- One size cushion — Larger fitting range.

MiniMe 2 Pediatric Nasal Mask



Vented

Size	Item#
Extra Small	CM60216
Small	CM60217
Medium	CM60218
Large	CM60221

Non-Vented

Size	Item#
Extra Small	CM60256
Small	CM60257
Medium	CM60258
Large	CM60261



Vented

Size	Item#
Extra Small	VCM055XS
Small	VCM055S
Medium	VCM055M
Large	VCM055L

Non-Vented

Size	Item#
Extra Small	VCM056XS
Small	VCM056S
Medium	VCM056M
Large	VCM056L



- Maximum comfort, minimum weight. Patented AIR°gel makes MiniMe masks some of the most gentle fitting masks on the market.
- The shapeable shell mask can be easily molded to fit individual needs.
- 100% designed and manufactured in the USA

Pediatric SleepWeaver® CPAP Mask

Item#
CM101325



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.

- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age




MiniMe Pediatric Nasal Mask



Headgear Size	Item#	OEM#
Small	CM004S	60213
Medium	CM004M	60214
Large	CM004L	60215

HCPCS# A7034



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

SleepWeaver® Advance Nasal Mask

Item#
CM017B

Mask with Adjustable Headgear

- Easy to store and pack for travel

OEM# 100274
HCPCS# A7034



NIV Tubing

Item#

TUB06 Each
SUH0610 10/pack

19mm with 22mm cuffs



- Sunset's most popular item
- ISO standard 22mm cuffs
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

Standard Elbow, Non-Vented

Item# **CAP009**
10/case



Single Limb Elbow, Non-Vented, AAV

Item# **CAP008**
10/case



CPAP Exhalation Port

Item# **VEN002-10**
10/case



Single Use Exhalation Port

Item# **VEN003**
10/case



Standard 4-Point Headgear

Item#

HG003S Small
HG003M Medium
HG003L Large



Cap Style Headgear

Item# **HG002**



Universal Neoprene Chinstrap

Item#

CS001
 Universal Chinstrap

HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows



Deluxe Chinstrap

Size	Item#
26in	CS004S
28in	CS004L

Compare to
 Respironics Deluxe

OEM# 302425

HCPCS# A7036

- Latex-free
- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort



Comfort Chinstrap

Item# **CS005**

Compare to
 ResMed Sullivan

OEM# 16015

HCPCS# A7036

- Latex-free
- One size fits most
- Easy to use —attaches to itself at the top of the head
- Soft cloth chin cup adds comfort



Premium Chinstrap

Item# **CS006**

Compare to
 Respironics Premium

OEM# 1012911

HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.



Ruby Style Adjustable Chinstrap




Size	Item#
Regular	CS007-ADJ
Extra Large	CS007-ADJXL

HCPCS# A7036

- Latex-free
- Fully adjustable
- Super soft materials for extra comfort



Ruby Style Chinstrap



For fitting, measure completely around the head, as indicated by the dashed line above.

Size	Item#
21.5-22.5in	CS007S
22-23in	CS007M
23-24in	CS007L

Compare to Puresom Ruby

OEM# TMS-09 • HCPCS# A7036

- Latex-free
- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

Halo Style Chinstrap

Item# **CS025**

Compare to BreatheWear Halo Chinstrap

OEM# H3407

HCPCS# A7036



Lynx Venturi Oxygen Delivery


Low Concentration

Item# **RES2400**

Each

Adapter Only **RES2401C**

- Design replaces multiple venturi adapters
- Connects to oxygen tubing and regulates flow
- Provides precision FiO₂ delivery from 24–60% with less than 2% variance during mild, moderate, and severe respiratory distress
- Excellent CO₂ washout



Unprecedented simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

Do no harm

Quicker de-escalation of O₂ levels and prevention of excessive CO₂ in the patient.

Reduction in waste

Remains with patient throughout EMS and hospital journey, replacing multiple adapters which may get lost or contaminated during transition.

High Concentration

Item# **RES2401**


Each

Adapter Only **RES241HC**

- Allows rapid de-escalation of FiO₂ in the EMS ambulance or ED to prevent CO₂ buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia (PaO₂/FiO₂) allows early intervention and better outcomes
- Excellent CO₂ washout

Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Capabilities could include replacing non-rebreathers and/or partial rebreather products



Unmatched precision

Provides specific levels of oxygen with less than 2% variance. Comparable products vary by as much as 30%.

Medium Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector
RES2100NT	Adult	None	Standard
RES2100U	Adult	7 ft	Universal
RES2100	Adult	7 ft	Standard
RES2203P	Pediatric	None	Standard
RES2100P	Pediatric	7 ft	Standard
RES2203A	Adult	None	Standard
RES2203A7	Adult	7 ft	Standard
RES2203I	Infant	None	Standard
RES2203I7	Infant	7 ft	Standard
RES2203P7	Pediatric	7 ft	Standard

- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



High Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length
RES2202A	Adult	None
RES2202A7	Adult	7 ft
RES2202I	Infant	None
RES2202I7	Infant	7 ft
RES2202P	Pediatric	None
RES2202P7	Pediatric	7 ft

Non-Rebreather Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	Safety Vent
RES2102I	Infant	7 ft	Standard	No
RES2102IM	Infant	None	Standard	No
RES2102P	Pediatric	7 ft	Standard	No
RES2102PS	Pediatric	7 ft	Standard	Yes
RES2102PSM	Pediatric	None	Standard	Yes
RES2102PSU	Pediatric	7ft	Universal	Yes
RES2102PU	Pediatric	7 ft	Universal	No
RES2102S	Adult	7 ft	Standard	Yes
RES2102SM	Adult	None	Standard	Yes
RES2102U	Adult	7 ft	Universal	No
RES2102	Adult	7 ft	Standard	No



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape

Humidoflo High Flow Nasal Cannula


10/case

Item#	Size	F&P Adapter
HFN001AS	Small	No
HFN001AM	Medium	No
HFN001AL	Large	No
HFN001ASA	Small	Yes
HFN001AMA	Medium	Yes
HFN001ALA	Large	Yes

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respironics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)



With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerable due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O₂ concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.



Adult Oxygen Cannula

Item#	Tubing Length	Connector	Qty
RES1100	None	Standard	50
RES1104	4 ft	Standard	50
RES1107N	7 ft	Standard	50
RES1107U	7 ft	Universal	50
RES1107NU	7 ft	Universal	50
RES1107	7 ft	Standard	50
RES1115U	15 ft	Universal	50
RES1115	15 ft	Standard	50
RES1125	25 ft	Standard	25
RES1135	35 ft	Standard	20
RES1150U	50 ft	Universal	20
RES1150Z	50 ft	Standard	20



Soft Cannula

Item#	Tubing Length	Connector	Qty
RES1100S	None	Standard	50
RES1107SU	7 ft	Universal	50
RES1107S	7 ft	Standard	50
RES1114S	14 ft	Standard	50
RES1114SU	14 ft	Universal	50
RES1115S	15 ft	Standard	50
RES1125S	25 ft	Standard	25
RES1125SU	25 ft	Universal	25

Flared Cannula

Item#	Tubing Length	Connector	Qty
RES1100F	None	Standard	50
RES1107F	7 ft	Standard	50
RES1107FU	7 ft	Universal	50
RES1114F	14 ft	Standard	50
RES1125F	25 ft	Standard	25
RES1150F	50 ft	Standard	20

Adult Cannula with Ear Cushions

Standard Touch, Standard Tip

Item#	Tubing Length	Connector	Qty
RES1100EC	None	Standard	50
RES1104EC	None	Standard	50
RES1107ECU	7 ft	Universal	50
RES1107EC	7 ft	Standard	50
RES1114EC	14 ft	Standard	50
RES1115EC	15 ft	Standard	50
RES1125EC	25 ft	Standard	25



Cannula Ear Cushions

Item# **RES020**
50 Pairs/case

- Add cushions to existing cannula



High Flow Soft Cannula

25/case

Length	Item#
4ft	RES1104SHF
7ft	RES1107SHF
15ft	RES1115SHF
25ft	RES1125SHF




- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

7 ft Cannula With CO₂ Sampling Line

50/case

Item#	Size	Connector	Luer	Filter
RES8007	Adult	Standard	Male	No
RES8007F	Adult	Standard	Male	Yes
RES8007I	Infant	Standard	Male	No
RES8007IF	Infant	Standard	Male	Yes
RES8007IU	Infant	Universal	Male	No
RES8007P	Pediatric	Standard	Male	No
RES8007PF	Pediatric	Standard	Male	Yes
RES8007PU	Pediatric	Universal	Male	No
RES8007U	Adult	Universal	Male	No
RES8107	Adult	Standard	Female	No
RES8107F	Adult	Standard	Female	Yes
RES8107I	Infant	Standard	Female	No
RES8107IF	Infant	Standard	Female	Yes
RES8107P	Pediatric	Standard	Female	No
RES8107PF	Pediatric	Standard	Female	Yes



Pediatric, Infant, and Neonate Cannula

50/case


Item#	Size	Tubing Length	Connector	Touch	Ear Cushions
RES1100P	Pediatric	None	Standard	Standard	No
RES1107I	Infant	7 ft	Standard	Standard	No
RES1107NE	Neonate	7 ft	Standard	Standard	No
RES1107P	Pediatric	7 ft	Standard	Standard	No
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes
RES1107SI	Infant	7 ft	Standard	Soft	No
RES1107SIU	Infant	7 ft	Universal	Soft	No
RES1107SNE	Neonate	7 ft	Standard	Soft	No
RES1107SP	Pediatric	7 ft	Standard	Soft	No
RES1107SPU	Pediatric	7 ft	Universal	Soft	No
RES1114SP	Pediatric	14 ft	Standard	Soft	No
RES1114SPU	Pediatric	14 ft	Universal	Soft	No

Secure Cannula Stickers

Item# **RES080**

25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



Oxygen Supply Tubing

Item#	Length	Connector	Color	Bore	Qty
RES3007SB	7 ft	Standard	None	Smooth	50
RES3007U	7 ft	Universal	None	Channeled	50
RES3007	7 ft	Standard	None	Channeled	50
RES3007G	7 ft	Standard	Green	Channeled	50
RES3015G	15 ft	Standard	Green	Channeled	50
RES3015	15 ft	Standard	None	Channeled	50
RES3025G	25 ft	Standard	Green	Channeled	25
RES3025U	25 ft	Universal	None	Channeled	25
RES3025	25 ft	Standard	None	Channeled	25
RES3030	30 ft	Standard	None	Channeled	20
RES3030G	30 ft	Standard	Green	Channeled	20
RES3040G	40 ft	Standard	Green	Channeled	20
RES3050G	50 ft	Standard	Green	Channeled	20
RES3050	50 ft	Standard	None	Channeled	20

Universal Connector



Finger Pulse Oximeter

Item# **RES5102**
Each





- Oversized 1.3" OLED display face
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- Includes lanyard and manual

Firesafe Thermal Shut-Off Valve

Item#
RES010-10
Bidirectional
10/case




- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire
- Smaller
- Laser etched
- Expiration date
- Human readable UDI
- Designed to last 5 years



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

Fatalities per 100,000 patients



FireSafe Nozzle

Item#
RES010FND
20/pack



FireSafe Tubing Connector

Item#
RES010TC 100/case



Disposable Bubble Humidifier 500 ml

Item#

RES050

Up to 15 LPM • 6 psi • 500 ml

50/case

HCPCS# E0560

OEM# HUM001



- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient

Bacterial/Viral HME Filter

Item#

BF2200**5/pack**

Oxygen Tubing Connector

Item# **RES012**

Fits 3/16in and 1/4in ID Tubing

10/pack

- Barbed ends ensure secure, leak-free connections

Christmas Tree Connector

Item# **RES001**

Green without Swivel

10/pack**RES001-30 30/pack**

- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Item# **RES002**

Green with Swivel

10/pack

Circuit Connectors

Item# **RES014**

22mm OD x 22mm OD

With Lip

10/pack

Humidifier Tubing Adapter

Item# **RES016**

22mm ID x 22mm OD

with Lip

10/pack

Bacterial/Viral HME Filter - Adult

Item#

BF104**Each**

- Bacterial efficiency 99.9999%
- Viral efficiency 99.999%
- Humidification efficiency 37.5 mg/l H₂O (500 ml tidal volume)
- Coupled with a hydrophobic bacterial viral filter membrane



Wrenches

Item# **RES031-1**

With Chain

1/pack

Tubing Reducer Connector

Item# **RES011**

Fits 3/16in, 1/4in,
5/16in, 3/8in, and 7/16in
ID Tubing

10/pack

- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing



Circuit Connectors

Item# **RES013**

22mm OD x 15mm ID/22mm OD

10/pack



Oxygen Bleed-In Adapter

Item# **RES019**

22mm with 6mm
Rotating Stem

10/pack

RES019-1 1/Pack



Oxygen Tubing Y-Connector

Item# **RES017**

Fits 3/16in and 1/4in ID Tubing
(All 3 sides are the same)

10/pack

- Connects 2 oxygen sources for higher flows



Connector with 6mm Oxygen Stem

Item# **RES058**

22mm with 6mm Stem

10/pack



Humidifier Connectors

Item# **HUM008**

L-Connector 12in

10/pk

HUM008-1 1/Pack

Item# **HUM008L**

L-Connector 21in

10/pk

Item# **HUM008S**

Straight Connector 9in

10/pk



Oxygen Tubing Connector — Swivel

Item# **RES018**

Male/Male

10/pack

RES018-1 1/Pack

- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing




Swivel Connector Silencer for Ventilator Circuits and CPAP

Item# **VEN002**

10/pack



Suction Machine Inline Filters



Item# **BF1412**
Final Bacteria Filter
1/4in to 1/2in Stepped Barb
on Both Sides
10/pack



Item# **BF3812**
3/8in to 1/2in
10/pack



Item# **BF1312**
1/3in to 1/2in
1/pack



Item# **BF516TS**
5/16in Barbed with
Threaded NPT
10/pack

HEPA Filter for 3B Stratus

Item# **BF6020**
HEPA Filter
Each
OEM# STR1020



6in HEPA Capsule with Removable Sticker

Item# **BFH6S**
Threaded
Each



3B RESmart Style Filter

Item# **CF8004**
Disposable Filter
3/pack • 3⁵/₁₆" x 1⁵/₈"
HCPCS# A7038



Oxygen Filters



Item# **BF316**
Final Bacteria Filter • 3/16 Stepped
Barb and 1/4in Straight Barb on
Both Sides
10/pack
BF316-1 1/pack



Item# **BF516**
1/4in Barbed on Both Sides
10/pack



Item# **BF1438**
1/4in to 3/8in
10/pack

Invacare Platinum 5 and 10 Style Filter

Item# **OF9002**
Foam Cabinet Filter
9" x 27⁷/₈" x 1¹/₂"
10/pack
OEM#1107412



Invacare Perfecto2 Style Filter

Item# **OF9004**
Foam Cabinet Filter
6⁷/₈" x 33³/₄" x 1¹/₂"
10/pack
OEM# 1143492



Extended Life Intake Bacteria Filter

Item# **OF1025D**
Each



DeVilbiss/Drive SolAiris III, 303 Series Style Filter

Item# **OF5004**
Foam Cabinet Filter
33³/₄" x 33³/₄" x 1¹/₂"
10/pack
OEM# 30302605



Filter for ResMed Astral Ventilator

Item# **VF003-4**
4/pack



Item# **VF003-H**
Filter Housing
Each



Item# **VF003-KIT**
4-Filter Kit
w/Filter Housing
OEM# 27939



Adult Air Entrainment Mask

50/case

Item#	Size	Kit	Tubing Length	Connector	Flex Tube
RES2400A	Adult	Mask Only	7 ft	Standard	Yes
RES2400AK	Adult	Yes	7 ft	Standard	Yes
RES2400AKU	Adult	Yes	7 ft	Universal	Yes
RES2400AU	Adult	Mask Only	7 ft	Universal	Yes
RES2400P	Pediatric	Mask Only	7 ft	Standard	Yes
RES2400PK	Pediatric	Yes	7 ft	Standard	Yes
RES2400PKU	Pediatric	Yes	7 ft	Universal	Yes
RES2400PU	Pediatric	Mask Only	7 ft	Universal	Yes



90-Degree CPAP Tubing Elbow

Item# **CAP1027**
Each

The Sunset 90-Degree CPAP Tubing Elbow Adapter connects standard 22mm or slim 15mm CPAP tubes to most CPAP machines at a 90-degree angle, relieving stress on the tubing connection and creating a more convenient pathway to the bed.

Comparable to OEM# 37394




Passive Non-Heated Breathing Circuits

Each

For use with Respiration Trilogy

Exhalation Port	Item#	OEM#
Swivel	VC1000	1073228 & 1073236
Fixed	VC1100	1069467 & 1069468



Non-Heated Breathing Circuit

For use with Respiration Trilogy


Item# **VC1001**

Active Adult, No Water Trap
22mm cuffs

30/case

OEM#

1073223 & 1073232



Adult Passive Heated Circuit with Exhalation Port

For use with Philips
Respiration Trilogy

Item#

VC1003

VC1003P

Pediatric Size


30/case

Comparable to OEM# 1076046,
1073227, 1073235



Includes:

- 1 VC1114 Passive Heated Breathing Circuit
- 1 VEN002 Swivel Exhalation Port



Passive Circuit w/Exhalation Port and Bacterial/Viral Filter

With 22mm cuffs

Item#

VC1005

Each

OEM# 1147378, 1132340, 7615





Non-Heated Breathing Circuits

For use with Respiration Trilogy
22mm cuffs


30/case

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
VC1001	1073223 & 1073232	Active	No	Adult	22mm
VC1004	1069443 & 1069444	Passive	Yes	Adult	22mm
VC1004P	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#
VC1001
Active Adult
OEM#
1073223
1073232




To use active circuits with Respiration Trilogy, simply cut ends off pressure and exhalation lines.



Item#
VC1004
Passive Adult
OEM#
1069443
1069444



Item#
VC1004P
Passive Pediatric
OEM#
1073224
1073233



Passive Heated Breathing Circuit

10/case


For use with MR850 Humidifier

Item#

VC1114

VC1114P

Pediatric Size



Adult NIV Patient Circuit


Item# **VC1160**

Compatible with V60, No Filters

10/case

Comparable to

OEM# 1069210



Adult NIV Single Limb Patient Circuit with Main Flow Filter

Item# **VC1161**

Includes 1 BF2200 22mm Bacteria Disk Filter

10/case




Adult NIV Single Limb Patient Circuit with Main Flow Filter and Exhalation Port Filter

Item# **VC1162**

Includes 2 BF2200 22mm Bacteria Disk Filters

10/case



Non-Heated Dual Limb Breathing Circuits

30/case

Item#	Size	Humidifier Limb	Water Trap
VC2001P	Pediatric	No	No
VC2001	Adult	No	No
VC2002P	Pediatric	Yes	No
VC2002	Adult	Yes	No
VC2003	Adult	Yes	Yes



Slim CPAP Tube

Item#

TUB006SS

6ft x 15mm with
22mm cuffs

Each

HCPCS# A7037



Highly durable

Strain relief cuff decreases separation and tearing

Virtually no memory

Designed to be the most flexible tube available

Minimal weight

One of the lightest tubes on the market at 2.6 ounces

Great value

Affordable and sustainable

CPAP Tubing

Item#

TUB008 8ft

TUB010 10ft

19mm with 22mm cuffs

Each

HCPCS# A7037



- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

Heated Tube for Respirationics

Item#


TUB06HPR

Heated Tube for Philips Respirationics
Dreamstation Series and System One

Each

OEM# HT15

HCPCS# A4604



- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off

ResMed S9 and AirSense 10 Series Filters

Item# CF2107

Disposable Filter

3/pack • 2 1/8" x 1 3/8"

HCPCS# A7038

CF2107-1 1/pack OEM# 36850

CF2107-2 2/pack OEM# 36851

CF2107-50 50/pack OEM# 36853

Item# CF2107HA

Hypoallergenic Filter

3/pack • 2 1/8" x 1 3/8"

HCPCS# A7038

CF2107HA-1 1/pack OEM# 36855

CF2107HA-2 2/pack OEM# 36856

CF2107HA-50 50/pack
OEM# 36858



ResMed AirSense 11 Series Filters

Item# CF2111

Disposable Filter

3/pack

HCPCS# A7038

CF2111-1 1/pack OEM# 39300

CF2111-2 2/pack OEM# 39301

CF2111-50 50/pack OEM# 39303

Item# CF2111HA-3

Hypoallergenic Filter

3/pack

1 1/4" x 1 3/16"

HCPCS# A7038

CF2111HA-1 1/pack

CF2111HA-2 2/pack



3B Luna G3 Style Filters

Item# CF8006-6

Disposable Filter

6/pack

HCPCS# A7038

CF8006-1 1/pack

CF8006-2 2/pack

Item# CF8006F-1

Reusable Filter

1/pack

HCPCS# A7039



3B Luna and Luna II Style Filters

Item# CF8005

Disposable Filter

12/pack • 2" x 1 3/16"

HCPCS# A7038

Item# CF8005F-1

Reusable Filter

1/pack • 2" x 1 3/16"

HCPCS# A7039

CF8005F-2 2/pack




3B Stratus 5 Style Filter

Item# OF6030

Foam Cabinet Filter

10/pack

OEM# STR1030



Respironics DreamStation 2 Series Filters

Item# **CF1012**

Disposable Filter

6/pack

OEM# 1142831

HCPSC# A7038

CF1012-1 1/pack OEM# 1142828

CF1012-2 2/pack OEM# 1142829

CF1012-3 3/pack

CF1012-50 50/pack



Item# **CF1012F**

Reusable Filter

2/pack

OEM# 1142687

HCPSC# A7039

CF1012F-1 1/pack



Respironics DreamStation Series Filters

Item# **CF1010**

Disposable Filter

6/pack

OEM# 1122519

HCPSC# A7038

CF1010-1 1/pack

CF1010-2 2/pack



Item# **CF1010F**

Reusable Filter

2/pack

OEM# 112446 (1-pack)

HCPSC# A7039

CF1010F-1 1/pack

CF1010F-6 6/pack



Flex HME

Item# **RES027FLEX**
40/case



TV Range (mL)	300-1500
Moisture Output (mg H2O/L)	33 mg/L
Weight (g)	17
Resistance (cm H2O)	0.5 @ 30 LPM 1.34 @ 60 LPM
Dead Space (mL)	75
Connection (mm)	22M,15F

Heat Moisture Exchangers

30/case

Item#					
Item#	RES027 T-HME	RES0270P T-HME with Oxygen Port	RES0271L Inline HME — Large Volume	RES027NE Inline HME — Neonatal	RES027HV2 Foam HME — Adult/Child
TV Range (mL)	250-1000	150-1250	50-250	15-50	250-1000
Moisture Output (mg H2O/L)	30, Vt = 0.1L	36.3 mg/L	31, Vt = 0.1L	30	29
Weight (g)	4.5	7.4	10	4.5	21
Resistance (cm H2O)	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.2 @ 30 LPM	0.4 @ 10 LPM	0.5 @ 10 LPM	0.52 @ 30 LPM 1.32 @ 60 LPM 2.54 @ 90 LPM
Dead Space (mL)	12.5	7	13	2.5	29
Connection (mm)	15F	15F 5.5mm O2 port	15M,22M/15F	15M,15F	22F/15M, 22M/15F

Adapter for T-HME

Item# **RES027A**
Each



Flexible Connector

Item# **TRAA004**
15mm F x 15mm M
Length 4.9-7.5cm
50/case
OEM# 3215



Tracheostomy Mask

50/case

Item# **RES2130**
Adult Regular Size

RES2130L
Adult Large Size

RES2230
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O2 supply tubing

Trach Ties

20/case


One-piece

- Easy on

Item#

RES029-17 2.5x43cm (17in)

RES029P-11 2.5x30cm (11.5in)



- Designed to provide secure stabilization of the tracheostomy tube

- Made from soft 60% Polyester, 20% Cotton, 20% Nylon

Two-piece

- Easy off


Item#

RES030P-10 2x25.5cm (10in)

RES030P-14 2x35.5cm (14in)

RES030P-18 2x45.5cm (18in)

RES030-20 2.5x51cm (20in)



Color Coded Dual Channel Berman Oral Airway

50/pack

Item#	Size	Color
AOA40BC-50	40mm	Pink
AOA50BC-50	50mm	Blue
AOA60BC-50	60mm	Black
AOA70BC-50	70mm	White
AOA80BC-50	80mm	Green
AOA90BC-50	90mm	Yellow
AOA100BC-50	100mm	Red
AOA110BC-50	110mm	Orange



- Smooth surface & rounded tip for improved patient comfort

- 100% Polyethylene (PE)
- Latex free

Oropharyngeal Guedel Airway

10/pack

Item#	Size	Color
AOP500GC	50mm	Blue
AOP601GC	60mm	Black
AOP702GC	70mm	White
AOP803GC	80mm	Green
AOP904GC	90mm	Yellow
AOP1005GC	100mm	Red



- Designed to overcome associated risks of loose or detached bite blocks
- 100% Polyethylene (PE)

- Smooth surface & rounded tip for improved patient comfort
- Latex free

Volumetric Incentive Spirometer

12/case

Item# **RES055**
2500ml

RES056
5000ml

RES055P
Pediatric

- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise.
- Dual-sided, clear and legible identification of volumes.
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration.
- Compact and ergonomic design.



SCAN FOR VIDEO



Three Ball Spirometer

Item# **RES057**
12/case

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient



- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

Valved T-Adapter, Aerosol 22mm ID x 22mm OD

Item# **NEB9001T**
30/case



Aerosol Mask

50/case

Item# **RES2110**
Adult Size

Item# **RES2110F**
Adult Size
w/ Flex Tube

Item# **RES2210**
Pediatric Size



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

SperaChamber

Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler

Item#
RES195-10
10/case

RES195-50
50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



Adult Mask



SperaChamber with Mask
10/case

Item#	Size
RES195MS	Small (Infant)
RES195MM	Medium (Pediatric)
RES195ML	Large (Adult)

Mask Only

Item#	Size	Qty
RES195MASK-S	Small (Infant)	Each
RES195MASK-S-10	Small (Infant)	10/box
RES195MASK-M	Medium (Pediatric)	Each
RES195MASK-M-10	Medium (Pediatric)	10/box
RES195MASK-L	Large (Adult)	Each
RES195MASK-L-10	Large (Adult)	10/box

Nebulizer Compressor

Each

Item#
NEB100

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Foam Filters (BF018 x5)

NEB100-R

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

Features

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

SCAN FOR VIDEO




Peak Flow Meter

Item# **RES052**
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**
Pediatric Mouthpiece for Peak Flow Meter
10/case

Plastic Intake Filter with Head For Respiration Inspiration 626 & 929

Item# **BF016**
10/pack



Valved T-Adapter

Item# **VEN001**
22F-18F-22M/15F
10/pack



Disposable Nebulizer Kits with Mouthpiece

50/case

All include 6cc jet nebulizer medicine cup

Accessory options

- T-piece
- Mouthpiece
- Corrugated tube
- 7ft supply tube
- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal



Item#	Size	Accessories
RES091	Adult	All
RES091NA	Adult	No corrugated tube
RES091U	Adult	All, Plus Universal Connector

Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**

Disposable
Single Patient
50/case

RES096U With Universal Connector

Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube



Disposable Nebulizer Kits with Mask

50/case

All include 6cc jet nebulizer medicine cup

Accessory options

- Mask
- 7ft supply tube

Item#	Size	Accessories
RES090	Adult	With Mask
RES092	Pediatric	All



- Disposable single patient use for long term delivery procedures
- 6cc bottle is easy to assemble and disassemble


Medicine Cup
Item# **RES090A**
50/case

Suction Catheters - Calibrated

50/case

• Sterile

Item#	Size	Tip	Valve	Type
TRASC106C	6fr	DeLee	Chimney	Coiled
TRASC108C	8fr	DeLee	Chimney	Coiled
TRASC108S	8fr	DeLee	Chimney	Straight
TRASC110C	10fr	Whistle	Y	Coiled
TRASC110S	10fr	Whistle	Y	Straight
TRASC112C	12fr	Whistle	Y	Coiled
TRASC112S	12fr	Whistle	Y	Straight
TRASC114C	14fr	Whistle	Y	Coiled
TRASC114CVS	14fr	Cloude	Chimney	Straight
TRASC114S	14fr	Whistle	Y	Straight
TRASC116C	16fr	Whistle	Y	Coiled
TRASC118C	18fr	Whistle	Y	Coiled
TRASC118S	18fr	DeLee	Chimney	Straight



Trach Care Kits

50/case

Sunset Trach Care Kit

Item# **TRATK003**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- Gauze
- 2 Cotton Tips
- Brush,
- Twill Tape
- Trach Dressing

Sunset Mini Trach Kit

Item# **TRATK004**

Includes:

- 3 Pipe Cleaners
- 4 Gauze
- 2 Cotton Tips
- Brush
- Twill Tape
- Dressing

Sunset Trach Care Kit with Forceps

Item# **TRATK005**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- 3 Gauze
- 3 Cotton Tips
- Brush
- Trach Dressing

Suction Catheter Kits with Cup and Gloves

50/case

• Sterile

Item#	Catheter Size	Tip	Valve
TRATK206	6fr	DeLee	Chimney
TRATK208	8fr	Whistle	Y
TRATK210	10fr	Whistle	Y
TRATK212	12fr	Whistle	Y
TRATK214	14fr	Whistle	Y
TRATK218	18fr	Whistle	Y



Pediatric Suction Catheter Kits

50/case

Item#	Catheter Size	Cup	
TRATK306	6fr	No	• Sterile • Y-Valve DeLee Catheter • Includes Gloves
TRATK308	8fr	No	
TRATK408	8fr	Yes	

Mini Catheter Tray with Gloves


50/case

Item#	Catheter Size	
TRATK310	10fr	• Coiled Catheter • Calibrated
TRATK312	12fr	
TRATK312VCVN	18fr, Vinyl Gloves	
TRATK314	14fr	
TRATK316	16fr	
TRATK318	18fr	

Sleeved Suction Coiled Catheter Kits with Cup

50/case

Item#	Catheter Size
TRATK410	10fr
TRATK412	12fr
TRATK414	14fr
TRATK418	18fr




Closed Ventilation Suction System

Standard Switch

Sterile and latex-free

ISO color coded for instant recognition of catheter sizes.



Closed Ventilation Suction System

Standard Switch

Sterile and latex-free

24 Hour Adult

Endotracheal 530mm With Swivel Connector

Item#	Description
VEN72453010M	10Fr, MDI Port
VEN72453012M	12Fr, MDI Port
VEN72453016M	16Fr, MDI Port

Endotracheal 530mm With T-Piece

Item#	Description
VEN82453010M	10Fr, MDI Port
VEN82453012M	12Fr, MDI Port
VEN82453014M	14Fr, MDI Port
VEN82453016M	16Fr, MDI Port

Tracheostomy With Swivel Connector

Item#	Description
VEN72430510M	305mm, 10Fr, MDI Port
VEN72430512M	305mm, 12Fr, MDI Port
VEN72430514M	305mm, 14Fr, MDI Port
VEN72430516M	305mm, 16Fr, MDI Port

Tracheostomy With T-Piece

Item#	Description
VEN82430510M	305mm, 10Fr, MDI Port
VEN82430512	305mm, 12Fr
VEN82430512M	305mm, 12Fr, MDI Port
VEN82430514	305mm, 14Fr
VEN82430514M	305mm, 14Fr, MDI Port
VEN82430516	305mm, 16Fr
VEN82430516M	305mm, 16Fr, MDI Port
VEN82440512M	405mm, 12Fr, MDI Port
VEN82440514M	405mm, 14Fr, MDI Port
VEN82440516M	405mm, 16Fr, MDI Port

72 Hour Adult

Endotracheal 530mm With Swivel Connector

Item#	Description
VEN77253010M	10Fr, MDI Port
VEN77253012M	12Fr, MDI Port
VEN77253014M	14Fr, MDI Port
VEN77253016M	16Fr, MDI Port

Tracheostomy 305mm With Swivel Connector

Item#	Description
VEN77230510M	10Fr, MDI Port
VEN77230512M	12Fr, MDI Port
VEN77230514M	14Fr, MDI Port
VEN77230516M	16Fr, MDI Port

72 Hour Neonatal/Pediatric

With Swivel Connector

Item#	Description
VEN77230505	305mm, 5Fr
VEN77230506	305mm, 6Fr
VEN77230507	305mm, 7Fr
VEN77230508	305mm, 8Fr



24 Hour Neonatal/Pediatric

With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600	305mm, 6Fr
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700	305mm, 7Fr
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800	305mm, 8Fr
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000	405mm, 10Fr
VEN66001000M	405mm, 10Fr, MDI Port

With Swivel Connector

Item#	Description
VEN72430505	305mm, 5Fr
VEN72430506	305mm, 6Fr
VEN72430507	305mm, 7Fr
VEN72430508	305mm, 8Fr

Closed Ventilation Suction System

Deluxe Switch (one-handed operation)

Sterile and latex-free

Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.


Isolation Valve*

Seals the irrigation chamber to provide a better cleaning of the catheter. Allows for continuous ventilation and prevents backflow of saline during cleaning to ensure patient comfort and safety.

*Only available with 72 Hour catheter

Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.



Closed Ventilation Suction System

Deluxe Switch (One-handed operation)

Sterile and latex-free

24 Hour Adult

Endotracheal 540mm

Item#	Description
VEN160100A	10Fr
VEN160100AM	10Fr, MDI Port
VEN162100A	10Fr, Wedge + Adapter
VEN162100AM	10Fr, MDI Port, Wedge + Adapter
VEN160120	12Fr
VEN160120M	12Fr, MDI Port
VEN162120	12Fr, Wedge + Adapter
VEN162120M	12Fr, MDI Port, Wedge + Adapter
VEN160140	14Fr
VEN160140M	14Fr, MDI Port
VEN162140	14Fr, Wedge + Adapter
VEN162140M	14Fr, MDI Port, Wedge + Adapter
VEN160160	16Fr
VEN160160M	16Fr, MDI Port
VEN162160	16Fr, Wedge + Adapter
VEN162160M	16Fr, MDI Port, Wedge + Adapter

Tracheostomy 305mm

Item#	Description
VEN161100A	10Fr
VEN161100AM	10Fr, MDI Port
VEN163100A	10Fr, Wedge + Adapter
VEN163100AM	10Fr, MDI Port, Wedge + Adapter
VEN161120	12Fr
VEN163120	12Fr, Wedge + Adapter
VEN161140	14Fr
VEN163140	14Fr, Wedge + Adapter
VEN161160	16Fr
VEN161160M	16Fr, MDI Port
VEN163160	16Fr, Wedge + Adapter

72 Hour Adult

Endotracheal 540mm

Item#	Description
VEN365100A	10Fr
VEN365100AM	10Fr, MDI Port
VEN367100A	10Fr, Wedge + Adapter
VEN367100AM	10Fr, MDI Port, Wedge + Adapter
VEN365120	12Fr
VEN365120M	12Fr, MDI Port
VEN367120	12Fr, Wedge + Adapter
VEN367120M	12Fr, MDI Port, Wedge + Adapter
VEN365140	14Fr
VEN365140M	14Fr, MDI Port
VEN367140M	14Fr, Wedge + Adapter
VEN365160	16Fr
VEN365160M	16Fr, MDI Port
VEN367160	16Fr, Wedge + Adapter
VEN367160M	16Fr, MDI Port, Wedge + Adapter

Tracheostomy 305mm

Item#	Description
VEN366100A	10Fr
VEN366100AM	10Fr, MDI Port
VEN368100A	10Fr, Wedge + Adapter
VEN368100AM	10Fr, MDI Port, Wedge + Adapter
VEN366120	12Fr
VEN366120M	12Fr, MDI Port
VEN368120	12Fr, Wedge + Adapter
VEN368120M	12Fr, MDI Port, Wedge + Adapter
VEN366140	14Fr
VEN366140M	14Fr, MDI Port
VEN368140	14Fr, Wedge + Adapter
VEN368140M	14Fr, MDI Port, Wedge + Adapter
VEN366160	16Fr
VEN366160M	16Fr, MDI Port
VEN368160	16Fr, Wedge + Adapter
VEN368160M	16Fr, MDI Port, Wedge + Adapter

Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port



Swivel Straight Connector

Item#
VEN005-20

Disposable
22mm ID, 15mm OD
20/pack



PeDIA - Pediatric Device for Induction of Anesthesia

Item#
ANEP021Y

With Mask Adapter Elbow

12/pack



Turn panic into play

For children, undergoing surgery is already unsettling, and the process of induction is often the most traumatizing part. Often kids will kick and scream, believing they're being trapped by the gas or suffocated by the mask. The experiences create negative emotional scarring that can last weeks, months or even years.

Specially designed for children undergoing surgery, the Pediatric Device for Induction of Anesthesia device, or "PeDIA," for short, includes a colorful whistling balloon, providing kids with something fun to focus on and turning a stressful moment into a lighter one.

Whether they're squeezing the balloon and focusing on its whistle or directly putting themselves to sleep by breathing the gas through balloon delivery, the PeDIA balloon captivates children's attention, allowing their active participation in the induction process.


It's a simple way to add a brighter touch to the medical experience, ease anxiety and bring smiles to patients, parents and staff.

This unique approach—presenting an interactive, familiar toy-like device—not only diverts children's focus from the impending procedure but also significantly reduces the anxiety and resistance commonly associated with mask induction. The aim is to decrease the anxiety felt by over 50% of pediatrics as well as parents and staff.

The PeDIA can be used alone as a mask-free device or, for younger children, with an anesthesia mask.

Use with
mask

Or use
without
mask



SCAN FOR VIDEOS

