



Motif Roam™

Breast Pump With Breast Milk Collection Cups

Item#
MATAAA0026



The Motif Roam is the durable and flexible breast pump designed to keep up with busy moms. It supports routine with ease and freedom, thanks to mom-favorite Breast Milk Collection Cups that eliminate the need for bulky bottles. Designed for freedom and flexibility, the Roam is ready to support moms' routines wherever life goes.

The Motif Roam provides true value with additional flange sizes and a price you won't be able to find anywhere else—especially not with a two-year warranty. The Roam also stands apart from competition with hospital-strength suction.



Features & Benefits

- **Hospital-Strength Suction:** Up to 280 mmHg suction strength, giving moms the efficiency they need for every pumping session
- **Portable:** Equipped with Breast Milk Collection Cups to give moms the freedom to pump hands-free
- **Discreet:** Collection Cups fit in nursing bra and hold up to 5 oz. of milk
- **More Options:** Resupply Kits available in both flange and collection cup styles
- **Variety of Flange Sizes:** Silicone inserts included in sizes 16 mm, 19 mm, 21 mm and 24mm so mom can find the perfect fit
- **Cluster Feeding and General Massage Mode:** 2 NEW Custom Modes designed by Lactation Professionals
- **Suction levels:** Works to express milk from the breast as efficiently as possible with 9 suction levels
- **Closed System:** Backflow prevention keeps the motor clean, preventing contamination

What's Included

- | | |
|-------------------------------|---------------------|
| 1 - Motif Roam Pump Motor | 2 - Diaphragm |
| 2 - Collection Cup | 2 - 16/19/21/24 mm |
| 2 - 28 mm Breast Shield | Flange Insert Pairs |
| 2 - Tubing w/ Roam Connectors | 1 - Lanyard |
| 2 - Tubing Connector | 1 - USB-C Cable |
| 4 - Duckbill Valve | 1 - Power Adapter |

Comparison Chart

				
BRAND	Motif	Zomee	Momcozy	Baby Buddha
MODEL	Roam	Z2	V2	2.0
WEIGHT	0.6 lbs	0.5 lbs	0.54 lbs	0.43 lbs
SUCTION STRENGTH	280 mmHg	260 mmHg	288 mmHg	320 mmHg
WARRANTY	2 years	2 year	1 year	1 year
NOISE LEVEL	<55 dB	N/A	<48 dB	<60 dB
MILK COLLECTION SIZE	5 oz.	8 oz.	7.7 oz.	N/A
INSURANCE UPGRADE NEEDED	✗	N/A	N/A	N/A
FLANGE SIZES	16 mm, 19 mm, 21 mm, 24 mm, 28 mm	24 mm and 28 mm	16 mm, 19 mm, 21 mm, 24 mm, 27 mm	24 mm
MODES	4 (Massage, Expression, General Pumping, Cluster Feeding)	3 (2 -Phase, Massage, Expression)	3 (Stimulation, Expression, Mixed)	3 (Soft Stimulation, Stimulation, Expression)
AUTO SHUT OFF	✓	✓	✓	✓
BPA FREE	✓	✓	✓	✓

*Designed to be fully covered by most insurance plans.

To order or request information

P 312.997.9980 E sales@sunsethcs.com W SUNSETHCS.COM

