

Ventilator Circuits Single Limb

Non-Heated

Active

Adult

VC1001

30/case

OEM#

1073223

1073232



Pediatric



VC1001P

30/case

OEM#

1073220

1073229



Passive

Adult


VC1004

30/case

OEM#

1069443

1069444



Pediatric


VC1004P

30/case

OEM#

1073224

1073233



VC1000

Each

OEM#

1073228

1073236



Swivel
exhalation port

VC1006

Each

OEM#

007936



Bacterial/Viral
Filter


VC1100

Each

OEM#

1069467

1069468



Fixed exhalation port

VC1007

Each

OEM#

007615



Bacterial/Viral
Filter

Fixed
exhalation port

VC1005

Adult

Each

OEM#

1147378

1132340

7615

VC1005P

Pediatric

Each

OEM#

1132339



Swivel exhalation port

Bacterial/Viral
Filter

Continued on next page

Heated

For use with F&P MR850

Passive

Adult

VC1114

10/case

22mm circuit
with optional connector

OEM# RT114



VC1003

10/case

Includes:



- 22mm circuit
- Swivel exhalation port
- Optional connector

OEM#

1076046

1073227

1073235




Pediatric

VC1114P

10/case

15mm circuit



VC1003P

10/case



Includes:

- 15mm circuit
- Swivel exhalation port
- 15mm/22mm connector
- 22mm/22mm connector

OEM#

1073225

1073234



Active

VC1002

30/case



Includes:

- 22mm circuit
- Active exhalation valve adapter
- Optional connector

OEM#

1073222

1073231



VC1002P

30/case


Includes:

- 15mm circuit
- Active exhalation valve adapter
- 15mm/22mm connector

OEM#

1073221

1073230



Adapters/Connectors

VEN001

Valved T-Adapter

Inline adapter for nebulizer treatment



VEN002

Swivel Exhalation Port

OEM# 332113



VEN003

Fixed Exhalation Port

OEM# 312149



VEN004

Active Exhalation

Valve Adapter

OEM# 1065659



VEN005

Straight Connector

22mm ID, 15mm OD



VEN006

Straight Connector

22mm ID, 22mm OD

