




Acute Care

Product Catalog



Contracted Supplier



vizient.
Awarded Supplier





Sales Representatives



Egle Sventnickaite
872-258-1654
esventnickaite@sunsethcs.com



Mason Moss
262-788-1862
mmoss@sunsethcs.com



Owen Ruth
312-533-2451
owenr@sunsethcs.com



Cindy Geighes
312-533-2478
cgeighes@sunsethcs.com



Hunter Baudino
262-788-1855
hbaudino@sunsethcs.com



Wyatt Kanemoto
262-788-8941
wkanemoto@sunsethcs.com



Director, Strategic Accounts
Lisa Iannazzone RN, BSN
636-222-3134
lisa@sunsethcs.com




Manager, Product & Business Development
Tom Munar
312-533-2476
tmunar@sunsethcs.com



Contract Manager, Strategic Accounts
Gina Rosenberg
914-837-5215
grosenberg@sunsethcs.com



RRT Clinical Specialist
Krystle Bodine BS, RRT
kbodine@sunsethcs.com



Corporate Office

Phone 872-275-3210

Fax 312-997-9985

141 W Jackson Blvd Ste 1950
Chicago IL 60604

Central Returns/Distribution Center

2725 Alft Ln
Elgin IL 60124

West Returns/Distribution Center

6750 W Van Buren St
Phoenix AZ 85043

Content

NIV Page

Adult NIV Masks	4-9
Pediatric NIV Masks	6-9
CPAP Tubing	10
Connectors and Ports	10
Headgear	11
Chinstraps	11-12

Microdial Flowmeter Page

Microflow	13
-----------	----

Respiratory Disposables Page

Vent Circuits	14-17
Lynx Venturi Oxygen Delivery	18
Oxygen Masks	19-20
Air Entrainment Masks	20
Oxygen Cannula and Accessories	21-23
Oxygen Supply and CPAP Tubing	24
Firesafe	25
Oxygen Connectors and Adapters	26
Oxygen Accessories	27
SpO2 Sensors	28
Bacteria Filters and Humidifiers	29-30
Oxygen Filters	30-31
Spirometers	32
Drainage Bags	32
Concentrator Filters	33-34
Berman Airways	35
Guedel and Nasopharyngeal Airways	36
Tracheostomy Kits and Accessories	36
HMEs and HME Accessories	37

Open Suction Page

Suction Catheter Kits	38
Suction Catheters	39
Trach Ties and Connectors	39
Suction Canisters and Liners	40
Suction Canister/Liner Kits	41
Suction Canister Brackets and Wall Mounts	42
Suction Tubing, Yankauers, and Bulk Roll	43
Solidifiers	44



Medication Delivery Page

Inhaler Chambers	45
Peak Flow Meters	45
Intake Filters and T-Adapters	45
Compressor Nebulizers	46
Handheld and Mini Nebulizers	46
Nebulizer Kits	47
Aerosol Masks	47

Closed Suction Catheters Page

Standard Switch	48-49
Deluxe Switch	50-51

Resuscitation Devices Page

Resuscitators	52
Manometers	53
PEEP Valves	53
Reservoir Bags	54
Pocket Resuscitator Set	54
Breathing Masks	54

Laryngeal Masks Page

Silicone Disposable	55
Silicone Reusable	55
Flexible Silicone Disposable	55
Flexible Silicone Reusable	55
PVC Disposable	56
Flexible PVC Disposable	56

Anesthesia Page

PeDIA	57
-------	----



CONTRACT

MS1087	Active Humidification + Respiratory Disposables
MS1097	Medication Delivery
MS1049	Non-Invasive Ventilation
MS7709	PeDIA Anesthesia Induction
MS1486	Open Suction
MS1488	Closed Suction Catheters



Contracted Supplier

CONTRACT

PP-NS-1741	Respiratory Disposables
PP-NS-1753	Medication Delivery
PP-NS-1762	Resuscitation Devices
PP-OR-2217	Laryngeal Masks
PP-FA-2138	Medical Gas Pipeline Equipment
BT-OR-019	Disposable Anesthesia

Full Face Mask

With Standard Headgear

Standard Elbow, Non-Vented

Size	Item#
Small	CM041S
Medium	CM041M
Large	CM041L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM041SEE
Medium	CM041MEE
Large	CM041LEE

Single Limb Elbow, Vented, AAV

Size	Item#
Small	CM041SEEV
Medium	CM041MEEV
Large	CM041LEEV

novaplus⁺

Standard Elbow, Non-Vented

Size	Item#
Small	VCM041S
Medium	VCM041M
Large	VCM041L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	VCM041SEE
Medium	VCM041MEE
Large	VCM041LEE

Single Limb Elbow, Vented, AAV

Size	Item#
Small	VCM041SEEV
Medium	VCM041MEEV
Large	VCM041LEEV

Replacement Cushions

Size	Item#
Small	CUCM005S-10
Medium	CUCM005M-10
Large	CUCM005L-10



With Cap Style Headgear

Standard Elbow, Non-Vented

Size	Item#
Small	CM041SC
Medium	CM041MC
Large	CM041LC

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM041SEEC
Medium	CM041MEEC
Large	CM041LEEC

Single Limb Elbow, Vented, AAV

Size	Item#
Small	CM041SEEV
Medium	CM041MEEV
Large	CM041LEEV

- Unobtrusive design offers a clear line of sight
- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone cushion increases comfort while maintaining optimal seal

Sunset Deluxe Nasal CPAP Masks

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM006S
Medium	CM006M
Large	CM006L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

CMSIZE01
 Deluxe Sizing Guide

HCPCS#
 MASK - A7034
 HEADGEAR - A7035
 CUSHION - A7032



Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.

Mask without Headgear

Size	Item#
Small	CM006SNH
Medium	CM006MNH
Large	CM006LNH

Replacement Cushions

Size	Item#
Small	CUCM006S
Medium	CUCM006M
Large	CUCM006L

Nasal Mask V

With Standard Headgear

Size	Item#
Small	CM008S
Medium	CM008M
Large	CM008L

novaplus+

Size	Item#
Small	VCM008S
Medium	VCM008M
Large	VCM008L



Standard Elbow, Non-Vented

Size	Item#
Small	CM008SVNV-10
Medium	CM008MVNV-10
Large	CM008LVNV-10

- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

Total Face Mask V

Standard Elbow

Size	Item#
Small	CM040S
Medium	CM040M
Large	CM040L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM040SEE
Medium	CM040MEE
Large	CM040LEE



Sunset Deluxe Nasal Pillow CPAP Mask V

Mask with Headgear and Small, Medium, and Large Pillows

Item# CM036

Headgear Only HG036

HCPCS#
MASK - A7034
HEADGEAR - A7035
PILLOWS - A7033

Replacement Pillows

Size	Item#
Small	CM036SP
Medium	CM036MP
Large	CM036LP



- Easy to fit, assemble, and maintain
- Three different pillow sizes to fit more people more easily
- 360° rotational elbow with flexible silicone tube
- Quiet vent design

Mask without Headgear CM036NH

Sunset Clearsight Deluxe Nasal Mask V

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM110S
Medium	CM110M
Large	CM110L

Replacement Cushions

Size	Item#
Small	CUCM110S
Medium	CUCM110M
Large	CUCM110L

HCPCS#
MASK - A7034
HEADGEAR - A7035



Clearsight design provides unobstructed view

Sol Full Face CPAP Mask

Mask with Headgear and Replaceable Cushion


10/case

Size	Item#
Small	CM205S
Medium	CM205M
Large	CM205L

Replacement Cushions

Size	Item#
Small	CUCM205S-10
Medium	CUCM205M-10
Large	CUCM205L-10

Easy Pop Cushions



A premium mask without the premium price

Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

Quiet vent

Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

Quick release connector

Rapidly and easily separates or connects mask and the tube

Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment

Ray Pediatric Masks Comparable to Nonny masks

Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	CM061S
Medium	CM061M
Large	CM061L
Extra Large	CM061XL

Replacement Cushions

Size	Item#
Small	CUCM061S
Medium	CUCM061M
Large	CUCM061L
Extra Large	CUCM061XL

Headgear Only

Size	Item#
Small	HG061S
Large	HG061L

SCAN CODE to see Ray Full Face Mask on YouTube



Nasal Mask with Headgear, Tube, and Replaceable Cushions


Comes with	Item#
S&M Cushions	CM060-KITS
L&XL Cushions	CM060-KITL

Replacement Cushions

Size	Item#
Small	CUCM060S
Medium	CUCM060M
Large	CUCM060L
Extra Large	CUCM060XL

Headgear Only

Size	Item#
Small	HG060S
Large	HG060L



Available Size Guide



SCAN CODE to see Ray Nasal Mask on YouTube



Sleepnet Innova V

Nasal Mask with Headgear

Vented

Size	Item#
Small/Medium	CM50180
Medium/Large	CM50181

novaplus+

Size	Item#
Small/Medium	VCM034SM
Medium/Large	VCM034ML

Non-Vented

Size	Item#
Small/Medium	CM50594
Medium/Large	CM50595

novaplus+

Size	Item#
Small/Medium	VCM034NVSM
Medium/Large	VCM034NVML

Innova® CPAP masks feature patented AIR®gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA



Sleepnet Ascend Mask Kit

Full Face System

Item# **CM037F**

Each

Vented

Mask Only

Size	Item#
Small	CM037S V
Medium	CM037M V
Large	CM037L V

Replacement Cushions

Size	Item#
Small	CM037CS V
Medium	CM037CM V
Large	CM037CL V



Nasal System

Item# **CM037N**

Each

Vented

Mask Only

Size	Item#
Small	CM037NS V
Medium	CM037NM V
Large	CM037NL V

Replacement Cushions

Size	Item#
Small	CM037NCS
Medium	CM037NCM
Large	CM037NCL



Sleepnet Phantom®2 and iQ®2 V

Phantom®2

Item# **CM050**

Replacement Cushion
CUCM050

Sleepnet's patented AIR®gel® technology provides maximum comfort while minimizing skin irritation and pressure points on the nasal bridge. Unlike many other gel masks, the unique qualities of AIR®gel® allow the mask to seal with minimal tension. By design AIR®gel® is twice as soft as competing gel products.



Replacement Frame/Tubing Assembly
CM051

iQ®2

Item# **CM049**

Replacement Cushion
CUCM049

- Sleepnet has the only mask that you can shape to each individual face. Flex, Fit, Forget!
- Cushions may be removed and replaced without purchasing a new mask.



Ask us about our 30-day mask guarantee

- Magnetic Clips — Easy on and off!
- 360° Swivel — Less tube drag for active sleepers.
- One size cushion — Larger fitting range.

Veraseal 2 Disposable Gel Mask with Headgear V



**Standard Elbow,
Non-Vented**



Size	Item#
Small	CM50933
Medium	CM50934
Large	CM50935
Extra Large	CM50936



Size	Item#
Small	VCM053S
Medium	VCM053M
Large	VCM053L
Extra Large	VCM053XL



**Single Limb Elbow,
Non-Vented, AAV**



Size	Item#
Small	CM50937
Medium	CM50938
Large	CM50939
Extra Large	CM50940



Size	Item#
Small	VCM054S
Medium	VCM054M
Large	VCM054L
Extra Large	VCM054XL



**Single Limb Elbow,
Vented, AAV**



Size	Item#
Small	CM50827
Medium	CM50828
Large	CM50829
Extra Large	CM50830



Size	Item#
Small	VCM052S
Medium	VCM052M
Large	VCM052L
Extra Large	VCM052XL

- The Veraseal® 2 disposable full face mask brings AIR°gel® Technology to the acute care environment for a giant step forward in patient comfort. Three color coded elbow designs make choosing the right mask for CPAP, single, or dual limb circuit NIV systems simple. The Veraseal® 2 fits quickly and easily while enhancing patient comfort. Both vented (with leak ports) and non-vented (without leak ports) versions are available for CPAP, Bilevel, and ventilator therapy.
- The AIR°gel® cushion provides the utmost in gentle cushioning to protect from pressure points and irritation. AIR°gel® promotes advanced skin protection, superior performance, and unbeatable comfort, all in a lightweight mask.
- The active headgear connector moves as the patient moves for a reliable seal.
- Breathable headgear lets patients enjoy soft, lightweight stability they can count on for a comfortable fit.

MiniMe 2 Pediatric Nasal Mask V



Vented

Size	Item#
Extra Small	CM60216
Small	CM60217
Medium	CM60218
Large	CM60221

Non-Vented

Size	Item#
Extra Small	CM60256
Small	CM60257
Medium	CM60258
Large	CM60261



Vented

Size	Item#
Extra Small	VCM055XS
Small	VCM055S
Medium	VCM055M
Large	VCM055L

Non-Vented

Size	Item#
Extra Small	VCM056XS
Small	VCM056S
Medium	VCM056M
Large	VCM056L



- Maximum comfort, minimum weight. Patented AIR°gel makes MiniMe masks some of the most gentle fitting masks on the market.
- The shapeable shell mask can be easily molded to fit individual needs.
- 100% designed and manufactured in the USA.

Pediatric SleepWeaver® CPAP Mask **V**

Item#
CM101325



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.


- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

MiniMe Pediatric Nasal Mask **V**



Headgear Size	Item#	OEM#
Small	CM004S	60213
Medium	CM004M	60214 P
Large	CM004L	60215

HCPCS# A7034



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

SleepWeaver® Advance Nasal Mask **V**

Item#
CM017B
Mask with Adjustable Headgear

- Easy to store and pack for travel

OEM# 100274
HCPCS# A7034



NIV Tubing

19mm with 22mm cuffs
10/case*

Item#

TUB0018 18in

TUB02 2ft

TUB004 4ft

TUB06 6ft *Individually Packed **V**

SUH0610 6ft **V**

TUB08 8ft **V** **P**

TUB010 10ft **V** **P**



- Sunset's most popular item
- ISO standard 22mm cuffs
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

Slim CPAP Tube **V** **P**

Item#

TUB006SS

6ft x 15mm with
22mm cuffs

HCPCS# A7037



Highly durable

Strain relief cuff decreases separation and tearing

Virtually no memory

Designed to be the most flexible tube available

Minimal weight

One of the lightest tubes on the market at 2.6 ounces

Great value

Affordable and sustainable

Heated Tube for Respironics **P** **V** Vizient Contract MS1087


Item#

TUB06HPR

Heated Tube for Philips
Respironics Dreamstation
Series and System One

OEM# HT15

HCPCS# A4604



- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off

Standard Elbow, **V** Non-Vented

Item# **CAP009**
10/case



Single Limb Elbow, **V** Non-Vented, AAV

Item# **CAP008**
10/case



Swivel Connector Silencer

Item#

VEN002 Each **P**

VEN002-10 10/pack **V** Contract MS1087



Single Use Exhalation Port **V**

Item# **VEN003**
10/pack



Standard 4-Point Headgear **V**

Item#
HG003S Small
HG003M Medium
HG003L Large



Cap Style Headgear **V**

Item# **HG002**



Universal Neoprene Chinstrap **V**

Item#
CS001
 Universal Chinstrap
 HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows



Deluxe Chinstrap **V**

Size	Item#
26in	CS004S
28in	CS004L

Compare to
 Respironics Deluxe
 OEM# 302425
 HCPCS# A7036

- Latex-free
- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort



Comfort Chinstrap **V**

Item# **CS005**
 Compare to
 ResMed Sullivan
 OEM# 16015
 HCPCS# A7036

- Latex-free
- One size fits most
- Easy to use—attaches to itself at the top of the head
- Soft cloth chin cup adds comfort



Premium Chinstrap **V**

Item# **CS006**
 Compare to
 Respironics Premium
 OEM# 1012911
 HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.



Ruby Style Adjustable Chinstrap




Size	Item#
Regular	CS007-ADJ
Extra Large	CS007-ADJXL

HCPCS# A7036

- Latex-free
- Fully adjustable
- Super soft materials for extra comfort



Ruby Style Chinstrap



For fitting, measure completely around the head, as indicated by the dashed line above.

Size	Item#
21.5-22.5in	CS007S
22-23in	CS007M
23-24in	CS007L

Compare to Puresom Ruby

OEM# TMS-09 • HCPCS# A7036

- Latex-free
- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

Halo Style Chinstrap V

Item# **CS025**

Compare to BreatheWear Halo Chinstrap

OEM# H3407

HCPCS# A7036



Microflow Microdial Flowmeter P

2024 Innovation Celebration Winner



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure Range	Flow Range*
RES8160048	345 – 500 kPa	A
RES8160050	350 – 550 kPa	E



DISS Inlet • DISS Outlet

Item#	Inlet Pressure Range	Flow Range*
RES8160049	345 – 500 kPa	A
RES8160051	250 – 450 kPa	E

*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
Range A	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
Range E	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

Assured, repeatable accuracy to enable safe, effective weaning.

What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.

Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned



Non-Heated Breathing Circuits

For use with Respiration Trilogy
22mm cuffs


30/case

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
VC1001	1073223 & 1073232	Active	No	Adult	22mm
VC1001P	1073220 & 1073229	Active	Yes	Pediatric	15mm
VC1004	1069443 & 1069444	Passive	Yes	Adult	22mm
VC1004P	1073224 & 1073233	Passive	Yes	Pediatric	15mm

Item#
VC1001 V P

Active Adult

OEM#
1073223
1073232




Item#


VC1001P

Active Pediatric

OEM#
1073220
1073229




To use active circuits with
Respiration Trilogy, simply
cut ends off pressure and
exhalation lines.



Item#
VC1004 V P

Passive Adult

OEM#
1069443
1069444




Item#

VC1004P V P

Passive Pediatric

OEM#
1073224
1073233



Adult Passive Heated Circuit with Exhalation Port **V** **P**

30/case

For use with Philips
 Respironics Trilogy

Item#
VC1003


VC1003P
 Pediatric Size

Comparable to OEM# 1076046,
 1073227, 1073235



Includes:

- 1 VC1114 Passive Heated Breathing Circuit
- 1 VEN002 Swivel Exhalation Port



Passive Circuit w/Exhalation Port and Bacterial/Viral Filter **V** **P**

With 22mm cuffs

Item#
VC1005
Each

OEM# 1147378, 1132340, 7615



Adult Passive Circuit with Bacteria Filter **V**

Comparable to Breas 007936

Item#
VC1006
Each



Adult Passive Circuit with Bacteria Filter and T-Adapter **V**

Comparable to Breas 007615


Item#
VC1007
Each



Passive Non-Heated Breathing Circuits **V** **P**

Exhalation Port

Item#	OEM#
Swivel	VC1000 1073228 & 1073236
Fixed	VC1100 1069467 & 1069468



Passive Heated Breathing Circuit **V** **P**


Item#

VC1114

For use with
MR850 Humidifier
10/case

VC1114P

Pediatric



Non-Heated Dual Limb Breathing Circuits **V** **P**

30/case

Item#	Size	Humidifier Limb	Water Trap
VC2001P	Pediatric	No	No
VC2001	Adult	No	No
VC2002P	Pediatric	Yes	No
VC2002	Adult	Yes	No
VC2003	Adult	Yes	Yes



Adult NIV Patient Circuit **V** **P**


Item# **VC1160**

Compatible with V60, No Filters

10/case

Comparable to

OEM# 1069210



Adult NIV Single Limb Patient Circuit with Main Flow Filter **V** **P**

Item# **VC1161**

Includes 1 BF2200 22mm Bacteria Disk Filter

10/case




Adult NIV Single Limb Patient Circuit with Main Flow Filter and Exhalation Port Filter **V** **P**

Item# **VC1162**

Includes 2 BF2200 22mm Bacteria Disk Filters

10/case



Lynx Venturi Oxygen Delivery V P


Low Concentration

Item# **RES2400**

Each

Adapter Only **RES240LC**

- Design replaces multiple venturi adapters
- Connects to oxygen tubing and regulates flow
- Provides precision FiO₂ delivery from 24-60% with less than 2% variance during mild, moderate, and severe respiratory distress
- Excellent CO₂ washout



Unprecedented simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

Do no harm

Quicker de-escalation of O₂ levels and prevention of excessive CO₂ in the patient.

Reduction in waste

Remains with patient throughout EMS and hospital journey, replacing multiple adapters which may get lost or contaminated during transition.

Unmatched precision

Provides specific levels of oxygen with less than 2% variance. Comparable products vary by as much as 30%.

High Concentration

Item# **RES2401**


Each

Adapter Only **RES241HC**

- Allows rapid de-escalation of FiO₂ in the EMS ambulance or ED to prevent CO₂ buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia (PaO₂/FiO₂) allows early intervention and better outcomes
- Excellent CO₂ washout

Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Capabilities could include replacing non-rebreathers and/or partial rebreather products



Simple O₂ Mask V P

With 7ft Safety Tubing
Elongated

Item# **RES2100**
Adult Size
50/case




- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Medium Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	V	P
RES2100NT	Adult	None	Standard	V	P
RES2100P	Pediatric	7 ft	Standard	V	P
RES2100U	Adult	7 ft	Universal	V	P
RES2203I	Infant	None	Standard	V	
RES2203P	Pediatric	None	Standard	V	P
RES2203A	Adult	None	Standard	V	P
RES2203A7	Adult	7 ft	Standard	V	P
RES2203I7	Infant	7 ft	Standard	V	P
RES2203P7	Pediatric	7 ft	Standard	V	P

- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



High Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	V	P
RES2202A	Adult	None	V	P
RES2202A7	Adult	7 ft	V	
RES2202I	Infant	None	V	P
RES2202I7	Infant	7 ft	V	P
RES2202P	Pediatric	None	V	
RES2202P7	Pediatric	7 ft	V	



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Non-Rebreather Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	Safety Vent	
RES2102	Adult	7 ft	Standard	No	V P
RES2102I	Infant	7 ft	Standard	No	V P
RES2102IM	Infant	7 ft	Standard	No	V
RES2102P	Pediatric	7 ft	Standard	No	V P
RES2102PS	Pediatric	7 ft	Standard	Yes	V P
RES2102PSM	Pediatric	None	Standard	Yes	V P
RES2102PSU	Pediatric	7ft	Universal	Yes	V P
RES2102PU	Pediatric	7 ft	Universal	No	V P
RES2102S	Adult	7 ft	Standard	Yes	V P
RES2102SM	Adult	None	Standard	Yes	V P
RES2102U	Adult	7 ft	Universal	No	V P

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



Adult Air Entrainment Mask

With Flex Tube
50/case



Item#	Size	Kit	Tubing Length	Connector	
RES2400A	Adult	Mask Only	7 ft	Standard	V P
RES2400AK	Adult	Yes	7 ft	Standard	V P
RES2400AKU	Adult	Yes	7 ft	Universal	V P
RES2400AU	Adult	Mask Only	7 ft	Universal	V P
RES2400P	Pediatric	Mask Only	7 ft	Standard	V
RES2400PK	Pediatric	Yes	7 ft	Standard	V P
RES2400PKU	Pediatric	Yes	7 ft	Universal	V
RES2400PU	Pediatric	Mask Only	7 ft	Universal	V P



Humidoflo High Flow Nasal Cannula **V** **P**

10/case

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respironics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)



Size	Item#
Small	HFN001AS
Medium	HFN001AM
Large	HFN001AL

With F&P Adapter

Size	Item#
Small	HFN001ASA
Medium	HFN001AMA
Large	HFN001ALA

Adult Oxygen Cannula

Item#	Tubing Length	Connector	Qty	
RES1100	None	Standard	50	V P
RES1100N	None	Standard	50	V
RES1104	4 ft	Standard	50	V P
RES1107	7 ft	Standard	50	V P
RES1107N	7 ft	Standard	50	V P
RES1107NU	7 ft	Universal	50	V
RES1107U	7 ft	Universal	50	V P
RES1110	10 ft	Standard	50	V P
RES1114N	14 ft	Standard	50	V P
RES1114U	14 ft	Universal	50	V P
RES1115	15 ft	Standard	50	V P
RES1115U	15 ft	Universal	50	V P
RES1125	25 ft	Standard	25	V P
RES1135	35 ft	Standard	20	V
RES1150Z	50 ft	Standard	20	V
RES1150U	50 ft	Universal	20	V P



- 1-6 LPM
- 6-channel supply tubing for safety

Soft Cannula

Item#	Tubing Length	Connector	Qty	
RES1100S	None	Standard	50	V P
RES1107S	7 ft	Standard	50	V P
RES1107SU	7 ft	Universal	50	V P
RES1114S	14 ft	Standard	50	V P
RES1114SU	14 ft	Universal	50	V P
RES1115S	15 ft	Standard	50	V P
RES1125S	25 ft	Standard	25	V P
RES1125SU	25 ft	Universal	25	V



- 1-6 LPM
- 6-channel supply tubing for safety

Pediatric, Infant, and Neonate Cannula

50/case

Item#	Size	Tubing Length	Connector	Touch	Ear Cushions
RES1100P	Pediatric	None	Standard	Standard	No
RES1107I	Infant	7 ft	Standard	Standard	No
RES1107NE	Neonate	7 ft	Standard	Standard	No
RES1107P	Pediatric	7 ft	Standard	Standard	No
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes
RES1107SI	Infant	7 ft	Standard	Soft	No
RES1107SIU	Infant	7 ft	Universal	Soft	No
RES1107SNE	Neonate	7 ft	Standard	Soft	No
RES1107SP	Pediatric	7 ft	Standard	Soft	No
RES1107SPU	Pediatric	7 ft	Universal	Soft	No
RES1114P	Pediatric	14 ft	Standard	Standard	No
RES1114PU	Pediatric	14 ft	Universal	Standard	No
RES1114SP	Pediatric	14 ft	Standard	Soft	No
RES1114SPU	Pediatric	14 ft	Universal	Soft	No

- 1-6 LPM
- 6-channel supply tubing for safety

Adult Cannula with Ear Cushions

Item#	Tubing Length	Connector	Qty
RES1100EC	None	Standard	50
RES1104EC	4 ft	Standard	50
RES1107EC	7 ft	Standard	50
RES1107ECU	7 ft	Universal	50
RES1114EC	14 ft	Standard	50
RES1115EC	15 ft	Standard	50
RES1125EC	25 ft	Standard	25

Flared Cannula V

Item#	Tubing Length	Connector	Qty
RES1100F	None	Standard	50
RES1107F	7 ft	Standard	50
RES1107FU	7 ft	Universal	50
RES1114F	14 ft	Standard	50
RES1114FU	14 ft	Universal	50
RES1125F	25 ft	Standard	25
RES1150F	50 ft	Standard	20

- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

High Flow Soft Cannula V

25/case


Item#	Tubing Length
RES1104SHF	4ft
RES1107SHF	7ft
RES1115SHF	15ft
RES1125SHF	25ft



7 ft Cannula With CO₂ Sampling Line **V**

50/case

Item#	Size	Connector	Luer	Filter
RES8007	Adult	Standard	Male	No
RES8007F	Adult	Standard	Male	Yes
RES8007I	Infant	Standard	Male	No
RES8007IF	Infant	Standard	Male	Yes
RES8007IU	Infant	Universal	Male	No
RES8007P	Pediatric	Standard	Male	No
RES8007PF	Pediatric	Standard	Male	Yes
RES8007PU	Pediatric	Universal	Male	No
RES8007U	Adult	Universal	Male	No
RES8107	Adult	Standard	Female	No
RES8107F	Adult	Standard	Female	Yes
RES8107I	Infant	Standard	Female	No
RES8107IF	Infant	Standard	Female	Yes
RES8107P	Pediatric	Standard	Female	No
RES8107PF	Pediatric	Standard	Female	Yes




Secure Cannula Stickers **V** **P**

Item# **RES080**

25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs




Ear Cushions for Cannula **V** **P**

Item# **RES020**

Cushions Only

50 Pairs/case



- Add cushions to existing cannula

Oxygen Supply Tubing

Item#	Length	Connector	Bore	Color	Qty	
RES3007	7 ft	Standard	Standard	None	50	V P
RES3007G	7 ft	Standard	Standard	Green	50	V P
RES3007SB	7 ft	Standard	Smooth	Green	50	V
RES3007U	7 ft	Universal	Standard	None	50	V P
RES3015	15 ft	Standard	Standard	None	50	V P
RES3015G	15 ft	Standard	Standard	Green	50	V
RES3025	25 ft	Standard	Standard	None	25	V P
RES3025G	25 ft	Standard	Standard	Green	25	V P
RES3025U	25 ft	Universal	Standard	None	25	V P
RES3030	30 ft	Standard	Standard	None	20	V
RES3030G	30 ft	Standard	Standard	Green	20	V
RES3040G	40 ft	Standard	Standard	Green	20	V P
RES3050	50 ft	Standard	Standard	None	20	V
RES3050G	50 ft	Standard	Standard	Green	20	V P

- 6-channel
- PHT/DEHP-free



Universal Connector

CPAP Tubing V P

Item#
TUB008 8ft
TUB010 10ft
19mm with 22mm cuffs
HPCPS# A7037





- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

Firesafe Thermal Shut-Off Valve V

Item#
RES010-10
Bidirectional
10/case




- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire
- Smaller
- Laser etched
- Expiration date
- Human readable UDI
- Designed to last 5 years



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

Fatalities per 100,000 patients



FireSafe Nozzle V P

Item#
RES010FND
10/pack



FireSafe Tubing Connector

Item#
RES010TC 100/case V
RES010TC-10 10/case



Christmas Tree Connector

Item# **RES001** V P
Green without Swivel
10/pack



RES001-30 P
30/pack

Item# **RES002** V P
Green with Swivel
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Oxygen Tubing Connector — Non-Swivel V P

Item# **RES012**
Fits 3/16in and 1/4in
ID Tubing
10/pack



- Barbed ends ensure secure, leak-free connections

Tubing Reducer Connector V P

Item# **RES011**
Fits 3/16in, 1/4in,
5/16in, 3/8in, and 7/16in
ID Tubing
10/pack



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

Humidifier Tubing Adapter V P

Item# **RES016**
22mm ID x
22mm OD
with Lip
10/pack



Circuit Connectors V P

Item# **RES013**
22mm OD x
15mm ID/
22mm OD
10/pack



Item# **RES014**
22mm OD x
22mm OD
With Lip
10/pack



Oxygen Tubing Connector — Swivel

Item# **RES018** V P
Male/Male
10/pack
RES018-1 1/pack



RES018-25 P
25/pack

Oxygen Tubing Y-Connector V P

Item# **RES017**
Fits 3/16in and 1/4in ID Tubing
(All 3 sides are the same)
10/pack



- Connects two oxygen sources for higher flows

Connector with 6mm Oxygen Stem V P

Item# **RES058**
22mm with 6mm Stem
10/pack



Oxygen Bleed-In Adapter V P

Item# **RES019**
22mm with 6mm
Rotating Stem
10/pack
RES019-1 1/pack



Finger Pulse Oximeter **P**

Item# **RES5102**
 Each



- Oversized 1.3" OLED display face
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- Includes lanyard and manual

Wrenches **V** **P**

Item# **RES031**
 With Chain
10/pack
RES031-1 1/pack




Humidifier Connectors **V** **P**

Item#
HUM008
 L-Connector 12in
10/pk
HUM008-1 1/pack

HUM008L
 L-Connector 21in
10/pk

HUM008S
 Straight Connector 9in
10/pk



Water Traps **V** **P**

Item# **RES037**
 Push On End Caps
 Each



- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

Oxygen Liter Meter **V** **P**

Item# **RES034**
 0-10 Liter/Minute
 Each

- Oxygen pen
- Works with tanks and concentrators



Covers for Concentrators **P**



Item# **EC251530C**
 25x15x30"
 1.5mil • Clear
250 bags/roll


SpO2 Sensors

NEW

Disposable SpO2 Sensor Neonatal/Adult

Item# **RES7100**


Similar to Nellcor MAX-N
20/box



Neo-M-1.5 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7200**

1.5ft cord • Non-replaceable
20/box



Neo-M-3 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7201**

3ft cord • Non-replaceable
20/box



Coban Tape Roll

Item# **RES7202CT**

Use for sensor life
extension or sensor
cord stability


Each roll 1" x 177"



Infant-M-1.5 Disposable SpO2 Sensor Infant (3-20 kg)

Item# **RES7203**

1.5ft cord
20/box




Pediatric-M-1.5 Disposable SpO2 Sensor Pediatric (10-50 kg)

Item# **RES7204**

1.5ft cord
20/box



Suction Machine Inline Filters V P



Item# **BF1412**
Final Bacteria Filter •
1/4in to 1/2in Stepped
Barb on Both Sides
10/pack



Item# **BF1312**
1/3in to 1/2in
1/pack



Item# **BF3812**
3/8in to 1/2in
10/pack



Item# **BF516TS**
5/16in Barbed with
Threaded NPT
10/pack

Oxygen Filters



Item# **BF316** V P
Final Bacteria Filter •
3/16 Stepped Barb and
1/4in Straight Barb on
Both Sides
10/pack
BF316-1 1/pack



Item# **BF516** V P
Final Bacteria Filter •
1/4in Barbed
10/pack
OEM# P0750X030



Item# **BF1438** P
1/4in to 3/8in
10/pack

HEPA Filter for 3B Stratus V P

Item# **BF6020**
HEPA Filter
Each
OEM# STR1020



6in HEPA Capsule with Removable Sticker V P

Item# **BFH6S**
Threaded
Each



Disposable Bubble Humidifier 500 ml P

Item#
RES050
Up to 15 LPM • 6 psi • 500 ml
50/case
HCPCS# E0560
OEM# HUM001

- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents

- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient



Extended Life Intake Bacteria Filter V P

Item# **OF1025D**
Each

**Foam Cabinet Filter** V P

Item# **OF9001**
6 3/4" x 33/8" x 1/2"
10/pack
OEM# 2000489



For Invacare
Mobilaire and
HomeFill oxygen
concentrators

Mobilaire 3 and 5 Kit P

Item# **OF9001SP**
2 Foam Cabinet Filters
1 HEPA Filter 6in Capsule
1 Final Bacteria Filter
Each

**Invacare Platinum 5 and 10 Style Filter** V P

Item# **OF9002**
Foam Cabinet Filter
9" x 27/8" x 1/2"
10/pack
OEM#1107412

**3B Stratus 5 Style Filter** P

Item# **OF6030**
Foam Cabinet Filter
10/pack
OEM# STR1030

**Invacare Perfecto2 Style Filter** V P

Item# **OF9004**
Foam Cabinet Filter
6 7/8" x 33/4" x 1/2"
10/pack
OEM# 1143492

**Perfecto2 Kit** P

Item# **OF9004SP**
2 Foam Cabinet Filters
1 HEPA Filter 6in Capsule
1 Final Bacteria Filter
Each

**Extended Life Bacteria Filter with Muffler** V P

Item# **OF6011**
22mm fittings
Each
Use if SolAiris was
manufactured
before Nov '03



For Drive/DeVilbiss
SolAiris and
525 Series oxygen
concentrators

EverFlo Kit P

Item# **OF1003SP**
1 Intake Filter
1 Final Bacteria Filter
Each

**DeVilbiss/Drive SolAiris III, 303 Series Style Filter** V P

Item# **OF5004**
Foam Cabinet Filter
33/4" x 33/4" x 1/2"
10/pack
OEM# 30302605

**SolAiris V Kit** P

Item# **OF5005SP**
2 Foam Cabinet Filters
1 Extended Life Bacteria Filter
w/Muffler
1 Final Bacteria Filter
Each



Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units **V** **P**

Item# **BF2200**
 22mm Bacteria
 Disk Filter
5/pack



Attach between the machine and tube for additional protection against bacteria and allergens

Filter Specifications

Bacterial Filtration Efficiency: (BFE) 99.99+%
 Viral Filtration Efficiency: (VFE) 99.9+%
 NaCl Particulate Filtration Efficiency: 96+%
 Flow Resistance: Approximately 8.8 mmH2O at 30 LPM
 Internal Volume: 66.3 mL without fittings / 78.3 mL with fittings

Connections:

- Patient End: 22 mm OD / 15 mm ID
- Machine End: 22 mm ID

Plastic Intake Filter with Head **V** **P**

Item# **BF017**
 For DeVilbiss Nebulizers
5/pack
 HCPCS# A7013



Foam Nebulizer Filter **V** **P**

Item# **BF018**
10/pack
 HCPCS# A7013



Plastic Intake Filter for Respireonics InnoSpire **V** **P**

Item# **BF0265**
10/pack
 HCPCS# A7013



- Mini
- MicroElite
- MiniElite

EverFlo Intake Filter **V** **P**

Item# **BF1001**
Each
 Respireonics OEM# 1038831



Bacterial/Viral HME Filter – Adult **V** **P**

Item# **BF104**
Each



- Bacterial efficiency 99.9999%
- Viral efficiency 99.999%
- Humidification efficiency 37.5 mg/l H2O (500 ml tidal volume)
- Coupled with a hydrophobic bacterial viral filter membrane

Filter for ResMed Astral Ventilator **V** **P**

Item# **VF003-4**
4/pack



Item# **VF003-H**
 Filter Housing
Each



Item# **VF003-KIT**
 4-Filter Kit
 w/Filter Housing
 OEM# 27939



Volumetric Incentive Spirometer V P

12/case

Item# **RES055**
2500ml

RES056
5000ml

RES055P
Pediatric



SCAN FOR VIDEO



- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise
- Dual-sided, clear and legible identification of volumes
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration
- Compact and ergonomic design

Three Ball Spirometer

Item# **RES057-10** P
10/pack

RES057 V P
12/case

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient



- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

Drainage Bag with Y-Adapter V P

Item# **RES071**
750cc
50/case



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

Drainage Bag with T-Adapter V P

Item# **RES070**
750cc
50/case



ResMed S9 and AirSense 10 Series Filters **V** **P**

Item# CF2107

Disposable Filter

3/pack • 2 1/8" x 1 3/8"

HCPCS# A7038

CF2107-1 1/pack OEM# 36850

CF2107-2 2/pack OEM# 36851

CF2107-50 50/pack OEM# 36853

Item# CF2107HA

Hypoallergenic Filter

3/pack • 2 1/8" x 1 3/8"

HCPCS# A7038

CF2107HA-1 1/pack OEM# 36855

CF2107HA-2 2/pack OEM# 36856

CF2107HA-50 50/pack

OEM# 36858



ResMed AirSense 11 Series Filters **P**

Item# CF2111

Disposable Filter

3/pack

HCPCS# A7038

CF2111-1 1/pack OEM# 39300

CF2111-2 2/pack OEM# 39301

CF2111-50 50/pack OEM# 39303

Item# CF2111HA-3

Hypoallergenic Filter

3/pack

1 1/4" x 1 3/16"

HCPCS# A7038

CF2111HA-1 1/pack

CF2111HA-2 2/pack



3B Luna G3 Style Filters **V** **P**

Item# CF8006-6

Disposable Filter

6/pack

HCPCS# A7038

CF8006-1 1/pack

CF8006-2 2/pack

Item# CF8006F-1

Reusable Filter

1/pack

HCPCS# A7039



3B Luna and Luna II Style Filters **P**

Item# CF8005

Disposable Filter

12/pack • 2" x 1 3/16"

HCPCS# A7038

Item# CF8005F-1

Reusable Filter

1/pack • 2" x 1 3/16"

HCPCS# A7039

CF8005F-2 2/pack



3B RESmart Style Filter **P**

Item#

CF8004

Disposable Filter

3/pack • 3 5/16" x 1 5/8"

HCPCS# A7038



Respironics DreamStation 2 Series Filters **V** **P**



Item# **CF1012**

Disposable Filter
6/pack

OEM# 1142831

HCPCS# A7038

CF1012-1 1/pack OEM# 1142828

CF1012-2 2/pack OEM# 1142829

CF1012-3 3/pack

CF1012-50 50/pack



Item# **CF1012F**

Reusable Filter
2/pack

OEM# 1142687

HCPCS# A7039

CF1012F-1 1/pack



Respironics DreamStation Series Filters **V** **P**



Item# **CF1010**

Disposable Filter
6/pack

OEM# 1122519

HCPCS# A7038

CF1010-1 1/pack

CF1010-2 2/pack



Item# **CF1010F**

Reusable Filter
2/pack

OEM# 112446 (1-pack)

HCPCS# A7039

CF1010F-1 1/pack

CF1010F-6 6/pack



90-Degree CPAP Tubing Elbow **V** **P**

Item# **CAP1027**

Each

The Sunset 90-Degree CPAP Tubing Elbow Adapter connects standard 22mm or slim 15mm CPAP tubes to most CPAP machines at a 90-degree angle, relieving stress on the tubing connection and creating a more convenient pathway to the bed.

Comparable to OEM# 37394



Color Coded Dual Channel Berman Oral Airway V P

50/pack

Item#	Size	Color
AOA40BC-50	40mm	Pink
AOA50BC-50	50mm	Blue
AOA60BC-50	60mm	Black
AOA70BC-50	70mm	White
AOA80BC-50	80mm	Green
AOA90BC-50	90mm	Yellow
AOA100BC-50	100mm	Red
AOA110BC-50	110mm	Orange



- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

White Dual Channel Berman Oral Airway V P

10/pack

Item#	Size
AOA40BD	40mm
AOA60BD	60mm
AOA80BD	80mm
AOA90BD	90mm
AOA100BD	100mm



Oropharyngeal Guedel Airway V P

10/pack

Item#	Size	Color
AOP500GC	50mm	Blue
AOP601GC	60mm	Black
AOP702GC	70mm	White
AOP803GC	80mm	Green
AOP904GC	90mm	Yellow
AOP1005GC	100mm	Red



- Designed to overcome associated risks of loose or detached bite blocks
- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

Nasopharyngeal Airway **V**

10/pack

Item#	Size	Length	PVC
ANP1472PNE	14fr	72mm	
ANP1685PNE	16fr	85mm	
ANP18105PNE	18fr	105mm	
ANP20115PNE	20fr	115mm	
ANP22125PNE	22fr	125mm	
ANP24130PNE	24fr	130mm	
ANP26140PNE	26fr	140mm	
ANP28155PNE	28fr	155mm	
ANP30165PNE	30fr	165mm	
ANP32170PNE	32fr	170mm	
ANP34175PNE	34fr	175mm	



Tracheostomy Swivel Elbow Connectors **V P**

Item# **RES021**

With suction port
15mm ID x 15mm OD
10/pack



Tracheostomy Ties with Fastener Tabs **V P**

20/case


One-piece

• Easy on

Item#

RES029 Adult

RES029P Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton

Two-piece

• Easy off

Item#

RES030 Adult

RES030P Pediatric



Endotracheal Tube Holders

20/case

Item#

TRAA005

Tube Holder

TRAA005B

Tube Holder with
Bite Block

Replaceable,
adjustable
neck strap

Moisture-resistant
upper lip retainer

Tube
runner

Soft protector
pads

Optional built-in
bite block

Tube
fastener

Sunset Trach Care Kits **V** Contract MS1486

50/case

Sunset Trach Care Kit

Item# **TRATK003**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- Gauze
- 2 Cotton Tips
- Brush,
- Twill Tape
- Trach Dressing

Sunset Trach Care Kit
with Forceps

Item# **TRATK005**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- 3 Gauze
- 3 Cotton Tips
- Brush
- Trach Dressing

Sunset Mini Trach Kit

Item# **TRATK004**

Includes:

- 3 Pipe Cleaners
- 4 Gauze
- 2 Cotton Tips
- Brush
- Twill Tape
- Dressing

Flex HME **V** **P**

Item# **RES027FLEX**
 40/case



TV Range (mL)	300-1500
Moisture Output (mg H2O/L)	33 mg/L
Weight (g)	17
Resistance (cm H2O)	0.5 @ 30 LPM 1.34 @ 60 LPM
Dead Space (mL)	75
Connection (mm)	22M,15F

Heat Moisture Exchangers **V** **P**

30/case

Item#					
Item#	RES027 T-HME	RES027OP T-HME with Oxygen Port	RES027IL Inline HME — Large Volume	RES027NE Inline HME — Neonatal	RES027HV2 Foam HME — Adult/Child
TV Range (mL)	250-1000	150-1250	50-250	15-50	250-1000
Moisture Output (mg H2O/L)	30, Vt = 0.1L	36.3 mg/L	31, Vt = 0.1L	30	29
Weight (g)	4.5	7.4	10	4.5	21
Resistance (cm H2O)	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.2 @ 30 LPM	0.4 @ 10 LPM	0.5 @ 10 LPM	0.52 @ 30 LPM 1.32 @ 60 LPM 2.54 @ 90 LPM
Dead Space (mL)	12.5	7	13	2.5	29
Connection (mm)	15F	15F 5.5mm O2 port	15M,22M/15F	15M,15F	22F/15M, 22M/15F

Tracheostomy Mask with O₂ Connector **V** **P**

50/case

Item#
RES2130
 Adult Regular Size

RES2130L
 Adult Large Size

RES2230
 Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O₂ supply tubing

Adapter for T-HME

Item#

RES027A **V** **P**
 Each

RES027A-30 30/pack **P**



Flexible Connector Length 4.9-7.5cm

50/case

Item#
TRAA004 **V** **P**

Adult — 22/15mm x 15mm

TRAA004P

Pediatric — 15mm F x 15mm M



Suction Catheter Kits P

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	TRATK305
6fr	Y	DeLee	Pediatric	TRATK306
8fr	Y	DeLee	Pediatric	TRATK308
8fr	Y	Whistle	Adult	TRATK108
8fr	Chimney	Whistle	Adult	TRATK108CV
10fr	Y	Whistle	Adult	TRATK110
10fr	Chimney	Whistle	Adult	TRATK110CV
12fr	Y	Whistle	Adult	TRATK112
12fr	Chimney	Whistle	Adult	TRATK112CV
14fr	Y	Whistle	Adult	TRATK114
14fr	Chimney	Whistle	Adult	TRATK114CV
16fr	Y	Whistle	Adult	TRATK116
16fr	Chimney	Whistle	Adult	TRATK116CV

- Sterile
- Includes Catheter and Gloves

Mini Catheter Tray with Gloves V

50/case


Item#	Catheter Size
TRATK310	10fr
TRATK312	12fr
TRATK312VCVN	18fr, Vinyl Gloves
TRATK314	14fr
TRATK316	16fr
TRATK318	18fr

- Coiled Catheter
- Calibrated

Sleeved Suction Coiled Catheter Kits with Cup V

50/case

Item#	Catheter Size
TRATK410	10fr
TRATK412	12fr
TRATK414	14fr
TRATK418	18fr



Suction Catheter Kits with Cup

50/case

Size	Valve Type	Tip	Intended Use	Item#	Contract
5fr	Y	DeLee	Pediatric	TRATK405	P
6fr	Y	DeLee	Pediatric	TRATK406	P
6fr	Chimney	DeLee	Pediatric	TRATK206	V
8fr	Y	DeLee	Pediatric	TRATK408	P
8fr	Y	Whistle	Adult	TRATK208	P V
8fr	Chimney	Whistle	Adult	TRATK208CV	P
10fr	Y	Whistle	Adult	TRATK210	P V
10fr	Chimney	Whistle	Adult	TRATK210CV	P
12fr	Y	Whistle	Adult	TRATK212	P V
12fr	Chimney	Whistle	Adult	TRATK212CV	P
14fr	Y	Whistle	Adult	TRATK214	P V
14fr	Chimney	Whistle	Adult	TRATK214CV	P
16fr	Y	Whistle	Adult	TRATK216	P
16fr	Chimney	Whistle	Adult	TRATK216CV	P
18fr	Y	Whistle	Adult	TRATK218	V



- Sterile
- Includes Catheter and Gloves

Suction Catheters

Size	Valve	Tip	Intended Use	Packaged	Item#	Qty	Contract
5fr	Y	DeLee	Pediatric	Coiled	TRASC305C	50	P
5fr	Y	DeLee	Pediatric	Straight	TRASC305S	100	P
6fr	Chimney	DeLee	Pediatric	Coiled	TRASC106C	50	V
6fr	Y	DeLee	Pediatric	Coiled	TRASC306C	50	P
6fr	Y	DeLee	Pediatric	Straight	TRASC306S	100	P
8fr	Chimney	DeLee	Pediatric	Coiled	TRASC108C	50	V
8fr	Chimney	DeLee	Pediatric	Straight	TRASC108S	50	V
8fr	Y	DeLee	Pediatric	Coiled	TRASC308C	50	P
8fr	Y	DeLee	Pediatric	Straight	TRASC308S	100	P
10fr	Y	Whistle	Adult	Coiled	TRASC110C	50	P
10fr	Y	Whistle	Adult	Straight	TRASC110S	100	P
12fr	Y	Whistle	Adult	Coiled	TRASC112C	50	P
12fr	Y	Whistle	Adult	Straight	TRASC112S	100	P
14fr	Y	Whistle	Adult	Coiled	TRASC114C	50	P
14fr	Chimney	Cloude	Adult	Straight	TRASC114CVS	50	V
14fr	Y	Whistle	Adult	Straight	TRASC114S	100	P
16fr	Y	Whistle	Adult	Coiled	TRASC116C	50	P V
16fr	Y	Whistle	Adult	Straight	TRASC116S	100	P
18fr	Y	Whistle	Adult	Coiled	TRASC118C	50	V
18fr	Chimney	DeLee	Adult	Straight	TRASC118S	50	V

- Sterile
- Single patient use
- Latex-free



Disposable Rigid Suction Canister

Item#	Size	Direct to Regulator	Qty
RES023S-24	800cc	No	24
RES025AW1200R	1200cc	Yes	40
RES025AW1200-40	1200cc	No	40
RES025AW2000	2000cc	No	40
RES025AW3000	3000cc	No	40



- Prevents backflow and premature shut off.
- Includes 20-degree patient port for 1200cc, or straight port for 800cc, to prevent tubing from kinking.
- Durability assurance for collecting and retaining aspirated fluids.
- Disposable.

Flexible Disposable Liner

50/case

Item#	Size	Includes Solidifier
LIN1000FD	1000cc	No
LIN1500FD	1500cc	No
LIN3000FD	3000cc	No
LIN1000FDS	1000cc	Yes
LIN1500FDS	1500cc	Yes
LIN3000FDS	3000cc	Yes



- Part of a complete fluid management system.
- These soft, plastic liners fit into the reusable cannisters, creating less plastic waste and greater material utilization.
- Not made with natural rubber latex.

Reusable Outer Canister for Disposable Liner with Built In On/Off Valve

10/case

For Flexible Disposable Liner

Item#	Size
LIN1000FDC	1000cc
LIN1500FDC	1500cc
LIN3000FDC	3000cc

For Semi-Rigid Disposable Liner

Item#	Size
LIN1000SRDC	1000cc
LIN1500SRDC	1500cc
LIN3000SRDC	3000cc



- Reusable outer cannisters provide a safe, convenient, and economical solution for fluid collection.
- Clear plastic cannisters are designed to hold the disposable liners.
- Reusable polycarbonate canister/bracket assembly has pre-attached tubing for connection to the liner.

Semi-Rigid Disposable Liner

Item#	Size
LIN1000SRD	1000cc
LIN1500SRD	1500cc
LIN3000SRD	3000cc

Canister/Liner Kits

Item#	Size and Type	Tubing Size and Length	Qty
LIN1000FD618	1000cc Flexible Liner	6mm x 1.8m	50
LIN1000FD627	1000cc Flexible Liner	6mm x 2.7m	50
LIN1500FD618	1500cc Flexible Liner	6mm x 1.8m	50
LIN1500FD627	1500cc Flexible Liner	6mm x 2.7m	50
LIN3000FD618	3000cc Flexible Liner	6mm x 1.8m	50
LIN1000SRD618	1000cc Semi-Rigid Liner Kit	6mm x 1.8m	72
LIN1000SRD627	1000cc Semi-Rigid Liner Kit	6mm x 2.7m	72
LIN1500SRD518	1500cc Semi-Rigid Liner Kit	5mm x 1.8m	72
LIN1500SRD618	1500cc Semi-Rigid Liner Kit	6mm x 1.8m	72
LIN1500SRD627	1500cc Semi-Rigid Liner Kit	6mm x 2.7m	72
LIN3000SRD618	3000cc Semi-Rigid Liner Kit	6mm x 1.8m	36
RES023S-T	800cc Disposable Suction Canister Kit	6mm x 1.8m	45
RES025AW1200-T518	1200cc Disposable Suction Canister Kit	5mm x 1.8m	30
RES025AW1200-T618	1200cc Disposable Suction Canister Kit	6mm x 1.8m	30
RES025AW1200-T518V	1200cc Disposable Suction Canister Kit	5mm x 1.8m 5mm x .46m Vacuum Tubing	30



Four-Canister Roll Stand

Item# **RES9100**
30"
Each



Wall Plate

Item# **RES9101**
10/case



Disposable Suction Canister Ring Bracket

20/case

Item#	Size
RES9102	1200cc
RES9103	2000-3000cc



Disposable Suction Canister Ring Bracket

20/case

Item#	Size
RES9102T	1200cc
RES9103T	2000cc



Disposable Suction Canister Ring Bracket Extended for Regulator with DISS

20/case

Item#	Size
RES9102D	1200cc
RES9103D	2000-3000cc



Disposable Suction Canister Ring Bracket with On/Off Valve

20/case

Item#	Size
RES9102V	1200cc
RES9103V	2000-3000cc



Semi-Rigid Suction Canister Ring Bracket Close to Wall

Each

Item#	Size
RES9104	1000cc
RES9105	1500cc



Semi-Rigid Suction Canister Ring Bracket Extended For Regulator

Item# **RES9105E**
1500cc
Each



Suction Tubing — Single Use Sterile

50/case


Item#	Diameter	Length
RES025-31618	5mm	46cm (18in)
RES025-31672	5mm	1.8m (6ft)
RES025-316120	5mm	3.0m (10ft)
RES025-316144	5mm	3.7m (12ft)
RES025-316240	5mm	6.1m (20ft)
RES025-14072	6mm	1.8m (6ft)
RES025-14120	6mm	3.0m (10ft)
RES025-14144	6mm	3.7m (12ft)
RES025-14240	6mm	6.1m (20ft)
RES025-93272	7mm	1.8m (6ft)
RES025-932120	7mm	3.0m (10ft)
RES025-932144	7mm	3.7m (12ft)
RES025-932240	7mm	6.1m (20ft)
RES025-38120	Large Bore 9.5mm	3.0m (10ft)



Yankauer — Single Use Sterile

100/case

Item#	Tip	Vented
TRAA002S	Bulbous	No
TRAA001S	Bulbous	Yes
TRAA004S	Straight	No
TRAA003S	Straight	Yes





- These lightweight yankauers are designed to provide optimal control and visibility while helping to minimize tissue trauma.

Bulk Roll — Non-Sterile

4/case

Item#	Diameter	Length
RES025B-530	5mm	30m
RES025B-630	6mm	30m
RES025B-730	7mm	30m


- Clear bulk tubing in rolls that can be cut into the exact length needed.
- Cut-to-size female connectors are located at 6ft intervals.
- Disposable accessories make it easy to add auxiliary items to your system.
- Flex Advantage liners, Guardian systems, Medi-Vac accessories are interchangeable to help you select the most effective setup for your needs.
- Allows you to adapt your system for various setups.
- Used to connect vacuum source to canisters and manifolds.



Solidifier — Dissolvable Sachet

100/case

Item#	Size
RES8001	150cc
RES8004	450cc
RES8015	1500cc



Economical

Pre-measured rupturing pouches eliminate the tearing and spillage associated with conventional solidifiers for medical waste.

Particulate-Free

Self-activating technology reduces particulates and bacterial exposure, and eliminates harmful dust.

Reduce Exposure

Helps to minimize the risk of cross contamination of infectious waste to medical staff.

Solidifier with Treatment — Dissolvable Sachet

100/case

Item#	Size
RES8002T	250cc
RES8005T	500cc
RES8007T	750cc
RES8015T	1500cc



Contain spills quickly with this astonishingly absorbent liquid solidifier packet

Super absorbent packets quickly absorb up to 400 times their own weight in aqueous solutions. These drop-in liquid absorbers offer the ideal solution for shipping, storing, and handling of aqueous solutions that must be contained quickly and effectively in the event of a leak.

Smart Design

Specially formulated pouch bursts open in the presence of a wide range of fluids to quickly immobilize them. Solidified material is easily contained by an outer vessel, shipping container, or other such device, resulting in safer handling.

Ideal for Medical Use

The most common applications for the use include medical shipping containers and bio-specimen transportation bags used for vials, test kits, and lab trays.

SperaChamber **V** **P**

Anti-Static Valved Inhaler Chamber
for Metered Dose Inhaler

Item# **RES195-10**
10/case

Item# **RES195-50**
50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



Adult Mask



SperaChamber with Mask 10/case

Item#	Size
RES195MS	Small (Infant)
RES195MM	Medium (Pediatric)
RES195ML	Large (Adult)

Mask Only

Item#	Size	Qty
RES195MASK-S	Small (Infant)	Each
RES195MASK-S-10	Small (Infant)	10/box
RES195MASK-M	Medium (Pediatric)	Each
RES195MASK-M-10	Medium (Pediatric)	10/box
RES195MASK-L	Large (Adult)	Each
RES195MASK-L-10	Large (Adult)	10/box

Peak Flow Meter **V** **P**

Item# **RES052**
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD** **V** **P**
Pediatric Mouthpiece
for Peak Flow Meter
10/case

Plastic Intake Filter with Head For Respiration Inspiration 626 & 929 **V** **P**

Item# **BF016**
10/pack



Valved T-Adapter, Aerosol **V** **P**

Item# **NEB9001T**
22mm ID x 22mm OD
10/case



Valved T-Adapter — Adult **V** **P**

Item# **VEN001**
22F-18F-22M/15F

Each Premier
10/pack Vizient



Compressor Nebulizers V P

Item# NEB100

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Foam Filters (BF018 x5)

NEB100-R

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)

SCAN FOR VIDEO



Sunny the Bear

Item# NEB300BEAR



Gordon the Fish

Item# NEB300FISH



NEB100-BAG V P
Carry Bag

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

Handheld Compressor Nebulizer V P

Item# NEB200

Included:

- Compressor
- Medication cup
- AC power adapter
- Air hose
- Elbow
- Mouthpiece
- Mask
- 4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (Approximately 3 microns)
- No replacement filter required



Mini Mesh Nebulizer with Rechargeable Battery V P

Item# NEB401

Included:

- Main unit
- Medicine cup
- Disposable mouthpiece
- Disposable adult mask
- Disposable pediatric mask
- USB cord
- Carrying bag
- Instruction manual



Small particle size
MMAD 2.2µm*

Efficient
Low Residual Liquid

Convenient
Small and Portable

Quiet
≤50dB (at one meter)

Modern
Rechargeable Lithium Battery

Disposable Nebulizer Kit

Item#

RES090 **V** **P**

Adult Size

50/case

- 6cc jet nebulizer bottle
- Mask
- 7ft supply tube

RES090A **V** **P**

No Accessories

50/case

RES090NM **P**

- 6cc jet nebulizer bottle
- 7ft supply tube

50/case

RES092 **V** **P**

Pediatric Size

50/case

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube



Disposable Nebulizer Kit with T-Piece **V** **P**

Item#

RES091

Disposable
Single Patient

50/case

RES091U

Disposable
Single Patient,
Universal Connector

50/case

RES091NA

Disposable
Single Patient,
No Aerosol

50/case

HPCPS# A7003

- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal



Includes:

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube

Reusable Nebulizer Kit with T-Piece **P**

Item# **RES093**

Reusable — Single Patient

10/case

HPCPS# A7005

Includes:

- 10cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube
- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°



Disposable Nebulizer Kit with Angled Mouthpiece **V** **P**

Item#

RES096

Disposable
Single Patient

50/case

HPCPS# A7003

RES096U

Disposable Single Patient,
Universal Connector

50/case

Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube



Aerosol Mask

50/case

Item#

RES2110 **V** **P**

Adult Size

RES2110F **P**

Adult Size
w/ Flex Tube

RES2210 **V** **P**

Pediatric Size

- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged




Closed Ventilation Suction System V

Standard Switch

Sterile and latex-free

ISO color coded for instant recognition of catheter sizes.



With the unique isolation valve, 72hr closed suction system features safe and simple access for mechanical ventilated patients for 72 hours of continuous use.

Once the airway is clear, withdraw the catheter completely into the plastic sheath. A black mark on the catheter indicates full withdrawal.

Smooth tubing surface and atraumatic rounded catheter tip for minimal tracheal wall injury.

Closed Ventilation Suction System V

Standard Switch

Sterile and latex-free

24 Hour Adult

Endotracheal 530mm With Swivel Connector

Item#	Description
VEN72453010M	10Fr, MDI Port
VEN72453012M	12Fr, MDI Port
VEN72453016M	16Fr, MDI Port

Endotracheal 530mm With T-Piece

Item#	Description
VEN82453010M	10Fr, MDI Port
VEN82453012M	12Fr, MDI Port
VEN82453014M	14Fr, MDI Port
VEN82453016M	16Fr, MDI Port

Tracheostomy With Swivel Connector

Item#	Description
VEN72430510M	305mm, 10Fr, MDI Port
VEN72430512M	305mm, 12Fr, MDI Port
VEN72430514M	305mm, 14Fr, MDI Port
VEN72430516M	305mm, 16Fr, MDI Port

Tracheostomy With T-Piece

Item#	Description
VEN82430510M	305mm, 10Fr, MDI Port
VEN82430512	305mm, 12Fr
VEN82430512M	305mm, 12Fr, MDI Port
VEN82430514	305mm, 14Fr
VEN82430514M	305mm, 14Fr, MDI Port
VEN82430516	305mm, 16Fr
VEN82430516M	305mm, 16Fr, MDI Port
VEN82440512M	405mm, 12Fr, MDI Port
VEN82440514M	405mm, 14Fr, MDI Port
VEN82440516M	405mm, 16Fr, MDI Port

72 Hour Adult

Endotracheal 530mm With Swivel Connector

Item#	Description
VEN77253010M	10Fr, MDI Port
VEN77253012M	12Fr, MDI Port
VEN77253014M	14Fr, MDI Port
VEN77253016M	16Fr, MDI Port

Tracheostomy 305mm With Swivel Connector

Item#	Description
VEN77230510M	10Fr, MDI Port
VEN77230512M	12Fr, MDI Port
VEN77230514M	14Fr, MDI Port
VEN77230516M	16Fr, MDI Port

72 Hour Neonatal/Pediatric

With Swivel Connector

Item#	Description
VEN77230505	305mm, 5Fr
VEN77230506	305mm, 6Fr
VEN77230507	305mm, 7Fr
VEN77230508	305mm, 8Fr



24 Hour Neonatal/Pediatric

With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600	305mm, 6Fr
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700	305mm, 7Fr
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800	305mm, 8Fr
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000	405mm, 10Fr
VEN66001000M	405mm, 10Fr, MDI Port

With Swivel Connector

Item#	Description
VEN72430505	305mm, 5Fr
VEN72430506	305mm, 6Fr
VEN72430507	305mm, 7Fr
VEN72430508	305mm, 8Fr

Closed Ventilation Suction System

Deluxe Switch (one-handed operation)

Sterile and latex-free

Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.


Isolation Valve*

Seals the irrigation chamber to provide a better cleaning of the catheter. Allows for continuous ventilation and prevents backflow of saline during cleaning to ensure patient comfort and safety.

*Only available with 72 Hour catheter

Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.



Closed Ventilation Suction System V

Deluxe Switch (One-handed operation)

Sterile and latex-free

24 Hour Adult

Endotracheal 540mm

Item#	Description
VEN160100A	10Fr
VEN160100AM	10Fr, MDI Port
VEN162100A	10Fr, Wedge + Adapter
VEN162100AM	10Fr, MDI Port, Wedge + Adapter
VEN160120	12Fr
VEN160120M	12Fr, MDI Port
VEN162120	12Fr, Wedge + Adapter
VEN162120M	12Fr, MDI Port, Wedge + Adapter
VEN160140	14Fr
VEN160140M	14Fr, MDI Port
VEN162140	14Fr, Wedge + Adapter
VEN162140M	14Fr, MDI Port, Wedge + Adapter
VEN160160	16Fr
VEN160160M	16Fr, MDI Port
VEN162160	16Fr, Wedge + Adapter
VEN162160M	16Fr, MDI Port, Wedge + Adapter

Tracheostomy 305mm

Item#	Description
VEN161100A	10Fr
VEN161100AM	10Fr, MDI Port
VEN163100A	10Fr, Wedge + Adapter
VEN163100AM	10Fr, MDI Port, Wedge + Adapter
VEN161120	12Fr
VEN163120	12Fr, Wedge + Adapter
VEN161140	14Fr
VEN163140	14Fr, Wedge + Adapter
VEN161160	16Fr
VEN161160M	16Fr, MDI Port
VEN163160	16Fr, Wedge + Adapter

72 Hour Adult

Endotracheal 540mm

Item#	Description
VEN365100A	10Fr
VEN365100AM	10Fr, MDI Port
VEN367100A	10Fr, Wedge + Adapter
VEN367100AM	10Fr, MDI Port, Wedge + Adapter
VEN365120	12Fr
VEN365120M	12Fr, MDI Port
VEN367120	12Fr, Wedge + Adapter
VEN367120M	12Fr, MDI Port, Wedge + Adapter
VEN365140	14Fr
VEN365140M	14Fr, MDI Port
VEN367140M	14Fr, Wedge + Adapter
VEN365160	16Fr
VEN365160M	16Fr, MDI Port
VEN367160	16Fr, Wedge + Adapter
VEN367160M	16Fr, MDI Port, Wedge + Adapter

Tracheostomy 305mm

Item#	Description
VEN366100A	10Fr
VEN366100AM	10Fr, MDI Port
VEN368100A	10Fr, Wedge + Adapter
VEN368100AM	10Fr, MDI Port, Wedge + Adapter
VEN366120	12Fr
VEN366120M	12Fr, MDI Port
VEN368120	12Fr, Wedge + Adapter
VEN368120M	12Fr, MDI Port, Wedge + Adapter
VEN366140	14Fr
VEN366140M	14Fr, MDI Port
VEN368140	14Fr, Wedge + Adapter
VEN368140M	14Fr, MDI Port, Wedge + Adapter
VEN366160	16Fr
VEN366160M	16Fr, MDI Port
VEN368160	16Fr, Wedge + Adapter
VEN368160M	16Fr, MDI Port, Wedge + Adapter

Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port



Swivel Straight Connector V

Item#
VEN005-20

Disposable
22mm ID, 15mm OD
20/pack



Adult Resuscitator P

12/case

Item#	Pop-Off	Mask
RES8100	No	None
RES8100H	Yes	None
RES8101	Yes	#4
RES8102	Yes	#5



Child Resuscitator P

12/case

Item#	Pop-Off	Mask
RES8110H	Yes	None
RES8111	Yes	#2
RES8112	Yes	#3



Infant Resuscitator P

12/case

Item#	Pop-Off	Mask
RES8120H	Yes	None
RES8121	Yes	#0
RES8122	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve

Disposable Pressure Manometer P

Item# **RES8201**
25/case



- Pressure range 0-60 cm H2O
- Indicator strip
- All models have the pressure range emboss printed on the housing
- Some users like to also have a color strip / indicator

Disposable Pressure Manometer, L-Shaped P

Item# **RES8201L**
25/case



- Pressure range 0-60 cm H2O
- Indicator strip
- All models have the pressure range emboss printed on the housing
- Some users like to also have a color strip / indicator

Reusable PEEP Valve, Orange P

2-10cm H2O
10/case

Item#
RES8203R
15mm ID



RES8203RL
30mm ID

- Available in orange color in 2-10cm H2O and blue color in 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Materials: Polycarbonate, silicone, stainless steel for spring

Reusable PEEP Valve, Blue P

5-20cm H2O
10/case

Item#
RES8204B
15mm ID



RES8204BL
30mm ID

- Available in orange color in 2-10cm H2O and blue color in 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Materials: Polycarbonate, silicone, stainless steel for spring

Disposable PEEP Valve, Orange P

2.5-10cm H2O
10/case

Item#
RES8205R
15mm ID



RES8205RL
30mm ID

- Available in 2.5-10cm H2O and 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Single patient use
- Materials: Polycarbonate, silicone, stainless steel

Disposable PEEP Valve, Blue P

5-20cm H2O
10/case

Item#
RES8205B
15mm ID



RES8205BL
30mm ID

- Available in 2.5-10cm H2O and 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Single patient use
- Materials: Polycarbonate, silicone, stainless steel

Reservoir Bag, PVC P

25/case

Item#	Size
RES8300P	0.5L
RES8301P	1L
RES8302P	2L



Reservoir Bag, Silicone P

25/case

Item#	Size
RES8300S	0.5L
RES8301S	1L
RES8302S	2L



Suitable for reusable resuscitator.

CPR Pocket Resuscitator Set P

Item# **RES8200**

12/case

- High transparency for better visualization
- Non-reusable, single use only
- 100% Medical level PVC material
- Latex-free components
- Economical choice
- Blow molding



Breathing Mask P

20/case

Item#	Size	Hooks
SBC01N	Neonatal (#0)	No
SBC01I	Infant (#1)	No
SBC01P	Pediatric (#2)	No
SBC01S	Adult Small (#3)	No
SBC01M	Adult Medium (#4)	No
SBC01L	Adult Large (#5)	No
SBC01XL	Adult Extra Large (#6)	No
SBC02N	Neonatal (#0)	Yes
SBC02I	Infant (#1)	Yes
SBC02P	Pediatric (#2)	Yes
SBC02S	Adult Small (#3)	Yes
SBC02M	Adult Medium (#4)	Yes
SBC02L	Adult Large (#5)	Yes
SBC02XL	Adult Extra Large (#6)	Yes



Silicone Disposable Laryngeal Masks P

Item#	Size
ALM1001	#1
ALM1015	#1.5
ALM1002	#2
ALM1025	#2.5
ALM1003	#3
ALM1004	#4
ALM1005	#5

- Medical grade silicone for super compatibility and patient comfort
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters
- Available with seven sizes for all patient use



Silicone Reusable Laryngeal Masks P

Item#	Size
ALM1001R	#1
ALM1015R	#1.5
ALM1002R	#2
ALM1025R	#2.5
ALM1003R	#3
ALM1004R	#4
ALM1005R	#5

- Reusable Silicone Laryngeal Masks are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Serial no. is printed on each connector for easy tracking
- Additional user card for record of device usage times and patient information is supplied in each pouch



Flexible Silicone Disposable Laryngeal Masks P

Item#	Size
ALM3001	#1
ALM3015	#1.5
ALM3002	#2
ALM3025	#2.5
ALM3003	#3
ALM3004	#4
ALM3005	#5

- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness



Flexible Silicone Reusable Laryngeal Masks P

Item#	Size
ALM3001R	#1
ALM3015R	#1.5
ALM3002R	#2
ALM3025R	#2.5
ALM3003R	#3
ALM3004R	#4
ALM3005R	#5

- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness



PVC Disposable Laryngeal Masks P

Item#	Size
ALM2001	#1
ALM2015	#1.5
ALM2002	#2
ALM2025	#2.5
ALM2003	#3
ALM2004	#4
ALM2005	#5

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



Flexible PVC Disposable Laryngeal Masks P

Item#	Size
ALM4001	#1
ALM4015	#1.5
ALM4002	#2
ALM4025	#2.5
ALM4003	#3
ALM4004	#4
ALM4005	#5

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



Size Patient Weight and Age

#1	Up to 5 kg	Pediatric
#1.5	5-10 kg	Pediatric
#2	10-20 kg	Pediatric
#2.5	20-30 kg	Pediatric
#3	30-50 kg	Pediatric
#4	50-70 kg	Adult
#5	70-100 kg	Adult

PeDIA - Pediatric Device for Induction of Anesthesia V P

Item#

ANEP021Y

With Mask Adapter Elbow


12/pack



2025 Innovation Celebration Winner

Use with
mask

Or use
without
mask



Turn panic into play

For children, undergoing surgery is already unsettling, and the process of induction is often the most traumatizing part. Often kids will kick and scream, believing they're being trapped by the gas or suffocated by the mask. The experiences create negative emotional scarring that can last weeks, months or even years.

Specially designed for children undergoing surgery, the Pediatric Device for Induction of Anesthesia device, or "PeDIA," for short, includes a colorful whistling balloon, providing kids with something fun to focus on and turning a stressful moment into a lighter one.

Whether they're squeezing the balloon and focusing on its whistle or directly putting themselves to sleep by breathing the gas through balloon delivery, the PeDIA balloon captivates children's attention, allowing their active participation in the induction process.

It's a simple way to add a brighter touch to the medical experience, ease anxiety and bring smiles to patients, parents and staff.

This unique approach—presenting an interactive, familiar toy-like device—not only diverts children's focus from the impending procedure but also significantly reduces the anxiety and resistance commonly associated with mask induction. The aim is to decrease the anxiety felt by over 50% of pediatrics as well as parents and staff.

The PeDIA can be used alone as a mask-free device or, for younger children, with an anesthesia mask.

SCAN FOR VIDEOS



