



Product Catalog

EUROPE



Policies

Ordering Options

Three Easy Ways To Order

Phone

(877) 5-SUNSET (877-578-6738)

Our customer service professionals are here to assist you
M-F 8:00 AM – 5:00 PM Central Time.

Email

intlsales@sunsethcs.com

customerservice@sunsethcs.com

*Contact us for more information
on fulfillment and compliance
solutions that best fit your needs.*

Payment Methods

Payment Portal and ACH

Visit www.sunsethcs.com/paymentportal to pay online with ACH or credit card.

Net 30-Day Payment Terms

For qualified companies. A 1.5% monthly finance charge will be assessed on all past due invoices.

COD

For non-qualified companies.

Credit Card

We accept all major credit cards. A processing fee may apply. After 10 days of invoice date, a surcharge of 2.5% will be applied.

VISA



DISCOVER

Customer Service

Returns

Contact us within 7 days of receiving your shipment.

Call us for a return authorization number. There is no restocking charge on orders returned within 20 days of invoice. Orders returned after 20 days may be assessed a restocking fee of 15%. No returns after 90 days.

(877) 5-SUNSET (877-578-6738)

Central Distribution Center

Sunset Healthcare Solutions
2725 Alft Ln
Elgin, IL 60124

Western Distribution Center

Sunset Healthcare Solutions
6750 W Van Buren St Ste 100
Phoenix, AZ 85043

Problems

Please call us immediately if you have any complaints (it is helpful to ask for the customer service representative with whom you placed your order).

Or email us at: customerservice@sunsethcs.com

Questions?

Our experienced customer service staff will be glad to provide you with information. Just give us a call.

Pricing

All prices are in U.S. Dollars.

Taxes are added where applicable.

We make every effort to hold pricing throughout the year.

However, prices are subject to change without notice.

Content

CPAP

	Page
Filters	5-9
Respironics	5
ResMed/Sullivan	6
ResMed	6
Fisher & Paykel	7
3B	7
Breas/GE	8
Breas	8
Covidien Nellcor Puritan Bennett	9
Beyond	9
Resvent	9
Löwenstein	9
Chinstraps	10-11

Oxygen

	Page
Connectors and Accessories	12-13

Heat Moisture Exchangers

	Page
Heat Moisture Exchangers	14

Respironics

DreamStation 2



Item# **CF1012**
 Disposable Filter
6/pack
 OEM# 1142831
CF1012-1 1/pack OEM# 1142828
CF1012-2 2/pack OEM# 1142829

Item# **CF1012F**
 Reusable Filter
2/pack
 OEM# 1142687
CF1012F-1 1/pack

DreamStation



Item# **CF1010**
 Disposable Filter
6/pack
 OEM# 1122519
CF1010-1 1/pack
CF1010-2 2/pack

Item# **CF1010F**
 Reusable Filter
2/pack
 OEM# 112446 (1-pack)
CF1010F-1 1/pack

M-Series



Item# **CF1007**
 Ultrafine Filter - Tabbed
6/pack • 1³/₄" x 7/₈"
CF1007-1 1/pack

Item# **CF1007F**
 Foam Filter
2/pack • 1⁵/₈" x 7/₈"
 OEM# 1029330
CF1007F-1 1/pack

PR System One



Item# **CF1009NT**
 Ultrafine Filter - No Tab
6/pack • 1³/₄" x 7/₈"
 OEM# 1029331
CF1009NT-1 1/pack

Item# **CF1009**
 Ultrafine Filter - Tabbed
6/pack • 1³/₄" x 7/₈"
CF1009-1 1/pack

Item# **CF1009F**
 Foam Filter
2/pack • 1⁵/₈" x 7/₈"
 OEM# 1029330
CF1009F-1 1/pack

REMstar Pro, REMstar Plus



Item# **CF1006**
 Ultrafine Filter
6/pack • 3¹/₂" x 1¹/₂"
CF1006-1 1/pack
 OEM# 1005945

Item# **CF1006F**
 Foam Filter
2/pack • 3⁵/₈" x 1⁵/₈"
CF1006F-1 1/pack
 OEM# 1005964

ResMed/Sullivan

AirSense 11



Item#
CF2111
Disposable Filter
3/pack
CF2111-1 1/pack OEM# 39300
CF2111-6 6/pack
CF2111-12 12/pack



S9 and AirSense 10 Series



Item# **CF2107**
Disposable Filter
3/pack • 2¹/₈" x 1³/₈"
CF2107-1 1/pack OEM# 36850
CF2107-2 2/pack OEM# 36851
CF2107-12 12/pack

Adapt SV



Item# **CF2106**
Disposable Filter
3/pack • 1¹/₄" x 1¹/₂"
OEM# 26926
CF2106-1 1/pack

Item# **CF2106HA**
Hypoallergenic Filter
3/pack • 1¹/₄" x 1¹/₂"
OEM# 26927
CF2106HA-1 1/pack

Item# **CF2107HA**
Hypoallergenic Filter
3/pack • 2¹/₈" x 1³/₈"
CF2107HA-1 1/pack OEM# 36855

AirMini



Item#
CF2108
Disposable Filter
2/pack • 1¹/₄" x 1³/₁₆"
OEM# 38815
CF2108-1 1/pack

ResMed

S8 Series



Item#
CF2105
Disposable Filter
3/pack • 1³/₈" x 1³/₈"
OEM# 33916
CF2105-1 1/pack OEM# 33915

Item#
CF2105HA
Hypoallergenic Filter
3/pack • 1³/₈" x 1³/₈"
OEM# 33920
CF2105HA-1 1/pack OEM# 33919

S7 Series



Item#
CF2103
Disposable Filter
3/pack • 1¹/₄" x 1¹/₂"
OEM# 30907
CF2103-1 1/pack OEM# 30919

Item#
CF2103HA
Hypoallergenic Filter
3/pack • 1¹/₄" x 1¹/₂"
CF2103HA-1 1/pack OEM# 30933

Fisher & Paykel

SleepStyle Series



Item# **CF4004-2**
Disposable Filter
2/pack
CF4004-1 1/pack



ICON Series



Item# **CF4003**
Disposable Filter
4/pack
2 1/16" x 5/8"
CF4003-1 1/pack OEM# 900ICON50



HC230 Series, HC600 Series



Item# **CF4002**
Disposable Filter
4/pack • 2 3/4" x 7/8"
OEM# 900HC240
CF4002-1 1/pack



HC200, HC220, HC221



Item# **CF4001**
Disposable Filter
4/pack • 5" x 3/4"
OEM# 900HC222
CF4001-1 1/pack



3B

Luna G3



Item# **CF8006-1**
Disposable Filter
1/pack



Item# **CF8006F-1**
Reusable Filter
1/pack



Luna and Luna II



Item# **CF8005-1**
Disposable Filter
1/pack • 2" x 1 3/16"



Item# **CF8005F-1**
Reusable Filter
1/pack • 2" x 1 3/16"



RESmart



Item# **CF8004-1**
Disposable Filter
1/pack • 3 5/16" x 1 5/8"



Breas/GE

iSleep 20 Series



Item# **CF7003-1**
Ultrafine Filter
1/pack • 2¹/₄" x 7⁷/₈"
OEM# 004153

Item# **CF7003F**
Foam Filter
1/pack • 2³/₄" x 1"
OEM# 004154

PV10



Item# **CF7002-1**
Ultrafine Filter
1/pack • 2³/₄" x 7⁷/₈"
OEM# 2B01975

Item# **CF7002F**
Foam Filter
Each • 2⁵/₈" x 7⁷/₈"
OEM# 2B01972

PV1




Item# **CF7001-1**
Ultrafine Filter
1/pack • 2¹/₂" x 1¹/₄"



Item# **CF7001F-1**
Foam Filter
1/pack • 2¹/₂" x 1¹/₄"

Breas

45LS Style Air Inlet Filter



Item# **VF010-5**
Disposable filter
5/pack
OEM# 007103



Item# **VF010F-5**
Reusable Filter
5/pack
OEM# 007104

Beyond

ResPlus



Item#
CF9005F-2
Reusable Filter
2/pack
CF9005F-1 1/pack

Resvent

iBreeze



Ultrafine Filter
CF9901-1 1/pack
CF9901-2 2/pack
CF9901-6 6/pack
CF9901-50 50/pack

Löwenstein

LUISA



Item# **VF008F-2**
Reusable Filter
2/pack
OEM# WM29928

Covidien Nellcor Puritan Bennett

GoodKnight 420/425



Item# **CF6003F**
Foam Filter
Each • 1³/₁₆" x 3³/₁₆"
OEM# M-413950-02

Sunset Chinstraps

Universal Neoprene Chinstrap



Item# **CS001**
Universal Chinstrap

- Latex-free
- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows

Adjustable Neoprene Chinstrap



Size	Item#
6in	CS003S
7in	CS003R
8in	CS003M
9in	CS003L
10in	CS003XL

Compare to Puresom Adjustable Ruby
OEM# TMS-09 ADJ

- Latex-free
- Thin adjustable strap attaches at the top of the head. Thicker bottom strap attaches to itself at the back of the head.
- Great for hard-to-fit patients

For fitting, measure along one side of the head, as indicated by the dashed line above.

Comfort Chinstrap



Item# **CS005**
Compare to ResMed Sullivan
OEM# 16015

- Latex-free
- One size fits most
- Easy to use —attaches to itself at the top of the head
- Soft cloth chin cup adds comfort

Basic Chinstrap



Item# **CS002**
Compare to Respirationics Basic
OEM# 302175

- Latex-free
- One size fits most
- Attaches to the top strap of headgear
- Cups the chin for comfort
- Works best with any nasal mask

Deluxe Chinstrap



Size	Item#
26in	CS004S
28in	CS004L

Compare to Respirationics Deluxe
OEM# 302425

- Latex-free
- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort

Premium Chinstrap




Item# **CS006**
Compare to Respirationics Premium
OEM# 1012911

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.

Sunset Chinstraps

Ruby Style Chinstrap



Size	Item#
21.5-22.5in	CS007S
22-23in	CS007M
23-24in	CS007L

Compare to Puresom Ruby
OEM# TMS-09

- Latex-free
- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

Halo Style Chinstrap



Item# **CS025**
Compare to BreatheWear
Halo Chinstrap
OEM# H3407

Ruby Style Adjustable Chinstrap



Size	Item#
Regular	CS007-ADJ
Extra Large	CS007-ADJXL

- Latex-free
- Fully adjustable
- Super soft materials for extra comfort



Oxygen Accessories

Wrenches and Chains

Item# **RES031**
With Chain **10/pack**

Item# **RES032**
Without Chain **10/pack**



- For use with post valve oxygen cylinders: A, B, C, D, JD or E tanks
- Made of tough, stress-resistant nylon

Item# **RES031M**
Metal Wrench with Chain for Post Valve Oxygen Cylinders and CGA 540 Style Regulators
5/pack



Humidifier Connectors

Item# **HUM008**
L-Connector 12in • **10/pk**

Item# **HUM008L**
L-Connector 21in • **10/pk**

Item# **HUM008S**
Straight Connector 11.625in • **10/pk**



Oxygen Liter Meter

Item# **RES034**
0-10 Liter/Minute
Each




- Oxygen pen
- Works with tanks and concentrators

Cannula Ear Cushions

Item# **RES020**
Cushions Only
50 Pairs/case

- Add cushions to existing cannula



Disposable Bubble Humidifier 500 ml

Item# **RES050**
Up to 15 LPM • 6 psi • 500 ml
50/case
OEM# HUM001



- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient

Water Traps

Item# **RES037**
Push On End Caps • **Each**



Item# **RES037T**
Threaded End Caps • **Each**

- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

Oxygen Connectors and Accessories

Christmas Tree Connector

Item# **RES001**
Green without Swivel
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections

Item# **RES002**
Green with Swivel
10/pack



- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Homefill Tubing Connector


Item# **RES039**
10/pack

- Compatible with Invacare, Drive, and Respirationics homefill machines
- Durable, coiled tubing



Homefill Racks

Item# **RES039-RACKB**
For Invacare Perfecto₂
4/case



RES039-RACKW
For Invacare Platinum
4/case



RES039-RACKD
For DeVilbiss 525/1025 Series
4/case



Easily attaches onto concentrator and is designed to conveniently hold Homefill and transfill oxygen systems.

Oxygen Tubing Connector — Non-Swivel

Item# **RES012**
Fits 3/16in and 1/4in ID Tubing
10/pack




- Barbed ends ensure secure, leak-free connections

Oxygen Tubing Y-Connector

Item# **RES017**
Fits 3/16in and 1/4in ID Tubing (All 3 sides are the same)
10/pack

- Connects 2 oxygen sources for higher flows



Oxygen Tubing Connector — Swivel

Item# **RES018**
Male/Male **10/pack**



RES018KIT-1
Female/Male **Each**



- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing

Oxygen Enrichment Adapter

Item# **RES019**
22mm with 6mm Rotating Stem **10/pack**



Washers with Rubber Ring

Item# **RES035**
Brass
25/pack



HMEs

Heat Moisture Exchangers

30/case*

*RES027FLEX 40/case



	RES027FLEX Flex HME — Foam	RES027 T-HME	RES027IL Inline HME with Gas Sampling Port	RES027NE Inline HME — Neonatal	RES027OP T-HME with O2 Port
TV Range (mL)	250-1500	150-1250	70-600	5-50	150-1250
Moisture Output (mg H₂O/L)	30 mg/L	36.4 at VT = 500 ml	34.5 at VT = 250 ml 32.5 at VT = 500 ml	32.8 at VT = 250 ml	32 mg/L
Weight (g)	15.5	6.6	10	1.9	7.4
Resistance (cm H₂O)	0.5 @ 30 LPM 1.4 @ 60 LPM	0.2 @ 30 LPM	0.8 @ 30 LPM 2.1 @ 60 LPM	0.8 @ 15 LPM 1.7 @ 30 LPM	0.2 @ 30 LPM
Dead Space (mL)	90	15	21	2.2	7
Connection (mm)	22F ,15F	15F	22M/15F-15M	15M,15F	15F 5.5mm O2 port

