

# Micro™ 510 Auto CPAP



Experience the freedom of sleep therapy without the bulk and hassle of a traditional CPAP machine. Designed with your active lifestyle in mind, the Transcend Micro 510 Auto CPAP is the epitome of convenience and portability, ensuring you receive the treatment you need, wherever life takes you.

## Features

- Ultra-Compact
- Whisper Quiet
- Universal Compatibility
- Travel-Friendly
- Advanced Comfort Features
- Easy Maintenance

### Compact and Lightweight Design

The Transcend Micro is the smallest and lightest travel CPAP on the market, weighing less than 0.5 pounds and measuring less than 4 inches wide. Its ultra-portable design makes it a breeze to pack and carry, offering a discreet and unobtrusive solution for on-the-go therapy.

### Universal Mask Compatibility

Enjoy the flexibility of using your preferred CPAP mask. The Transcend Micro is designed to work seamlessly with any CPAP mask with a standard 22mm connection.

### Versatile Power Options

Whether you're at home, traveling abroad, or venturing off the grid, the Transcend Micro adapts to your environment. It comes equipped with an AC power adapter for traditional use and offers the Transcend Micro PowerAway™ Battery as an ideal companion for camping, flights, remote locations or use every day as a backup power source.

#### Item#

- CAPSOM507000** White
- CAPSOM507001** Pink
- CAPSOM507002** Gray
- CAPSOM507003** Blue

## What's Included

- Transcend Micro CPAP and filter
- Transcend Micro AirFlex™ 6-foot hose
- Transcend WhisperSoft™ Muffler
- AC power adapter
- 2 Year manufacturer's warranty



### Whisper-Quiet Operation

The Transcend Micro operates quietly, ensuring that both you and your companions can enjoy a peaceful night's sleep. When used with the WhisperSoft muffler (included), the device's noise level is further minimized, allowing for peaceful rest.

### Maintenance Made Simple

With user-friendly features and easy maintenance, the Transcend Micro is designed to simplify your sleep therapy routine.

### Health Benefits of CPAP Therapy

Consistent use of a sleep apnea device can lead to numerous health benefits, including lower blood pressure, reduced risk of stroke, decreased daytime sleepiness, and improved mental health. By maintaining open airways during sleep, CPAP therapy can also reduce snoring and enhance overall sleep quality.

### Warranty and Support

Your investment in the Transcend Micro is protected by a 2-year warranty and backed by comprehensive support.

# PowerAway® Battery

Item# **CAPSOM505001**

## Extended Run Time

- **Long-Lasting Power:** Offers 2-3 nights or 17.7 hours of run time on a single charge at a pressure setting of 9cm, ensuring restful sleep without frequent recharging.
- **Quick Charge Time:** Recharges in less than 4 hours, so it's ready to go when you are.

## Reliability and Safety

- **FAA-Approved:** Cleared for in-flight use, making it an essential for travelers
- **Backup Power:** Provides a permanent backup power source during outages, ensuring uninterrupted sleep therapy.
- **Long Life Span:** Capable of 1,000+ full charge/discharge cycles for long-term reliability.



## Essential Pack

Item# **CAPSOM505053**

### Includes:

- AirMist HME Cartridges Pack (4 filters + 1 holder)
- Micro WhisperSoft Muffler
- (1) 6ft AirFlex Hose
- (1) 4ft AirFlex Hose
- (2) Hose Clips
- Travel Power Supply Plug Adapters
- Standard Woven filters (2 pack)
- Hypoallergenic filters (2 pack)
- Transcend Micro Travel Pouch



## Essential Bundle

Item# **CAPSOM505055**

### Includes:

- Essential Pack
- White Micro 510
- SleepPak Travel Bag



## Power Bundle

Item# **CAPSOM505056**

### Includes:

- Essential Pack
- White Micro 510
- PowerAway Battery
- SleepPak Travel Bag



## Resupply Pack

Item# **CAPSOM505054**

### Includes:

- Micro WhisperSoft Muffler
- 6ft Airflex Hose
- AirMist HME Cartridges Pack (4 filters and 1 holder)
- Standard Woven filters (2 pack)



SCAN FOR MORE:  
VIDEOS AND SPECS



To order or for more information

P 312.997.9980 E sales@sunsethcs.com W SUNSETHCS.COM

