



# Brighter

Pediatrics by Sunset



**Contracted Supplier**



**vizient.**  
Awarded Supplier





Sales Representatives



**Egle Sventnickaite**  
872-258-1654  
esventnickaite@sunsethcs.com



**Mason Moss**  
262-788-1862  
mmoss@sunsethcs.com



**Owen Ruth**  
312-533-2451  
owenr@sunsethcs.com



**Cindy Geighes**  
312-533-2478  
cgeighes@sunsethcs.com




**Hunter Baudino**  
262-788-1855  
hbaudino@sunsethcs.com



**Wyatt Kanemoto**  
262-788-8941  
wkanemoto@sunsethcs.com

*Director, Strategic Accounts*  
**Lisa Iannazzone** RN, BSN  
636-222-3134  
lisa@sunsethcs.com

*Manager, Product & Business Development*  
**Tom Munar**  
312-533-2476  
tmunar@sunsethcs.com



*Contract Manager, Strategic Accounts*  
**Gina Rosenberg**  
914-837-5215  
grosenberg@sunsethcs.com

*RRT Clinical Specialist*  
**Krystle Bodine** BS, RRT  
kbodine@sunsethcs.com

**Corporate Office**

Phone 872-275-3210  
Fax 312-997-9985

141 W Jackson Blvd Ste 1950  
Chicago IL 60604

**Central Returns/Distribution Center**

2725 Alft Ln  
Elgin IL 60124

**West Returns/Distribution Center**

6750 W Van Buren St  
Phoenix AZ 85043

# Content

<b>NIV</b>	<b>Page</b>
NIV Masks	4
CPAP/NIV Masks	5
NIV Tubing and Connectors	5

<b>Oxygen Regulators</b>	<b>Page</b>
Microdial Flowmeter	6
Regulators	7

<b>Respiratory Disposables</b>	<b>Page</b>
Spirometers	8
Oxygen Cannula and Accessories	8
Oxygen Masks	9
Ventilator Filters	9
Air Entrainment Masks	10
Oxygen Supply Tubing	10
Breathing Circuits	11
Trach Ties	12
HMEs	12
Trach Masks	12
Adapters and Connectors	12
SpO2 Sensors	13

<b>Medication Delivery</b>	<b>Page</b>
Nebulizers	14-15
Nebulizer Kits	15
Aerosol Masks	15
Inhaler Chamber	15

<b>Closed Suction Catheters</b>	<b>Page</b>
Deluxe Switch	16
Standard Switch	17

<b>Open Suction</b>	<b>Page</b>
Suction Catheters	18
Trach Care Kits	18
Suction Catheter Kits and Trays	18-19

<b>Laryngeal Masks</b>	<b>Page</b>
Silicone Disposable	20
Silicone Reusable	20
Flexible Silicone Disposable	20
Flexible Silicone Reusable	20
PVC Disposable	21
Flexible PVC Disposable	21
Silicone Disposable Dual Seal	21
PVC Disposable Dual Seal	21

<b>Resuscitation Devices</b>	<b>Page</b>
Resuscitators	22
Breathing Masks	22

<b>Anesthesia</b>	<b>Page</b>
PeDIA	23



## CONTRACT

PP-NS-1741	<b>Respiratory Disposables</b>
PP-NS-1753	<b>Medication Delivery</b>
PP-NS-1762	<b>Resuscitation Devices</b>
PP-OR-2217	<b>Laryngeal Masks</b>
PP-FA-2138	<b>Medical Gas Pipeline Equipment</b>
BT-OR-019	<b>Disposable Anesthesia</b>



## CONTRACT

MS1087	<b>Active Humidification + Respiratory Disposables</b>
MS1097	<b>Medication Delivery</b>
MS1049	<b>Non-Invasive Ventilation</b>
MS7709	<b>PeDIA Anesthesia Induction</b>
MS1486	<b>Open Suction</b>

## SleepNet MiniMe 2 Pediatric Nasal Mask V

10/pack

### Vented Mask with Headgear

Size	Item#
Extra Small	<b>CM60216</b>
Small	<b>CM60217</b>
Medium	<b>CM60218</b>
Large	<b>CM60221</b>

### Non-Vented Mask with Headgear

Size	Item#
Extra Small	<b>CM60256</b>
Small	<b>CM60257</b>
Medium	<b>CM60258</b>
Large	<b>CM60261</b>

novaplus<sup>+</sup>

### Vented Mask with Headgear

Size	Item#
Extra Small	<b>VCM055XS</b>
Small	<b>VCM055S</b>
Medium	<b>VCM055M</b>
Large	<b>VCM055L</b>

### Non-Vented Mask with Headgear

Size	Item#
Extra Small	<b>VCM056XS</b>
Small	<b>VCM056S</b>
Medium	<b>VCM056M</b>
Large	<b>VCM056L</b>



HCPCS#  
MASK - A7034  
HEADGEAR - A7035

### Replacement Headgear

Size	Item#
Extra Small/Small	<b>HG027S</b>
Medium/Large	<b>HG027L</b>

## Ray Pediatric Masks Comparable to Nonny masks V

### Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	<b>CM061S</b>
Medium	<b>CM061M</b>
Large	<b>CM061L</b>
Extra Large	<b>CM061XL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM061S</b>
Medium	<b>CUCM061M</b>
Large	<b>CUCM061L</b>
Extra Large	<b>CUCM061XL</b>

### Headgear Only

Size	Item#
Small	<b>HG061S</b>
Large	<b>HG061L</b>



SCAN CODE to see  
Ray Full Face Mask  
on YouTube

### Nasal Mask with Headgear, Tube, and Replaceable Cushions

Comes with	Item#
S&M Cushions	<b>CM060-KITS</b>
L&XL Cushions	<b>CM060-KITL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM060S</b>
Medium	<b>CUCM060M</b>
Large	<b>CUCM060L</b>
Extra Large	<b>CUCM060XL</b>

### Headgear Only

Size	Item#
Small	<b>HG060S</b>
Large	<b>HG060L</b>





SCAN CODE to see  
Ray Nasal Mask  
on YouTube

### Non-Vented Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Large	<b>CM061L-NV</b>
Extra Large	<b>CM061XL-NV</b>

Large and Extra Large sizes help fill a critical pediatric sizing gap, with more potential to avoid intubation.



Available Size Guide

**FinTek allows more adjustability to conform the mask to the face**

**One-piece molded silicone mask is small, light and comfortable**

**Headgear with 6-point velcro provides maximum comfort and versatility**

**Quick release clips**

**Dual 360 degree swivel elbow**

- New non-vented configuration supports use in closed-circuit ventilation systems, including ICU settings
- Designed specifically for pediatric patients, reducing reliance on ill-fitting small adult masks
- Compatible with dual-limb circuits and adaptable to single-limb systems with an exhalation valve
- May support improved therapy delivery through better fit and system flexibility
- SKU reduction by reducing the need for multiple mask types per patient

## NIV Tubing V P

19mm with 22mm cuffs  
10/case\*

Item#

**TUB06** 6ft \*Individually Packed

**SUH0610** 6ft

**TUB010** 10ft

- Sunset's most popular item
- ISO standard 22mm cuffs
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free




## MiniMe Pediatric Nasal Mask



Headgear Size	Item#	OEM#
Small	<b>CM004S</b>	60213 <span>V</span>
Medium	<b>CM004M</b>	60214 <span>V</span> <span>P</span>
Large	<b>CM004L</b>	60215 <span>V</span>

HCPCS# A7034



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

## Pediatric SleepWeaver® CPAP Mask V

Item#

**CM101325**



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.

- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

## Standard Elbow, Non-Vented V



Item# **CAP009**  
10/case

## Single Limb Elbow, Non-Vented, AAV V



Item# **CAP008**  
10/case

## Swivel Connector Silencer

Item#

**VEN002** Each P

**VEN002-10** 10/pack V Contract MS1087



## Single Use Exhalation Port V

Item# **VEN003**  
10/pack



## Microflow Microdial Flowmeter P

2024 Innovation Celebration Winner



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure Range	Flow Range*
<b>RES8160048</b>	345 – 500 kPa	<b>A</b>
<b>RES8160050</b>	350 – 550 kPa	<b>E</b>



DISS Inlet • DISS Outlet

Item#	Inlet Pressure Range	Flow Range*
<b>RES8160049</b>	345 – 500 kPa	<b>A</b>
<b>RES8160051</b>	250 – 450 kPa	<b>E</b>

\*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
<b>Range A</b>	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
<b>Range E</b>	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

## Assured, repeatable accuracy to enable safe, effective weaning.

### What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.

### Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned



## Mini Oxygen Regulator CGA 870 Style

Item# **RES046**

Pediatric • Barbed

0-4 LPM • **Each**

0,  $1/32$ ,  $1/16$ ,  $1/8$ ,  $1/4$ ,  $1/2$ ,  $3/4$ , 1, 1  $1/2$ , 2, 3, 4 LPM

HCPCS# E1353

For use with E size or smaller cylinders: A, B, C, D, JD, E, M2, M4, M6, ML6, M7, M9, M15, M22, or M24.



Barbed hose outlet  
for easy connection to  
the oxygen tubing

"T" Handle for easy  
mounting — no tools  
needed

### Both Regulators Feature

- Pressure gauge reads to 3000 PSI
- Easy-click flow selector knob
- 2-year warranty

## Mini Oxygen Regulator CGA 540 Style

Item# **RES047**

Pediatric • DISS

0-4 LPM • **Each**

0,  $1/32$ ,  $1/16$ ,  $1/8$ ,  $1/4$ ,  $1/2$ ,  $3/4$ , 1, 1  $1/2$ , 2, 3, 4 LPM

HCPCS# E1353

For use with M60 and larger cylinders: M60, M90, M, MM, M122, M250, H, J, T or other cylinders with CGA 540 style valve.



Threaded  
DISS outlet

## Volumetric Incentive Spirometer P

12/case

Item#  
**RES055P**

Pediatric

- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise
- Dual-sided, clear and legible identification of volumes
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration
- Compact and ergonomic design



## Pediatric, Infant, and Neonate Cannula V

50/case

Item #	Size	Tubing Length	Connector	Touch	Ear Cushions
<b>RES1100P</b>	Pediatric	None	Standard	Standard	No
<b>RES1107I</b>	Infant	7 ft	Standard	Standard	No
<b>RES1107P</b>	Pediatric	7 ft	Standard	Standard	No
<b>RES1107PEC</b>	Pediatric	7 ft	Standard	Standard	Yes
<b>RES1107PECU</b>	Pediatric	7 ft	Universal	Standard	Yes
<b>RES1107SI</b>	Infant	7 ft	Standard	Soft	No
<b>RES1107SIU</b>	Infant	7 ft	Universal	Soft	No
<b>RES1107NE</b>	Neonate	7 ft	Standard	Standard	No
<b>RES1107SNE</b>	Neonate	7 ft	Standard	Soft	No
<b>RES1107SP</b>	Pediatric	7 ft	Standard	Soft	No
<b>RES1107SPU</b>	Pediatric	7 ft	Universal	Soft	No
<b>RES1114SP</b>	Pediatric	14 ft	Standard	Soft	No
<b>RES1114SPU</b>	Pediatric	14 ft	Universal	Soft	No



## Ear Cushions for Cannula V P

Item# **RES020**  
Cushions Only  
50 Pairs/case



- Add cushions to existing cannula



## Secure Cannula Stickers V P

Item# **RES080**  
25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



### High Concentration Oxygen Mask V

50/case

Item #	Size	Tubing Length
<b>RES2202I</b>	Infant	None
<b>RES2202I7</b>	Infant	7 ft
<b>RES2202P</b>	Pediatric	None
<b>RES2202P7</b>	Pediatric	7 ft



### Medium Concentration Oxygen Mask V

50/case

Item #	Size	Tubing Length	Connector
<b>RES2203P</b>	Pediatric	None	Standard
<b>RES2100P</b>	Pediatric	7 ft	Standard
<b>RES2203I</b>	Infant	None	Standard
<b>RES2203I7</b>	Infant	7 ft	Standard
<b>RES2203P7</b>	Pediatric	7 ft	Standard



- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

### Non-Rebreather Oxygen Mask V

50/case

Item #	Size	Tubing Length	Connector	Safety Vent
<b>RES2102I</b>	Infant	7 ft	Standard	No
<b>RES2102IM</b>	Infant	None	Standard	No
<b>RES2102P</b>	Pediatric	7 ft	Standard	No
<b>RES2102PS</b>	Pediatric	7 ft	Standard	Yes
<b>RES2102PSM</b>	Pediatric	None	Standard	Yes
<b>RES2102PSU</b>	Pediatric	7ft	Universal	Yes
<b>RES2102PU</b>	Pediatric	7 ft	Universal	No



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape

### Filter for ResMed Astral Ventilator V P

Item# **VF003-4**  
4/pack



Item# **VF003-H**  
Filter Housing  
Each



Item# **VF003-KIT**  
4-Filter Kit  
w/Filter Housing  
OEM# 27939



## Air Entrainment Mask

With Flex Tube

50/case

Item#	Kit	Tubing Length	Connector	
<b>RES2400P</b>	Mask Only	7 ft	Standard	V
<b>RES2400PK</b>	Yes	7 ft	Standard	V P
<b>RES2400PKU</b>	Yes	7 ft	Universal	V
<b>RES2400PU</b>	Mask Only	7 ft	Universal	V P



## Oxygen Supply Tubing V

Item #	Length	Connector	Bore	Color	Qty
<b>RES3007</b>	7ft	Standard	Channelled	None	50
<b>RES3007SB</b>	7 ft	Standard	Smooth	None	50
<b>RES3007U</b>	7 ft	Universal	Channelled	None	50
<b>RES3015</b>	15 ft	Standard	Channelled	None	50
<b>RES3015G</b>	15 ft	Standard	Channelled	Green	50
<b>RES3025</b>	25 ft	Standard	Channelled	None	25
<b>RES3025G</b>	25 ft	Standard	Channelled	Green	25
<b>RES3025U</b>	25 ft	Universal	Channelled	None	25
<b>RES3030</b>	30 ft	Standard	Channelled	None	20
<b>RES3030G</b>	30 ft	Standard	Channelled	Green	20
<b>RES3040G</b>	40 ft	Standard	Channelled	Green	20
<b>RES3050</b>	50 ft	Standard	Channelled	None	20
<b>RES3050G</b>	50 ft	Standard	Channelled	Green	20

- 6-channel
- PHT/DEHP-free





Universal Connector

### Non-Heated Breathing Circuits


For use with Respiroics Trilogy  
22mm cuffs  
**30/case**

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
<b>VC1001P</b>	1073220 & 1073229	Active	Yes	Pediatric	15mm
<b>VC1004P</b>	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#  
**VC1001P**  
Active Pediatric  
OEM#  
1073220  
1073229



To use active circuits with Respiroics Trilogy, simply cut ends off pressure and exhalation lines.



Item#  
**VC1004P**  
Passive Pediatric  
OEM#  
1073224  
1073233



### Passive Heated Breathing Circuit P V

Item#  
**VC1114P**  
For use with  
MR850 Humidifier  
**10/case**



### Passive Heated Breathing Circuit With Swivel Exhalation Port P V

Item# **VC1003P**  
For use with Philips  
Respiroics Trilogy  
**Each**

- Includes:**
- 1 VC1114P Passive Heated Breathing Circuit
  - 1 VEN002 Swivel Exhalation Port

OEM# 1076046,  
1073227, 1073235



### Non-Heated Dual Limb Breathing Circuits V

**30/case**

REF	Humidifier Limb	Water Trap
<b>VC2001P</b>	No	No
<b>VC2002P</b>	Yes	No




## Tracheostomy Ties with Fastener Tabs P

**20/case**

### One-piece

- Easy on

Item#  
**RES029P** Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton

### Two-piece

- Easy off

Item#  
**RES030P** Pediatric



## Heat Moisture Exchanger V

**30/case**

Item#  
**RES027NE-10**  
Inline HME — Neonatal



<b>TV Range (mL)</b>	15-50
<b>Moisture Output (mg H<sub>2</sub>O/L)</b>	30
<b>Weight (g)</b>	4.5
<b>Resistance (cm H<sub>2</sub>O)</b>	0.5 @ 10 LPM
<b>Dead Space (mL)</b>	2.5
<b>Connection (mm)</b>	15M,15F

## Tracheostomy Mask with O<sub>2</sub> Connector V

**50/case**

Item#  
**RES2230**  
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O<sub>2</sub> supply tubing

## Oxygen Tubing Connector — Swivel P V

Item# **RES018-1**  
Male/Male  
**1/pack**



## Flexible Connector P V

**50/case**

Item# **TRAA004P**  
15mm F x 15mm M  
Length 4.9-7.5cm  
OEM# 3215



## Oxygen Bleed-In Adapter P

Item# **RES019**  
22mm with 6mm  
Rotating Stem  
**10/pack**

**RES019-1** 1/pack P V



## Connector with 6mm Oxygen Stem V P

Item# **RES058**  
22mm with 6mm Stem  
**10/pack**




## SpO2 Sensors

NEW


### Disposable SpO2 Sensor Neonatal/Adult

Item# **RES7100**  
Similar to Nellcor MAX-N  
20/box




### Neo-M-1.5 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7200**  
1.5ft cord • Non-replaceable  
20/box



### Neo-M-3 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7201**  
3ft cord • Non-replaceable  
20/box




#### Coban Tape Roll

Item# **RES7202CT**  
Use for sensor life extension  
or sensor cord stability  
Each roll 1" x 177"

### Infant-M-1.5 Disposable SpO2 Sensor Infant (3-20 kg)

Item# **RES7203**  
1.5ft cord  
20/box



### Pediatric-M-1.5 Disposable SpO2 Sensor Pediatric (10-50 kg)

Item# **RES7204**  
1.5ft cord  
20/box



## Compressor Nebulizers

Item#

**NEB100**   
Each

### Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

### NEB100-R

- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

### Sunny the Bear

#### NEB300BEAR

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

#### NEB300BEAR-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



### Gordon the Fish

#### NEB300FISH

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

#### NEB300FISH-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



## Handheld Compressor Nebulizer

Item#

**NEB200**  
Each

### Included

Compressor  
Medication cup  
AC power adapter  
Air hose  
Elbow  
Mouthpiece  
Mask  
4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (approximately 3 microns)
- No replacement filter required



## Handheld Mesh Nebulizer

Item#

**NEB401**  
Each

### With internal rechargeable battery

### Included

Main unit  
Medicine cup  
Disposable mouthpiece  
Disposable adult mask  
Disposable pediatric mask  
USB cord  
Carrying bag  
Instruction manual

- Small particle size: MMAD 2.2µm\*
- Efficient: Low Residual Liquid
- Convenient: Small and Portable
- Quiet: ≤50dB (at one meter)
- Modern: Rechargeable Lithium Battery



## Handheld Mesh Nebulizer

Item#  
**NEB402**

With internal rechargeable battery

**Particle Size**  
MMAD < 5µm

**Nebulization Rate**  
0.25ml/min-0.9ml/min

**Medication Cup Capacity**  
10ml (Max)

**Device Weight**  
Approx. 1/4 lb

**Device Size**  
Approx. 1 1/2 x 2 x 3 7/8 inches



**Included**  
Medication cup  
USB cord  
Mask  
Mouthpiece  
Portable bag  
Manual

## Disposable Nebulizer Kit with T-Piece V P

Item#  
**RES091U**  
Disposable  
Single Patient,  
Universal Connector  
**50/case**



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

**Includes:**

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube

## Aerosol Mask V P

**50/case**

Item#  
**RES2210**  
Pediatric Size



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

## Disposable Nebulizer Kit V

Item# **RES092**  
Pediatric Size  
**50/case**

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube



## SperaChamber V P

Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler

Item#  
**RES195-10**  
**10/case**

**RES195-50**  
**50/case**

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



**SperaChamber with Mask**  
**10/case**

Item#	Size
<b>RES195MS</b>	Small (Infant)
<b>RES195MM</b>	Medium (Pediatric)



**Mask Only**

Item#	Size	Qty
<b>RES195MASK-S</b>	Small (Infant)	<b>Each</b>
<b>RES195MASK-S-10</b>	Small (Infant)	<b>10/box</b>
<b>RES195MASK-M</b>	Medium (Pediatric)	<b>Each</b>
<b>RES195MASK-M-10</b>	Medium (Pediatric)	<b>10/box</b>

## Closed Ventilation Suction System

# Deluxe Switch (One-handed operation)

Sterile and latex-free

### Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port

#### Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

#### MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

#### Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

#### Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.

#### Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.



## Closed Ventilation Suction System V

# Standard Switch

Sterile and latex-free

### 72 Hour Neonatal/Pediatric

#### With Swivel Connector

Item#	Description
VEN77230505	305mm, 5Fr
VEN77230506	305mm, 6Fr
VEN77230507	305mm, 7Fr
VEN77230508	305mm, 8Fr



### 24 Hour Neonatal/Pediatric

#### With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600	305mm, 6Fr
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700	305mm, 7Fr
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800	305mm, 8Fr
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000	405mm, 10Fr
VEN66001000M	405mm, 10Fr, MDI Port

#### With Swivel Connector

Item#	Description
VEN72430505	305mm, 5Fr
VEN72430506	305mm, 6Fr
VEN72430507	305mm, 7Fr
VEN72430508	305mm, 8Fr

#### Swivel Straight Connector V


Item#  
**VEN005-20**  
Disposable  
22mm ID, 15mm OD  
**20/pack**



## Suction Catheters - Calibrated

50/case

Item#	Size	Tip	Valve	Type	Sterile
<b>TRASC106C</b>	6fr	DeLee	Chimney	Coiled	No
<b>TRASC108C</b>	8fr	DeLee	Chimney	Coiled	No
<b>TRASC108S</b>	8fr	DeLee	Chimney	Straight	No
<b>TRASC110C</b>	10fr	Whistle	Y	Coiled	No
<b>TRASC110S</b>	10fr	Whistle	Y	Straight	Yes



## Trach Care Kits

50/case

### Sunset Trach Care Kit

Item# **TRATK003**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- Gauze
- 2 Cotton Tips
- Brush,
- Twill Tape
- Trach Dressing

### Sunset Mini Trach Kit

Item# **TRATK004**

Includes:

- 3 Pipe Cleaners
- 4 Gauze
- 2 Cotton Tips
- Brush
- Twill Tape
- Dressing

### Sunset Trach Care Kit with Forceps

Item# **TRATK005**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- 3 Gauze
- 3 Cotton Tips
- Brush
- Trach Dressing

## Suction Catheter Kits with Cup and Gloves

50/case

Item#	Catheter Size	Tip	Valve	Sterile
<b>TRATK206</b>	6fr	DeLee	Chimney	Yes
<b>TRATK208</b>	8fr	Whistle	Y	No
<b>TRATK210</b>	10fr	Whistle	Y	No



### Pediatric Suction Catheter Kits **V**

**50/case**

Item#

**TRATK306**

Kit with 6fr Catheter

- Sterile
- Y-Valve DeLee Catheter
- Includes Gloves

**TRATK308**

Kit with 8fr Catheter

**TRATK408**

Kit with 8fr Catheter and Cup

### Mini Catheter Tray with Gloves **V**

**50/case**

Item#

**TRATK310**

Kit with 10fr Catheter

- Coiled Catheter
- Calibrated


### Sleeved Suction Coiled Catheter Kits with Cup **V**

**50/case**

Item#

**TRATK410**

Kit with 10fr Catheter



## Silicone Disposable Laryngeal Masks P

Item #	Size
<b>ALM1001</b>	#1
<b>ALM1015</b>	#1.5
<b>ALM1002</b>	#2
<b>ALM1025</b>	#2.5
<b>ALM1003</b>	#3

- Medical grade silicone for super compatibility and patient comfort
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters
- Available with seven sizes for all patient use



## Silicone Reusable Laryngeal Masks P

Item #	Size
<b>ALM1001R</b>	#1
<b>ALM1015R</b>	#1.5
<b>ALM1002R</b>	#2
<b>ALM1025R</b>	#2.5
<b>ALM1003R</b>	#3

- Reusable Silicone Laryngeal Masks are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Serial no. is printed on each connector for easy tracking
- Additional user card for record of device usage times and patient information is supplied in each pouch



## Flexible Silicone Disposable Laryngeal Masks P

Item #	Size
<b>ALM3001</b>	#1
<b>ALM3015</b>	#1.5
<b>ALM3002</b>	#2
<b>ALM3025</b>	#2.5
<b>ALM3003</b>	#3



## Flexible Silicone Reusable Laryngeal Masks P

Item #	Size
<b>ALM3001R</b>	#1
<b>ALM3015R</b>	#1.5
<b>ALM3002R</b>	#2
<b>ALM3025R</b>	#2.5
<b>ALM3003R</b>	#3




- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness

Size	Patient Weight
#1	Up to 5 kg
#1.5	5-10 kg
#2	10-20 kg
#2.5	20-30 kg
#3	30-50 kg

### PVC Disposable Laryngeal Masks **P**

Item #	Size
<b>ALM2001</b>	#1
<b>ALM2015</b>	#1.5
<b>ALM2002</b>	#2
<b>ALM2025</b>	#2.5
<b>ALM2003</b>	#3

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



### Flexible PVC Disposable Laryngeal Masks **P**

Item #	Size
<b>ALM4001</b>	#1
<b>ALM4015</b>	#1.5
<b>ALM4002</b>	#2
<b>ALM4025</b>	#2.5
<b>ALM4003</b>	#3



### Silicone Disposable Dual Seal Laryngeal Masks **P**

Item #	Size
<b>ALM5015</b>	#1.5
<b>ALM5002</b>	#2
<b>ALM5025</b>	#2.5
<b>ALM5003</b>	#3



### PVC Disposable Dual Seal Laryngeal Masks **P**

Item #	Size
<b>ALM6003</b>	#3



Size	Patient Weight
#1	Up to 5 kg
#1.5	5-10 kg
#2	10-20 kg
#2.5	20-30 kg
#3	30-50 kg

## Child Resuscitator P

12/case

Item #	Pop-Off	Mask
<b>RES8110H</b>	Yes	None
<b>RES8111</b>	Yes	#2
<b>RES8112</b>	Yes	#3



## Infant Resuscitator P

12/case

Item #	Pop-Off	Mask
<b>RES8120H</b>	Yes	None
<b>RES8121</b>	Yes	#0
<b>RES8122</b>	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve

## Breathing Mask P

20/case

Item #	Size	Hooks
<b>SBC01N</b>	Neonatal (#0)	No
<b>SBC01I</b>	Infant (#1)	No
<b>SBC01P</b>	Pediatric (#2)	No
<b>SBC02N</b>	Neonatal (#0)	Yes
<b>SBC02I</b>	Infant (#1)	Yes
<b>SBC02P</b>	Pediatric (#2)	Yes



# PeDIA - Pediatric Device for Induction of Anesthesia V P

Item#

**ANEP021Y**

With Mask Adapter Elbow

**12/pack**



*2025 Innovation Celebration Winner*

Use with mask



Or use without mask



## Turn panic into play

For children, undergoing surgery is already unsettling, and the process of induction is often the most traumatizing part. Often kids will kick and scream, believing they're being trapped by the gas or suffocated by the mask. The experiences create negative emotional scarring that can last weeks, months or even years.

Specially designed for children undergoing surgery, the Pediatric Device for Induction of Anesthesia device, or "PeDIA," for short, includes a colorful whistling balloon, providing kids with something fun to focus on and turning a stressful moment into a lighter one.

Whether they're squeezing the balloon and focusing on its whistle or directly putting themselves to sleep by breathing the gas through balloon delivery, the PeDIA balloon captivates children's attention, allowing their active participation in the induction process.

It's a simple way to add a brighter touch to the medical experience, ease anxiety and bring smiles to patients, parents and staff.

This unique approach—presenting an interactive, familiar toy-like device—not only diverts children's focus from the impending procedure but also significantly reduces the anxiety and resistance commonly associated with mask induction. The aim is to decrease the anxiety felt by over 50% of pediatrics as well as parents and staff.

The PeDIA can be used alone as a mask-free device or, for younger children, with an anesthesia mask.

### SCAN FOR VIDEOS



