

# Ray Pediatric Non-Vented NIV Mask

- FinTek allows more adjustability to conform the mask to the face
- One-piece molded silicone mask is small, light and comfortable
- Headgear with 6-point velcro provides maximum comfort and versatility
- Quick release clips
- Dual 360 degree swivel elbow

Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Large	<b>CM061L-NV</b>
Extra Large	<b>CM061XL-NV</b>

Replacement Cushions


Size	Item#
Large	<b>CUCM061L</b>
Extra Large	<b>CUCM061XL</b>

Replacement Headgear

Size	Item#
Small	<b>HG061S</b>
Large	<b>HG061L</b>

**Large and Extra Large sizes help fill a critical pediatric sizing gap, with more potential to avoid intubation.**

- New non-vented configuration supports use in closed-circuit ventilation systems, including ICU settings
- Designed specifically for pediatric patients, reducing reliance on ill-fitting small adult masks
- Compatible with dual-limb circuits and adaptable to single-limb systems with an exhalation valve
- May support improved therapy delivery through better fit and system flexibility
- SKU reduction by reducing the need for multiple mask types per patient



Available size guide