

Active Heated Ventilator Circuits



Active Heated Circuit

Includes VEN004 Active Exhalation Valve Adapter

VC1002 Adult

VC1002P Pediatric

Designed to provide heated active alternatives to the Respiration circuits 1073222 & 1073231 & 1073221 & 1073230. VC1002 and VC1002P can be used on a wide variety of units (excluding VOCSN and Resmed Astral).



VC1002 and VC1002P are compatible with a wide range of ventilator units.




Scan or click for video

VEN004

Active Exhalation Valve Adapter

Designed for seamless integration with heated circuits from other manufacturers to create an active circuit, can be used with MR810, HC550, and other humidifiers.




Active Heated Circuit for Resmed Units Only

Includes VEN004-R Active Exhalation Valve Adapter

VC1002-R Adult

VC1002P-R Pediatric

Designed to provide heated active alternative to the Respiration circuits 1073222 & 1073231 & 1073221 & 1073230 for use exclusively with Resmed Astral with the single limb valve adapter OEM 27941.



VEN004-R

Active Exhalation Valve Adapter for Resmed

For use with Resmed units only. Can be combined with heated circuits from other manufacturers to create an active circuit, to be used with MR810, HC550, and other humidifiers.

