



# Homecare Product Catalog



# Policies

## Ordering Options

---

### Two Easy Ways To Order

#### Phone

(877) 5-SUNSET (877-578-6738)

Our customer service professionals are here to assist you  
M-F 8:00 AM – 5:00 PM Central Time.

#### Email

sales@sunsethcs.com

*Contact us for more information  
on fulfillment and compliance  
solutions that best fit your needs.*

## Shipping

---

Ask your rep about shipping to various locations and our drop shipping partners.

#### UPS

Whenever possible, orders ship via UPS. If another carrier is desired, please contact your sales rep.

All orders are packed in the smallest boxes available and combined with other products on the order to reduce shipping costs.

All products ship prepaid and are added to the invoice.

All products ship from Elgin, IL or Phoenix, AZ. Please refer to UPS.com for transit times.

#### Damaged Merchandise

Call us within 10 days. Keep original cartons for UPS inspection.

*We proudly partner with*



**MCKESSON**



## Payment Methods

---

#### Payment Portal and ACH

Visit [www.sunsethcs.com/paymentportal](http://www.sunsethcs.com/paymentportal) to pay online with ACH or credit card.

#### Net 30-Day Payment Terms

For qualified companies. A 1.5% monthly finance charge will be assessed on all past due invoices.

#### COD

For non-qualified companies.

#### Credit Card

We accept all major credit cards. A processing fee may apply. After 10 days of invoice date, a surcharge of 2.5% will be applied.

**VISA**



**DISCOVER**

## Customer Service

---

#### Returns

Contact us within 7 days of receiving your shipment.

Call us for a return authorization number. There is no restocking charge on orders returned within 20 days of invoice. Orders returned after 20 days may be assessed a restocking fee of 15%. No returns after 90 days.

(877) 5-SUNSET (877-578-6738)

#### Central Distribution Center

Sunset Healthcare Solutions  
2725 Alft Ln  
Elgin, IL 60124

#### Western Distribution Center

Sunset Healthcare Solutions  
6750 W Van Buren St Ste 100  
Phoenix, AZ 85043

#### Problems

Please call us immediately if you have any complaints (it is helpful to ask for the customer service representative with whom you placed your order).

Or email us at: [customerservice@sunsethcs.com](mailto:customerservice@sunsethcs.com)

#### Questions?

Our experienced customer service staff will be glad to provide you with information. Just give us a call.

#### Pricing

All prices are in U.S. Dollars.

Taxes are added where applicable.

We make every effort to hold pricing throughout the year.

However, prices are subject to change without notice.

# Content

## CPAP

	Page
Masks and Replacement Parts	4-14
Sunset Healthcare Solutions	4-7
Sleepnet	8-10
Circadiance	11-12
Airway Management	13
Bleep	14
Sleep Apnea Appliances	15-16
Tubing	17-18
Filters	19-27
Respironics	19-20
ResMed	21-22
Fisher & Paykel	23
React Health/3B	23
AEIOMed	24
Beyond	24
Breas/GE	24
CareFusion/Viasys	25
Covidien Nellcor Puritan Bennett	25
DeVilbiss	26
Evo Medical	26
Healthdyne	27
Invacare	27
Löwenstein	27
Resvent	27
ProBasics	27
Humidifier Chambers	28
Chinstraps	29-30
Headgear	30
Cash Items	31-39
Cleaners and Accessories	31-32
Battery/Power Inverter	33
Transcend Micro™ Auto CPAP	34
Mask Liners, Soft Headgear Covers	35
Aromatherapy	36
Connectors, Adapters, Inline Filter	36
Pillows	37-38
Tubing Wraps	38
Storage	39

## Oxygen

	Page
Pulse Oximeter	40
SpO2 Sensors	41
Firesafe	42
Portable Concentrator	43
Concentrators	44
Masks	45
Cannula and Cannula Accessories	46-47
Tubing	47
Microdial Flowmeter	48
Regulators and Conserver	49
Oxygen HCPCS Chart	50
Cylinder Carriers	51-52
Connectors	53
Accessories	53-55
Cleaners	55
Analyzers	56
Celon Wet Bands	56
Bacteria and HEPA Filters	57
Concentrator Cabinet Filters	58-64
Caire/AirSep/SeQual	58
CHAD	59
Covidien Nellcor Puritan Bennett	59
Dynarex	59
OxyGo	59
DeVilbiss/Drive	60
Invacare	61
Rhythm, Jumao, Attendant, Lifestyle, Medacure	62
Respironics	63
Healthdyne	64
React Health/3B	64
Nidek	64
Repair Parts	65-69

## Suction

	Page
Suction Machine	70
Suction Canister and Kits	70-71
Open Suction Catheters	72
Open Suction Catheter Kits	73
Suction Filters	74
Suction Adapters and Connectors	74
Yankauers	74
Tracheostomy Supplies	75
Heat and Moisture Exchangers	76
Closed Suction Catheters	77-78

## Ventilator

	Page
LUISA Ventilator	79
Ventilator Filters	80-81
Adapters and Connectors	81
Ventilator Mask	81
Breathing Circuits	82-85
Heated Single-Limb	82
Non-Heated Single-Limb	83-84
Dual-Limb	85

## Other

	Page
Disinfectants	86
Cough Assist Kits	87
Breathing Masks	87
Aerosol	88
Tabletop Nebulizers	89
Handheld Nebulizers	90
Nebulizer Kits	91
Nebulizer and Anesthesia Filters	92
StrongArm Mobility and Comfort	93-94
Bath Safety	95
Maternity	96-97
Equipment Covers	98-99
Blood Testing	99

## Sunset Healthcare Solutions

### Vesper CPAP Masks

**NEW**

#### Full Face Mask with Headgear and Cushion

Size	Item#
Small	<b>CM206S</b>
Medium	<b>CM206M</b>
Large	<b>CM206L</b>

**CM206-FP**  
**Full Face Fit Pack**  
 (Mask and S, M, L cushions)

#### Replacement Cushions

Size	Item#
Small	<b>CUCM206S</b>
Medium	<b>CUCM206M</b>
Large	<b>CUCM206L</b>

#### CM206-CLIPS

Replacement Headgear Clips

#### IFU-00133

Sizing Guide

#### Features

**Soft and skin-friendly headgear fabric** helps patients comfortably adapt to treatment

**Quiet honeycomb vent** reduces airflow erosion intensity for peaceful sleep

**Easy pop cushions** attach or detach from mask in an easy pop; no need to painstakingly align pieces

**Quick release connector** rapidly and easily separates or connects mask and the tube



HCPCS#  
 MASK - A7030  
 HEADGEAR - A7035  
 CUSHION - A7031

#### Nasal Mask with Headgear and Cushion

Size	Item#
Small	<b>CM207S</b>
Medium	<b>CM207M</b>
Large	<b>CM207L</b>

**CM207-FP**  
**Nasal Fit Pack**  
 (Mask and S, M, L cushions)

#### Replacement Cushions

Size	Item#
Small	<b>CUCM207S</b>
Medium	<b>CUCM207M</b>
Large	<b>CUCM207L</b>

#### CM207-CLIPS

Replacement Headgear Clips

#### IFU-00136


Sizing Guide



HCPCS#  
 MASK - A7034  
 HEADGEAR - A7035  
 CUSHION - A7032



Available in  
**SleepGlad**



*Sunset's Vesper Full Face and Nasal Masks present the latest in CPAP mask design to provide patients with fresh options for maximum comfort.*

Vesper masks introduce a pared back fit, with up-to-date design details like the ergonomic stabilizing arms that curve up to the temples for patients who prefer to forego a forehead piece.

Vesper masks provide more visibility, a feeling of openness and a lighter fit. The flexible interface moves with side or stomach sleepers. The open structure eliminates forehead marks. Vesper masks also feature a 360-degree swivel elbow to prevent awkward tube tangles.

Vesper masks bring balance and calm to the nighttime, softening the edges of sleep.

Vesper comes in attractive retail packaging



SCAN FOR VESPER  
 FULL FACE VIDEO



SCAN FOR VESPER  
 NASAL VIDEO



# Sunset Healthcare Solutions

## Sol Full Face CPAP Mask

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM205S</b>
Medium	<b>CM205M</b>
Large	<b>CM205L</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM205S</b>
Medium	<b>CUCM205M</b>
Large	<b>CUCM205L</b>

**HG205** Replacement Headgear - Fits Most  
**CM205-CLIPS** Replacement Headgear Clips  
**IFU-00057** Sizing Guide

### Sol without Headgear

Size	Item#
Small	<b>CM205SNH</b>
Medium	<b>CM205MNH</b>
Large	<b>CM205LNH</b>

HCPCS#  
 MASK - A7030  
 HEADGEAR - A7035  
 CUSHION - A7031



*A premium mask without the premium price*

### Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

### Quiet vent

Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

### Quick release connector

Rapidly and easily separates or connects mask and the tube


### Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment

SCAN TO WATCH SOL VIDEO



Easy Pop Cushions



Available to add to your formulary in SleepGlad

Sol comes in attractive retail packaging



Ask us about our 30-day mask guarantee

## Ray Pediatric CPAP/NIV Masks Comparable to Nonny masks

### Full Face Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#	OEM#
Small	<b>CM061S</b>	AG-PEDFF-S
Medium	<b>CM061M</b>	AG-PEDFF-M
Large	<b>CM061L</b>	AG-PEDFF-L
Extra Large	<b>CM061XL</b>	AG-PEDFF-XL

### Replacement Cushions

Size	Item#	OEM#
Small	<b>CUCM061S</b>	AG-PEDFF-CS
Medium	<b>CUCM061M</b>	AG-PEDFF-CM
Large	<b>CUCM061L</b>	AG-PEDFF-CL
Extra Large	<b>CUCM061XL</b>	AG-PEDFF-CXL

### Replacement Headgear

Size	Item#	OEM#
Sm/Med	<b>HG061S</b>	AG-PEDFF-HGS
Lg/XL	<b>HG061L</b>	AG-PEDFF-HGL

HCPCS#  
 MASK - A7030 • HEADGEAR - A7035 • CUSHION - A7031



### Nasal Mask with Headgear, Tube, and Replaceable Cushions

Size	Item#	OEM#
S&M Cushions	<b>CM060-KITS</b>	AG-PEDKIT-S
L&XL Cushions	<b>CM060-KITL</b>	AG-PEDKIT-L

### Replacement Cushions

Size	Item#	OEM#
Small	<b>CUCM060S</b>	AG-PEDKIT-CS
Medium	<b>CUCM060M</b>	AG-PEDKIT-CM
Large	<b>CUCM060L</b>	AG-PEDKIT-CL
Extra Large	<b>CUCM060XL</b>	AG-PEDKIT-CXL

### Replacement Headgear

Size	Item#	OEM#
Small	<b>HG060S</b>	AG-PEDKIT-HGS
Large	<b>HG060L</b>	AG-PEDKIT-HGL

HCPCS#  
 MASK - A7034 • HEADGEAR - A7035 • CUSHION - A7032



## Sunset Healthcare Solutions

### Sunset Deluxe Full Face CPAP Mask

#### Mask with Flex Tube, Headgear, and Replaceable Cushion

Size	Item#
Small	<b>CM005S</b>
Medium	<b>CM005M</b>
Large	<b>CM005L</b>

**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

**CMSIZE01**  
Deluxe Sizing Guide

HCPCS#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031



Unobtrusive design offers a clear line of sight

- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone replaceable cushion increases comfort while maintaining optimal seal

#### Replacement Cushions

Size	Item#
Small	<b>CUCM005S</b>
Medium	<b>CUCM005M</b>
Large	<b>CUCM005L</b>

#### Mask without Headgear

Size	Item#
Small	<b>CM005SNH</b>
Medium	<b>CM005MNH</b>
Large	<b>CM005LNH</b>

### Sunset Deluxe Nasal CPAP Masks

#### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM006S</b>
Medium	<b>CM006M</b>
Large	<b>CM006L</b>

**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

**CMSIZE01**  
Deluxe Sizing Guide

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
CUSHION - A7032

#### Mask without Headgear

Size	Item#
Small	<b>CM006SNH</b>
Medium	<b>CM006MNH</b>
Large	<b>CM006LNH</b>

Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.



#### Mask with Headgear, Adjustable Forehead T, and Replaceable Cushion

Size	Item#
Small	<b>CM106S</b>
Medium	<b>CM106M</b>
Large	<b>CM106L</b>



#### Replacement Cushions

Size	Item#
Small	<b>CUCM006S</b>
Medium	<b>CUCM006M</b>
Large	<b>CUCM006L</b>

# Sunset Healthcare Solutions

## Sunset Deluxe Nasal Pillow CPAP Mask

Mask with Headgear and Small, Medium, and Large Pillows

Item# **CM036**

Headgear Only  
**HG036**

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
PILLOWS - A7033

### Replacement Pillows

Size	Item#
Small	<b>CM036SP</b>
Medium	<b>CM036MP</b>
Large	<b>CM036LP</b>



- Easy to fit, assemble, and maintain
- Three different pillow sizes to fit more people more easily
- 360° rotational elbow with flexible silicone tube
- Quiet vent design

Mask without Headgear  
**CM036NH**

## Sunset Clearsight Deluxe Nasal Mask

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM110S</b>
Medium	<b>CM110M</b>
Large	<b>CM110L</b>

**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

### Replacement Cushions

Size	Item#
Small	<b>CUCM110S</b>
Medium	<b>CUCM110M</b>
Large	<b>CUCM110L</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035



*Clearsight design provides unobstructed view*

Mask Series	Mask Type	Replaceable Cushion	Flex Tube Included	Headgear Included
<b>CM005</b>	Full Face	✓	✓	✓
<b>CM006</b>	Nasal - Forehead Support	✓	✓	✓
<b>CM106</b>	Nasal - Adjustable	✓		✓
<b>CM007</b>	Full Face			✓
<b>CM008</b>	Nasal			✓
<b>CM110</b>	Nasal	✓		✓
<b>CM036</b>	Nasal Pillow	✓	✓	✓

## Sunset Classic CPAP Masks

Full Face Mask with Headgear

Size	Item#
Small	<b>CM007S</b>
Medium	<b>CM007M</b>
Large	<b>CM007L</b>
Extra Large	<b>CM007XL</b>

**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

HCPCS#  
MASK - A7030  
HEADGEAR - A7035



- Conventional forehead pad removed to create a clear line of sight
- Minimalist design and ultra-lightweight silicone
- Quick-release clips make it easy to remove
- 360° dual-swivel elbow port allows for free movement
- One-piece mask and cushion design is durable and easy

Nasal Mask with Headgear

Size	Item#
Medium	<b>CM008M</b>
Large	<b>CM008L</b>

**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

HCPCS#  
MASK - A7034  
HEADGEAR - A7035



- Great value for cash patients or those with lower insurance reimbursement
- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

# Sleepnet

## Phantom®2 and iQ®2

### Phantom®2

Item#  
**CM050**

Replacement  
Cushion  
**CUCM050**



### iQ®2

Item#  
**CM049**

Replacement  
Cushion  
**CUCM049**



Replacement Frame/Tubing Assembly  
**CM051**

*Ask us about  
our 30-day  
mask guarantee*

- Sleepnet's patented AIR°gel® technology provides maximum comfort while minimizing skin irritation and pressure points on the nasal bridge. Unlike many other gel masks, the unique qualities of AIR°gel allow the mask to seal with minimal tension. By design, AIR°gel is twice as soft as competing gel products.
- Sleepnet has the only mask that you can shape to each individual face. Flex, Fit, Forget!

- Cushions may be removed and replaced without purchasing a new mask.
- Magnetic Clips — Easy on and off!
- 360° Swivel — Less tube drag for active sleepers.
- One size cushion — Larger fitting range.



# Sleepnet

## Innova

### Nasal Mask with Headgear

Size	Item#
Small/Medium	<b>CM034SM</b>
Medium/Large	<b>CM034ML</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035



### Full Face Mask with Headgear

Size	Item#
Small	<b>CM035S</b>
Medium	<b>CM035M</b>
Large	<b>CM035L</b>

HCPCS#  
MASK - A7030  
HEADGEAR - A7035



Innova® CPAP masks feature patented AIR<sup>o</sup>gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.


- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA

## MiniMe Pediatric Nasal Mask



### Mask with Headgear

Headgear Size	Item#	OEM#
Small	<b>CM004S</b>	60213
Medium	<b>CM004M</b>	60214
Large	<b>CM004L</b>	60215



### Replacement Headgear

Size	Item#	OEM#
Small	<b>HG010S</b>	55241
Medium	<b>HG010M</b>	55242
Large	<b>HG010L</b>	55243

HCPCS#  
MASK - A7034  
HEADGEAR - A7035

- Non-vented version available

MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

## MiniMe 2 Pediatric Nasal Mask

### Mask with Headgear

Size	Item#	OEM#
Extra Small	<b>CM027XS</b>	60216
Small	<b>CM027S</b>	60217
Medium	<b>CM027M</b>	60218
Large	<b>CM027L</b>	60221

### Replacement Headgear

Size	Item#
Extra Small/Small	<b>HG027S</b>
Medium/Large	<b>HG027L</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035

- Non-vented version available



## Sleepnet

### Ascend

#### Nasal Mask with Headgear and All Cushions

Item# **CM037N**


#### Replacement Cushions

Size	Item#
Small	<b>CM037NCS</b>
Medium	<b>CM037NCM</b>
Large	<b>CM037NCL</b>

#### Nasal Mask with Headgear and 1 Cushion

Size	Item#
Small	<b>CM037NS</b>
Medium	<b>CM037NM</b>
Large	<b>CM037NL</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
CUSHION - A7032



- One frame with three sizes of replaceable cushions: S, M, L
- Unique flexible spacebar design provides superior stability, holding a seal without contacting the skin
- Sleepnet's patented AIR<sup>o</sup>gel minimizes pressure points and skin irritation while offering maximum comfort with minimum weight
- Lightweight, breathable quick-release headgear slides off easily with one hand
- 360-degree rotating elbow for optimal tube positioning and easier mobility
- Active Headgear Connector moves as the patient moves for a reliable seal
- Made in the USA

#### Full Face Mask with Headgear and All Cushions

Item# **CM037F**

#### Replacement Cushions

Size	Item#
Small	<b>CM037CS</b>
Medium	<b>CM037CM</b>
Large	<b>CM037CL</b>



#### Full Face Mask with Headgear and 1 Cushion

Size	Item#
Small	<b>CM037S</b>
Medium	<b>CM037M</b>
Large	<b>CM037L</b>

HCPCS#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031



### iQ Nasal Mask

Item#  
**CM001**

iQ Nasal Mask with Headgear

OEM# TMS-889

HCPCS#  
MASK - A7034  
HEADGEAR - A7035



### Mojo Full Face Mask

#### Mask with Headgear

Size	Item#	OEM#
Small	<b>CM003SH</b>	TMS-835
Medium	<b>CM003MH</b>	TMS-836
Large	<b>CM003LH</b>	TMS-837

#### Replacement Headgear

**HG003**  
OEM# TMS-89

HCPCS#  
MASK - A7030  
HEADGEAR - A7035



### Mojo 2 Full Face Mask

#### Mask with Headgear

Size	Item#	OEM#
Small	<b>CM50817</b>	50817
Medium	<b>CM50818</b>	50818
Large	<b>CM50819</b>	50819
Extra Large	<b>CM50820</b>	50820

- Maximum comfort, minimum weight. Patented AIR<sup>o</sup>gel minimizes pressure points and skin irritation
- Forehead adjustment screw removes pressure from the bridge of the nose for personalized fit
- The shapeable shell of each mask can be easily molded to fit individual needs: bend it, mold it, shape it!
- Lightweight, breathable headgear disconnects from mask easily with one hand thanks to the quick-release magnetic headgear clips
- Oxygen port allows easy delivery of supplemental oxygen
- 360° rotating elbow for optimal tube positioning and easier mobility



# Circadiance

## SleepWeaver® Advance Nasal Mask



- One size fits most (Small Tan available for smaller faces)
- Includes Advance headgear that can be adjusted to fit a broader range of patients



### Mask with Adjustable Headgear

Mask Color	Item#	OEM#
Blue	<b>CM017B</b>	100274
Camo	<b>CM017C</b>	100281
Leopard	<b>CM017L</b>	100285
Pink	<b>CM017P</b>	100277
Tan	<b>CM017T</b>	100332

**CM017T-S**  
Small Tan  
OEM# 101224

HCPCS#  
MASK - A7034  
HEADGEAR - A7035

**HG017**  
Headgear Only  
OEM# 100435

**CM017TNH**  
Tan — No Headgear  
OEM# 100336

- Easy to store and pack for travel

## SleepWeaver® Élan Nasal Mask

### Mask Kit

Mask Color	Item#	OEM#
Tan	<b>CM022T-KIT</b>	100563
Blue	<b>CM022B-KIT</b>	100782

### Kits Include:

- Mask
- Small Cushion
- Regular Cushion
- Large Cushion
- Headgear
- Feather Weight Tube®
- Overhead Tube Tether Strap
- Quick Start Illustration Guide
- Instructional DVD

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
TUBE - A7037

### Replacement Cushions

Size and Color	Item#	OEM#
Small Blue	<b>CUCM022B-S</b>	100743
Regular Blue	<b>CUCM022B-R</b>	100654
Large Blue	<b>CUCM022B-L</b>	100664
Small Tan	<b>CUCM022T-S</b>	100742
Regular Tan	<b>CUCM022T-R</b>	100653
Large Tan	<b>CUCM022T-L</b>	100663
Small Leopard	<b>CUCM022L-S</b>	100746
Regular Leopard	<b>CUCM022L-R</b>	100657
Large Leopard	<b>CUCM022L-L</b>	100667

**TUB0018F**  
Replacement Feather Weight Tube®  
OEM# 100477

**HG022**  
Headgear Only  
OEM# 100630

### Mask with Headgear only

Size and Color	Item#	OEM#
Small Blue	<b>CM022B-S</b>	100728
Regular Blue	<b>CM022B-R</b>	100519
Large Blue	<b>CM022B-L</b>	100482
Small Tan	<b>CM022T-S</b>	100727
Regular Tan	<b>CM022T-R</b>	100518
Large Tan	<b>CM022T-L</b>	100481



## All SleepWeaver Masks

- Made of soft cloth
- Allow sleeping on side, stomach, and even wearing glasses
- Easy to store and pack for travel
- Hand or machine washable gentle cycle, mild detergent, hang to dry



- Detachable 90° elbow with 360° rotation allows for over-the-head hose placement

# Circadiance

## SleepWeaver® Anew™ Full Face Mask

### Blue Mask with Headgear

Size	Item#	OEM#
Small	<b>CM026B-S</b>	100955
Regular	<b>CM026B-R</b>	100960
Large	<b>CM026B-L</b>	100965

### Headgear

Size	Item#	OEM#
Small	<b>HGCM026-S</b>	101145
Reg/Lg	<b>HGCM026-R</b>	101084

### TUB0018F

Feather Weight Tube®  
OEM# 100477

HCPCS#  
MASK - A7027  
HEADGEAR - A7035  
TUBE - A7037



- Soft material wicks away moisture and leaves no pressure marks on the face
- Loose fit provides ultimate comfort and eliminates strap marks
- Open line of sight accommodates eyeglasses and is ideal for people who like to fall asleep while reading or watching TV
- 90° elbow swivel connector rotates 360° allowing users to sleep in more natural positions; easy to remove and clean

## Pediatric SleepWeaver® CPAP Mask

Item#  
**CM101325**



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.

- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

## SleepWeaver® 3D®

Item#  
**CM027**

### Blue Mask without Headgear

100/packs

Headgear	Item#
With Large Headgear	<b>CM027-LB</b>
With Regular Headgear	<b>CM027-RB</b>

Headgear	Item#
With Large Headgear	<b>CM027-L-100</b>
With Regular Headgear	<b>CM027-R-100</b>

- New wing-shaped stay routes the headgear under the ears, ensuring comfort and a great fit
- Permanent 360° swivel provides reliable performance in any position, including over the headboard
- Breathable cloth design eliminates fogging that is common in hard plastic masks - nothing breathes like cloth, keeping skin dry and comfortable
- For those who want to read or watch TV while in bed, SleepWeaver 3D can be worn with glasses
- Users prefer a soft cloth against the face over rubber, silicone and other gels. Cloth masks eliminate the need for a cushion or nasal pad



## Airway Management

### Tap Pap

#### No headgear needed!

Item#

**CM024NH**

Mask Kit with Small, Medium, and Large Pillows

OEM# PAP-NP1-000

#### Replacement Pillows

Size	Item#	OEM#
Small	<b>CUCM024S</b>	PAP-NP1-101
Medium	<b>CUCM024M</b>	PAP-NP1-102
Large	<b>CUCM024L</b>	PAP-NP1-103

**CM024MP** Replacement Mouthpiece

OEM# PAP-NP1-201

HCPCS#

MASK - A7034

PILLOWS - A7033

HEADGEAR - A7035



- Mouthpiece molds to the upper teeth to provide a custom and secure fit even if the bottom teeth are not in contact with the mouthpiece
- Unique nut and bolt joint allows the nasal pillows to be set to the desired angle
- Provides a stable fit regardless of the sleep position
- Patients report not being able to feel the mask on their face
- 3 sizes of nasal pillows included (Small, Medium, and Large)
- 90-day compliance guarantee for new users of the mask, 30-day guarantee for experienced users

## Bleep

### Bleep Eclipse Starter Kit

#### Eclipse Frame and Two 32/Packs of Halos

Item#  
**CAP100589**

The BLEEP Eclipse™ Starter Kit includes the headgear-free Eclipse™ mask frame and two 32-packs of Halos™ (32-day total supply) for immediate use. This clinically-proven NO LEAKS PAP solution is the only one of its kind on the market. No headgear, no mask lines and no broken hair. The comfortable Eclipse mask frame includes a diffuser and flexible tube and connects to Halos in lieu of headgear. Halos affix to the nasal area without inserting into the nostrils, significantly reducing irritation and claustrophobia, driving patient comfort and health, and preventing face marks. Halos are hypoallergenic, with gentle, surgical grade adhesive, and free from BPA, corn or latex. Halos utilize innovative MagSeal™ technology for rapid connection and disconnection. BLEEP sleep solutions are covered by Medicare and insurance, and the only CPAP solutions fully made in the USA.

HPCS# A7034



Item# **CAP100590**  
Bleep Halos for  
Eclipse Mask Frame  
**32/Pack**  
HPCS# A7033

### Bleep DreamPort Sleep Solution

#### DreamWay Interface with 32/Pack of DreamPorts

Item#  
**CMB100381**

Make the mask a thing of the past.

Unlike traditional CPAP masks the DreamPort Sleep Solution from Bleep is headgear free, lightweight, flexible and easy to use. The patented solution uses medical grade adhesive strips to hold the DreamPorts (seals) in place eliminating headgear. And unlike traditional seals with the DreamPorts you will not wake up with marks on your face.

#### Includes:

- Bleep DreamWay Interface with 10 inch hose and connector
- Bleep DreamPorts (32-Pack, 16-Night Supply)
- Instruction Guide



Item# **CUB100379**  
Bleep DreamPorts for  
DreamWay Interface  
**32/Pack**

### Bleep DreamPreps Witch Hazel Pads

Item#  
**CAPB100433**  
**100 Pads**



### Bleep Vent Diffuser Kit

Item#  
**CAP100604**  
**1 Hook and Loop  
Black Strap and  
3 White Filters**



Item# **CAP100662**  
Bleep Vent Filter  
Replacement Pack  
**12/pack**

## Sleep Apnea Appliances

### Hush Strips®

Item#

**CAP1025**



**Box of 32 Strips**

**Hush Strips® are designed to improve sleep quality by reducing mouth breathing and snoring during sleep**

Hush Strips® sleep strips are made from medical grade nonwoven fabric and gentle easy-release adhesive, which moves with the skin flexibly, making them easy to reposition and remove.

**Easy to apply residue-free material**

The simple, non-invasive and easy to use anti-snoring strips are made of medical-grade hypoallergenic material which doesn't leave residue and is safe for both kids and adults!



#### How to use

1. Wash and dry lips
2. Remove strip from the liner
3. Close your mouth and turn lips slightly inward
4. Apply strip on lips and press to seal

## Airway Management

### myTAP



Item# **CM039**

#### Precision-fit, immediate, simple and affordable solution for snoring and mild to moderate sleep apnea.

- Offers immediate treatment
- Can be fitted in 15 minutes in a sleep lab, physician's office or DME
- Clinically proven in over 40 independent peer reviewed studies to treat snoring and obstructive sleep apnea (OSA)
- Patients control and remold trays with Precision-Fit, enabling overall comfort, compliance and success: 90% of patients achieve ideal fit upon first try, no bite registration needed
- Thin, sleek trays create a comfortable, low-profile appliance
- Simple design makes it an affordable and convenient option
- Includes an AM Aligner to help offset potential bite changes
- Includes a mouth shield accessory to prevent dry mouth or excess saliva. A secondary benefit is the promotion of nasal breathing, the best natural way to breathe

## AirAvant

### Bongo Rx EPAP

Item#  
**CAPBNG500**

Bongo Rx all sizes starter kit

#### Includes:

- All four sizes: small, medium, large, and extra large
- Travel case
- Drying stand
- Headgear
- No machine, tube, electricity, or batteries needed
- Use at home, traveling, camping... use anytime, anywhere
- Small enough to fit in a shirt pocket
- Reusable and easy to clean with simply soap and water
- Designed, molded, and assembled in the USA
- Clinically tested to demonstrate safety and effectiveness
- Supplemental lifestyle option for your compliant patients
- Patient friendly option for your non-compliant patients



The Bongo Rx is a simple device with soft nasal seals and is available in various sizes (S, M, L, and XL). The Bongo Rx is an EPAP (expiratory positive airway pressure) device and works as follows:

- During inhalation, small valves open to allow normal breathing through the device.
- During exhalation, the valves close and direct exhaled air through small vents to generate EPAP, which keeps the airway open.

#### Annual Replenishment Pack

Includes 4 pairs, travel case, drying stand, and headgear

Size	Item#
Small	<b>CAPBNG502</b>
Medium	<b>CAPBNG503</b>
Large	<b>CAPBNG504</b>
Extra Large	<b>CAPBNG505</b>



#### Single Replenishment Unit

Includes 1 pair and headgear

Size	Item#
Small	<b>CAPBNG512</b>
Medium	<b>CAPBNG513</b>
Large	<b>CAPBNG514</b>
Extra Large	<b>CAPBNG515</b>



## CPAP Tubing

### Durable CPAP Tubing — Gray



Item#  
**TUB06**  
6ft  
19mm with 22mm cuffs

HCPCS# A7037


- Sunset's most popular item
- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

#### All Lengths

<b>TUB0018</b>	18in
<b>TUB02</b>	2ft
<b>TUB004</b>	4ft
<b>TUB06</b>	6ft
<b>TUB006.5</b>	6ft 6in — Compatible with ResMed tube wraps
<b>TUB007</b>	7ft
<b>TUB08</b>	8ft
<b>TUB10</b>	10ft

Compatible with Sunset's new Quick Release Connector!

Item#  
**CAP1016**



### Slim CPAP Tube

Item#  
**TUB006SS**  
6ft x 15mm with  
22mm cuffs  
HCPCS# A7037



#### Highly durable

Strain relief cuff decreases separation and tearing

#### Virtually no memory

Designed to be the most flexible tube available

#### Minimal weight

One of the lightest tubes on the market at 2.6 ounces

#### Great value

Affordable and sustainable

### Pressure Port Tubing

Item# **TUB006PP**  
6ft with  
Pressure Sensor Lines  
22mm Cuffs  
HCPCS# A7037



- Built-in sensor line for use with the Puritan Bennett: GoodKnight 420S, 420E and 425 and KnightStar 330

## Heated CPAP Tubing

### Heated Tube for Respironics

---

Item#


**TUB06HPR**

Heated Tube for  
Philips Respironics  
Dreamstation Series  
and System One

**Each**

OEM# HT15  
HCPCS# A4604

- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off



### Heated Tube for Fisher & Paykel ICON

---

Item#

**TUBHC900ICON208**

OEM# 900ICON208



## Respironics

### DreamStation 2



**Item# CF1012**

Disposable Filter  
**6/pack**

OEM# 1142831

HCPCS# A7038

**CF1012-1** 1/pack OEM# 1142828

**CF1012-2** 2/pack OEM# 1142829

**CF1012-3** 3/pack

**CF1012-50** 50/pack

**Item# CF1012F**

Reusable Filter  
**2/pack**

OEM# 1142687

HCPCS# A7039

**CF1012F-1** 1/pack

### DreamStation



**Item# CF1010**

Disposable Filter  
**6/pack**

OEM# 1122519

HCPCS# A7038

**CF1010-1** 1/pack

**CF1010-2** 2/pack

**Item# CF1010F**

Reusable Filter  
**2/pack**

OEM# 112446 (1-pack)

HCPCS# A7039

**CF1010F-1** 1/pack

**CF1010F-6** 6/pack

### M-Series



**Item# CF1007**

Ultrafine Filter - Tabbed  
**6/pack** • 1<sup>3</sup>/<sub>4</sub>" x 7<sup>7</sup>/<sub>8</sub>"

HCPCS# A7038

**CF1007-1** 1/pack

**CF1007-2** 2/pack

**Item# CF1007F**

Foam Filter  
**2/pack** • 1<sup>5</sup>/<sub>8</sub>" x 7<sup>7</sup>/<sub>8</sub>"

OEM# 1029330

HCPCS# A7039

**CF1007F-1** 1/pack

### PR System One



**Item# CF1009NT**

Ultrafine Filter - No Tab  
**6/pack** • 1<sup>3</sup>/<sub>4</sub>" x 7<sup>7</sup>/<sub>8</sub>"

OEM# 1029331

HCPCS# A7038

**CF1009NT-1** 1/pack

**CF1009NT-2** 2/pack

**Item# CF1009**

Ultrafine Filter - Tabbed  
**6/pack** • 1<sup>3</sup>/<sub>4</sub>" x 7<sup>7</sup>/<sub>8</sub>"

HCPCS# A7038

**CF1009-1** 1/pack

**CF1009-2** 2/pack

**Item# CF1009F**

Foam Filter  
**2/pack** • 1<sup>5</sup>/<sub>8</sub>" x 7<sup>7</sup>/<sub>8</sub>"

OEM# 1029330

HCPCS# A7039

**CF1009F-1** 1/pack

*All filters can be ordered in any pack quantity.  
Contact your sales rep for details.*

## Respironics

### REMstar Pro, REMstar Plus



Item# **CF1006**  
Ultrafine Filter  
**2/pack** • 3<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
OEM# 1005945  
HCPCS# A7038

Item# **CF1006F**  
Foam Filter  
**2/pack** • 3<sup>5</sup>/<sub>8</sub>" x 1<sup>5</sup>/<sub>8</sub>"  
OEM# 1005964  
HCPCS# A7039

### BiPAP Pro, Plus, Duet LX, Synchrony



Item# **CF1005**  
Ultrafine Filter  
**6/pack** • 4<sup>3</sup>/<sub>4</sub>" x 3<sup>1</sup>/<sub>4</sub>"  
OEM# 622219  
HCPCS# A7038

Item# **CF1005F**  
Foam Filter  
**2/pack** • 4<sup>3</sup>/<sub>4</sub>" x 2<sup>3</sup>/<sub>8</sub>"  
OEM# 622220  
HCPCS# A7039

### Aria, Virtuoso, Duet, Quartet



Item# **CF1004**  
Ultrafine Filter  
**6/pack** • 3<sup>1</sup>/<sub>2</sub>" x 2<sup>3</sup>/<sub>16</sub>"  
OEM# 532311  
HCPCS# A7038

Item# **CF1004F**  
Foam Filter  
**2/pack** • 3<sup>1</sup>/<sub>4</sub>" x 2<sup>1</sup>/<sub>8</sub>"  
OEM# 532247  
HCPCS# A7039

### REMstar



Item# **CF1003**  
Ultrafine Filter  
**6/pack** • 5<sup>1</sup>/<sub>2</sub>" x 4<sup>3</sup>/<sub>4</sub>"  
OEM# 362522  
HCPCS# A7038

Item# **CF1003F**  
Foam Filter  
**1/pack** • 5<sup>3</sup>/<sub>4</sub>" x 5<sup>1</sup>/<sub>4</sub>"  
OEM# 362521  
HCPCS# A7039

### Solo, Aria LX, Virtuoso LX, REMstar LX



Item# **CF1002**  
Ultrafine Filter  
**6/pack** • 3<sup>3</sup>/<sub>8</sub>" x 2<sup>3</sup>/<sub>4</sub>"  
OEM# 622017  
HCPCS# A7038

Item# **CF1002F**  
Foam Filter  
**2/pack** • 3<sup>5</sup>/<sub>16</sub>" x 2<sup>5</sup>/<sub>16</sub>"  
OEM# 622018  
HCPCS# A7039

### SleepEasy/BiPAP



Item# **CF1001**  
Ultrafine Filter  
**6/pack** • 7<sup>3</sup>/<sub>4</sub>" x 4<sup>7</sup>/<sub>8</sub>"  
OEM# 302064  
HCPCS# A7038

## ResMed

### AirSense 11



Item#  
**CF2111**  
Disposable Filter  
**3/pack**  
HCPCS# A7038  
**CF2111-1** 1/pack OEM# 39300  
**CF2111-2** 2/pack OEM# 39301  
**CF2111-50** 50/pack OEM# 39303

Item#  
**CF2111HA-3**  
Hypoallergenic Filter  
**3/pack**  
1<sup>1</sup>/<sub>4</sub>" x 1<sup>3</sup>/<sub>16</sub>"  
HCPCS# A7038  
**CF2111HA-1** 1/pack  
**CF2111HA-2** 2/pack

AirSense 11 Carry Bag  
**CAP1028-S11**



### S9 and AirSense 10 Series



Item#  
**CF2107**  
Disposable Filter  
**3/pack** • 2<sup>1</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>8</sub>"  
HCPCS# A7038  
**CF2107-1** 1/pack OEM# 36850  
**CF2107-2** 2/pack OEM# 36851  
**CF2107-50** 50/pack OEM# 36853

Item#  
**CF2107HA**  
Hypoallergenic Filter  
**3/pack** • 2<sup>1</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>8</sub>"  
HCPCS# A7038  
**CF2107HA-1** 1/pack OEM# 36855  
**CF2107HA-2** 2/pack OEM# 36856  
**CF2107HA-50** 50/pack  
OEM# 36858

### AirMini



Item#  
**CF2108**  
Disposable Filter  
**2/pack** • 1<sup>1</sup>/<sub>4</sub>" x 1<sup>3</sup>/<sub>16</sub>"  
OEM# 38815  
HCPCS# A7038

### Adapt SV



Item#  
**CF2106**  
Disposable Filter  
**3/pack** • 1<sup>1</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
OEM# 26926  
HCPCS# A7038

Item#  
**CF2106HA**  
Hypoallergenic Filter  
**3/pack** • 1<sup>1</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
OEM# 26927  
HCPCS# A7038

## ResMed

### S8 Series



Item#  
**CF2105**  
Disposable Filter  
**3/pack** • 1<sup>3</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>8</sub>"  
OEM# 33916  
HCPCS# A7038

**CF2105-1** 1/pack OEM# 33915  
**CF2105-2** 2/pack OEM# 33973  
**CF2105-50** 50/pack OEM# 33918



Item#  
**CF2105HA**  
Hypoallergenic Filter  
**3/pack** • 1<sup>3</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>8</sub>"  
OEM# 33920  
HCPCS# A7038

**CF2105HA-1** 1/pack OEM# 33919  
**CF2105HA-2** 2/pack OEM# 33975  
**CF2105HA-50** 50/pack OEM# 33922

### S7 Series



Item#  
**CF2103**  
Disposable Filter  
**3/pack** • 1<sup>1</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
OEM# 30907  
HCPCS# A7038

**CF2103-1** 1/pack OEM# 30919  
**CF2103-2** 2/pack  
**CF2103-50** 50/pack OEM# 30908



Item#  
**CF2103HA**  
Hypoallergenic Filter  
**3/pack** • 1<sup>1</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
HCPCS# A7038

**CF2103HA-1** 1/pack OEM# 30933

### S6 Series



Item#  
**CF2102**  
Disposable Filter  
**3/pack** • 5<sup>3</sup>/<sub>4</sub>" x 5<sup>5</sup>/<sub>8</sub>"  
OEM# 21935  
HCPCS# A7038

**CF2102-50** 50/pack  
OEM# 21936

## Fisher & Paykel

### SleepStyle Series



Item#  
**CF4004-2**  
Disposable Filter  
**2/pack**  
HCPCS# A7038



### ICON Series



Item#  
**CF4003**  
Disposable Filter  
**4/pack**  
2 1/16" x 5/8"  
HCPCS# A7038



**CF4003-1** 1/pack  
OEM# 900ICON503

**CF4003-6** 6/pack  
OEM# 900ICON212

### HC230 Series, HC600 Series



Item#  
**CF4002**  
Disposable Filter  
**4/pack** • 2 3/4" x 7/8"  
OEM# 900HC240  
HCPCS# A7038  
**CF4002-2** 2/pack



### HC200, HC220, HC221



Item#  
**CF4001**  
Disposable Filter  
**4/pack** • 5" x 3/4"  
OEM# 900HC222  
HCPCS# A7038  
**CF4001-2** 2/pack



## React Health/3B

### 3B Luna G3 and React G3X



Item# **CF8006-6**  
Disposable Filter  
**6/pack**  
HCPCS# A7038  
**CF8006-1** 1/pack  
**CF8006-2** 2/pack



Item# **CF8006F-1**  
Reusable Filter  
**1/pack**  
HCPCS# A7039



### Luna and Luna II



Item#  
**CF8005**  
Disposable Filter  
**12/pack** • 2" x 1 3/16"  
HCPCS# A7038



Item#  
**CF8005F-1**  
Reusable Filter  
**1/pack** • 2" x 1 3/16"  
HCPCS# A7039  
**CF8005F-2** 2/pack



### RESmart




Item#  
**CF8004**  
Disposable Filter  
**3/pack** • 3 5/16" x 1 5/8"  
HCPCS# A7038



## AEIOMed

### Everest Aura, Everest III and Integrated



Item#  
**CF8001**  
Foam Filter  
**2/pack** • 2" x 1<sup>9</sup>/<sub>16</sub>"  
OEM# 100446  
HCPCS# A7039

## Beyond

### ResPlus



Item#  
**CF9005F-2**  
Reusable Filter  
**2/pack**  
HCPCS# A7039  
**CF9005F-1** 1/pack

## Breas/GE

### iSleep 20 Series



Item# **CF7003**  
Ultrafine Filter  
**2/pack** • 2<sup>1</sup>/<sub>4</sub>" x 7/8"  
OEM# 004153  
HCPCS# A7038

Item# **CF7003F**  
Foam Filter  
**Each** • 2<sup>3</sup>/<sub>4</sub>" x 1"  
OEM# 004154  
HCPCS# A7039

### PV10



Item# **CF7002**  
Ultrafine Filter  
**2/pack** • 2<sup>3</sup>/<sub>4</sub>" x 7/8"  
OEM# 2B01975  
HCPCS# A7038

Item# **CF7002F**  
Foam Filter  
**Each** • 2<sup>5</sup>/<sub>8</sub>" x 7/8"  
OEM# 2B01972  
HCPCS# A7039

### PV1



Item# **CF7001**  
Ultrafine Filter  
**2/pack** • 2<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
HCPCS# A7038



Item# **CF7001F**  
Foam Filter  
**3/pack** • 2<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
HCPCS# A7039

## CareFusion/Viasys

### PureSom



Item#  
**CF9003**  
Ultrafine Filter  
**3/pack** • 1<sup>3</sup>/<sub>16</sub>" x 1<sup>7</sup>/<sub>16</sub>"  
OEM# 7511  
HCPCS# A7038

Item#  
**CF9003F**  
Foam Filter  
**3/pack** • 1<sup>3</sup>/<sub>16</sub>" x 1<sup>7</sup>/<sub>16</sub>"  
OEM# 7510  
HCPCS# A7039

### Orion and Pegasus



Item#  
**CF8002**  
Ultrafine Filter  
**6/pack** • 2<sup>7</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>4</sub>"  
OEM# 774688  
HCPCS# A7038



Item#  
**CF8002F**  
Foam Filter  
**Each** • 2<sup>7</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>4</sub>"  
OEM# 774689  
HCPCS# A7039

## Covidien Nellcor Puritan Bennett

### GoodKnight 314/318



Item#  
**CF6005F**  
Foam Filter  
**2/pack** • 6<sup>1</sup>/<sub>4</sub>" x 1"  
OEM# S-290609-00  
HCPCS# A7039

### GoodKnight 420/425



Item#  
**CF6003**  
Ultrafine Filter  
**2/pack** • 1<sup>3</sup>/<sub>16</sub>" x 3<sup>1</sup>/<sub>16</sub>"  
OEM# M-413950-04  
HCPCS# A7038

Item#  
**CF6003F**  
Foam Filter  
**Each** • 1<sup>3</sup>/<sub>16</sub>" x 3<sup>1</sup>/<sub>16</sub>"  
OEM# M-413950-02  
HCPCS# A7039

### GoodKnight 418 A/G/P/S



Item#  
**CF6002**  
Ultrafine Filter  
**4/pack** • 5<sup>3</sup>/<sub>16</sub>" x 3<sup>3</sup>/<sub>4</sub>"  
OEM# Y-413560-01  
HCPCS# A7038

### GoodKnight 418



Item#  
**CF6001**  
Ultrafine Filter  
**4/pack** • 2<sup>1</sup>/<sub>4</sub>" x 1<sup>2</sup>/<sub>5</sub>"  
OEM# S-645030-00  
HCPCS# A7038

Item#  
**CF6001F**  
Foam Filter  
**6/pack** • 2<sup>7</sup>/<sub>16</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
OEM# S-645029-00  
HCPCS# A7039

## DeVilbiss

### IntelliPAP



Item#  
**CF5004**  
Ultrafine Filter  
**4/pack** • 1<sup>3</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
OEM# DV51D-603  
HCPCS# A7038

Item#  
**CF5004F**  
Foam Filter  
**4/pack** • 1<sup>5</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>16</sub>"  
OEM# DV51D-602  
HCPCS# A7039

### Horizon LT



Item#  
**CF5003**  
Ultrafine Filter  
**4/pack** • 4<sup>1</sup>/<sub>4</sub>" x 1<sup>5</sup>/<sub>16</sub>"  
OEM# 8000D-603  
HCPCS# A7038

Item#  
**CF5003F**  
Foam Filter  
**4/pack** • 4<sup>1</sup>/<sub>4</sub>" x 1<sup>5</sup>/<sub>16</sub>"  
OEM# 8000D-602  
HCPCS# A7039

## Evo Medical

### RemRest



Item#  
**CF8003F**  
Foam Filter  
**3/pack** • 2<sup>3</sup>/<sub>16</sub>" x 1<sup>1</sup>/<sub>16</sub>"  
OEM# 21117  
HCPCS# A7039

### Comfort PAP



Item#  
**CF9003**  
Ultrafine Filter  
**3/pack** • 1<sup>3</sup>/<sub>16</sub>" x 1<sup>7</sup>/<sub>16</sub>"  
OEM# 7511  
HCPCS# A7038

Item#  
**CF9003F**  
Foam Filter  
**3/pack** • 1<sup>3</sup>/<sub>8</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
OEM# 7510  
HCPCS# A7039

Enhance your branding!

## FREE PRIVATE LABELING

Sunset's free service allows you to put your logo and contact information on our products.

Increase brand awareness and allow patients to easily contact you.

Work with our marketing team to customize the layout with your logo, phone number, website address, and reorder information.

Available on all CPAP tubes and filters.

Minimum orders apply.

Certain restrictions may apply.

Contact your sales rep for details.



## Healthyne

### Quest



Item# **CF3003**  
Ultrafine Filter  
with Curves  
**5/pack** • 4<sup>1</sup>/<sub>2</sub>" x 4<sup>3</sup>/<sub>8</sub>"  
OEM# 7302  
HCPCS# A7038

Item# **CF3003F**  
Foam Filter  
with Curves  
**5/pack** • 4<sup>1</sup>/<sub>2</sub>" x 4<sup>1</sup>/<sub>4</sub>"  
OEM# 7301  
HCPCS# A7039

### Tranquility



Item#  
**CF3001**  
Ultrafine Filter  
**5/pack** • 5<sup>1</sup>/<sub>8</sub>" x 1<sup>7</sup>/<sub>8</sub>"  
OEM# 7170  
HCPCS# A7038

## Löwenstein

### prisma



Item# **VF008-12**  
Disposable Filter  
**12/pack**  
OEM# WM29652

Item# **VF008F-2**  
Reusable Filter  
**2/pack**  
OEM# WM29928

## Resvent

### iBreeze



Ultrafine Filter  
**CF9901-1** 1/pack  
**CF9901-2** 2/pack  
**CF9901-6** 6/pack  
**CF9901-50** 50/pack  
HCPCS# A7038

## Invacare

### Polaris LT



Item# **CF9002**  
Ultrafine Filter  
**3/pack** • 5" x 2"  
OEM# 1082428  
HCPCS# A7038

Item# **CF9002F**  
Foam Filter  
**3/pack** • 5" x 2"  
OEM# 1082995  
HCPCS# A7039

### Polaris EX



Item# **CF9001**  
Ultrafine Filter  
**3/pack** • 2<sup>3</sup>/<sub>16</sub>" x 1"  
HCPCS# A7038

Item# **CF9001F**  
Foam Filter  
**3/pack** • 2<sup>1</sup>/<sub>4</sub>" x 1<sup>11</sup>/<sub>16</sub>"  
OEM# 1124211  
HCPCS# A7039

### Polaris TR



Item# **CF8001**  
Foam Filter  
**2/pack** • 2" x 1<sup>9</sup>/<sub>16</sub>"  
OEM# 100446  
HCPCS# A7039

## ProBasics

### Zzz-PAP



Item# **CF9003**  
Ultrafine Filter  
**3/pack** • 1<sup>3</sup>/<sub>16</sub>" x 1<sup>7</sup>/<sub>16</sub>"  
OEM# 7511  
HCPCS# A7038

Item# **CF9003F**  
Foam Filter  
**3/pack** • 1<sup>3</sup>/<sub>16</sub>" x 1<sup>7</sup>/<sub>16</sub>"  
OEM# 7510  
HCPCS# A7039

## Respironics

### Humidifier Chamber for PR System One



Item# **HUM1136502**  
Humidifier Chamber for PR System  
One 60 Series  
Each HCPCS# A7046

### Humidifier Chamber for DreamStation



Item# **HUM1122520**  
Humidifier Chamber  
for DreamStation  
Each HCPCS# A7046

## Fisher & Paykel

### Humidifier Chambers for HC100 Series



Item# **HUMHC300**  
Extended Life Humidifier Chamber  
for the Fisher & Paykel HC100 Series  
Each OEM# HC 300

Item# **HUMHC325S**  
Replacement Chamber for the  
Fisher & Paykel HC100 Series  
Each OEM# HC 325S  
HCPCS# A7046

Compatible with:

- HC100
- HC105
- HC150



### Humidifier Chambers for HC200 Series



Item# **HUMHC355**  
Extended Life Humidifier Chamber for  
the Fisher & Paykel HC200 Series  
Each OEM# HC 355



Item# **HUMHC385**  
Replacement Chamber for the Fisher  
& Paykel SleepStyle HC200 Series  
Each OEM# HC 385S  
HCPCS# A7046



Compatible with:

- HC200
- HC201
- HC202
- HC203
- HC204
- HC205
- HC210
- HC211
- HC212
- HC213
- HC214
- HC215
- HC220
- HC221
- HC222
- HC223
- HC224
- HC225
- HC230
- HC231
- HC232
- HC233
- HC234
- HC235
- HC240
- HC241
- HC242
- HC243
- HC244
- HC245
- HC250
- HC251
- HC252
- HC253
- HC254
- HC255

### Humidifier Chamber for HC600 Series



Item# **HUMHC365**  
Replacement Chamber for the  
Fisher & Paykel SleepStyle 600 Series  
Each OEM# HC 365S · HCPCS# A7046

Compatible with:

- HC604
- HC608



### Humidifier Chamber for ICON Series



Item# **HUMHC900ICON200**  
Replacement Chamber for the  
Fisher & Paykel ICON Series  
Each OEM# 900ICON200 · HCPCS# A7046

Compatible with:

- ICON Premo
- ICON Novo
- ICON Auto



## Sunset Chinstraps

### Universal Neoprene Chinstrap



**Item# CS001**  
Universal Chinstrap  
HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows

### Basic Chinstrap



**Item# CS002**  
Compare to  
Respironics Basic  
OEM# 302175  
HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to the top strap of headgear
- Cups the chin for comfort
- Works best with any nasal mask

### Adjustable Neoprene Chinstrap



Size	Item#
6in	<b>CS003S</b>
7in	<b>CS003R</b>
8in	<b>CS003M</b>
9in	<b>CS003L</b>
10in	<b>CS003XL</b>

Compare to Puresom  
Adjustable Ruby  
OEM# TMS-09 ADJ  
HCPCS# A7036

- Latex-free
- Thin adjustable strap attaches at the top of the head. Thicker bottom strap attaches to itself at the back of the head
- Great for hard-to-fit patients

For fitting, measure along one side of the head, as indicated by the dashed line above.

### Deluxe Chinstrap



Size	Item#
26in	<b>CS004S</b>
28in	<b>CS004L</b>

Compare to  
Respironics Deluxe  
OEM# 302425  
HCPCS# A7036

- Latex-free
- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort

### Comfort Chinstrap



**Item# CS005**  
Compare to  
ResMed Sullivan  
OEM# 16015  
HCPCS# A7036

- Latex-free
- One size fits most
- Easy to use —attaches to itself at the top of the head
- Soft cloth chin cup adds comfort

### Premium Chinstrap




**Item# CS006**  
Compare to  
Respironics Premium  
OEM# 1012911  
HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted
- Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.

## Sunset Chinstraps and Headgear

### Ruby Style Chinstrap



For fitting, measure completely around the head, as indicated by the dashed line above.

Size	Item#
21.5-22.5in	<b>CS007S</b>
22-23in	<b>CS007M</b>
23-24in	<b>CS007L</b>

Compare to Puresom Ruby

OEM# TMS-09 · HCPCS# A7036

- Latex-free
- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

### Halo Style Chinstrap



Item# **CS025**

Compare to BreatheWear Halo Chinstrap

OEM# H3407

HCPCS# A7036

### Ruby Style Adjustable Chinstrap



Size	Item#
Regular	<b>CS007-ADJ</b>
Extra Large	<b>CS007-ADJXL</b>

HCPCS# A7036

- Latex-free
- Fully adjustable
- Super soft material for extra comfort



### Universal Mesh Headgear



Item# **HG001**

Compare to Respironics SimpleStrap Headgear

OEM# 302328

HCPCS# A7035

- Latex-free
- One size fits most
- Can be used with most 4-point nasal or full face masks

### Universal Headgear



Size	Length	Item#
Regular	22.5in	<b>HG003</b>
Extra Large	24.5in	<b>HG003XL</b>

HCPCS# A7035

- Latex-free
- Can be used with most 4-point nasal or full face masks

## Cleaners

### Sunset Cleaning Wipes — Flowpack

Item#

**CAP1003S**

Flowpack of 64 wipes  
**12/case**

- Ingredients are safe and gentle for daily cleaning of all types of CPAP and oxygen masks and equipment



### Private Label CPAP Wipes

Order the flowpack wipes with your logo and contact info on a lid sticker, with a minimum quantity commitment. Talk to your sales rep for more info.

### Citrus II Cleaning Wipes

Item# **CAP1003**

Canister of 62 wipes  
**12/case**



### Cleaning Wipes — Single Packs

Item# **CAP1005**

Citrus II Box of 12 single-pack wipes  
**12/case**

- Convenient individually wrapped wipes are great for patient samples or travel use



### CleanSmart CPAP Disinfectant Spray

Item#

**CAP1018-L**

16oz • **30/case**

**CAP1018-S**

2oz • **180/case**

- Kills 99.9% of viruses & bacteria
- Included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2
- Active ingredient Hypochlorous Acid does not leave a chemical residue
- EPA Registration Number 89896-2
- Formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment



### Citrus II Mask Cleaning Spray

Item# **CAP1001**

1.5oz • **24/case**

Item# **CAP1002**

8oz • **12/case**

- All-natural cleanser for daily mask cleaning
- Available in 1.5oz travel size or 8oz non-aerosol spray cans




### Filter and Check Valve Kit for SoClean 2 CPAP Cleaning Device

Item# **CAP1007S-KIT**

**Each**

- With filter, check valve, retail packaging and instruction manual
- Replace every 6 months



**SUNSET**  
HEALTHCARE SOLUTIONS

## CPAP Cash Items

### CPAP Tube Cleaning Brush

Item#  
**CAP1013**  
Standard 6ft  
Each

- Brushes all the way through a 6ft CPAP tube without leaving debris in the middle
- Stainless steel wire and soft nylon bristles



### Handheld CPAP Pressure Meter

Item#  
**CAP003**  
Each

- Measures flow 0-18 cm/H2O
- Connects directly to the CPAP tubing without an adapter
- Lightweight device used to verify the pressure of a CPAP unit



### SleepRes V-Com™

Item#  
**CAP1026**  
Each

#### Helping patients adjust to CPAP on their own

Adding V-Com to a patient's therapy is like Training Wheels for CPAP™. This simple device can mean the difference between your patient giving up on PAP or continuing their lifesaving treatment.

Within minutes, a respiratory therapist can show a patient what the PAP pressure feels like with and without the V-Com device inline. This feedback mechanism allows them to better understand how their airflow should feel.

As you know, comfort equals adherence and adherence equals better treatment outcomes.



#### The V-Com Difference

Improve early PAP tolerance with V-Com by providing your patients with initial comfort that doesn't compromise their treatment.

Traditionally, when a patient struggles with the constant flow of air pressure, pressure settings are reduced.

Reducing pressure settings to increase a patient's initial tolerance improves comfort, but at a cost.

Decreasing pressure settings also means decreasing expiratory positive airway pressure (EPAP) which is principally responsible for therapy.

The V-Com effect is flow-dependent with maximum comfort during the high inspiratory flows. As a result, V-Com has minimal effect on expiratory positive airway pressure.

### 90-Degree CPAP Tubing Elbow

Item#  
**CAP1027**  
Each

The Sunset 90-Degree CPAP Tubing Elbow Adapter connects standard 22mm or slim 15mm CPAP tubes to most CPAP machines at a 90°, relieving stress on the tubing connection and creating a more convenient pathway to the bed.

OEM# 37394



### Nasal Pad

Item#  
**CAP1029**  
5/pack



**NEW**

- Compatible with most CPAP masks
- Soft and gentle, sensitive skin-approved, latex- and silicone-free
- Protects against skin breakdown, redness, irritation and pressure marks
- Conveniently reusable and simple to apply, remove and clean
- Gel can improve mask stability, leak-prevention and seal for comfortable, restful sleep

## CPAP Battery and Power Inverter/Charger

### Portable Outlet UPS CPAP Battery

Item#  
**CAP1600X9**

Includes battery, wall charger, and travel bag

- 159 Watt, 14.4 Ah UPS and Portable CPAP Battery
- Works with all models of CPAP and BiPAP machines
- Does not require additional adapters or connectors
- Power a CPAP and 2 USB devices at the same time
- Can power laptops and other small electronics
- Pure Sine Wave technology adds battery life
- Color display & power percentage indicator
- Meets TSA & FAA requirements for air travel
- Works directly with international voltages
- Fast charging USB-A and USB-C, 5V/3.0 A
- 1 year replacement warranty



#### Specs


<b>Battery Capacity</b>	159 Wh/14.4Ah
<b>AC Input</b>	100-240V/50-60Hz
<b>DC Input</b>	12-19V, 2-3.2A
<b>AC Output</b>	110V/60Hz (will match AC input)
<b>USB Output</b>	5V, 3.0A
<b>Rated/Peak Power Output</b>	100W/145W
<b>Optimal Operating Temp</b>	40°F to 105°F
<b>Dimensions</b>	W 6" X D 9.25" X H 1.5"
<b>Weight</b>	3lbs



#### When power goes out, your CPAP stays ON!®



Plug your CPAP into the Portable Outlet. The battery stays fully charged.



When there is an outage, the **battery provides power**. Your therapy continues without any interruption.



### DC to AC Power Inverter, Car Charger, Lighter Plug-in

Item# **CAPDC2AC150**

Recharge your Portable Outlet CPAP Battery using a car's 12V lighter outlet

- Can also power most CPAPs and other small electronics
- 150 Watts
- 3-prong outlet
- 2 USB ports



## CPAP Cash Items

### Transcend Micro™ 510 Auto CPAP

#### Each

##### Item#

**CAPSOM507000** White

**CAPSOM507001** Pink

**CAPSOM507002** Gray

**CAPSOM507003** Blue

- **Small:** About the size of a baseball
- **Lightweight:** Weighs less than a mobile phone
- **Works with any CPAP mask:** Use your existing mask
- **Quiet operation:** The Transcend WhisperSoft™ muffler is included

- Smallest and lightest travel CPAP device
- Compatible with any CPAP mask
- Unplug and sleep anywhere with more power options
- Adjustable GentleRise and AirRelief technology
- Stay dry, all the time
- AirMist HME – waterless humidification

#### Includes:

- WhisperSoft™ Muffler
- Hard Shell Case
- Power Supply
- TransFlex 6ft Hose
- AirMist HME (1 Adapter +1 Cartridge)
- Woven Filter (Inside the unit)
- Extra Hypoallergenic Filter
- 3 Year manufacturer's warranty

★★★  
MADE IN  
USA



#### SCAN FOR MORE: VIDEOS AND SPECS



#### Transcend Micro™ Accessories

##### Item#

**CAPSOM505001** PowerAway® Battery

**CAPSOM505049** Sun Watt Solar Panel Charger

**CAPSOM505035** Plug Adapters (AU/EU/UK)

**CAPSOM505010** Power Supply

**CAPSOM505030** AirMist Starter Pack (2 filters and 1 holder)

**CAPSOM505009** AirMist HME Pack (4 filters and 1 holder)

**CAPSOM505014** AirMist HME Cartridges (4 Pack)

**CAPSOM505015** AirMist HME Cartridges (8 Pack)

**CAPSOM505016** AirMist HME Cartridges (12 Pack)

**CAPSOM505008** SleepPak Travel Bag

**CAPSOM505038** WhisperSoft Muffler Kit: 6' Hose and Muffler

**CAPSOM505051** Standard Woven Filters - 2/PK

**CAPSOM505052** Hypoallergenic Filters - 2/PK

#### Item# **CAPSOM505001**

PowerAway® Battery



#### Extended Run Time

- **Long-Lasting Power:** Offers 2-3 nights or 17.7 hours of run time on a single charge at a pressure setting of 9cm, ensuring restful sleep without frequent recharging.
- **Quick Charge Time:** Recharges in less than 4 hours, so it's ready to go when you are.

#### Reliability and Safety

- **FAA-Approved:** Cleared for in-flight use, making it an essential for travelers
- **Backup Power:** Provides a permanent backup power source during outages, ensuring uninterrupted sleep therapy.
- **Long Life Span:** Capable of 1,000+ full charge/discharge cycles for long-term reliability.

#### Item# **CAPSOM505054**

Resupply Pack

#### Includes:

- Micro WhisperSoft Muffler
- 6ft AirFlex Hose
- AirMist HME Cartridges Pack (4 filters and 1 holder)
- Standard Woven filters (2 pack)



#### Item# **CAPSOM505053**

Essential Pack

#### Includes:

- AirMist HME Cartridges Pack (4 filters + 1 holder)
- Micro WhisperSoft Muffler
- (1) 6ft AirFlex Hose
- (1) 4ft AirFlex Hose
- (2) Hose Clips
- Travel Power Supply Plug Adapters
- Standard Woven filters (2 pack)
- Hypoallergenic filters (2 pack)
- Transcend Micro Travel Pouch



#### Item# **CAPSOM505055**

Essential Bundle

#### Includes:

- Essential Pack
- White Micro 510
- SleepPak Travel Bag



#### Item# **CAPSOM505056**

Power Bundle

#### Includes:

- Essential Pack
- White Micro 510
- PowerAway Battery
- SleepPak Travel Bag



## CPAP Cash Items

### Snugz CPAP Mask Liners

Item#

#### **CAP7001**

Machine-Washable Liners  
for Most Full Face Masks  
**2/pack**



#### **CAP7002**

Machine-Washable Liners  
for Most Nasal Masks  
**2/pack**

- 2 per pack, 90-day supply
- Machine washable
- Fits MOST nasal or full face CPAP masks (Not recommended for ResMed Mirage Quattro or Respirationics Amara View and DreamWear)
- Blend of spandex, moisture-wicking polyester and lycra fibers
- Latex-free

Snugz is a machine washable, one-size-fits-most liner for CPAP masks that is designed to help reduce noisy leaks, prevent redness and irritation, and feel great against the skin. The liner helps to grip the skin and keep the mask in place around the nose and mouth, reducing the need for a chinstrap. The improved seal is especially beneficial for high-pressure/high-flow PAP patients and can reduce nightly sleep disturbances and AHI.


### Soft Headgear Strap Covers

Item#

#### **CAP2004**



- Reduces strap marks
- No bulky and uncomfortable velcro
- Reduces strap irritation on sensitive skin
- Thick padding for extra comfort
- No need for clips – just slip wraps around headgear straps
- Easy and helpful for patients with arthritis in hands or fingers
- Can be machine washed and dried
- Wraps securely grip headgear and slide to adjust
- One size fits all
- Made from 230g Coral Velvet



### RemZzzs CPAP Mask Liners

Increases compliance by adding a comfortable and absorbent pad to CPAP masks. Reduces air leaks and protects skin from red marks.

Item#

#### **CAP3001(S,M,L)**

Liners for Full Face Resmed,  
Respirationics, DeVilbiss, and  
Fisher & Paykel

OEM# RZO-FFK1-SM32 (SMALL)  
RZO-FFK2-MD32 (MEDIUM)  
RZO-FFK3-LG32 (LARGE)

#### **CAP3002(S,M,L)**

Liners for Full Face Sunset, Fisher &  
Paykel FlexiFit, Forma, and Hans Rudolph

OEM# RZO-FFK4-SM32 (SMALL)  
RZO-FFK5-MD32 (MEDIUM)  
RZO-FFK6-LG-32 (LARGE)

#### **CAP3003(S,M,L)**

Liners for Nasal Respirationics Comfort  
and DeVilbiss Easy Fit

OEM# RZO-NAK7-SM32 (SMALL)  
RZO-NAK8-MD32 (MEDIUM)  
RZO-NAK9-LG32 (LARGE)

#### **CAP3004(XS, S,M,L)**

Liners for Nasal ResMed, Respirationics,  
Fisher & Paykel

OEM# RZO-NAK10A-XS32 (EXTRA SMALL)  
RZO-NAK10-SM32 (SMALL)  
RZO-NAK11-MD32 (MEDIUM)  
RZO-NAK12-LG32 (LARGE)

#### **CAP3005(L, XL)**

Liners for Respirationics  
FitLife Whole Face Masks

OEM# RZO-TF6A-SM32 (SMALL AND LARGE)  
RZO-TF6B-LG32 (EXTRA LARGE)

#### **CAP3006**

Liners for Nasal Pillows, Most  
Respirationics, ResMed, and TapPac

OEM# RZO-NP7A-UN32

#### **CAP3007**

Liners for Hybrid ResMed AirFit/  
AirTouch, Respirationics Amara View

OEM# RZO-MF6C-UN32



Sold in  
attractive  
packaging

## CPAP Cash Items

### Pur-Sleep CPAP Aromatherapy

#### Item# **CAP5001**

VaporClear Aromatherapy  
Starter Kit  
OEM# VAPOR



#### 30ml Aromatic Refills

Scent/OEM#	Item#
Black Cherry	<b>CAP5004-BCH</b>
Calm	<b>CAP5004-CLM</b>
Clear	<b>CAP5004-CLR</b>
Dark Chocolate	<b>CAP5004-DCH</b>
Deep	<b>CAP5004-DEE</b>
French Vanilla	<b>CAP5004-FRV</b>
Fresh	<b>CAP5004-FSH</b>
Peace	<b>CAP5004-PCE</b>
Spice	<b>CAP5004-SPC</b>
Strawberry	<b>CAP5004-STR</b>

#### Kit Includes:

- 1 aluminum diffusion stand
- 4 reusable diffusion pads
- 1.5ml bottle of Clear (peppermint, lime)
- 1.5ml bottle of Calm (eucalyptus, lavender, rosemary)


#### Item# **CAP5002**

Reusable Diffusion Pads  
**10/pack** • OEM# DP10

#### Item# **CAP5003**


Replacement Aluminum  
Diffusion Tray  
**Each** • OEM# TRAY

- Aromatherapy products specifically for use with CPAP
- Set up the diffuser right below the CPAP filter intake, add a drop of essential oil to the pad and the fragrance is drawn into the CPAP machine



### Quick Tube Release Connector

Item# **CAP1016**  
**Each**



### Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item# **BF2200**  
22mm  
**5/pack**

Attach between the machine and tube for additional protection against bacteria and allergens



### Circuit Connectors

Item# **RES013**  
22mm OD x  
15mm ID/  
22mm OD  
**10/pack**



Item# **RES014**  
22mm OD x  
22mm OD  
With Lip  
**10/pack**



### Oxygen Bleed-In Adapter

Item# **RES019**  
22mm with 6mm  
Rotating Stem  
**10/pack**



**RES019-CAP**  
With Stem Cap  
**10/pack**

### Swivel Exhalation Port

Item# **VEN002**  
For CPAP and  
Ventilator Circuits



### Humidifier Tubing Adapter

Item# **RES016**  
22mm ID x 22mm OD • with Lip  
**10/pack**



# CPAP Cash Items

## CPAPMax Pillow

Item# **CAP4003**  
Each



OEM# 15-501R-1-920

Traditional pillows interfere with mask performance and can cause a stiff neck, shifting, leaks, and facial discomfort.



The CPAPMax Pillow offers superb mask accommodation, excellent support, and unparalleled comfort.



**Cool 3D Mesh Air-Flow**  
Reduces heat and perspiration



**Fiberfill**  
Soft microfiber pillow top with hypoallergenic fiberfill

**CPAPMax Replacement Cover**  
**CAP4003RC**

Each  
OEM# 1-630-101-1-933

**Replacement Pillow Case**

**CAP4003C**  
Each  
OEM# 1-626-900-1-02

**2 Comfortable Sides:**



**Cozy and Soft Super Plush Side**  
Pillow top with hypoallergenic fiberfill



**NEW! Cooling Gel Side**  
Reduces heat and perspiration

**Adjust Height to Your Comfort**



## 10-in-1 Flip Pillow

Item# **CAP4006**  
Each



**The most versatile wedge and bed pillow available!**

10 support positions with 1 pillow! Comfort and align your neck, shoulder and lumbar region when used as back support cushion. Use it as a knee or leg pillow wedge or even in an office chair for added support. Filled with a soft hypo-allergenic fiber-fill.

**10-in-1 Flip Pillowcase Cover**  
**CAP4006CN** Navy  
**CAP4006CW** White  
Each



## Leg Pillow

Item# **CAP29450R**  
Each



**New & Improved Patented Shape**

Stays put - even when you bend your knees!

**Breathable, Washable cover**

Plush, soft cover easily removes for washing.



- Restores Your Natural Alignment for a Healthier Night's Sleep
- Provides perfect alignment of back and hips, relieves discomfort in hips and joints

**Doctor Recommended for Back, Hip & Knee Comfort**

- Helps Align Lower Back
- Stabilizes Hips
- Cushions Knees
- 25% More Leg Support
- Ventilated Cushion & Cover
- Cool Dry Comfort



**New Packaging!**  
Clamshell design conveys high quality of product & efficient use of shelf space

## CPAP Cash Items

### Folding Wedge Pillow

Each



Height	Item#
7 in	<b>CAP3092407R</b>
10 in	<b>CAP3092410R</b>
12 in	<b>CAP3092412R</b>

The Folding Wedge will help your customers sleep comfortably elevated, help aid digestion, and ease breathing and sinus pressure. The upright, angled position will provide comfort while supporting the head and chest. Plus, they'll sleep with their spine properly aligned which helps relieve muscle tension and back pain. The Folding Wedge easily unfolds for nighttime use or folds to store during the day.

#### Features & Benefits:

- Helps aid digestion
- Also helps with other conditions often associated with Acid Reflux such as respiratory issues, nasal and sinus congestion, and snoring
- Soothes muscle tension and supports the spine
- Includes a soft, velour zippered cover for easy removal
- Folds for convenient storage and fits under most beds
- Cover is machine washable





**New Packaging!**  
Merchandise Bulky Wedges Neatly on a Shelf

### CoolPAP Pillow

Item# **CAP4007**

Each

**Works with any style and brand of mask!**



Rolls up for travel or storage



New retail packaging!

**CoolPAP Pillow Cover**  
**CAP4007C**  
Each



### Deluxe Tubing Wrap

Item# **CAP2002**  
6ft Nylon Tubing Wrap  
Each  
Compare to ResMed Tubing Wrap



- 6ft nylon wrap with fleece lining maximizes heat retention and reduces noise
- Zipper closure for ease of use



Both Tubing Wraps come in attractive retail packaging

### Comfort Tubing Wrap

Item# **CAP2001**  
6ft Velour Tubing Wrap  
Each  
Compare to Respironics Tubing Wrap



- 6ft plush velour material adds comfort while preventing rainout
- Zipper closure for ease of use

## CPAP Cash Items

### CPAP Bedside Storage Bag

Item#

**CAP1015**

- Well-made and stylish
- Designed to hold CPAP supplies such as a mask and tube
- Weighted lid allows bag to hang off top of nightstand or tabletop



### Carry Bag for AirSense 11

Item#

**CAP1028-S11**

**NEW**



## CPAP HCPCS and Replacement Schedule

HCPCS	Item	Recommended Replacement Schedule
A4604	Heated Circuit Tubing	1 per 3 months
A7027	Combination Oral/Nasal Mask	1 per 3 months
A7028	Replacement Oral Cushion for Combination Oral/Nasal Mask	2 per month
A7029	Replacement Nasal Pillows for Combination Oral/Nasal Mask — Pair	2 per month
A7030	Full Face Mask	1 per 3 months
A7031	Replacement Full Face Mask Cushion or Liner	1 per month
A7032	Replacement Nasal Cushion or Liner	2 per month
A7033	Replacement Pillows — Pair	2 per month
A7034	Nasal or Pillows Mask	1 per 3 months
A7035	Headgear	1 per 6 months
A7036	Chinstrap	1 per 6 months
A7037	Tubing	1 per 3 months
A7038	Filter — Disposable	2 per month
A7039	Filter — Nondisposable	1 per 6 months
A7044	Oral Interface	
A7045	Replacement Exhalation Port	
A7046	Replacement Water Chamber	1 per 6 months
A9279	Monitoring Feature/Device	
E0470	BiLevel Device without Backup Rate	
E0471	BiLevel Device with Backup Rate	
E0561	Humidifier - Non-Heated	
E0562	Humidifier - Heated	
E0601	CPAP Device	

*Due to the rapid changes in reimbursement, this information is for reference only. Sunset Healthcare Solutions cannot guarantee its comprehensiveness or accuracy. Sunset urges you to contact a reimbursement specialist for up-to-date information and advice.*

# Oxygen Products

Oxygen	Page
Pulse Oximeter	40
SpO2 Sensors	41
Firesafe	42
Portable Concentrator	43
Concentrators	44
Masks	45
Cannula and Cannula Accessories	46-47
Tubing	47
Microdial Flowmeter	48
Regulators and Conserver	49
Oxygen HCPCS Chart	50
Cylinder Carriers	51-52
Connectors	53
Accessories	53-55
Cleaners	55
Analyzers	56
Celon Wet Bands	56
Bacteria and HEPA Filters	57
Concentrator Cabinet Filters	58-64
Caire/AirSep/SeQual	58
CHAD	59
Covidien Nellcor Puritan Bennett	59
Dynarex	59
OxyGo	59
DeVilbiss/Drive	60
Invacare	61
Rhythm, Jumao, Attendant, Lifestyle, Medacure	62
Respironics	63
Healthdyne	64
React Health/3B	64
Nidek	64
Repair Parts	65-69

## Pulse Oximeter

### Finger Pulse Oximeter

#### Item# **RES5102**

- Oversized 1.3" OLED display face is the largest of any finger pulse oximeter
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- CE, ISO and FDA approved
- Includes lanyard and manual



# SpO2 Sensors

## Disposable SpO2 Sensor Neonatal/Adult

Item# **RES7100**  
Similar to Nellcor MAX-N  
**20/box**




## Neo-M-1.5 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7200**  
1.5ft cord • Non-replaceable  
**20/box**



## Neo-M-3 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7201**  
3ft cord • Non-replaceable  
**20/box**




### Coban Tape Roll

Item# **RES7202CT**  
Use for sensor life extension or sensor cord stability  
**Each roll 1" x 177"**



## Infant-M-1.5 Disposable SpO2 Sensor Infant (3-20 kg)

Item# **RES7203**  
1.5ft cord  
**20/box**



## Pediatric-M-1.5 Disposable SpO2 Sensor Pediatric (10-50 kg)

Item# **RES7204**  
1.5ft cord  
**20/box**



## SpO2 Adapter Cable

Item#  
**RES7200R97-10** 10ft  
**RES7200R97-14** 14ft  
**Each**



**NEW**

## SpO2 Adapter Cable

Item#  
**RES7200R57-4** 4ft  
**RES7200R8-10** 10ft  
**Each**



**NEW**

# Firesafe

## Firesafe Thermal Shut-Off Valve


Item#  
**RES010**  
Each

Bidirectional



- Now smaller
- Laser etched expiration date
- Human readable UDI
- Designed to last 5 years

- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

### Fatalities per 100,000 patients



\* 1 death was reported among 73% of the patient population (2013-17)

## FireSafe Kits

Item#  
**RES010KIT**  
2 Firesafe Valves  
1 Tubing Connector



Item#  
**RES010CANNULAKIT**  
2 Firesafe Valves  
1 Tubing Connector  
1 Adult Soft Cannula



Item#  
**RES010VAKIT**  
2 Firesafe Cannula Valves  
1 Tubing Connector  
1 Cannula with Female Port  
1 7ft Oxygen Supply Tube



## FireSafe Tubing Connector

Item#  
**RES010TC**  
100/case

**RES010TC-10**  
10/case



## FireSafe Nozzle

Item#  
**RES010FND-20**  
20/pack



# Belluscura Oxygen Concentrators

## X-PLOR® Portable Oxygen Concentrator

Item#

**RES52051001**

### X-PLOR

Lightweight and compact, the X-PLOR® is ideal for patients who want to replace their heavy oxygen tanks without missing a moment. The X-PLOR is also Bluetooth®-capable and can be connected to mobile devices through Belluscura's newly developed Nomad Biometrics™ App. With a rechargeable battery, easy-to-use interface, and high oxygen output, supplementary oxygen has never been more approachable or affordable.

### What's Included

- X-PLOR® Unit
- 8-Cell Battery
- User-Replaceable Sieve Bed Cartridge
- Air Intake Filter
- AC Power Adaptor
- DC Power Adaptor
- Carry Bag With Handle & Shoulder Strap

### Accessories

Item#

<b>RES50073</b>	User-Replaceable Sieve Bed Cartridge
<b>RES50188</b>	8-Cell Battery
<b>RES53012</b>	DC Power Adaptor
<b>RES53011</b>	AC Power Adaptor
<b>RES52191</b>	Carry Bag
<b>RES52291</b>	Backpack



### Features

- Oxygen output 82%-92% purity
- User-Replaceable Sieve Bed Cartridge
- Airgonomically™ Designed to be lightweight and compact for easy travel
- Long-Lasting, Direct Charge Battery: Up to 4.5 hours of use when fully charged
- Pulse flow settings 1 through 5 delivering 215 to 900 ml/min
- No Breath Detected (NBD) mode is a battery-saving mode that runs at Auto-pulse15 BPM at the units current setting
- FAA Compliant

## Oxygen Concentrators

### Drive DeVilbiss Oxygen Concentrators

Item#  
**RES525DS**  
**5 Liter**  
OEM# 525DS



Item#  
**RES1025DS**  
**10 Liter**  
OEM# 1025DS

- Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter and 10 Liter Oxygen Concentrators were designed with enhanced quality, durability and simplicity in mind
- Patented DeVilbiss Turn-Down Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator
- Exclusive DeVilbiss OSD (Oxygen Sensing Device), standard on every unit, ensures patient safety and reliability for longer service intervals
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device

### Drive DeVilbiss iG02 Portable Oxygen Concentrator

Item#  
**RES125D**  
**iG02**

Item#  
**RES125D-XB**  
**iG02 with**  
**Extra Battery**



- Settings 1-5 Pulse Dose
- SmartDose Auto-Adjusting Conserver Technology included
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive conserver trigger on the POC market

### Covers for Concentrators



25×15×30"  
1.5mil  
**250 bags/roll**

Item#	Color
<b>EC251530C</b>	Clear
<b>EC251530T</b>	Tan Tint

## Oxygen Masks

### Simple O2 Mask

With 7ft Safety Tubing

Item# **RES2100**  
Adult Size  
**50/case**

Item# **RES2200V**  
Pediatric Size  
**50/case**

HCPCS# A4620



- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

### Partial Rebreather / Non-Rebreather O2 Mask with 7ft Tubing

Oxygen High Concentration  
Mask with 7ft Safety Tubing  
and Reservoir Bag

Item# **RES2101V**  
Adult Size  
**50/case**

OEM# 8120-7-50  
HCPCS# A4620



- 5-15 LPM
- Safety vent
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

### Non-Rebreather O2 Mask

With 7ft Safety Tubing  
& Reservoir Bag

Item# **RES2102**  
Adult Size  
**50/case**

Item# **RES2202V**  
Pediatric Size  
**50/case**

HCPCS# A4620



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape

# Oxygen Cannula

## Cannula with Tubing

### Adult Cannula with Tubing

Length	Material	Item#	Qty
4ft	Standard	<b>RES1104</b>	<b>50/case</b>
4ft	Soft	<b>RES1104S</b>	<b>50/case</b>
7ft	Standard	<b>RES1107</b>	<b>50/case</b>
7ft	Soft	<b>RES1107S</b>	<b>50/case</b>
15ft	Standard	<b>RES1115</b>	<b>50/case</b>
25ft	Standard	<b>RES1125</b>	<b>25/case</b>
25ft	Soft	<b>RES1125S</b>	<b>25/case</b>

HCPCS# A4615

- Clean, fresh look
- 1-6 LPM
- Curved nasal prongs for added comfort
- 6-channel safety tubing designed to prevent kinks
- Soft, lightweight material is memory free for long-term comfort



Item#

**RES1207V**

Pediatric Cannula  
with 7ft Supply Tube  
**50/case**



Item#

**RES1307V**

Infant Cannula  
with 7ft Supply Tube  
**50/case**




## Cannula without Tubing

### Adult Cannula without Tubing

Port	Material	Item#	Qty
Male	Standard	<b>RES1100</b>	<b>50/case</b>
Male	Soft	<b>RES1100S</b>	<b>50/case</b>
Female	Soft	<b>RES1100FS</b>	<b>50/case</b>

HCPCS# A4615



Male port

Female port

## Dual-Lumen (Demand) Cannula

Item# **RES4147V**

7ft

**25/case**

HCPCS# A4615



- 1-6 LPM
- Designed for use with dual port oxygen conservation delivery systems
- Dual tubing set provides independent tubes for oxygen delivery and sensing of patient breathing
- Smoothbore tubing

## High Flow Soft Cannula

Length	Item#	Qty
4ft	<b>RES1104SHF</b>	<b>25/case</b>
7ft	<b>RES1107SHF</b>	<b>25/case</b>
15ft	<b>RES1115SHF</b>	<b>25/case</b>
25ft	<b>RES1125SHF</b>	<b>25/case</b>

HCPCS# A4615



- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

# Oxygen Cannula

## Humidiflo High Flow Nasal Cannula

10/case

Size	Item#
Small	<b>HFN001AS</b>
Medium	<b>HFN001AM</b>
Large	<b>HFN001AL</b>

With F&P Adapter

Size	Item#
Small	<b>HFN001ASA</b>
Medium	<b>HFN001AMA</b>
Large	<b>HFN001ALA</b>

**HFN001-FPA** F&P Adapter



With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerated due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 - 60 L/min and O2 concentrations range from 21 - 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respironics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)


See page 85 for **LUISA Ventilator**



## Secure Cannula Stickers

Item# **RES080**  
25 pairs/case

- Hold cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



## Comfort Cannula with Ear Cushions

Item# **RES1107EC**  
7ft Cannula with Ear Cushions  
25/case  
HCPCS# A4615



Item# **RES020**  
Cushions Only  
50 Pairs/case

- Add cushions to existing cannula

# Oxygen Tubing

## Oxygen Supply Tube



Length	Item#	Qty
7ft	<b>RES3007</b>	<b>50/case</b>
15ft	<b>RES3015</b>	<b>50/case</b>
25ft	<b>RES3025</b>	<b>25/case</b>
40ft	<b>RES3040</b>	<b>20/case</b>
50ft	<b>RES3050</b>	<b>20/case</b>

HCPCS# A4616

- Crush resistant with low memory to prevent kinking
- Great durability and comfort level

## Oxygen Supply Tube — Green



Length	Item#	Qty
7ft	<b>RES3007G</b>	<b>50/case</b>
15ft	<b>RES3015G</b>	<b>50/case</b>
25ft	<b>RES3025G</b>	<b>25/case</b>
40ft	<b>RES3040G</b>	<b>20/case</b>
50ft	<b>RES3050G</b>	<b>20/case</b>

HCPCS# A4616

- Latex-free
- 1-6 LPM
- 6 channels

# Microdial Flowmeter

## Microflow Microdial Flowmeter

2024 Innovation Celebration Winner



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure	Inlet Pressure Range	Flow Range*
<b>RES8160048</b>			<b>A</b>
<b>RES8160050</b>	400 kPa (58 psi)	350 – 550 kPa	<b>E</b>



DISS Inlet • DISS Outlet

Item#	Inlet Pressure	Inlet Pressure Range	Flow Range*
<b>RES8160049</b>	400 kPa (58 psi)	345 – 500 kPa	<b>A</b>
<b>RES8160051</b>	350 kPa (51 psi)	250 – 450 kPa	<b>E</b>

\*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
<b>Range A</b>	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
<b>Range E</b>	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

**Assured, repeatable accuracy to enable safe, effective weaning.**

### What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.

### Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned



# Oxygen Regulators and Conserver

## Mini Oxygen Regulator CGA 870 Style



"T" Handle for easy mounting: no tools needed

Barbed hose outlet for easy connection to the oxygen tubing

**Item# RES042**  
 Adult • Barbed  
 0-8 LPM • **Each**  
 0, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8 LPM

**Item# RES043**  
 Adult • Barbed  
 0-15 LPM • **Each**  
 0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

**Item# RES046**  
 Pediatric • Barbed  
 0-4 LPM • **Each**  
 0, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1 1/2, 2, 3, 4 LPM

For use with E size or smaller cylinders: A, B, C, D, JD, E, M2, M4, M6, ML6, M7, M9, M15, M22, or M24.



HCPCS# E1353



Pediatric models have a green dial cap

## Mini Oxygen Regulator CGA 540 Style



Threaded DISS outlet

**Item# RES044**  
 Adult • DISS  
 0-8 LPM • **Each**  
 0, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8 LPM

**Item# RES045**  
 Adult • DISS  
 0-15 LPM • **Each**  
 0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

**Item# RES047**  
 Pediatric • DISS  
 0-4 LPM • **Each**  
 0, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1 1/2, 2, 3, 4 LPM



HCPCS# E1353

For use with M60 and larger cylinders: M60, M90, M, MM, M122, M250, H, J, T or other cylinders with CGA 540 style valve.

### Both Regulators Feature

- Pressure gauge reads to 3000 PSI
- Easy-click flow selector knob
- 2-year warranty



**Private label laser etching now available!**

Ask your rep for details.

## Oxygen Regulator CGA 540 Style Long Stem



**Item# RES045L**  
 Adult • DISS  
 0-15 LPM • **Each**  
 0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

## Bonsai Pneumatic Conserver

**Item# RES049**  
 Single Lumen • Barbed Outlet  
**Each**



- Delivers a uniform oxygen pulse in increments: 1, 2, 3, 4, 5, and 6
- Two continuous flow (CF) settings preset at 2 LPM and 4 LPM
- Delivers up to a 6:1 savings ratio from 14-40 breaths per minute
- Senses breath/delivers oxygen in first half of inspiratory cycle
- Operating pressure range: 200-3000 PSI
- Light weight and quiet operation
- Easy to operate single selector knob
- Compact design fits easily in virtually any oxygen carry bag
- Limited 2-year warranty

## Bags for Regulators and Conserver



**Item# EC1012C**  
 10×12"  
 1.5mil, Clear  
**1000 bags/box**

## Sample Oxygen Cylinder Specifications<sup>\*\*</sup>

Use the handy chart below to determine which accessories to purchase for your oxygen cylinder.

\* This information is intended to be used as a guide. Dimensions may vary by manufacturer. Sunset Healthcare Solutions does not sell oxygen tanks.



Older Name	A	B				C	D	JD	E		
<b>Newer Name</b>	M2	M4	M6	ML6	M7	M9	M15	M22	M24	M60	M/MM/M122
<b>Diameter (in.)</b>	2.5	3.2	3.2	4.3	4.3	4.3	4.3	5.3	4.3	7.3	8
<b>Height (in.)</b>	5.3	8.5	11.5	7.6	9.1	11	16.5	16.5	25.5	23	36
<b>Empty Weight (lb)</b>	0.7	1.6	2.2	2.8	3.3	3.7	5.3	8	7.9	22.3	39.5
<b>Capacity (L)</b>	42	113	164	165	198	255	425	640	680	1738	3455
<b>Transport Method</b>	CSBM6, CSBM6BP	CSBM6, CSBM6BP	CSBM6, CSBM6BP, CSBM6Flex, CSBM6-Carrier	CSBML6M7, CSBML6-Hardshell	CSBML6M7	CSBM9, CSBM6BP, CSBM6Flex, RES038B	CSBD, CSBDBP, CSBD1Roller, CSBD2Roller, RES038B, RES060		CSBE, CSBD1Roller, CSBD2Roller, RES038B, RES060	Not Portable	Not Portable
<b>Regulator Type</b>	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES044, RES045	RES044, RES045

## Oxygen HCPCS

HCPCS	Item
A4615	Cannula, nasal
A4616	Oxygen tubing, per foot
A4620	Oxygen mask, variable concentration
E0424	Stationary compressed gaseous oxygen system, rental
E0431	Portable gaseous oxygen system, rental
E0443	Portable oxygen contents, gaseous (for use only with portable gaseous systems when no stationary gas or liquid system is used), 1 month supply = 1 unit
E0445	Oximeter device for measuring blood oxygen levels non-invasively
E0560	Humidifier, durable for supplemental humidification during IPPB treatment or oxygen delivery
E1353	Oxygen supply regulator
E1354	Oxygen supply wheeled cart
E1355	Oxygen supply stand/rack
E1390	Oxygen concentrator, single delivery port, capable of delivering 85% or greater oxygen concentration at the prescribed flow rate
E1391	Oxygen concentrator, dual delivery port, capable of delivering 85% or greater oxygen concentration at the prescribed flow rate
E1392	Portable oxygen concentrator, rental
K0738	Portable gaseous oxygen system, rental; home compressor used to fill portable oxygen cylinders; includes portable containers, regulators, flow meter, humidifier, cannula or mask and tubing

# Cylinder Bags

## Comfort Shoulder Bags

Cylinder	Item#
ML6/M7	<b>CSBML6M7</b>
C/M9	<b>CSBM9</b>
B/M6	<b>CSBM6</b>
D	<b>CSBD</b>



**CASH ITEM**

- Comfort flex shoulder strap
- Double zipper for easy access
- Machine washable



CSBML6M7

CSBM9

CSBM6

CSBD

## Flex O2 Bag

Item# **CSBM6Flex**  
 Reconfigurable Bag for M6 and M9 Cylinders  
**Each**

- Adjusts to 5 carrying modes



**CASH ITEM**

## Backpack

Cylinder	Item#
B/M6/C/M9	<b>CSBM6BP</b>
D	<b>CSBDBP</b>

- Fully adjustable and comfortable backpack straps
- Machine washable
- Converts to either backpack or shoulder bag



**CASH ITEM**

## Wheelchair/Scooter Carrier

Item# **CSBE**  
 For E Cylinders  
**Each**

- Multi-part stabilizer straps for easy chair attachment
- Machine washable



**CASH ITEM**

## Cylinder Carriers

### ML6 Hardshell Tote

Item#  
**CSBML6-Hardshell**  
For Homefill Cylinders  
Each



**CASH ITEM**



### Double Tank Rolling/Shoulder Bag

Item#  
**CSBD2Roller**  
For D/E Cylinders  
Each

- May also be used with shoulder strap



**CASH ITEM**



### Oxygen Cylinder Carrier

Item#  
**CSBM6-Carrier**  
Soft Polyester Tote  
for M6 Cylinders  
Each

- Tote handles
- Padded shoulder strap
- Fastening strap for mobility
- Holds up to 8 M6 tanks



### Cylinder Rack

Item# **RES038B**  
For D, E, or M9 Cylinders —  
Lightweight, Durable  
Plastic Composite Rack  
Each



### D/E Cart

Item#  
**RES060**  
4/case

- Accommodates one D or E tank
- Sturdy yet lightweight (4.18 LBS welded metal construction)
- Adjustable height handle with PVC grip
- Large EVA Foam wheels have just the right amount of resilience for good traction and maneuverability
- Elegant, grey satin finish
- Bottom fork rest with rubber caps



## Oxygen Connectors and Accessories

### Christmas Tree Connector

Item# **RES001**  
Green without Swivel  
**10/pack**



Item# **RES002**  
Green with Swivel  
**10/pack**



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

### Oxygen Tubing Y-Connector

Item# **RES017**  
Fits 3/16in and 1/4in ID Tubing  
(All 3 sides are the same)  
**10/pack**



- Connects 2 oxygen sources for higher flows

### Oxygen In Use Signs

Item# **RES003**  
6in x 9in  
With String  
**50/pack**



- Double-sided (English and Spanish)

### Washer with Rubber Ring

Item# **RES035**  
Brass  
**25/pack**



### Oxygen Accessory Kit

Item# **RES1000**  
RES002 Christmas Tree Connector - Swivel  
RES012 Tubing Connector - Non-Swivel  
RES017 Tubing Y-Connector  
RES018 Tubing Connector - Swivel  
RES032 Wrench  
RES035 Brass Washer

Each



**CUSTOMIZABLE KIT**



- Order all common oxygen accessories in one easy-to-handle kit
- Contents can be customized to fit your needs

### Oxygen Tubing Connector — Non-Swivel

Item# **RES012**  
Fits 3/16in and 1/4in  
ID Tubing  
**10/pack**



- Barbed ends ensure secure, leak-free connections

### Oxygen Tubing Connector — Swivel

Item# **RES018**  
Male/Male  
**10/pack**  
**RES018-25 25/pack**



**RES018KIT-1**  
Female/Male  
**Each**



- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing

## Oxygen Accessories

### O2 Starter Kit

**Item# RES02KIT**

RES1107 Cannula with 7ft Tube x2  
 RES3015 15ft Supply Tube x2  
 RES002 Connector — Swivel  
 RES003 Oxygen In Use Sign  
 RES012 Tubing Connector — Non-Swivel  
 RES017 Tubing Y-Connector  
 RES018 Tubing Connector — Swivel  
 RES032 Wrench  
 RES035 Brass Washer



**CUSTOMIZABLE KIT**

### Wrenches and Chains

**Item# RES031**

With Chain **10/pack**

**Item# RES032**

Without Chain **10/pack**



- For use with post valve oxygen cylinders: A, B, C, D, JD or E tanks
- Made of tough, stress-resistant nylon

**Item# RES031M**

Metal Wrench with Chain for Post Valve Oxygen Cylinders and CGA 540 Style Regulators **5/pack**



**Item# RES031-CHAIN**

12in Ball Chain with Coupling Connector **50/pack**

- Add chain to an oxygen valve wrench



**Item# RES031-BW**

Cylinder Wrench with Two-Ring Bungee and Two Yoke Washers **10/pack**



### Disposable Bubble Humidifier 500 ml

**Item#**

**RES050**

Up to 15 LPM • 6 psi • 500 ml

**50/case**

HCPCS# E0560

OEM# HUM001



- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient

### Humidifier Connectors

**Item# HUM008**

L-Connector 12in • **10/pk**

**Item# HUM008L**

L-Connector 21in • **10/pk**

**Item# HUM008S**

Straight Connector 11.625in • **10/pk**



### Water Traps

**Item# RES037**

Push On End Caps • **Each**

**Item# RES037T**

Threaded End Caps • **Each**



- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

## Accessories and Cleaners

### Homefill Tubing Connector

Item# **RES039**  
10/pack

- Compatible with Invacare, Drive, and Respiration Homefill machines
- Durable, coiled tubing



### Homefill Racks

Item# **RES039-RACKB**  
For Invacare Perfecto2  
4/case




**RES039-RACKW**  
For Invacare Platinum  
4/case



**RES039-RACKD**  
For DeVilbiss  
525/1025 Series  
4/case

Easily attaches onto concentrator and is designed to conveniently hold Homefill and Transfill oxygen systems.



### Oxygen Liter Meter

Item# **RES034**  
0-10 Liter/Minute  
Each



- Oxygen pen
- Works with tanks and concentrators

### Germicidal Deodorizing Cleaner

Item# **RES6100**  
22oz spray bottles  
12/case

- Ready-to-use spray cleans, disinfects and deodorizes
- Formulated for use on hard, non-porous surfaces
- Bactericidal, Tuberculocidal, Fungicidal, Mildewstat
- Does not contain abrasives, bleach, phosphates or phosphorous
- Effective against SARS-CoV-2 (COVID 19 virus) on hard non-porous surfaces



### Germicidal Deodorizing Cleaner — Gallon

Item# **RES6102**  
1 Gallon

- No-scrub, hospital strength germicidal cleaner
- Cleans, deodorizes, and disinfects surfaces, equipment and non-critical instruments



# Oxygen Accessories

## Smart Check Analyzer

Item# **RES7003** **NEW**  
Each

Portable, handheld oxygen analyzer designed specifically for testing and maintaining oxygen concentrators.

The ultimate oxygen + flow + pressure analyzer.

The Smart Check is an all-in-one oxygen analyzer for testing continuous and pulse (POC) concentrators. Measure oxygen, flow, pressure and pulse volume with a single device. Use the Smart Check to also test CPAP, BIPAP and nebulizers.

The Smart Check is the only device that also offers a built-in flashlight, backlight, and integrated tube management system.

With its fast temperature response, SmartCheck quickly adapts to temperature changes, ensuring accurate results in less time.

Offering all these features at a very competitive price, SmartCheck delivers exceptional performance at an incredible value.



	Smart Check	UltraMax O2	Pro2 Check Elite
Measure Oxygen & Flow	✓	✓	✓
Test Pressure	✓	✓	✓
Quick Touch Pressure Check	✓	✓	✗
Calibrated for Life	✓	✓	✗
Test Pulse Concentrators	✓	✗	✗
Test CPAP + Nebulizer	✓	✗	✗
Flashlight	✓	✗	✗
Backlight	✓	✗	✗
Tube Management	✓	✗	✗
Fast Environmental Response	✓	✗	✗
Cost	\$	\$\$\$	\$\$

**Features:**

- LED Display
- Calibration Check Button
- Sample Inlet
- Flashlight
- Sample Outlet
- Battery Cover

## Celon Wet Bands

Create a perfect seal to protect oxygen tank valves

**Easy to apply, safe and effective. Biodegradable and sustainable material.**

- Eliminates dust, dirt & dampness
- Hand applied
- Heat-free solution
- Creates a secure and wrinkle-free seal

Item# **WB66120**  
66 x 120mm for CGA 540 Valves  
**Bucket of 1,500**



- Professional finish
- Environmentally friendly
- Biodegradable

Item# **WB30540**  
30.5 x 40mm for CGA 870 Valves  
**Bucket of 11,500**



Item# **WB30525**  
30.5 x 25mm for CGA 870 Valves  
**Bucket of 18,000**



# Oxygen Bacteria And HEPA Filters

## 6in HEPA Capsule with Removable Sticker

Item# **BFH6S**

Threaded

Each



## EverFlo Intake Filter

Item# **BF1001**

Each

Respironics OEM# 1038831



## Extended Life Bacteria Filter

Item# **BF5000**

Nonthreaded

Each



## Extended Life Intake Bacteria Filter

Item# **BF1103**

Each



## Final Bacteria Filter

Item# **BF516**

1/4in Barbed on Both Sides

10/pack



## Respironics SimplyGo Style Inlet Filter

Item#

**BF8000**

Each



## HEPA Filter for React Health/3B Stratus

Item# **BF6020**

HEPA Filter

Each

OEM# STR1020



## Intake Bacteria Filter

Item# **BF2222**

22mm fittings

5/pack



## Extended Life Bacteria Filter with Muffler

Item#

**OF6011**

Each



## Extended Life Intake Bacteria Filter

Item#

**OF1025D**

Each



## Intake Bacteria Filter

Item# **BF7000**

Jumao JMC5A, Attendant 8X899, Lifestyle LM5A, Medacure OCSL, Rhythm LM5CA, Drive DECO2 5L Concentrators

Each




## Final Bacteria Filter — 1/4in to 3/8in

Item# **BF1438**

1/4in to 3/8in Stepped Barb on Both Sides

10/pack



## Final Bacteria Filter — 3/16in

Item# **BF316**

3/16 Stepped Barb and 1/4in Straight Barb on Both Sides

10/pack



## Millennium M10 Style HEPA Filter

Item# **BF1002**

Long Life Filter

Each



## Caire/AirSep/SeQual

### SeQual Integra



Item# **OF11005**  
Foam Cabinet Filter  
**10/pack**  
75/8" x 31/2"-43/4 x 3/4"  
OEM# 2185



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**  
OEM# 2033



Item# **BF516**  
Final Bacteria Filter  
1/4in Barbed  
**10/pack**

### AirSep Newlife, QuietLife, Intensity



Item# **OF8002**  
Foam Cabinet Filter  
**10/pack**  
51/4" x 37/8" x 7/8"  
OEM# F1-50-08-05



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in  
Barbed  
**10/pack**  
OEM# F1-30-03



Item# **OF8001**  
Felt Filter  
**10/pack**  
2" OD x 5/8" ID x 1"  
OEM# M1 161-1



**KIT**

Item# **OF8002SP**  
Newlife Kit  
2 - OF8002  
1 - BF1438  
1 - OF8001  
Each

### CAIRE Companion 5



Item# **BF1103**  
Extended Life Intake  
Bacteria Filter  
**Each**  
OEM# 2185

### SeQual Eclipse



Item# **OF11006**  
Foam Cabinet Filter  
**10/pack**  
21/2" x 21/2" x 3/4"

## CHAD

Total O2



Item# **BF2222**  
Intake Bacteria Filter  
22mm fittings  
**5/pack**  
Healthdyne OEM# 100-6631-00  
Respironics OEM# H620



Item# **OF3003**  
Foam Intake Filter  
for use with BF2222  
25/8" OD x 3/4" ID x 23/8"  
**10/pack**  
OEM# H613

## Covidien Nellcor Puritan Bennett

590 and 592A



Item# **OF6001**  
Foam Cabinet Filter  
63/4" x 43/4" x 3/8"  
**10/pack**  
OEM# 492672



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**



**KIT**



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb and  
1/4in Straight Barb  
**10/pack**  
OEM# 492141

Item#  
**OF6001SP**  
590 and 592 Kit  
2 - OF6001  
1 - BF5000  
1 - BF316  
Each

## Dynarex and React Health

Resp-02 and PHOENIX 5L



Item# **OF6013-10**  
Foam Cabinet Filter  
5.75" x 3.625"  
**10/pack**



Item# **OF6011**  
Air Intake Filter  
**Each**  
OEM# RPOC-14AIF

## OxyGo

OxyHome



Item# **OF8006-10**  
Foam Cabinet Filter  
63/4" x 41/4" x 3/8"  
**10/pack**  
OEM# 1170-5401



Item# **BFH6S**  
HEPA Filter  
6in Capsule  
**Each**  
OEM# 2000806

## DeVilbiss/Drive

### SolAiris III, 303 Series



Item# **OF5004**  
Foam Cabinet Filter  
33/4" x 33/4" x 1/2"  
**10/pack**  
OEM# 30302605



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651



Item# **BF5000**  
Extended Life Bacteria Filter  
**Each**

### 1025 Series



Item# **OF1025D**  
Extended Life Bacteria Filter  
**Each**



Item# **OF5004**  
Foam Cabinet Filter  
33/4" x 33/4" x 1/2"  
**10/pack**  
OEM# 303DZ-605



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651

### SolAiris V, 505, 515 Series



Item# **OF5005**  
Foam Cabinet Filter  
5 1/2" x 35/8" x 5/8"  
**10/pack**  
OEM# 50502694



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651



Item# **BF5000**  
Extended Life Bacteria Filter  
**Each**



Item# **OF6011**  
Extended Life Bacteria Filter—22mm fittings  
**Each**  
Use if SolAiris was manufactured before Nov '03

Item# **OF5005SP**  
SolAiris V Kit  
2 - OF5005  
1 - OF6011  
1 - BF1438  
Each


### 525 DS Series



Item# **OF5004**  
Foam Cabinet Filter  
33/4" x 33/4" x 1/2"  
**10/pack**  
OEM# 303DZ-605



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651



Item# **OF6011**  
Extended Life Bacteria Filter—22mm fittings  
**Each**

### All MC Units



Item# **OF5001**  
Foam Cabinet Filter  
**10/pack**  
6 1/8" x 3 1/2" x 1/4"  
OEM# 444-514



Item# **BF5000**  
Extended Life Bacteria Filter • **Each**



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651

# Invacare

## Mobilaire 3 and 5



Item# **OF9001SP**  
 Mobilaire Kit  
 2 - OF9001  
 1 - BFH6S  
 1 - BF516  
 Each



Item# **OF9001**  
 Foam Cabinet Filter  
 6 1/2" x 33/8" x 1/2"  
**10/pack**  
 OEM# 2000489



Item# **BFH6S**  
 HEPA Filter  
 6in Capsule  
**Each**  
 OEM# 2000806



Item# **BF516**  
 Final Bacteria Filter  
 1/4in Barbed  
**10/pack**  
 OEM# P0750X030

## Perfecto2



Item# **OF9004SP**  
 Perfecto2 Kit  
 2 - OF9004  
 1 - BFH6S  
 1 - BF516  
 Each



Item# **OF9004**  
 Foam Cabinet Filter  
 7" x 33/4" x 1/2"  
**10/pack**  
 OEM# 1143492



Item# **BFH6S**  
 HEPA Filter 6in Capsule  
**Each**  
 OEM# 2000806



Item# **BF516**  
 Final Bacteria Filter  
 1/4in Barbed  
**10/pack**  
 OEM# P0750X030

## Platinum 5 and 10



Item# **OF9002SP**  
 Platinum Kit  
 2 - OF9002  
 1 - BFH6S  
 1 - BF516  
 Each



Item# **OF9002**  
 Foam Cabinet Filter 9" x 27/8" x 1/2"  
**10/pack**  
 OEM#1107412



Item# **BFH6S**  
 HEPA Filter 6in Capsule with Sticker  
**Each**  
 OEM# 2000806



Item# **BF516**  
 Final Bacteria Filter 1/4" Barbed  
**10/pack**  
 OEM# P0750X030



Item# **HUM008S**  
 Straight Humidifier Connector  
 w/ 9in Tube  
**10/pack**

## HomeFill



Item# **OF9001**  
 Foam Cabinet Filter  
 6 1/2" x 33/8" x 1/2"  
**10/pack**  
 OEM# 2000489



## Rhythm HC, Jumao, Attendant, Lifestyle, Medacure

### Rhythm LM5CA, Jumao JMC5A, Attendant 8X899, Lifestyle LM5A, Medacure OC5L

Compatible with

- Rhythm LM5CA



Item# **OF7001-10**  
Foam Cabinet Filter  
**10/pack**  
OEM# LM209069



Compatible with

- Jumao JMC5A
- Attendant 8X899
- Lifestyle LM5A
- Medacure OC5L
- Drive DECO2



Item# **OF7000**  
Foam Cabinet Filter  
**10/pack**

Compatible with

- Rhythm LM5CA
- Jumao JMC5A
- Attendant 8X899
- Lifestyle LM5A
- Medacure OC5L



Item# **BF7000**  
Intake Filter  
**Each**  
OEM# LM209066



Item# **BF516**  
Final Bacteria Filter,  
1/4in Barbed  
**10/pack**  
OEM# LM209068

## Rhythm

### Rhythm HC LM5BA



Item# **OF8000**  
Foam Cabinet Filter  
**10/pack**  
OEM# LM5BA-ITF



Item# **OF6011**  
Extended Life  
Bacteria Filter —  
22mm fittings  
**Each**  
OEM# LM5BA-ILF



Item# **BF516**  
Final Bacteria Filter,  
1/4in Barbed  
**10/pack**  
OEM# LM5BA-BF

### Rhythm LM10A



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**

OEM# LM10A-ILF

### Rhythm P2, P2-E6, and P2-E7



Item# **OF7002-5**  
Cotton Intake Filter  
**5/pack**  
OEM# P2FC-1

# Respironics

## EverGo



Item# **OF1005SP**  
EverGo Kit  
2 - OF1005  
1 - BF316  
Each



Item# **OF1005**  
Foam Cabinet Filter  
3" x 2" x 3/8"  
**10/pack**  
OEM# 900-150



Item# **BF316**  
Final Bacteria Filter  
3/16in Stepped Barb  
and 1/4in Straight Barb  
**10/pack**  
OEM# H621

## SimplyGo and SimplyGo Mini



Item# **BF8000**  
Inlet Filter for SimplyGo  
**Each**  
OEM# 1095248



Item# **RES1008**  
Compressor Filter for  
SimplyGo Mini  
**Each**

## Millennium



Item# **OF1001SP**  
Millennium Kit  
2 - OF1001  
1 - BF2222  
1 - OF3003  
1 - BF316  
Each



Item# **OF1001**  
Foam Cabinet Filter  
8<sup>3</sup>/<sub>16</sub>" x 4" x 3/8"  
(1/4" tabs)  
**10/pack**  
OEM# H619



Item# **BF2222**  
Intake Bacteria Filter  
22mm fittings  
**5/pack**  
OEM# H620



Item# **OF3003**  
Foam Intake Filter  
for use with BF2222  
2<sup>5</sup>/<sub>8</sub>" OD x 3/4" ID x 2<sup>3</sup>/<sub>8</sub>"  
**10/pack**  
OEM# H618



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb,  
1/4in Straight Barb  
**10/pack**  
OEM# H621



Item# **OF1002**  
Felt Filter  
2<sup>1</sup>/<sub>8</sub>" OD x 1/2"  
**6/pack**



Item# **BF1002**  
Millennium M10 Style  
Long Life Filter  
**Each**

## EverFlo



Item# **OF1003SP**  
EverFlo Kit  
1 - BF1001  
1 - BF316  
Each



Item# **OF1006SP**  
Cabinet Filter Kit  
(No Cover)  
**2 Filters, 1 Casing**  
OEM# 1055625



Item# **OF1006**  
Cabinet Filter  
5<sup>11</sup>/<sub>16</sub>" x 4<sup>1</sup>/<sub>2</sub>" x 7/16"  
**10/pk**



Item# **OF1006-CASE**  
Foam Cabinet Filter  
Casing  
**10/pack**



Item# **BF1001**  
EverFlo Intake Filter  
**Each**  
OEM# 1038831

Item# **HUM008**  
Humidifier  
L-Connector  
With 12in  
Tube  
**10/pk**



Item# **BF316**  
Final Bacteria Filter  
3/16in Stepped Barb  
and 1/4in Straight Barb  
**10/pack**  
OEM# H621

## Healthdyne

### Alliance



Item# **OF3002**  
Foam Cabinet Filter  
8<sup>1</sup>/<sub>4</sub>" x 3<sup>3</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
**10/pack**  
OEM# 100-0078-10



Item# **BF2222**  
Intake Bacteria Filter  
22mm fittings  
**5/pack**  
OEM# 100-6631-00



Item# **OF3003**  
Foam Intake Filter  
for use with BF2222  
2<sup>5</sup>/<sub>8</sub>" OD x 3<sup>3</sup>/<sub>4</sub>" ID x 2<sup>3</sup>/<sub>8</sub>"  
**10/pack**  
OEM# H618



Item# **BF316**  
Final Bacteria Filter  
3/16in Stepped Barb  
and 1/4in Straight Barb  
**10/pack**  
OEM# H621

## React Health/3B

### Stratus 5



Item# **OF6030**  
Foam Cabinet Filter  
**10/pack**  
OEM# STR1030



Item# **BF6020**  
HEPA Filter  
**Each**  
OEM# STR1020



Item# **BF516**  
Final Bacteria Filter,  
1/4in Barbed  
**10/pack**  
OEM# P0750X030

## Nidek

### Nuvo



Item# **OF11003**  
Foam Cabinet Filter  
9<sup>3</sup>/<sub>8</sub>" x 3<sup>15</sup>/<sub>16</sub>" x 5<sup>1</sup>/<sub>16</sub>"  
**10/pack**  
OEM# 9250-1025



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# 7631-1053

Item#  
**OF11003SP**  
Nuvo Kit  
2 - OF11003  
1 - BF1438  
1 - BF5000  
Each



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**  
OEM# 9250-1028

### Mark IV, V, V+



Item# **OF11004**  
Foam Cabinet Filter  
10" x 4<sup>1</sup>/<sub>8</sub>" x 5<sup>1</sup>/<sub>8</sub>"  
**10/pack**  
OEM# 9007-6G37



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in  
Barbed  
**10/pack**  
OEM# 7631-1053

# Repair Parts NEW



## Invacare Platinum 5 / Platinum 10

**Item# RP001**  
Pilot Valve and  
1 O-Ring  
**Each**  
OEM# 2003052



**Item# RP006**  
O2 Conc Hose Barb  
Fitting 1/8" Conc Port  
**Each**  
OEM# 1104151



**Item# RP019**  
Pressure Regulator  
**Each**



**Item# RP002**  
O Ring for Pilot Valve  
**25/pack**  
OEM# 1085112



**Item# RP009**  
Power Cord  
**Each**  
OEM# 1108764



**Item# RP031**  
Cabinet Clips  
**2/pack**



**Item# RP003**  
PE Valve  
**Each**  
OEM# 1143578



**Item# RP010**  
4-Way Valve  
Rebuild Gasket Kit  
**Each**



**Item# RP044**  
Filter Access Door  
**Each**  
OEM# 110766



**Item# RP045**  
Sieve Bed Filter Grid  
**Each**  
OEM# 1109520



**Item# RP004**  
On/Off Switch  
**Each**  
OEM# 1139950



**Item# RP014**  
Concentrator Capacitor  
15MFD  
**Each**  
OEM# 2003019



**Item# RP081**  
Piston Cup Seals for  
Thomas Compressor  
2650/2660  
**50/pack**



**Item# RP005**  
Hour Meter  
**Each**  
OEM# 2002939



**Item# RP016**  
Concentrator Capacitor  
30MFD \*For Platinum 10  
**Each**  
OEM# 1128740



**Item# RP083**  
Brass O2 Outlet  
**10/pack**



## Invacare HomeFill

**Item# RP004**  
Power Switch  
**Each**  
OEM# 2001134



**Item# RP013**  
Concentrator Capacitor  
10MFD  
**Each**



**Item# RP028**  
O-Ring for  
Homefill Tank  
**10/pack**



**Item# RP007**  
O2 Conc Fitting Receptacle Port  
**Each**  
OEM# 1104152



**Item# RP026**  
Cylinder Dust Cap  
**10/pack**



**Item# RP031**  
Cabinet Clips  
**2/pack**



**Item# RP008**  
2000psi Pressure Switch  
**Each**  
OEM# 1089467



**Item# RP027**  
Compressor Dust Cap  
**10/pack**  
OEM# 1156644



**Item# RP042**  
Male, QDV,  
High Pressure, INTL  
with Dust Cover  
**Each**



## Repair Parts NEW



### Invacare Perfecto2 / Perfecto2V

**Item# RP001**  
Pilot Valve and  
1 O-Ring  
**Each**  
OEM# 2003052



**Item# RP010**  
4-Way Valve Rebuild  
Gasket Kit  
**Each**



**Item# RP043**  
Cabinet Access Door  
\*for Perfecto2  
**Each**  
OEM# 1161831



**Item# RP002**  
O Ring for Pilot Valve  
**25/pack**  
OEM# 1085112



**Item# RP013**  
Concentrator Capacitor  
10MFD  
\*for Perfecto2  
**Each**



**Item# RP045**  
Sieve Bed Filter Grid  
\*for Perfecto2  
**Each**  
OEM# 1109520



**Item# RP003**  
PE Valve  
**Each**  
OEM# 1143578



**Item# RP015**  
Concentrator Capacitor  
22MFD  
\*for Perfecto2V  
**Each**



**Item# RP046**  
PVC F-Tube  
**Each**  
OEM# 1143446



**Item# RP004**  
On/Off Switch  
**Each**  
OEM# 1139950



**Item# RP018**  
Compressor  
Rebuild Kit  
\*for Perfecto2V  
**Each**  
OEM# 1191916



**Item# RP047**  
Control Panel  
**Each**  
OEM# 1143469



**Item# RP005**  
Hour Meter  
**Each**  
OEM# 2002939



**Item# RP019**  
Pressure Regulator  
**Each**



**Item# RP081**  
Piston Cup Seals for  
Thomas Compressor  
2650/2660  
**50/pack**



**Item# RP006**  
O2 Conc Hose Barb  
Fitting 1/8"  
Conc Port  
**Each**  
OEM# 1104151



**Item# RP031**  
Cabinet Clips  
**2/pack**



**Item# RP083**  
Brass O2 Outlet  
**10/pack**



**Item# RP009**  
Power Cord  
**Each**  
OEM# 1108764



**Item# RP036**  
Complete  
Cabinet Set  
**3/case**  
OEM# 1153341



**Item# RP101**  
Sticker Labels  
for Perfecto2 Shroud  
**10/pack**



**Item# RP102**  
Sticker Labels  
for Perfecto2V Shroud  
**10/pack**



# Repair Parts NEW



## Drive/DeVilbiss 525 / 1025

**Item# RP004**  
On/Off Switch  
**Each**  
OEM# 1139950



**Item# RP022**  
Exhaust Muffler  
**Each**  
OEM# 515A-705



**Item# RP048**  
Power Cord  
**Each**



**Item# RP018**  
Compressor  
Rebuild Kit  
*\*also compatible  
with 303DZ*  
**Each**  
OEM# 1191916



**Item# RP023**  
Circuit Breaker  
**Each**  
OEM# 525D-613



**Item# RP049**  
DS SSP Capacitor  
*\*for 1025DS*  
**Each**



**Item# RP020**  
Flow Meter  
**Each**  
OEM# 505DZ-607



**Item# RP024**  
Oxygen Outlet Port  
**5/case**  
OEM# 525DD-606



**Item# RP093**  
Rebuilt Drive  
Rotary Valve  
**12/case**



**Item# RP021**  
Cooling Fan  
**Each**  
OEM# 525DD-634



**Item# RP032**  
Pressure Regulator  
8.5 PSI  
**Each**



**Item# RP094**  
Rebuild Kit for  
Rotary Valve  
**Each**



**Item# RP095**  
Spacers for Rebuilt  
Rotary Valves  
**10/pack**



## Respironics EverFlo (Post-2013)

**Item# RP006**  
O2 Conc Hose  
Barb Fitting 1/8"  
Conc Port  
**Each**  
OEM# 1104151



**Item# RP035**  
Pressure Regulator  
Adapter Kit  
**Each**  
OEM# 1038824



**Item# RP012**  
Flow Meter  
**Each**  
OEM# 0209184



**Item# RP082**  
SMC Pilot and  
Body Valve  
**Each**



## Respironics Trilogy

**Item# RP025**  
O2 Quick Connect  
L-Shape Adapter  
**2/pack**



**Item# RP025S**  
O2 Quick Connect  
Straight Adapter  
**2/pack**



**RP025S is Also Compatible with:**

- Löwenstein **LUISA**
- ResMed **Astral**
- React Health **V\*Home**

# Repair Parts NEW



## Rhythm LM5BA & LM5CA, Lifestyle LM5A, Jumao JMC5A, Direct Supply Attendant 8X899, Medacure OC5L

**Item# RP004**  
On/Off Switch  
**Each**  
OEM# 1139950



**Item# RP040**  
Power Cord  
**Each**



**Item# RP085**  
Power Switch  
*\*for Rhythm LM5BA*  
**Each**



**Item# RP005**  
Hour Meter  
**Each**  
OEM# 2002939



**Item# RP041**  
Casters  
**4/pack**



**Item# RP086**  
Cooling Fan  
*\*for Rhythm LM5BA*  
**Each**



**Item# RP019**  
Pressure Regulator  
*\*For Rhythm LM5BA*  
**Each**



**Item# RP075**  
Cup Seal for Compressor Rebuild Kit  
**Each**



**Item# RP087**  
Solenoid Valve  
*\*for Rhythm LM5BA*  
**Each**



**Item# RP023**  
Circuit Breaker  
*\*for Rhythm LM5BA*  
**Each**



**Item# RP079**  
Pressure Regulator, 7 PSI  
**Each**



**Item# RP088**  
Oxygen Sensor  
*\*for Rhythm LM5BA*  
**Each**



**Item# RP037**  
Muffler  
**Each**



**Item# RP080**  
O-Rings for Pressure Regulator  
**50/pack**



**Item# RP089**  
Circuit Board  
*\*for Rhythm LM5BA*  
**Each**



**Item# RP038**  
PE Valve  
**Each**



**Item# RP084**  
Power Cord  
*\*for Rhythm LM5BA*  
**Each**



**Item# RP090**  
Display Screen  
*\*for Rhythm LM5BA*  
**Each**



**Item# RP039**  
Control Valve  
**Each**



**Item# RP096**  
Cabinet with Door for Direct Supply Attendant  
**2/box**



## React Stratus 5

**Item# RP071**  
Compressor  
**Each**



**Item# RP074**  
Cooling Fan  
**Each**



**Item# RP072**  
Compressor Seal - Piston Ring  
**Each**



**Item# RP073**  
Regulator Valve  
**Each**



**Item# RP076**  
Power Switch  
**25/pack**



**Item# RP077**  
Pneumatic Valve  
**Each**



# Liquid Oxygen Repair Parts NEW



## Caire Companion C1000 / C1000T

**Item# RP050**  
Cover Side B  
**Each**



**Item# RP054**  
Top Bezel  
\*for C1000  
**Each**  
OEM# B-775003-00



**Item# RP062**  
Manifold O-Ring  
**Each**  
OEM# 775037



**Item# RP051**  
Side Cover Decal for Companion T  
**10/box**  
OEM# B-775499-00



**Item# RP055**  
Cover Side A  
\*for C1000  
**Each**  
OEM# B-775001-00



**Item# RP063**  
Vent / Fill Assembly  
\*for C1000T  
**Each**  
OEM# 775037



**Item# RP053**  
Moisture Control Pad  
\*for C1000  
**Each**  
OEM# B-775092-00



**Item# RP058**  
Tube Fill  
**Each**



**Item# RP065**  
Lens Cap  
**Each**



**Item# RP052**  
Carrying Strap  
**Each**  
OEM# B-775097-00



**Item# RP060**  
Quick Connect Adapter  
**Each**  
OEM# B-775263



**Item# RP067**  
Portable Rubber Foot  
**10 Sheets of 56**



## Caire HELIOS

**Item# RP056M**  
Lip Seal - Metal  
**Each**  
OEM# B-775260-00



**Item# RP056R**  
Lip Seal - Rubber  
**Each**  
OEM# B-775260-00



## Puritan Bennett 41L Base Liquid Oxygen Vessel

**Item# RP057**  
Rebuilt Contents Gauge  
**Each**  
OEM# B-775847-SV



**Item# RP066**  
Contents Gauge Tubing  
**2/pack**  
OEM# B-775856-SV, PB21, PB41, H46



*We're always adding new repair parts!*

[Click here to see our full lineup](#)

# Suction and Tracheostomy

Suction	Page
Suction Machine	70
Suction Canister and Kits	70-71
Open Suction Catheters	72
Open Suction Catheter Kits	73
Suction Filters	74
Suction Adapters and Connectors	74
Yankauers	74
Tracheostomy Supplies	75
Heat and Moisture Exchangers	76
Closed Suction Catheters	77-78

## Suction

### Suction Machine

Item#  
**SU100DC**  
With Internal  
Rechargeable  
Battery

**SU100AC**  
Without Battery

**SU100DPA**  
Includes DC Mobile  
Power Adapter  
and Internal  
Rechargeable Battery



- Effective removal of respiratory fluids
- Durable yet lightweight
- Variable vacuum control
- Quiet and easy to use
- Two-year warranty

**Comes with**  
Machine  
Carry Bag  
800cc Suction Canister  
AC Adapter & Power Cord  
6ft Suction Tubing  
3.5in Suction Tube Connector  
4.5in Suction Tube Connector  
BF1438S Bacteria Filter

#### Accessories and replacement parts

<b>SU100DC-PA</b>	DC Mobile Power Adapter (for car power outlet)
<b>SU100DC-Adapt</b>	Replacement Power Adapter (without cord)
<b>SU100DC-Cord</b>	Replacement Power Cord (without adapter)
<b>SU100DC-BAG</b>	Carrying/Storage Bag
<b>SU100DC-BATTERY</b>	Internal rechargeable battery replacement

### Suction Canister

Item#  
**RES023AW**  
800cc • Single Use  
**48/case**  
HCPCS# A7000

- Float valve shutoff



### Covers for Suction Machines



20×24"  
2mil  
**500 bags/roll**

Item#	Color
<b>EC2024C</b>	Clear
<b>EC2024T</b>	Tan Tint

# Suction

## Suction Kit for SU100DC

**Item# RES026S-SS**

- RES023AW Suction Canister
- RES024 Suction Tubing Connector 3.5in
- RES024M Suction Tubing Connector 4.5in
- BF1438S Bacteria Filter
- RES025 Suction Tubing 1/4in x 6ft

HCPCS#  
CANISTER - A7000  
TUBING - A7002



## Suction Filter Component Kit for SU100DC

**Item# RES024-MKIT**

- RES024 Suction Tubing Connector 3.5in
- RES024M Suction Tubing Connector 4.5in
- BF1438S Bacteria Filter



## Suction Kit

**Item# RES026S-MULTI**

- RES023AW Suction Canister 800cc
- RES024XL Suction Tube Connector 16in
- RES025 Suction Tubing 6ft x 1/4in
- BF516TS Bacteria Filter
- RES024S Suction Tube Connector 1.5in



## Suction Component Kit

**Item# RES024-MULTI**

- RES024S Connector Tube 1.5in
- RES024XL Connector Tube 16in
- BF516TS Bacteria Filter



**RES026S-MULTI and RES024-MULTI are compatible with:**

**Precision Medical**

Easy Vac	PM60
Power Vac	PM61
Easy Go	PM66 series, PM65 PM65E, PM65HG

**Compass/Roscoe**

Roscoe Portable Suction	50004 (AC model) 50006 (DC model)
Cliq Single Cylinder Aspirator	ASP-ROS

**Medline**

Heavy Duty	HCS7000
------------	---------

**Invacare**

Aspirator	IRC1135
Durable Suction Unit	IRC1136

**Drive/DeVilbiss**

Heavy-Duty Suction Machine	18600
VacuMax™ Go Suction Machine	18615

**Allied Healthcare**

Schuco Portable Aspirator with Sturdy Legs	S430A
Schuco Portable Suction Unit with Metal Base	S130A
Schuco Portable Aspirator with Molded Plastic Base	S330A
Gomco OptiVac® AC/DC Portable Suction Unit & 305 Aspirator	G180 G180CE G180CE-JAR G180-JAR 01-22-3005 01-22-0305
Gomco® Tabletop Aspirator Model 300	01-22-0300
Gomco® General Use Aspirators Models 3020 & 3040	01-62-3020 01-62-3040
Gomco® 405 Aspirator	01-12-0405 01-12-4005 01-23-0405
Gomco® Stand-Mounted Models 4032 & 4042	01-52-4032 01-52-4042
Gomco Table Top Gastric Aspirator 270	01-22-0270

**RES026S-SS and RES024-MKIT are compatible with:**

**Drive/DeVilbiss**

Vacu-Aide® QSU Quiet Suction Unit	7314P-D 7314P-D-EXF 7314D-D
Vacu-Aide® Suction Unit	7305D-D 7305P-D-EXF 7305D-D-EXF 7305P-D

**Sunset**

Suction Machine (Portable Unit)	SU100DC
Suction Machine (Stationary Unit)	SU100AC

## Open Suction Catheters

50/case

Size	Length	Valve	Tip	Intended Use	Packaged	Item#
6fr	44 cm	Chimney	DeLee	Pediatric	Coiled	<b>TRASC6-CC</b>
6fr	44 cm	Chimney	DeLee	Pediatric	Straight	<b>TRASC6-CS</b>
8fr	44 cm	Chimney	DeLee	Pediatric	Coiled	<b>TRASC8-CC</b>
8fr	44 cm	Chimney	DeLee	Pediatric	Straight	<b>TRASC8-CS</b>
8fr	44 cm	Directional	DeLee	Pediatric	Straight	<b>TRASC8-DS</b>
10fr	58 cm	Chimney	Whistle	Adult	Coiled	<b>TRASC10-CC</b>
10fr	58 cm	Chimney	Whistle	Adult	Straight	<b>TRASC10-CS</b>
10fr	58 cm	Directional	Whistle	Adult	Straight	<b>TRASC10-DS</b>
12fr	58 cm	Chimney	Whistle	Adult	Coiled	<b>TRASC12-CC</b>
12fr	58 cm	Chimney	Whistle	Adult	Straight	<b>TRASC12-CS</b>
12fr	58 cm	Directional	Whistle	Adult	Straight	<b>TRASC12-DS</b>
14fr	58 cm	Chimney	Whistle	Adult	Coiled	<b>TRASC14-CC</b>
14fr	58 cm	Chimney	Whistle	Adult	Straight	<b>TRASC14-CS</b>
16fr	58 cm	Chimney	Whistle	Adult	Straight	<b>TRASC16-CS</b>
18fr	58 cm	Chimney	Whistle	Adult	Straight	<b>TRASC18-CS</b>

- Sterile
- Single patient use
- Latex-free

**NEW**




Coiled Packed

Straight Packed




Chimney Valve



Directional Valve

Whistle Tip



DeLee Tip

### Open Suction Catheter Kits with Cup

50/case

Size	Length	Valve Type	Tip	Intended Use	Item#
5fr	44 cm	Chimney	DeLee	Pediatric	<b>TRAKIT5-C</b>
6fr	44 cm	Chimney	DeLee	Pediatric	<b>TRAKIT6-C</b>
8fr	44 cm	Directional	DeLee	Pediatric	<b>TRAKIT8-D</b>
8fr	44 cm	Chimney	DeLee	Pediatric	<b>TRAKIT8-C</b>
10fr	58 cm	Directional	Whistle	Adult	<b>TRAKIT10-D</b>
10fr	58 cm	Chimney	Whistle	Adult	<b>TRAKIT10-C</b>
12fr	58 cm	Directional	Whistle	Adult	<b>TRAKIT12-D</b>
12fr	58 cm	Chimney	Whistle	Adult	<b>TRAKIT12-C</b>
14fr	58 cm	Directional	Whistle	Adult	<b>TRAKIT14-D</b>
14fr	58 cm	Chimney	Whistle	Adult	<b>TRAKIT14-C</b>
16fr	58 cm	Chimney	Whistle	Adult	<b>TRAKIT16-C</b>



- Sterile
- Includes Catheter and Gloves

### Sleeved Coiled Open Suction Catheter Kits with Cup

50/case


Size	Length	Valve Type	Tip	Intended Use	Item#
8fr	44 cm	Directional	DeLee	Pediatric	<b>TRAKIT8-SL-NG</b>
10fr	58 cm	Directional	Whistle	Adult	<b>TRAKIT10-SL-NG</b>
12fr	58 cm	Directional	Whistle	Adult	<b>TRAKIT12-SL-NG</b>
14fr	58 cm	Directional	Whistle	Adult	<b>TRAKIT14-SL-NG</b>



## Suction

### Suction Machine Inline Filters

Item# **BF1412**  
Final Bacteria Filter  
1/4in to 1/2in  
Stepped Barb on  
Both Sides  
**10/pack**



Item# **BF1438S**  
1/4in to 3/8in  
**10/pack**



Item# **BF3812**  
3/8in to 1/2in  
**10/pack**



Item# **BF1312**  
1/3in to 1/2in  
**1/pack**



Item# **BF516TS**  
5/16in Barbed with  
Threaded NPT  
**10/pack**



### Suction Tubing

#### 1/4in

Length	Item#	Qty
18in	<b>RES025S</b>	<b>50/case</b>
6ft	<b>RES025</b>	<b>10/case</b>



#### 3/16in

Length	Item#	Qty
18in	<b>RES025-31618</b>	<b>100/case</b>
20in	<b>RES025-31620</b>	<b>100/case</b>
6ft	<b>RES025-31672</b>	<b>50/case</b>

### Suction Tube Connectors

#### 10/pack

Size	Item#
1.5in	<b>RES024S</b>
3.5in	<b>RES024</b>
4.5in	<b>RES024M</b>
6in	<b>RES024L</b>
10in	<b>RES024L10</b>
16in	<b>RES024XL</b>



### Rubber Suction Elbow Connector Drive DeVilbiss Style

Item#  
**RES028**  
**6/pack**



### Suction Elbow Connector Precision Medical Style

Item#  
**RES033**  
**10/pack**



### Bulb Tip Yankauers

Item# **TRAA001S**  
Vented • **50/case**

Item# **TRAA002S**  
Non-Vented • **50/case**

HCPCS# A4628



- Sterile
- Single patient use
- Not made with natural latex rubber

# Tracheostomy

## Tracheostomy Mask

Item#  
**RES2130**  
Adult Size  
**50/case**  
HCPCS# A7525



**RES2230**  
Pediatric Size  
**50/case**



## Flexible Connector

**50/case**

Item#  
**TRAA004** Adult  
22/15mm x 15mm  
Length 4.9-7.5cm  
OEM# 3222



**TRAA004P** Pediatric  
15mm F x 15mm M  
Length 4.9-7.5cm  
OEM# 3215



## Adapter for T-HME

Item# **RES027A**  
**Each**



## Flexible Trach Adapter

Item#  
**TRAA003**  
22mm Cuff  
**Each**

OEM# 1073902, 1147898



**TRAA003P** Pediatric  
15mm Cuff  
OEM# 1073863



## Bacterial/Viral HME Filter

Item#  
**BF104**  
**Each**

With hydrophobic filter media

**Bacterial Filtration Efficiency (BFE)**  
99.99+%

**Viral Filtration Efficiency (VFE)**  
99.99+%

**Flow Resistance** 0.09 (30L/min kpa)

**Dead Space** 38ml

**Weight** 30.6g

**Tidal Volume** 150-1500ml

**Moisture Output** 35.1 mg/L




## Tracheostomy Ties with Fastener Tabs

**20/case**

One-piece • Easy on

Item#  
**RES029P-11**  
1 in x 11.5 in  
**RES029-17**  
1 in x 17 in



**Softer, more  
breathable  
100% cotton  
patient-facing  
material in 2026**

Two-piece • Easy off

Item#  
**RES030P-10**  
0.78 in x 10 in  
**RES030P-14**  
0.78 in x 14 in  
**RES030P-18**  
0.78 in x 18 in  
**RES030-20**  
1 in x 20 in



Sunset's trach ties deliver value and comfort, with soft straps and light Velcro closures for customization. Secure, but flexible for patient movement, the ties help prevent sliding, skin irritation and accidental extubation. 60% Polyester, 20% Cotton, 20% Nylon.

# Heat and Moisture Exchangers

## Tracheostomy Heat Moisture Exchangers

30/case\*

\*RES027UOP 50/case



**T-HME  
RES027**



**T-HME with O2 Port  
RES027OP**



**HME with O2 Port  
RES027UOP**

<b>TV Range (mL)</b>	150-1250	150-1250	50-1000
<b>Moisture Output (mg H<sub>2</sub>O/L)</b>	36.4 at VT = 500 ml	32 mg/L	13.2 mg/L
<b>Weight (g)</b>	6.6	7.4	5
<b>Resistance (cm H<sub>2</sub>O)</b>	0.2 @ 30 LPM	0.2 @ 30 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM
<b>Dead Space (mL)</b>	15	7	17
<b>Connection (mm)</b>	15F	15F 5.5mm O2 port	15F

## Inline Heat Moisture Exchangers

30/case

\*RES027FLEX 40/case



**Flex HME — Foam  
RES027FLEX**



**Inline Small Volume  
HME with Gas Sampling  
Port  
RES027IL**




**Inline Large Volume HME  
with Gas Sampling Port  
RES027ILLV**



**Inline HME —  
Neonatal  
RES027NE**

<b>TV Range (mL)</b>	250-1500	70-600	150-1500	5-50
<b>Moisture Output (mg H<sub>2</sub>O/L)</b>	30 mg/L	34.5 at VT = 250 ml 32.5 at VT = 500 ml	37.7 at VT = 250 ml 35.5 at VT = 500 ml 33.5 at VT = 1000 ml	32.8 at VT = 250 ml
<b>Weight (g)</b>	15.5	10	15	1.9
<b>Resistance (cm H<sub>2</sub>O)</b>	0.5 @ 30 LPM 1.4 @ 60 LPM	0.8 @ 30 LPM 2.1 @ 60 LPM	1.0 @ 30 LPM 2.6 @ 60 LPM	0.8 @ 15 LPM 1.7 @ 30 LPM
<b>Dead Space (mL)</b>	90	21	34	2.2
<b>Connection (mm)</b>	22F, 15F	22M/15F-15M	22M/15F-15M	15M, 15F

## Standard Closed Ventilation Suction System



# Standard Closed Ventilation Suction System

## 24 Hour Adult

### Endotracheal 540mm

Item#	Description
VEN46001000A	10Fr, Flex Tube
VEN46001000AM	10Fr, MDI Port
VEN46001200	12Fr, Flex Tube
VEN46001200M	12Fr, MDI Port
VEN46001400	14Fr, Flex Tube
VEN46001400M	14Fr, MDI Port
VEN46001600	16Fr
VEN46001600M	16Fr, MDI Port

### Tracheostomy 305mm

Item#	Description
VEN46011000A	10Fr
VEN46011000AM	10Fr, MDI Port
VEN46011200	12Fr, Flex Tube
VEN46011200M	12Fr, MDI Port
VEN46011400	14Fr, Flex Tube
VEN46011400M	14Fr, MDI Port, Flex Tube
VEN46011600	16Fr
VEN46011600M	16Fr, MDI Port

## 72 Hour Adult

### Endotracheal 540mm

Item#	Description
VEN56051000A	10Fr
VEN56051000AM	10Fr, MDI Port
VEN56051200	12Fr
VEN56051200M	12Fr, MDI Port
VEN56051400	14Fr
VEN56051400M	14Fr, MDI Port
VEN56051600	16Fr
VEN56051600M	16Fr, MDI Port

### Tracheostomy 305mm

Item#	Description
VEN56061000A	10Fr, Flex Tube
VEN56061000AM	10Fr, MDI Port
VEN56061200	12Fr, Flex Tube
VEN56061200M	12Fr, MDI Port
VEN56061400	14Fr, Flex Tube
VEN56061400M	14Fr, MDI Port
VEN56061600	16Fr
VEN56061600M	16Fr, MDI Port



## Neonatal/Pediatric

### With DSE Connection

Item#	Description
VENCSS6-305	6Fr, 305mm
VENCSS8-203	8Fr, 203mm
VENCSS8-355	8Fr, 355mm
VENCSS10-203	10Fr, 203mm
VENCSS10-405	10Fr, 405mm



# Ventilator

## Ventilator

## Page

LUISA Ventilator	79
Ventilator Filters	80-81
Adapters and Connectors	81
Ventilator Mask	81
Breathing Circuits	82-85
Heated Single-Limb	82
Non-Heated Single-Limb	83-84
Dual-Limb	85

## Ventilator

### Movair / Löwenstein LUISA®

Item# **VNT31390**  
Each

**NEW**

LUISA with High Flow Therapy provides long term care with easy serviceability. Routine patient maintenance consisting of only filter changes and a device wipe down provides uncomplicated upkeep.

With simple, routine maintenance and no adapters required between circuits, the LUISA provides clinicians and patients' peace of mind with an all-encompassing system in a compact, home ventilator with the added benefit of High Flow Therapy.

#### Includes:

- Exhalation Module
- Leakage Circuit
- Power Supply/Power Cord
- O2 Nozzles
- Fine and Course Filters
- Protection Bag
- Accessories Bag
- USB Drive



- High Flow Therapy: all circuit systems
- Connectivity: customized patient download reports
- Uncomplicated: no adapters, intuitive operation
- Battery run time: up to 18 hours
- Monitoring: FiO2 and SpO2 measurement
- Safety: alarm displayed in two languages simultaneously
- Display: 10" angled touchscreen
- Functionality: 100 - 3,000 mL Vt without additional flow sensor
- Responsive: adjustable target volume speed
- Lightweight: ~ 8 lbs.
- Positioning: vertical or horizontal; Circuit orientation options (left or right side)

#### Accessories and Replacement Parts

<b>VNT31404</b>	Exhalation Module - Disposable
<b>VNT31569</b>	Power Supply
<b>VNT27027</b>	Power Cord
<b>VNT30669</b>	Oxygen Connection Nozzle
<b>VF29652</b>	Set of 12 Pollen/Fine Filters - 12/PK
<b>VF29928</b>	Set of 2 Air/Coarse Dust Filters - 2/PK
<b>VNT31417</b>	Protective Bag
<b>VNT31440</b>	Accessories Bag
<b>VNT31414</b>	USB Flash Drive

## ResMed

### Filter for Astral Ventilator

Item# **VF003-4**  
4/pack

OEM# 27939



NEW

4-Filter Kit with  
Filter Housing

Item# **VF003-KIT**



Item# **VF003-H**  
Filter Housing  
Each



### 02 Quick Connect Straight Adapter 1/4" Barb

Item# **RP025S**

1/4" Barb, for use with Astral  
2/pack



NEW

## Respironics

### Air Inlet Vent Filter for V60

Item# **VF002**

For use with Philips Respironics V60  
5/pack



### Filters for Trilogy/Trilogy EV300 Ventilators

Item# **VF005-10**  
Particulate Filter for  
Trilogy Ventilator  
10/pack



NEW

Item# **VF007-1**  
Reusable Filter for Trilogy EV300  
1/pack

**VF007-10** 10/pack  
OEM# 1134591



Item# **CF1009F**  
Reusable Filter for Trilogy  
2/pack

OEM# 1029330



### 02 Quick Connect Straight Adapter 1/4" Barb

Item# **RP025S**

1/4" Barb, for use with Trilogy  
2/pack



NEW

## Breas

### Vivo 45/45LS Style Air Inlet Filters

Item# **VF010-5**

Disposable filter for Vivo 45LS  
5/pack

OEM# 007103



NEW

Item# **VF010F-5**

Reusable Filter for Vivo 45LS  
5/pack

OEM# 007104



NEW

Item# **VF011F-5**

Reusable filter for Vivo 45  
5/pack

OEM# 007105



NEW

## Movair / Löwenstein

### Filters for LUISA

Item# **VF008-12**  
Disposable filter  
12/pack

OEM# WM29652



**3-Filter Kit**

1 VF008 Disposable Filter  
1 VF008F Reusable Filter  
1 VF014F Gross Particle Filter  
Item# **VF008-KIT**

Item# **VF008F-2**

Reusable Filter  
2/pack

OEM# WM29928



Item# **VF014F**  
Gross Particle Filter  
Each



NEW

### 02 Quick Connect Straight Adapter 1/4" Barb

Item# **RP025S**

1/4" Barb, for use with LUISA  
2/pack



NEW

## Coviden/Nellcor

### Ventilator Felt Filter

Item# **VF006**

For use with NPB LP6 & LP10  
10/pack

OEM# L-002917-000



## Newport Medical

### Proximal Inline Filter for Newport HT50

Item# **VF009-5**  
5/pack



## Universal filter

### Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item# **BF2200**  
22mm  
5/pack



Attach between the machine and tube for additional protection against bacteria and allergens

## Ventec VOSCN

### Filters for Ventec VOCSN

Item# **VF015-10**  
Internal Bacteria Filter  
10/pack  
OEM# PRT-00976-000



Item# **VF016-10**  
External Bacteria Filter  
10/pack  
OEM# PRT-00850-000



## Ventilator Adapters, Connectors, and Mask

### Valved T-Adapter — Adult

Item# **VEN001**  
Inline adapter for nebulizer treatment  
22F-18F-22M/15F  
Each



### Circuit Connectors

Item# **RES013**  
15mm ID x 22mm OD  
10/pack



Item# **RES014**  
22mm ID x 22mm OD  
10/pack



### Exhalation Ports For Ventilator Circuits & CPAP

Item# **VEN002**  
Swivel  
Each



Item# **VEN003**  
Fixed  
22 mm OD • 15 mm ID  
10/pack  
OEM# 312149



### Straight Connectors

Item# **VEN005**  
22mm ID, 15mm OD  
Each



Item# **VEN006**  
22mm ID, 22mm OD  
Each



### Active Exhalation Valve Adapter

Item# **VEN004**  
Each  
OEM# 1065659



**VEN004-R**  
Compatible with ResMed

### NIV Mask with Headgear — Non-vented

Each

Item#	Size
<b>CM038M</b>	Medium
<b>CM038L</b>	Large



## Heated Single Limb Ventilator Circuits without Humidifier Chamber

### Active Heated Breathing Circuits

#### 30/case

For use with MR850 Humidifier

Item#

**VC1002** Active Adult

**With Exhalation Valve Adapter**

NEW

#### Includes:

- 22mm circuit
- Active exhalation valve adapter
- Optional connector

OEM# 1073222, 1073231



**VC1002-R**  
Compatible with ResMed

Item#

**VC1002P** Active Pediatric

**With Exhalation Valve Adapter**



NEW

#### Includes:

- 15mm circuit
- Active exhalation valve adapter
- 15mm/22mm connector

OEM# 1073221, 1073230



**VC1002P-R**  
Compatible with ResMed

### Passive Heated Breathing Circuits

#### 30/case

For use with MR850 Humidifier

Item#

**VC1003** Adult

**With Swivel Exhalation Port**

#### Includes:

- 22mm circuit
- Swivel exhalation port
- Optional connector

OEM# 1076046,  
1073227, 1073235



### Adult Single Limb Heated Circuit

#### 10/case

For use with MR850 Humidifier

Item#

**VC1114** Adult

OEM# RT114



Item#

**VC1003P** Pediatric

**With Swivel Exhalation Port**



NEW

#### Includes:

- 15mm circuit
- Swivel exhalation port
- 15mm/22mm connector
- 22mm/22mm connector

OEM# 1073225, 1073234



## Passive Ventilator Circuits available from Sunset

Item#	Heating	Exhalation Port	Filter	Water Trap	Size	Tube Diameter	OEM#
<b>VC1000</b>	Non-Heated	Swivel	None	No	Adult	22mm	1073228 & 1073236
<b>VC1003</b>	Heated	Swivel	None	No	Adult	22mm	1076046, 1073227, 1073235
<b>VC1003P</b>	Heated	Swivel	None	No	Pediatric	15mm	1073225, 1073234
<b>VC1004</b>	Non-Heated	Fixed	None	Yes	Adult	22mm	1069443 & 1069444
<b>VC1004P</b>	Non-Heated	Swivel	None	Yes	Pediatric	15mm	1073224 & 1073233
<b>VC1005</b>	Non-Heated	Swivel	Bacterial/Viral	No	Adult	19mm	1147378, 1132340, 7615
<b>VC1005P</b>	Non-Heated	Swivel	Bacterial/Viral	No	Pediatric	15mm	1132339
<b>VC1006</b>	Non-Heated	None	Bacterial/Viral	No	Adult	22mm	
<b>VC1007</b>	Non-Heated	Fixed	Bacterial/Viral	No	Adult	22mm	
<b>VC1114</b>	Heated	None	None	No	Adult	22mm	RT114
<b>VC1100</b>	Non-Heated	Fixed	None	No	Adult	22mm	1069467 & 1069468
<b>VC1160</b>	Non-Heated	Fixed	None	No	Adult	22mm	1069210

## Non-Heated Single Limb Ventilator Circuits without Humidifier Chamber

### Passive Non-Heated Breathing Circuit

Item#  
**VC1160**  
Compatible with V60,  
No Filters  
OEM# 1069210



### Active Non-Heated Breathing Circuits

30/case  
22mm cuffs



Item#  
**VC1001** Adult  
OEM# 1073223, 1073232

**VC1001-R**  
Compatible with ResMed




Item#  
**VC1001P** Pediatric  
OEM# 1073220, 1073229

**VC1001P-R**  
Compatible with ResMed




To use active circuits with Respiration Trilogy, simply cut ends off pressure and exhalation lines.




# Non-Heated Single Limb Ventilator Circuits without Humidifier Chamber

## Passive Non-Heated Breathing Circuits with Water Trap

Item# **VC1004** Adult  
 OEM# 1069443, 1069444



Item# **VC1004P** Pediatric  
 OEM# 1073224, 1073233



## Passive Non-Heated Breathing Circuit

For use with Respironics Trilogy  
**Each**

Item# **VC1000**  
 with Swivel Exhalation Port  
 OEM# 1073228 & 1073236

Item# **VC1100**  
 with Fixed Exhalation Port  
 OEM# 1069467 & 1069468



## Passive Circuit with Exhalation Port and Bacterial/Viral Filter

With 22mm cuffs  
**Each**

Item# **VC1005**  
 Adult • 19mm Tube Diameter  
 OEM# 1147378, 1132340, 7615

Item# **VC1005P**  
 Pediatric • 15mm Tube Diameter  
 OEM# 1132339



## Adult Passive Circuit with Bacteria Filter

Item# **VC1006**  
 Comparable to Breas 007936  
**Each**



## Adult Passive Circuit with Bacteria Filter and Fixed Exhalation Port

Item# **VC1007**  
 Comparable to Breas 007615  
**Each**



## Dual Limb Ventilator Circuits

**NEW**


### Adult Dual Limb Active Circuit 6ft

Item#

**VC1009**

For use with Vyair LTV

OEM# 29657-001



### Pediatric Dual Limb Active Circuit 6ft

Item#

**VC1010**

For use with Vyair LTV

OEM# 29658-001




### Pediatric Dual Limb Active Circuit Water Trap

Item#

**VC1011**

For use with Vyair LTV

OEM# 29662-001




### Adult Dual Limb Active Circuit with Water Trap

Item#

**VC1012**

For use with Vyair LTV

OEM# 29663-001



# Other

Other	Page
Disinfectants	86
Cough Assist Kits	87
Breathing Masks	87
Aerosol	88
Tabletop Nebulizers	89
Handheld Nebulizers	90
Nebulizer Kits	91
Nebulizer and Anesthesia Filters	92
StrongArm Mobility and Comfort	93-94
Bath Safety	95
Maternity	96-97
Equipment Covers	98-99
Blood Testing	99

## Other HCPCS

A6402	Gauze - sterile, pad size 16 sq. in. or less, without adhesive border, each dressing
A4483	Heat and Moisture Exchanger — Disposable
A4626	Tracheostomy Cleaning Brush
A4628	Oropharyngeal Suction Catheter (Yankauer)
A4629	Tracheostomy Care Kit - for established tracheostomy
A7000	Suction Canister — Disposable
A7002	Tubing for Suction Pump
A7003	Nebulizer Set - Small Volume—Disposable
A7005	Nebulizer Set - Small Volume—Nondisposable

A7013	Nebulizer Compressor Filter— Disposable
A7015	Aerosol Mask used with Nebulizer
A7525	Tracheostomy Mask
E0570	Nebulizer with Compressor
E0600	Respiratory Suction Pump - home model, portable or stationary, electronic
E2601	General Use Wheelchair Seat Cushion - width less than 22 inches, any depth
K1005	Disposable collection and storage bag for breast milk, any size, any type, each

## Disinfectants

### Simpli Germicidal Disposable Wipes

Item#  
**RES6108W**  
Canister of 160 wipes

- Bactericidal, Tuberculocidal, Fungicidal, Virucidal\*
- Premoistened wipes containing a quaternary/alcohol-based solution.
- Made in the USA
- Kills \*SARS-CoV-2, Covid 19, Mycobacterium bovis - BCG (Tuberculosis—TB) in 1 Minute
- Kills Candida Auris in 2 Minutes



### CleanSmart Disinfectant

Item#  
**CAP1019-C**  
23oz Spray Bottle  
6/case

- CleanSmart is EPA listed as a broad-spectrum and hospital disinfectant
- EPA Registration Number 89896-2
- CleanSmart is formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment
- CleanSmart is included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2, the virus that causes COVID-19
- Active ingredient Hypochlorous Acid does not leave a chemical residue



## Cough Assist Kits

**Cough Assist Kits** Compatible with ABM Respiratory Care BiWaze Cough System and Philips Respironics Cough Assist T70

**All kits come with 1 BF2200 Bacterial/Viral Filter and 1 RES013 Circuit Connector**

Item#

**CAKIT001**

With 1 TUB06 6ft Tube  
**Each**

OEM# 1090835



Item#

**CAKIT002**

With 1 TUB06 6ft Tube and  
1 RES019 Oxygen Enrichment Adapter  
**Each**



**CAKIT003**

With 1 TUB06 6ft Tube  
and 1 TRAA003 Flexible  
Trach Adapter  
**Each**



### ABM Respiratory Care BiWaze Cough System



Item# **OF8003-3**  
Air Inlet Foam Filter  
**3/pack**

### Kits With Tube and SBC01 Series Breathing Mask

With 1 TUB06 6ft Tube

Mask Size	Item#	Comparable to OEM#
Infant	<b>CAKIT004I</b>	1090830
Toddler	<b>CAKIT004T</b>	1090831
Adult Small	<b>CAKIT004S</b>	1090832
Adult Medium	<b>CAKIT004M</b>	1090833
Adult Large	<b>CAKIT004L</b>	1090834



With 1 TUB04 4ft Tube

Mask Size	Item#	Comparable to OEM#
Infant	<b>CAKIT005I</b>	325-9237
Toddler	<b>CAKIT005T</b>	325-9236
Adult Small	<b>CAKIT005S</b>	325-9235
Adult Medium	<b>CAKIT005M</b>	325-9234
Adult Large	<b>CAKIT005L</b>	325-9217



**CAKIT006**

with 1 CAMP001 Mouthpiece  
and 1 TUB06 6ft Tube  
**Each**

OEM# 1090836



## Breathing Mask

### Air Cushion Breathing Mask with Valve - No Hooks

**20/case**

Mask Size	Item#
Infant	<b>SBC01I</b>
Toddler	<b>SBC01T</b>
Adult Small	<b>SBC01S</b>
Adult Medium	<b>SBC01M</b>
Adult Large	<b>SBC01L</b>



- Ergonomic cushion design provides an optimal seal with minimal pressure
- Clear mask permits easy patient observation
- Available in a variety of sizes for the best possible fit
- Inflation port for optimized seal

## Aerosol

### Aerosol Mask

Item# **RES2110**  
Adult Size  
50/case

Item# **RES2210**  
Pediatric Size  
50/case

HCPCS# A7015

- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



### Peak Flow Meter

Item# **RES052**  
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**  
Pediatric Mouthpiece  
for Peak Flow Meter  
10/case



### Drainage Bag with T-Adapter

Item# **RES070**  
750cc • 50/case



### Drainage Bag with Y-Adapter

Item# **RES071**  
750cc • 50/case



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

### SperaChamber

Anti-Static Valved Inhaler Chamber  
for Metered Dose Inhaler

Chamber Only

Item#  
**RES195-10** 10/pack  
**RES195-50** 50/pack

Chamber with Mask  
10/pack

Item#  
**RES195MS** Small (Infant)  
**RES195MM** Medium (Pediatric)

**RES195ML** Large (Adult)

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free
- Dog mask appeals to kids



Mask Only

Size	Item#
Small	<b>RES195S</b>
Medium	<b>RES195M</b>
Large	<b>RES195L</b>



# Tabletop Compressor Nebulizers

## Compressor Nebulizers

Item#

### NEB100

- Disposable Nebulizer Kit
- Foam Filters

### NEB100-R

- Disposable Nebulizer Kit
- Reusable Nebulizer Kit
- Foam Filters x 5



**NEB100-BAG**  
Carry Bag

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

### Sunny the Bear

#### NEB300BEAR

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Foam Filters x 5

#### NEB300BEAR-R

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Reusable T-Piece Kit
- Foam Filters x 5



**NEB300-BAG**  
**Bag for Nebulizer Compressor**  
To include bag in configuration, just add "B " to Item Number:

- NEB300FISH-B**
- NEB300FISH-RB**
- NEB300BEAR-B**
- NEB300BEAR-RB**



### Gordon the Fish

#### NEB300FISH

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Foam Filters x 5

#### NEB300FISH-R

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Reusable T-Piece Kit
- Foam Filters x 5



## Covers for Nebulizers



20x24"  
2mil  
**500 bags/roll**

Item#	Color
<b>EC2024C</b>	Clear
<b>EC2024T</b>	Tan Tint



## Handheld Nebulizers

### Handheld Mesh Nebulizer

Item#  
**NEB402**

**With internal rechargeable battery**

**Included**

Medication cup  
USB cord  
Mask  
Mouthpiece  
Portable bag  
Manual



**Particle Size** MMAD < 5µm  
**Nebulization Rate** 0.25ml/min-0.9ml/min  
**Medication Cup Capacity** 10ml (Max)  
**Device Weight** Approx. 1/4 lb  
**Device Size** Approx. 1<sup>1</sup>/<sub>2</sub> x 2 x 3<sup>7</sup>/<sub>8</sub> inches

### Handheld Compressor Nebulizer

Item#  
**NEB200**

**Included**

Compressor  
Medication cup  
AC power adapter  
Air hose  
Elbow  
Mouthpiece  
Mask  
4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (Approximately 3 microns)
- No replacement filter required



## Nebulizer Kits

### Disposable Nebulizer Kit with T-Piece

Item# **RES091**  
 Disposable  
 Single Patient  
**50/case**  
 HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
  - T-piece with mouthpiece and corrugated reservoir tube
  - 7ft supply tube



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

### Reusable Nebulizer Kit with T-Piece

Item# **RES093**  
 Reusable  
 Single Patient  
**10/case**  
 HCPCS# A7005

- Includes:**
- 10cc jet nebulizer bottle
  - T-piece with mouthpiece and corrugated reservoir tube
  - 7ft supply tube



- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°

### Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**  
 Disposable  
 Single Patient  
**50/case**  
 HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
  - Angled mouthpiece
  - 7ft supply tube



### Disposable Nebulizer Kit with Mask

Item# **RES090**  
 Adult Size  
**50/case**

Item# **RES092**  
 Pediatric Size  
**50/case**



HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
  - Mask
  - 7ft supply tube

Disposable single patient use for long term delivery procedures

6cc bottle is easy to assemble and disassemble



## Nebulizer Filters

### Plastic Intake Filter with Head

Item# **BF016**

**10/pack**

HCPCS# A7013



### Plastic Intake Filter with Head

Item# **BF017**

**5/pack**

HCPCS# A7013



### Foam Filter

Item# **BF018**

**10/pack**

HCPCS# A7013



### Intake Filter for 50 PSI Compressor

Item#

**BF020**

**3/pack**

Compatible with CHAD Drive  
50 PSI Compressor  
OEM# 18450-FILTER



### Plastic Intake Filter for Respironics InnoSpire

Item# **BF026S**

**10/pack**

HCPCS# A7013

- Mini
- MicroElite
- MiniElite



### Foam Inlet Filter for Invacare 50 PSI

Item#

**BF027-10**

**10/pack**



## Nebulizer Filter Compatibility

### Caire

Medi-Mist	BF018
-----------	-------

### CHAD

Drive 50psi	BF020
-------------	-------

### Drive/DeVilbiss

PulmoNeb	BF017
----------	-------

PulmoNeb LT	BF017
-------------	-------

Sunrise Compact	BF018
-----------------	-------

Pulmo-Aide	BF017
------------	-------

Pulmo-Aide Compact	BF017
--------------------	-------

PulmoMate	BF017
-----------	-------

Power Neb Ultra 18080	BF017
-----------------------	-------

CHAD 50psi	BF020
------------	-------

### Evo

AeroMax	BF018
---------	-------

SportsNeb	BF018
-----------	-------

Super SportsNeb	BF018
-----------------	-------

### Invacare

Passport	BF018
----------	-------

Envoy	BF018
-------	-------

Envoy Jr.	BF018
-----------	-------

Select	BF018
--------	-------

### John Bunn

Nebulite EV	BF018
-------------	-------

LT2	BF018
-----	-------

LX	BF018
----	-------

### Medline

SportsMist	BF018
------------	-------

### Respironics

Inspiration 626 & 929	BF016
-----------------------	-------

MisterNeb	BF016
-----------	-------

Inspiration Elite	BF016
-------------------	-------

MiniElite	BF018
-----------	-------

InnoSpire	BF026S
-----------	--------

### Salter

AIRE Plus	BF018
-----------	-------

### Schuco

2000	BF018
------	-------

3000	BF018
------	-------

### Sunset

NEB100 Series	BF018
---------------	-------

NEB300 Series	BF018
---------------	-------

# StrongArm®

## StrongArm® Comfort Cane

12/case

Item#	Color
<b>BM001</b>	Black
<b>BM001-BL</b>	Blue
<b>BM001-PR</b>	Purple

**Color pack (one of each)**

**BM001-3MC**



- Ergonomic handgrip supports forearm and relieves pressure from wrist
- Offset bend centers grip over tip for maximum stability when cane is at the user's side
- Lightweight and strong
- 12 height adjustments
- Comfort Fit Handle
- Supports up to 500 lbs
- For left and right arm
- Self-standing base: Larger 'footprint' greatly improves traction and overall stability
- Safely supports user walking on grass, gravel, sand, carpet, hard surfaces, and stairs



**Now comes in retail packaging**

## StrongArm® Grabber Reacher Tool

Item#  
**BM002**  
12/case

The ultimate companion for effortless grabbing and reaching! Engineered with precision and designed for practicality, this innovative tool is perfect for those hard-to-reach items and provides extra support for heavier items. Grab with Ease and Confidence in Each Squeeze.



- Transfers weight for improved support
- Rotating jaw for vertical or horizontal use
- Integrated hook and magnetic retriever
- Durable, lightweight aluminum construction
- Folds for easy storage and portability
- 32" of added reach

## StrongArm® Clear View Universal Walker Tray

12/case

Item#	Color
<b>BM003-BL</b>	Blue
<b>BM003-BK</b>	Grey
<b>BM003-PL</b>	Purple



**Unobstructed view for safer mobility**

The Clear View Walker Tray is designed to enhance the user's ability to see obstacles ahead and reduce the risk of falls. Its transparency helps maintain the user's line of sight, providing added safety and confidence while moving. Compatible with the leading brand walkers.

**Transparent**

See where you are walking. Avoid tripping hazards.

**Universal Fit**

Tray works with the leading brand walkers.

**Convenience**

Built-in cup holder and storage bag.

# StrongArm®

## StrongArm® Universal Clip-On Light

Item#  
**BM004**  
4/box

Featuring a super-bright LED light, this clip-on accessory significantly enhances visibility, providing added safety during evening strolls or walking throughout the house. Effortlessly clip onto canes, walkers, rollators, crutches, grabber reachers, knee scooters, wheelchairs, and more, ensuring powerful and reliable illumination on any mobility device.



## StrongArm® LED Magnifying Glass

Item#  
**BM005**  
10/box

- Handheld 4x LED Magnifier
- 6 LED Lights
- Compact, Lightweight, and Portable
- Lens Diameter: 2.5"
- Weighs Only 1.7 oz
- Uses 2 AAA Batteries (Included)



## StrongArm® Retail Starter Kit

Item#  
**BM006**  
Each

**NEW**

Save when you order the StrongArm Retail Starter Kit!

**Includes:**

- 2 BM001** Black StrongArm Comfort Canes
- 2 BM001-BL** Blue StrongArm Comfort Canes
- 2 BM001-PR** Purple StrongArm Comfort Canes
- 2 BM002** StrongArm Grabber Reacher Tools
- 2 BM004** StrongArm Universal Clip-On Lights
- 2 BM005** StrongArm LED Magnifying Glasses



## Bath Safety

### PreserveTech™ 360° Swivel Bath Chair

Item#

**BSE001**

Compact and customizable, the PreserveTech™ Swivel Bath Chair provides a safe transfer in and out of the tub for those with limited mobility. The seat rotates a full 360° in either direction and automatically locks at 90° intervals to help prevent skin irritation and friction wounds during transfers. The seat height can be adjusted from 15"-22", and the padded armrests can be raised out of the way.



### Splash Defense™ Transfer Bench

Item#

**BSE002**

Splash Defense™ Transfer Bench is the first of its kind to allow the shower curtain liner to remain inside the bath tub while a user is bathing. This unique product enhancement is essential to help prevent water from spilling out onto the bathroom floor, creating a danger zone for users and caretakers. The Splash Defense™ Transfer Bench provides security for the caretaker, while giving users back their independence and allowing them to bathe safely and privately again.



### PreserveTech™ 360° Universal Raised Toilet Seat

Item#

**BSE003**

Introducing the PreserveTech™ Universal Raised Toilet Seat. With this truly universal raised toilet seat users can set three different seat heights (2, 4 or 6") by maneuvering just two components. This seat can fit on most standard and elongated toilets and features two hygienic cutaways for easy cleaning. The pivoting armrests can be raised or removed, depending on user needs and bathroom space. This item also includes a lid with a recessed grip for effortless opening.



## Maternity

### Motif Roam Breast Pump with Breast Milk Collection Cups

Item#  
**MATAAA0026**

The Motif Roam is the durable and flexible breast pump designed to keep up with busy moms. It supports routine with ease and freedom, thanks to mom-favorite Breast Milk Collection Cups that eliminate the need for bulky bottles. Designed for freedom and flexibility, the Roam is ready to support moms' routines wherever life goes.

The Motif Roam™ comes equipped with everything you need for a seamless pumping experience; 28 mm breastshields, tubing with connectors, duckbill valves, diaphragms, and flange inserts in sizes 16 mm, 19 mm, 21 mm, and 24 mm for a perfect fit. A lanyard adds portability, while the USB-C cable and power adapter ensure quick, easy charging. Designed for freedom and flexibility, the Roam is ready to support your routine wherever life takes you.

- **Hospital-Strength Suction:** Up to 280 mmHg suction strength, giving moms the efficiency they need for every pumping session
- **Portable:** Equipped with Breast Milk Collection Cups to give moms the freedom to pump hands-free
- **Discreet:** Collection Cups fit in nursing bra and hold up to 5 oz. of milk
- **More Options:** Resupply Kits available in both flange and collection cup styles
- **Variety of Flange Sizes:** Silicone inserts included in sizes 16 mm, 19 mm, 21 mm and 24mm so mom can find the perfect fit
- **Cluster Feeding and General Massage Mode:** 2 NEW Custom Modes designed by Lactation Professionals
- **Suction levels:** Works to express milk from the breast as efficiently as possible with 9 suction levels
- **Closed System:** Backflow prevention keeps the motor clean, preventing contamination

**NEW**



#### What's Included

- |                              |  |
|------------------------------|--|
| 1 - Motif Roam pump motor    | 4 - Pairs of flange inserts (16 mm, 19 mm, 21 mm, 24 mm) |
| 2 - Collection cups          | 1 - Lanyard  |
| 2 - 28 mm breast shield      | 1 - USB-C cable  |
| 2 - Tubing w/Roam connectors | 1 - Power adapter  |
| 2 - Tubing connector         | 1 - Nipple Ruler   |
| 4 - Duckbill valve           |  |
| 2 - Diaphragm                |  |

## Maternity

### Motif Luna Double Electric Breast Pump

Item#  
**MATAAA0013**

Sophisticated, strong, and soothing, the Motif Luna is Motif's most powerful breast pump. With a battery life of 2+ hours between charges, the Luna was built for modern motherhood. The Luna is available in both battery and non-battery versions. Lightweight and portable, the Luna weighs less than 2 pounds. With a suction strength of 280 mmHg, the Luna is proven to provide More Milk in Less Time®.

- Very quiet
- Customizable pumping sessions
- Compact size and light weight
- Less guilt, pain, & anxiety
- Convenient nightlight & digital screen



### Motif Milk Storage Bags

Designed to make breastmilk storage as simple as possible, these bags feature an easy-pour spout for a no-mess experience. Easily transfer milk from the bottle to the bag, seal with the leak-proof double zipper lock, label, and store in your fridge or freezer. No spills. No hassle. Made with your busy day in mind. Comes in a box of 40, 90, 100 or 300 count.

HCPCS: K1005

Item#  
**MATAAC000806**  
**100/pack**

Other pack counts are available.  
Ask your rep for details.



- Easy-pour spout design
- BPA-free
- Leak-proof with double zipper lock
- Stores up to 8 oz. of milk
- Write-on label
- Freezer safe
- Lays flat for compact freezer storage
- Self stand for refrigerator storage

## Equipment Covers

### Bags for Regulators, Conservers, and Small Items



Item#  
**EC1012C**  
10×12"  
1.5mil, Clear  
**1000 bags/box**

### Covers for Nebulizers and Suction Machines



20×24"  
2mil  
**500 bags/roll**

Item#	Color
<b>EC2024C</b>	Clear
<b>EC2024T</b>	Tan Tint

### Respiratory Set Up Bags with Draw Tape



Item#  
**EC1216RB**  
12×16"  
1mil  
**500 bags/box**

### Covers for Concentrators



25×15×30"  
1.5mil  
**250 bags/roll**

Item#	Color
<b>EC251530C</b>	Clear
<b>EC251530T</b>	Tan Tint

### Covers for CPAP/Vents and Walkers



Item#  
**EC2430C**  
24×30"  
1.5mil, Clear  
**500 bags/roll**

## Equipment Covers

### Covers for Wheelchairs, Walkers, and Commodes



Item#	Thickness	Size	Color	Quantity
<b>EC302035C</b>	1.5mil	30x20x35"	Clear	<b>200 bags/roll</b>
<b>EC302035T</b>	1.5mil	30x20x35"	Tan Tint	<b>200 bags/roll</b>
<b>EC5045</b>	1mil	50x45"	Clear	<b>250 bags/roll</b>

- Protect equipment and materials from dust and moisture
- Packed on convenient rolls
- Low Density Polyethylene
- Packaged in dispenser carton

### Mattress Covers



Item#	Vented	Thickness	Size	Color	Quantity
<b>EC39990C</b>	Yes	1.5mil	39x9x90"	Clear	<b>100 bags/box</b>
<b>EC7252C</b>	No	1mil	72x52"	Clear	<b>200 bags/roll</b>
<b>EC7252T</b>	No	1mil	72x52"	Tan Tint	<b>100 bags/roll</b>


## Lancets

### Pressure Activated Safety Lancet

Item# **DIALAN001**  
21G, 2.0mm Depth  
**10/pack**



**NEW**



**Connect** with us



## Your sales representatives



**Sean O'Donnell**

Direct 312-637-6558  
Email sodonnell@sunsethcs.com



**Madison Gwilliam**

Direct 312-533-2459  
Email mgwilliam@sunsethcs.com



**Evan Berger**

Direct 262-788-8948  
Email eberger@sunsethcs.com



**John Bisceglie**

Direct 312-533-2462  
Email jbisceglie@sunsethcs.com



**Jon Kmiec**

Direct 312-533-2482  
Email jkmiec@sunsethcs.com



**Phil Cosens**

Direct 312-533-2470  
Email pcosens@sunsethcs.com



**Brian Lipson**

Direct 312-533-2463  
Email blipson@sunsethcs.com



**Jonny Deutsch**

Direct 872-259-3691  
Email jdeutsch@sunsethcs.com



**Michelle Chojnacki**

Direct 312-533-2486  
Email mchojnacki@sunsethcs.com



**Chris Crawford**

Direct 872-258-1650  
Email ccrawford@sunsethcs.com



**Pavle Pantovic**

Direct 872-282-0751  
Email ppantovic@sunsethcs.com



**Danielle Dubowski**

Direct 312-533-2479  
Email ddubowski@sunsethcs.com



Phone 877-578-6738  
Email sales@sunsethcs.com  
Visit www.sunsethcs.com

Corporate Headquarters  
141 W Jackson Blvd Ste 1950  
Chicago IL 60604

Central Distribution Center  
2725 Alft Ln  
Elgin IL 60124

Western Distribution Center  
6750 W Van Buren St Ste 100  
Phoenix AZ 85043

## National Accounts



**Gabrielle Murphy**

Direct 312-638-8680  
Email gmurphy@sunsethcs.com



**Shannon St. Clair**

Direct 312-533-2449  
Email sstclair@sunsethcs.com

## International Accounts



**Bryan Vigna**

Direct 312-533-2455  
Email bvigna@sunsethcs.com