




# Active Heated Ventilator Circuits

## Active Heated Circuit

Includes VEN004 Active Exhalation Valve Adapter



**VC1002** Adult

**VC1002P** Pediatric



Designed to provide heated active alternatives to the Resprionics circuits 1073222 & 1073231 & 1073221 & 1073230. VC1002 and VC1002P can be used on a wide variety of units (excluding VOC SN and Resmed Astral).


**VC1002 and VC1002P are compatible with a wide range of ventilator units.**



## Active Exhalation Valve Adapter VEN004



Designed for seamless integration with heated circuits from other manufacturers to create an active circuit, can be used with MR810, HC550, and other humidifiers.




## Active Heated Circuit for Resmed Units Only

Includes VEN004-R Active Exhalation Valve Adapter

**VC1002-R** Adult

**VC1002P-R** Pediatric

Designed to provide heated active alternative to the Resprionics circuits 1073222 & 1073231 & 1073221 & 1073230 for use exclusively with Resmed Astral with the single limb valve adapter OEM 27941.



## Active Exhalation Valve Adapter for Resmed VEN004-R

**For use with Resmed units only.** Can be combined with heated circuits from other manufacturers to create an active circuit, to be used with MR810, HC550, and other humidifiers.

