

Available from



CAIRE®

HOMESTYLE™ Aspen

Stationary Oxygen Concentrator

Item#
RESSCE1000



5-Year Warranty*



Quiet Operation



Simple User Interface



Your O₂ Fleet Solution

Features

- Continuous flow 0.5-5 LPM
- Up to 93%±3% oxygen purity
- Sound less than or equal to 40 dB(A)*
- Power supply: 110V~,60Hz
- Power consumption: 320VA
- 32 lbs.
- 15 inches L x 9 inches W x 23 inches H
- Built-in cord carasel
- Elegant and compact frame



*Aspen Warranty Coverage

Proven performance, quality, and warranty, with Sunset's 5-year coverage for the main unit and sieve beds, and 1 year for accessories.

Serving Your Customers Across the Continuum Of Caire.

Non-Delivery Solution

Optimize field efficiency and productivity while improving customer care with the CAIRE non-delivery solution, which combines the industry-leading FreeStyle® Comfort® portable oxygen concentrator with the HOMESTYLE™ Aspen stationary oxygen concentrator.



FreeStyle® Comfort®



HOMESTYLE™ Aspen