

Available from



CAIRE®

NewLife® Intensity® 10

High Flow Stationary Concentrator

Item#

RESAS0994 Newlife Intensity 10L SOC Dual Flow 120V/60Hz

RESAS0991 Newlife Intensity 10L SOC Single Flow 120V/60Hz



Trusted therapeutic oxygen solution designed for medical equipment fleets world wide

Built-to-last Oxygen Delivery Innovation

Since being introduced in 2009, the New Life Intensity 10 has been a trusted high flow stationary concentrator addressing patients' oxygen delivery needs in over 100 countries including specialty applications in home and medical care facilities.



Time-proven design

- **Industry-leading reliability** — greater than 99% reliability rate out of the box and less than 2%* warranty repair rate
- **Comprehensive warranty** — 3-year warranty on the device
- **Durability** — engineered with an oversized compressor and sieve beds to minimize service



Versatile oxygen delivery

- **Oxygen delivery** — up to 95.5% concentration at up to 10 LPM, meeting a wide range of patient prescriptions
- **High outlet pressure** — the first oxygen concentrator offering 20 psi to support specialty applications and up to 200 feet of oxygen tubing
- **Energy efficient** — reduced power consumption, minimizing the burden on the patient's energy bill



Ease of service

- **Minimal maintenance** — proven quality design and performance to minimize required service
- **Modular design** — fewer components reduce time and effort to service the device
- **Enhanced diagnostics** — onboard alarm alerts enable quick troubleshooting

SPECIFICATIONS

Flow Rate 2.0 – 10.0 LPM

Dimensions

27.5 x 16.5 x 14.5 in.
(69.9 x 41.9 x 36.8 cm)

Weight 58 lbs (26.3 kg)

Sound Level 58 dB(a)

Outlet Pressure
20 psig (138 kPa)

Power 600 W

Maximum Tubing
200 ft (61 m)

Warranty 3 years

Maintenance Schedule

Felt filter: 1 year replacement
Intake filter: Clean weekly