



PORTABLE
OUTLET ALL NEW
X9

Portable Outlet CPAP Battery

Rechargeable Battery with
AC/USB Outlets & Autoswitch

Item# **CAP1600X9**

- International medical grade certified IEC60601-1-2



- Larger display, USB-C charging output, runs cooler

- Completely new design and electronics
- Less buttons, simplified features and functions

One battery. Every CPAP. One simple message for your customers.

The Portable Outlet X9 combines universal compatibility with true Uninterruptible Power Supply (UPS) technology, giving dealers a simpler product to stock and an easier solution to recommend. Instead of managing multiple batteries and proprietary adapters, the X9 works with every CPAP, BiPAP and home ventilator using the manufacturer's original power supply.

- ✓ **Universal Compatibility for every CPAP and BiPAP**
- ✓ **True UPS automatically keeps therapy running during outages**
- ✓ **One SKU replaces multiple proprietary battery solutions**



What's Included

- Portable Outlet X9
- AC Power Cord
- Protective Travel Bag
- Owner's Manual

When power goes out, your CPAP stays ON!®

Today's CPAP customers expect more than a travel battery. They want confidence that their therapy will continue whether they have an unexpected, overnight power outage at home or they are preparing for powerful storms in the forecast. Portable Outlet helps dealers meet that need with one universal solution that's easy to explain and simple to support.

Unlike machine-specific batteries, the Portable Outlet X9 uses a standard AC outlet, allowing patients to use their existing manufacturer power supply instead of proprietary adapters or cables. Because it functions as a true Uninterruptible Power Supply (UPS), customers simply leave their CPAP plugged into the battery. If utility power fails, therapy continues automatically.

For dealers, that means fewer compatibility questions, fewer SKUs to inventory and one straightforward recommendation for virtually every CPAP customer.



Key Features & Benefits Explained

True Uninterruptible Power Supply (UPS)

Unlike standard CPAP batteries that require users to reconnect their machine after an outage begins, the X9 automatically switches to battery power so therapy continues without interruption.

Universal Compatibility

The standard AC outlet works with every CPAP and BiPAP machine using the manufacturer's original power supply, eliminating proprietary adapters and reducing compatibility questions.

One Solution Instead of Multiple Batteries

One universal battery can replace multiple machine-specific solutions, simplifying inventory management and making it easier for staff to recommend the right product.

Engineered for Sensitive Medical Equipment

Pure Sine Wave technology and IEC-certified electronics deliver clean, stable power designed to help protect sensitive CPAP electronics.

Full-Night Battery Performance

Typical runtime provides 10-12 hours on most standard CPAP machines (without heated tubing) and over 16 hours on many travel CPAP machines.*

Travel Ready

TSA/FAA compliant for carry-on air travel**, with optional vehicle and solar charging accessories available for off-grid use.

Charge More Than Your CPAP

USB-A and USB-C ports let users conveniently charge phones, tablets and other small electronics while their CPAP remains connected.

Added Peace of Mind

Backed by the industry's first \$10,000 Safe CPAP Power Guarantee™, giving customers additional confidence in their purchase. See details here.

*Actual runtime varies based on machine type, pressure setting and heated humidity/tubing.

**Always check with your airline for any additional battery requirements before flying.

Why Dealers Prefer Portable Outlet

- ✓ One universal battery instead of multiple machine-specific batteries
- ✓ Future-proof compatibility
- ✓ Premium differentiation for safe power
- ✓ Easier inventory management
- ✓ Easier employee training
- ✓ Easier recommendation for customers



Compatibility

Works with:

- Resmed
- Philips Respironics
- React Health
- Fisher & Paykel
- DeVilbiss
- 3B Medical
- Most other CPAP/BiPAP machines



Q "What happens if my power goes out?"

A Leave your CPAP plugged into the Portable Outlet all the time. If your home's power goes out, your CPAP stays running automatically.

Q "Will this work with my machine?"

A Yes. It works with every CPAP because you use your manufacturer's original power supply.

Q "What if I buy a different CPAP later?"

A You can keep using the same Portable Outlet!

To order or request information

P 312.997.9980 E sales@sunsethcs.com W sunsethcs.com

