

PORTABLE OUTLET X9



Most Competitors




True Uninterruptible Power Supply (UPS)

- Powers your CPAP at all times while it charges.
- UPS and travel battery provides 10-12 hours of therapy.
- When power goes out, your CPAP stays ON!

Not a UPS (Stand alone battery)

- Cannot recharge while powering a CPAP.
- Backup only for travel or after an outage.
- When power goes out, find your battery.



Pure Sine Wave (Clean, Stable Power)

- Cleanest power minimizing risk to devices and patients.
- Recommended by ResMed and other manufacturers.*
- Protects motors, blower assemblies, humidifiers, and electronic control boards.

Modified Sine Wave (Dirty Power)

- Waveform distortion, voltage instability.
- May cause heating & performance issues.
- Users risk electrical noise, interference, or even long-term cumulative damage.




Use Manufacturer's Power Supply

- Accepts standard 110V AC plugs.
- Use the manufacturer's power supply (safest method)^
- Protects motors, blower assemblies, humidifiers, and electronic control boards.

Requires Aftermarket Adapters

- Cannot use genuine OEM power supplies.
- CPAP brands warn of potential damage.
- Risk of incorrect voltage delivery which could void manufacturer's warranty



Built-In Power Supply

- Uses a thin universal power cord.
- Inexpensive to replace and readily available.
- Fits with battery in the included durable travel bag.

Requires Bulky Transformer

- Large power supplies clutter nightstands.
- Expensive to replace, hard to find locally.
- Charger won't fit with battery in travel bag.



Universal 110V AC Outlet

- Simple, like plugging into your home's wall outlet.
- Powers all sleep therapy devices without adapters.
- Adapters won't be needed for future therapy devices.

Requires Adapter To Connect

- Confusing for patients and technicians.
- Adapters restrict CPAP compatibility.
- New adapters or batteries may be needed.



IEC 60601-1-2 Electronic Certification

- The global electromagnetic and medical authority.
- Certified to the same medical-device electronic standards used by ResMed and other manufacturers.°
- Tests include extreme heat, shock, vibration & overload.

Not Medically Certified



\$10,000 Safe CPAP Power Guarantee™

- If your CPAP, BiPAP or BiPAP ASV or ST is damaged by the X9 CPAP Battery, we will replace it up to \$10,000.
- Piece of mind: Covered by Safe CPAP Power Guarantee.

Not Guaranteed

* Published in ResMed 198103/7/1 2017-08

° ResMed AirSense & AirCurve devices are also IEC 60601-1-2 certified.

^ Published in ResMed 10113885/1 2023-10