



Homecare Product Catalog



Policies

Ordering Options

Two Easy Ways To Order

Phone

(877) 5-SUNSET (877-578-6738)

Our customer service professionals are here to assist you
M-F 8:00 AM – 5:00 PM Central Time.

Email

sales@sunsethcs.com

*Contact us for more information
on fulfillment and compliance
solutions that best fit your needs.*

Shipping

Ask your rep about shipping to various locations and
our drop shipping partners.

UPS

Whenever possible, orders ship via UPS. If another
carrier is desired, please contact your sales rep.

All orders are packed in the smallest boxes available
and combined with other products on the order to
reduce shipping costs.

All products ship prepaid and are added to the invoice.

All products ship from Elgin, IL or Phoenix, AZ.
Please refer to UPS.com for transit times.

Damaged Merchandise

Call us within 10 days. Keep original cartons for
UPS inspection.

We proudly partner with



McKESSON



Payment Methods

Payment Portal and ACH

Visit www.sunsethcs.com/paymentportal to pay online with ACH
or credit card.

Net 30-Day Payment Terms

For qualified companies. A 1.5% monthly finance charge will be
assessed on all past due invoices.

COD

For non-qualified companies.

Credit Card

We accept all major credit cards. A processing fee may apply.
After 10 days of invoice date, a surcharge of 2.5% will be applied.



DISCOVER

Customer Service

Returns

Contact us within 7 days of receiving your shipment.

Call us for a return authorization number. There is no
restocking charge on orders returned within 20 days of invoice.
Orders returned after 20 days may be assessed a restocking
fee of 15%. No returns after 90 days.
(877) 5-SUNSET (877-578-6738)

Central Distribution Center

Sunset Healthcare Solutions
2725 Alft Ln
Elgin, IL 60124

Western Distribution Center

Sunset Healthcare Solutions
6750 W Van Buren St Ste 100
Phoenix, AZ 85043

Problems

Please call us immediately if you have any complaints (it is
helpful to ask for the customer service representative with
whom you placed your order).
Or email us at: customerservice@sunsethcs.com

Questions?

Our experienced customer service staff will be glad to provide
you with information. Just give us a call.

Pricing

All prices are in U.S. Dollars.
Taxes are added where applicable.
We make every effort to hold pricing throughout the year.
However, prices are subject to change without notice.

Content

CPAP

	Page
Masks and Replacement Parts	4-14
Sunset Healthcare Solutions	4-7
Sleepnet	8-10
Circadiance	11-12
Airway Management	13
Bleep	14
Sleep Apnea Appliances	15-16
Tubing	17-18
Filters	19-27
Respironics	19-20
ResMed	21-22
Fisher & Paykel	23
React Health/3B	23
AEIOMed	24
Beyond	24
Breas/GE	24
CareFusion/Viasys	25
Covidien Nellcor Puritan Bennett	25
DeVilbiss	26
Evo Medical	26
Healthdyne	27
Invacare	27
Löwenstein	27
Resvent	27
ProBasics	27
Humidifier Chambers	28
Chinstraps	29-30
Headgear	30
Cash Items	31-39
Cleaners and Accessories	31-32
Battery/Power Inverter	33
Transcend Micro™ Auto CPAP	34
Mask Liners, Soft Headgear Covers	35
Aromatherapy	36
Connectors, Adapters, Inline Filter	36
Pillows	37-38
Tubing Wraps	38
Storage	39

Oxygen

	Page
Pulse Oximeter	40
SpO2 Sensors	41
Firesafe	42
Concentrators	43-44
Masks	45
Cannula and Cannula Accessories	46-47
Tubing	47
Microdial Flowmeter	48
Regulators and Conserver	49
Oxygen HCPCS Chart	50
Cylinder Carriers	51-52
Connectors	53
Accessories	53-55
Cleaners	55
Analyzers	56
Celon Wet Bands	56
Bacteria and HEPA Filters	57
Concentrator Cabinet Filters	58-64
Caire/AirSep/SeQual	58
CHAD	59
Covidien Nellcor Puritan Bennett	59
Dynarex	59
OxyGo	59
DeVilbiss/Drive	60
Invacare	61
Rhythm, Jumao, Attendant, Lifestyle, Medacure	62
Respironics	63
Healthdyne	64
React Health/3B	64
Nidek	64
Repair Parts	65-69

Suction

	Page
Suction Machine	70
Suction Canister and Kits	70-71
Open Suction Catheters	72
Open Suction Catheter Kits	73
Suction Filters	74
Suction Adapters and Connectors	74
Yankauers	74
Tracheostomy Supplies	75
Heat and Moisture Exchangers	76
Closed Suction Catheters	77-78

Ventilator

	Page
LUSA Ventilator	79
Ventilator Filters	80-81
Adapters and Connectors	81
Ventilator Mask	81
Breathing Circuits	82-85
Heated Single-Limb	82
Non-Heated Single-Limb	83-84
Dual-Limb	85

Other

	Page
Disinfectants	86
Cough Assist Kits	87
Breathing Masks	87
Aerosol	88
Tabletop Nebulizers	89
Handheld Nebulizers	90
Nebulizer Kits	91
Nebulizer and Anesthesia Filters	92
StrongArm Mobility and Comfort	93-94
Bath Safety	95
Maternity	96-97
Equipment Covers	98-99
Blood Testing	99

Sunset Healthcare Solutions

Vesper CPAP Masks

NEW

Full Face Mask with Headgear and Cushion

Size	Item#
Small	CM206S
Medium	CM206M
Large	CM206L

CM206-FP Full Face Fit Pack (Mask and S, M, L cushions)

Replacement Cushions

Size	Item#
Small	CUCM206S
Medium	CUCM206M
Large	CUCM206L

CM206-CLIPS

Replacement Headgear Clips

IFU-00133

Sizing Guide

Vesper Full Face without Headgear

Size	Item#
Small	CM206SNH
Medium	CM206MNH
Large	CM206LNH



HCPCS#
MASK - A7030
HEADGEAR - A7035
CUSHION - A7031

Vesper Full Face Headgear HG206

Nasal Mask with Headgear and Cushion

Size	Item#
Small	CM207S
Medium	CM207M
Large	CM207L

CM207-FP Nasal Fit Pack (Mask and S, M, L cushions)

Replacement Cushions

Size	Item#
Small	CUCM207S
Medium	CUCM207M
Large	CUCM207L

CM207-CLIPS

Replacement Headgear Clips

IFU-00136

Sizing Guide

Vesper Nasal without Headgear

Size	Item#
Small	CM207SNH
Medium	CM207MNH
Large	CM207LNH




HCPCS#
MASK - A7034
HEADGEAR - A7035
CUSHION - A7032

Vesper Nasal Headgear HG207

Vesper masks introduce a pared back fit, with up-to-date design details like the fabric-wrapped ergonomic stabilizing arms that curve up to the temples for patients who prefer to forego a forehead piece.

Vesper masks provide more visibility, a feeling of openness and a lighter fit. The flexible interface moves with side or stomach sleepers. The open structure eliminates forehead marks. Vesper masks also feature a 360-degree swivel elbow to prevent awkward tube tangles.



Features

Soft and skin-friendly headgear fabric helps patients comfortably adapt to treatment

Quiet honeycomb vent reduces airflow erosion intensity for peaceful sleep

Easy pop cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

Quick release connector rapidly and easily separates or connects mask and the tube

Vesper comes in attractive retail packaging



SCAN FOR VESPER
FULL FACE VIDEO



SCAN FOR VESPER
NASAL VIDEO



Sunset Healthcare Solutions

Sol Full Face CPAP Mask

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM205S
Medium	CM205M
Large	CM205L

Replacement Cushions


Size	Item#
Small	CUCM205S
Medium	CUCM205M
Large	CUCM205L

HG205 Replacement Headgear - Fits Most
CM205-CLIPS Replacement Headgear Clips
IFU-00057 Sizing Guide

SCAN TO WATCH
SOL VIDEO



Easy Pop Cushions



Sol without Headgear

Size	Item#
Small	CM205SNH
Medium	CM205MNH
Large	CM205LNH

HCPCS#
 MASK - A7030
 HEADGEAR - A7035
 CUSHION - A7031



Available to add
to your formulary
in SleepGlad

Ask us about our 30-day
mask guarantee

*A premium mask without
the premium price*

Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

Quiet vent

Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

Quick release connector

Rapidly and easily separates or connects mask and the tube

Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment

Sol comes in attractive retail packaging



Ray Pediatric CPAP/NIV Masks Comparable to Nonny masks

Full Face Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#	OEM#
Small	CM061S	AG-PEDFF-S
Medium	CM061M	AG-PEDFF-M
Large	CM061L	AG-PEDFF-L
Extra Large	CM061XL	AG-PEDFF-XL

Replacement Cushions

Size	Item#	OEM#
Small	CUCM061S	AG-PEDFF-CS
Medium	CUCM061M	AG-PEDFF-CM
Large	CUCM061L	AG-PEDFF-CL
Extra Large	CUCM061XL	AG-PEDFF-CXL

Replacement Headgear

Size	Item#	OEM#
Sm/Med	HG061S	AG-PEDFF-HGS
Lg/XL	HG061L	AG-PEDFF-HGL

HCPCS#
 MASK - A7030 • HEADGEAR - A7035 • CUSHION - A7031



Nasal Mask with Headgear, Tube, and Replaceable Cushions

Size	Item#	OEM#
S&M Cushions	CM060-KITS	AG-PEDKIT-S
L&XL Cushions	CM060-KITL	AG-PEDKIT-L

Replacement Cushions

Size	Item#	OEM#
Small	CUCM060S	AG-PEDKIT-CS
Medium	CUCM060M	AG-PEDKIT-CM
Large	CUCM060L	AG-PEDKIT-CL
Extra Large	CUCM060XL	AG-PEDKIT-CXL

Replacement Headgear

Size	Item#	OEM#
Small	HG060S	AG-PEDKIT-HGS
Large	HG060L	AG-PEDKIT-HGL

HCPCS#
 MASK - A7034 • HEADGEAR - A7035 • CUSHION - A7032



Sunset Healthcare Solutions

Sunset Deluxe Full Face CPAP Mask

Mask with Flex Tube, Headgear, and Replaceable Cushion

Size	Item#
Small	CM005S
Medium	CM005M
Large	CM005L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

CMSIZE01
 Deluxe Sizing Guide

HCPCS#
 MASK - A7030
 HEADGEAR - A7035
 CUSHION - A7031

Unobtrusive design
 offers a clear
 line of sight



- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone replaceable cushion increases comfort while maintaining optimal seal

Replacement Cushions

Size	Item#
Small	CUCM005S
Medium	CUCM005M
Large	CUCM005L

Mask without Headgear

Size	Item#
Small	CM005SNH
Medium	CM005MNH
Large	CM005LNH

Sunset Deluxe Nasal CPAP Masks

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM006S
Medium	CM006M
Large	CM006L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

CMSIZE01
 Deluxe Sizing Guide

HCPCS#
 MASK - A7034
 HEADGEAR - A7035
 CUSHION - A7032

Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.



Mask without Headgear

Size	Item#
Small	CM006SNH
Medium	CM006MNH
Large	CM006LNH

Mask with Headgear, Adjustable Forehead T, and Replaceable Cushion

Size	Item#
Small	CM106S
Medium	CM106M
Large	CM106L

Replacement Cushions

Size	Item#
Small	CUCM006S
Medium	CUCM006M
Large	CUCM006L



Sunset Healthcare Solutions

Sunset Deluxe Nasal Pillow CPAP Mask

Mask with Headgear and Small, Medium, and Large Pillows

Item# **CM036**

Headgear Only
HG036

HCPCS#
MASK - A7034
HEADGEAR - A7035
PILLOWS - A7033

Replacement Pillows

Size	Item#
Small	CM036SP
Medium	CM036MP
Large	CM036LP



- Easy to fit, assemble, and maintain
- Three different pillow sizes to fit more people more easily
- 360° rotational elbow with flexible silicone tube
- Quiet vent design

Mask without Headgear
CM036NH

Sunset Clearsight Deluxe Nasal Mask

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM110S
Medium	CM110M
Large	CM110L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

Replacement Cushions

Size	Item#
Small	CUCM110S
Medium	CUCM110M
Large	CUCM110L

HCPCS#
MASK - A7034
HEADGEAR - A7035



Clearsight design provides unobstructed view

Mask Series	Mask Type	Replaceable Cushion	Flex Tube Included	Headgear Included
CM005	Full Face	✓	✓	✓
CM006	Nasal - Forehead Support	✓	✓	✓
CM106	Nasal - Adjustable	✓		✓
CM007	Full Face			✓
CM008	Nasal			✓
CM110	Nasal	✓		✓
CM036	Nasal Pillow	✓	✓	✓

Sunset Classic CPAP Masks

Full Face Mask with Headgear

Size	Item#
Small	CM007S
Medium	CM007M
Large	CM007L
Extra Large	CM007XL

HG003 Headgear - Fits Most
HG003XL Headgear - XL

HCPCS#
MASK - A7030
HEADGEAR - A7035



- Conventional forehead pad removed to create a clear line of sight
- Minimalist design and ultra-lightweight silicone
- Quick-release clips make it easy to remove
- 360° dual-swivel elbow port allows for free movement
- One-piece mask and cushion design is durable and easy

Nasal Mask with Headgear

Size	Item#
Medium	CM008M
Large	CM008L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

HCPCS#
MASK - A7034
HEADGEAR - A7035



- Great value for cash patients or those with lower insurance reimbursement
- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

Sleepnet

Phantom®2 and iQ®2

Phantom®2

Item#
CM050

Replacement
Cushion
CUCM050



iQ®2

Item#
CM049

Replacement
Cushion
CUCM049



*Ask us about
our 30-day
mask guarantee*

Replacement Frame/Tubing Assembly
CM051

- Sleepnet's patented AIR°gel® technology provides maximum comfort while minimizing skin irritation and pressure points on the nasal bridge. Unlike many other gel masks, the unique qualities of AIR°gel allow the mask to seal with minimal tension. By design, AIR°gel is twice as soft as competing gel products.
- Sleepnet has the only mask that you can shape to each individual face. Flex, Fit, Forget!

- Cushions may be removed and replaced without purchasing a new mask.
- Magnetic Clips — Easy on and off!
- 360° Swivel — Less tube drag for active sleepers.
- One size cushion — Larger fitting range.



Sleepnet

Innova

Nasal Mask with Headgear

Size	Item#
Small/Medium	CM034SM
Medium/Large	CM034ML

HCPCS#
MASK - A7034
HEADGEAR - A7035



Full Face Mask with Headgear

Size	Item#
Small	CM035S
Medium	CM035M
Large	CM035L

HCPCS#
MASK - A7030
HEADGEAR - A7035



Innova® CPAP masks feature patented AIR°gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA

MiniMe Pediatric Nasal Mask



Mask with Headgear


Headgear Size	Item#	OEM#
Small	CM004S	60213
Medium	CM004M	60214
Large	CM004L	60215

Replacement Headgear

Size	Item#	OEM#
Small	HG010S	55241
Medium	HG010M	55242
Large	HG010L	55243

HCPCS#
MASK - A7034
HEADGEAR - A7035

- Non-vented version available



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

MiniMe 2 Pediatric Nasal Mask

Mask with Headgear

Size	Item#	OEM#
Extra Small	CM027XS	60216
Small	CM027S	60217
Medium	CM027M	60218
Large	CM027L	60221

Replacement Headgear

Size	Item#
Extra Small/Small	HG027S
Medium/Large	HG027L

HCPCS#
MASK - A7034
HEADGEAR - A7035

- Non-vented version available



Sleepnet

Ascend

Nasal Mask with Headgear and All Cushions

Item# **CM037N**

Replacement Cushions

Size	Item#
Small	CM037NCS
Medium	CM037NCM
Large	CM037NCL

Nasal Mask with Headgear and 1 Cushion

Size	Item#
Small	CM037NS
Medium	CM037NM
Large	CM037NL

HCPCS#
MASK - A7034
HEADGEAR - A7035
CUSHION - A7032

- One frame with three sizes of replaceable cushions: S, M, L
- Unique flexible spacebar design provides superior stability, holding a seal without contacting the skin
- Sleepnet's patented AIR^ogel minimizes pressure points and skin irritation while offering maximum comfort with minimum weight
- Lightweight, breathable quick-release headgear slides off easily with one hand
- 360-degree rotating elbow for optimal tube positioning and easier mobility
- Active Headgear Connector moves as the patient moves for a reliable seal
- Made in the USA



Full Face Mask with Headgear and All Cushions

Item# **CM037F**

Replacement Cushions

Size	Item#
Small	CM037CS
Medium	CM037CM
Large	CM037CL

Full Face Mask with Headgear and 1 Cushion

Size	Item#
Small	CM037S
Medium	CM037M
Large	CM037L

HCPCS#
MASK - A7030
HEADGEAR - A7035
CUSHION - A7031



iQ Nasal Mask

Item#
CM001

iQ Nasal Mask with Headgear

OEM# TMS-889

HCPCS#
MASK - A7034
HEADGEAR - A7035



Mojo Full Face Mask

Mask with Headgear

Size	Item#	OEM#
Small	CM003SH	TMS-835
Medium	CM003MH	TMS-836
Large	CM003LH	TMS-837

Replacement Headgear

HG003

OEM# TMS-89

HCPCS#
MASK - A7030
HEADGEAR - A7035



Mojo 2 Full Face Mask

Mask with Headgear

Size	Item#	OEM#
Small	CM50817	50817
Medium	CM50818	50818
Large	CM50819	50819
Extra Large	CM50820	50820

- Maximum comfort, minimum weight. Patented AIR^ogel minimizes pressure points and skin irritation
- Forehead adjustment screw removes pressure from the bridge of the nose for personalized fit
- The shapeable shell of each mask can be easily molded to fit individual needs: bend it, mold it, shape it!
- Lightweight, breathable headgear disconnects from mask easily with one hand thanks to the quick-release magnetic headgear clips
- Oxygen port allows easy delivery of supplemental oxygen
- 360° rotating elbow for optimal tube positioning and easier mobility



Circadiance

SleepWeaver® Advance Nasal Mask



- One size fits most (Small Tan available for smaller faces)
- Includes Advance headgear that can be adjusted to fit a broader range of patients



Mask with Adjustable Headgear

Mask Color	Item#	OEM#
Blue	CM017B	100274
Camo	CM017C	100281
Leopard	CM017L	100285
Pink	CM017P	100277
Tan	CM017T	100332

CM017T-S
Small Tan
OEM# 101224

HCPCS#
MASK - A7034
HEADGEAR - A7035

HG017
Headgear Only
OEM# 100435

CM017TNH
Tan — No Headgear
OEM# 100336

- Easy to store and pack for travel

SleepWeaver® Élan Nasal Mask

Mask Kit

Mask Color	Item#	OEM#
Tan	CM022T-KIT	100563
Blue	CM022B-KIT	100782

Kits Include:

- Mask
- Small Cushion
- Regular Cushion
- Large Cushion
- Headgear
- Feather Weight Tube®
- Overhead Tube Tether Strap
- Quick Start Illustration Guide
- Instructional DVD

HCPCS#
MASK - A7034
HEADGEAR - A7035
TUBE - A7037

Replacement Cushions

Size and Color	Item#	OEM#
Small Blue	CUCM022B-S	100743
Regular Blue	CUCM022B-R	100654
Large Blue	CUCM022B-L	100664
Small Tan	CUCM022T-S	100742
Regular Tan	CUCM022T-R	100653
Large Tan	CUCM022T-L	100663
Small Leopard	CUCM022L-S	100746
Regular Leopard	CUCM022L-R	100657
Large Leopard	CUCM022L-L	100667

TUB0018F
Replacement
Feather Weight Tube®
OEM# 100477

HG022
Headgear Only
OEM# 100630

Mask with Headgear only

Size and Color	Item#	OEM#
Small Blue	CM022B-S	100728
Regular Blue	CM022B-R	100519
Large Blue	CM022B-L	100482
Small Tan	CM022T-S	100727
Regular Tan	CM022T-R	100518
Large Tan	CM022T-L	100481



All SleepWeaver Masks

- Made of soft cloth
- Allow sleeping on side, stomach, and even wearing glasses
- Easy to store and pack for travel
- Hand or machine washable gentle cycle, mild detergent, hang to dry



- Detachable 90° elbow with 360° rotation allows for over-the-head hose placement

Circadiance

SleepWeaver® Anew™ Full Face Mask

Blue Mask with Headgear

Size	Item#	OEM#
Small	CM026B-S	100955
Regular	CM026B-R	100960
Large	CM026B-L	100965

Headgear

Size	Item#	OEM#
Small	HGCM026-S	101145
Reg/Lg	HGCM026-R	101084

TUB0018F

Feather Weight Tube®
OEM# 100477

HCPCS#
MASK - A7027
HEADGEAR - A7035
TUBE - A7037



- Soft material wicks away moisture and leaves no pressure marks on the face
- Loose fit provides ultimate comfort and eliminates strap marks
- Open line of sight accommodates eyeglasses and is ideal for people who like to fall asleep while reading or watching TV
- 90° elbow swivel connector rotates 360° allowing users to sleep in more natural positions; easy to remove and clean

Pediatric SleepWeaver® CPAP Mask

Item#
CM101325



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.

- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

SleepWeaver® 3D®

Item#
CM027

Blue Mask without Headgear

Headgear	Item#
With Large Headgear	CM027-LB
With Regular Headgear	CM027-RB

100/packs

Headgear	Item#
With Large Headgear	CM027-L-100
With Regular Headgear	CM027-R-100

- New wing-shaped stay routes the headgear under the ears, ensuring comfort and a great fit
- Permanent 360° swivel provides reliable performance in any position, including over the headboard
- Breathable cloth design eliminates fogging that is common in hard plastic masks - nothing breathes like cloth, keeping skin dry and comfortable
- For those who want to read or watch TV while in bed, SleepWeaver 3D can be worn with glasses
- Users prefer a soft cloth against the face over rubber, silicone and other gels. Cloth masks eliminate the need for a cushion or nasal pad



Airway Management

Tap Pap

No headgear needed!

Item#

CM024NH

Mask Kit with Small, Medium, and Large Pillows

OEM# PAP-NP1-000

Replacement Pillows

Size	Item#	OEM#
Small	CUCM024S	PAP-NP1-101
Medium	CUCM024M	PAP-NP1-102
Large	CUCM024L	PAP-NP1-103

CM024MP Replacement Mouthpiece

OEM# PAP-NP1-201

HCPCS#

MASK - A7034

PILLOWS - A7033

HEADGEAR - A7035



- Mouthpiece molds to the upper teeth to provide a custom and secure fit even if the bottom teeth are not in contact with the mouthpiece
- Unique nut and bolt joint allows the nasal pillows to be set to the desired angle
- Provides a stable fit regardless of the sleep position
- Patients report not being able to feel the mask on their face
- 3 sizes of nasal pillows included (Small, Medium, and Large)
- 90-day compliance guarantee for new users of the mask, 30-day guarantee for experienced users

Bleep

Bleep Eclipse Starter Kit

Eclipse Frame and Two 32/Packs of Halos

Item#
CAP100589

The BLEEP Eclipse™ Starter Kit includes the headgear-free Eclipse™ mask frame and two 32-packs of Halos™ (32-day total supply) for immediate use. This clinically-proven NO LEAKS PAP solution is the only one of its kind on the market. No headgear, no mask lines and no broken hair. The comfortable Eclipse mask frame includes a diffuser and flexible tube and connects to Halos in lieu of headgear. Halos affix to the nasal area without inserting into the nostrils, significantly reducing irritation and claustrophobia, driving patient comfort and health, and preventing face marks. Halos are hypoallergenic, with gentle, surgical grade adhesive, and free from BPA, corn or latex. Halos utilize innovative MagSeal™ technology for rapid connection and disconnection. BLEEP sleep solutions are covered by Medicare and insurance, and the only CPAP solutions fully made in the USA.

HPCS# A7034



Item# **CAP100590**
Bleep Halos for
Eclipse Mask Frame
32/Pack
HPCS# A7033

Bleep DreamPort Sleep Solution

DreamWay Interface with 32/Pack of DreamPorts

Item#
CMB100381

Make the mask a thing of the past.

Unlike traditional CPAP masks the DreamPort Sleep Solution from Bleep is headgear free, lightweight, flexible and easy to use. The patented solution uses medical grade adhesive strips to hold the DreamPorts (seals) in place eliminating headgear. And unlike traditional seals with the DreamPorts you will not wake up with marks on your face.

Includes:

- Bleep DreamWay Interface with 10 inch hose and connector
- Bleep DreamPorts (32-Pack, 16-Night Supply)
- Instruction Guide



Item# **CUB100379**
Bleep DreamPorts for
DreamWay Interface
32/Pack



Bleep DreamPreps Witch Hazel Pads

Item#
CAPB100433
100 Pads



Bleep Vent Diffuser Kit

Item#
CAP100604
**1 Hook and Loop
Black Strap and
3 White Filters**



Item# **CAP100662**
Bleep Vent Filter
Replacement Pack
12/pack



Sleep Apnea Appliances

Hush Strips®

Item#

CAP1025


Box of 32 Strips

Hush Strips® are designed to improve sleep quality by reducing mouth breathing and snoring during sleep

Hush Strips® sleep strips are made from medical grade nonwoven fabric and gentle easy-release adhesive, which moves with the skin flexibly, making them easy to reposition and remove.

Easy to apply residue-free material

The simple, non-invasive and easy to use anti-snoring strips are made of medical-grade hypoallergenic material which doesn't leave residue and is safe for both kids and adults!



How to use

1. Wash and dry lips
2. Remove strip from the liner
3. Close your mouth and turn lips slightly inward
4. Apply strip on lips and press to seal

Airway Management

myTAP



Precision-fit, immediate, simple and affordable solution for snoring and mild to moderate sleep apnea.

- Offers immediate treatment
- Can be fitted in 15 minutes in a sleep lab, physician's office or DME
- Clinically proven in over 40 independent peer reviewed studies to treat snoring and obstructive sleep apnea (OSA)
- Patients control and remold trays with Precision-Fit, enabling overall comfort, compliance and success: 90% of patients achieve ideal fit upon first try, no bite registration needed
- Thin, sleek trays create a comfortable, low-profile appliance
- Simple design makes it an affordable and convenient option
- Includes an AM Aligner to help offset potential bite changes
- Includes a mouth shield accessory to prevent dry mouth or excess saliva. A secondary benefit is the promotion of nasal breathing, the best natural way to breathe

AirAvant

Bongo Rx EPAP

Item# CAPBNG500

Bongo Rx all sizes starter kit

Includes:

- All four sizes: small, medium, large, and extra large
- Travel case
- Drying stand
- Headgear
- No machine, tube, electricity, or batteries needed
- Use at home, traveling, camping... use anytime, anywhere
- Small enough to fit in a shirt pocket
- Reusable and easy to clean with simply soap and water
- Designed, molded, and assembled in the USA
- Clinically tested to demonstrate safety and effectiveness
- Supplemental lifestyle option for your compliant patients
- Patient friendly option for your non-compliant patients



The Bongo Rx is a simple device with soft nasal seals and is available in various sizes (S, M, L, and XL). The Bongo Rx is an EPAP (expiratory positive airway pressure) device and works as follows:

- During inhalation, small valves open to allow normal breathing through the device.
- During exhalation, the valves close and direct exhaled air through small vents to generate EPAP, which keeps the airway open.

Annual Replenishment Pack

Includes 4 pairs, travel case, drying stand, and headgear

Size	Item#
Small	CAPBNG502
Medium	CAPBNG503
Large	CAPBNG504
Extra Large	CAPBNG505



Single Replenishment Unit

Includes 1 pair and headgear

Size	Item#
Small	CAPBNG512
Medium	CAPBNG513
Large	CAPBNG514
Extra Large	CAPBNG515



CPAP Tubing

Durable CPAP Tubing — Gray



Item#

TUB06

6ft

19mm with 22mm cuffs

HCPCS# A7037

- Sunset's most popular item
- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free


All Lengths

TUB0018	18in
TUB02	2ft
TUB004	4ft
TUB06	6ft
TUB006.5	6ft 6in — Compatible with ResMed tube wraps
TUB007	7ft
TUB08	8ft
TUB10	10ft

Compatible with Sunset's new Quick Release Connector!

Item#

CAP1016



Slim CPAP Tube

Item#

TUB006SS

6ft x 15mm with
22mm cuffs

HCPCS# A7037



Highly durable

Strain relief cuff decreases separation and tearing

Virtually no memory

Designed to be the most flexible tube available

Minimal weight

One of the lightest tubes on the market at 2.6 ounces

Great value

Affordable and sustainable

Pressure Port Tubing

Item# **TUB006PP**

6ft with
Pressure Sensor Lines
22mm Cuffs

HCPCS# A7037



- Built-in sensor line for use with the Puritan Bennett: GoodKnight 420S, 420E and 425 and KnightStar 330

Heated CPAP Tubing

Heated Tube for Respirationics

Item#

TUB06HPR


Heated Tube for
Philips Respirationics
Dreamstation Series
and System One

Each

OEM# HT15

HPCPS# A4604

- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off



Heated Tube for Fisher & Paykel ICON

Item#

TUBHC900ICON208

OEM# 900ICON208



Respironics

DreamStation 2



Item# CF1012

Disposable Filter

6/pack

OEM# 1142831

HCPCS# A7038

CF1012-1 1/pack OEM# 1142828

CF1012-2 2/pack OEM# 1142829

CF1012-3 3/pack

CF1012-50 50/pack

Item# CF1012F

Reusable Filter

2/pack

OEM# 1142687

HCPCS# A7039

CF1012F-1 1/pack

DreamStation



Item# CF1010

Disposable Filter

6/pack

OEM# 1122519

HCPCS# A7038

CF1010-1 1/pack

CF1010-2 2/pack

Item# CF1010F

Reusable Filter

2/pack

OEM# 112446 (1-pack)

HCPCS# A7039

CF1010F-1 1/pack

CF1010F-6 6/pack

M-Series



Item# CF1007

Ultrafine Filter - Tabbed

6/pack • 1³/₄" x 7⁷/₈"

HCPCS# A7038

CF1007-1 1/pack

CF1007-2 2/pack

Item# CF1007F

Foam Filter

2/pack • 1⁵/₈" x 7⁷/₈"

OEM# 1029330

HCPCS# A7039

CF1007F-1 1/pack

PR System One



Item# CF1009NT

Ultrafine Filter - No Tab

6/pack • 1³/₄" x 7⁷/₈"

OEM# 1029331

HCPCS# A7038

CF1009NT-1 1/pack

CF1009NT-2 2/pack

Item# CF1009

Ultrafine Filter - Tabbed

6/pack • 1³/₄" x 7⁷/₈"

HCPCS# A7038

CF1009-1 1/pack

CF1009-2 2/pack

Item# CF1009F

Foam Filter

2/pack • 1⁵/₈" x 7⁷/₈"

OEM# 1029330

HCPCS# A7039

CF1009F-1 1/pack

*All filters can be ordered in any pack quantity.
Contact your sales rep for details.*

Respironics

REMstar Pro, REMstar Plus



Item# CF1006
Ultrafine Filter
2/pack • 3¹/₂" x 1¹/₂"
OEM# 1005945
HCPCS# A7038

Item# CF1006F
Foam Filter
2/pack • 3⁵/₈" x 1⁵/₈"
OEM# 1005964
HCPCS# A7039

BiPAP Pro, Plus, Duet LX, Synchrony



Item# CF1005
Ultrafine Filter
6/pack • 4³/₄" x 3¹/₄"
OEM# 622219
HCPCS# A7038

Item# CF1005F
Foam Filter
2/pack • 4³/₄" x 2³/₈"
OEM# 622220
HCPCS# A7039

Aria, Virtuoso, Duet, Quartet



Item# CF1004
Ultrafine Filter
6/pack • 3¹/₂" x 2³/₁₆"
OEM# 532311
HCPCS# A7038

Item# CF1004F
Foam Filter
2/pack • 3¹/₄" x 2¹/₈"
OEM# 532247
HCPCS# A7039

REMstar



Item# CF1003
Ultrafine Filter
6/pack • 5¹/₂" x 4³/₄"
OEM# 362522
HCPCS# A7038

Item# CF1003F
Foam Filter
1/pack • 5³/₄" x 5¹/₄"
OEM# 362521
HCPCS# A7039

Solo, Aria LX, Virtuoso LX, REMstar LX



Item# CF1002
Ultrafine Filter
6/pack • 3³/₈" x 2³/₄"
OEM# 622017
HCPCS# A7038

Item# CF1002F
Foam Filter
2/pack • 3⁵/₁₆" x 2⁵/₁₆"
OEM# 622018
HCPCS# A7039

SleepEasy/BiPAP



Item# CF1001
Ultrafine Filter
6/pack • 7³/₄" x 4⁷/₈"
OEM# 302064
HCPCS# A7038

ResMed

AirSense 11



Item#
CF2111

Disposable Filter

3/pack

HCPCS# A7038

CF2111-1 1/pack OEM# 39300

CF2111-2 2/pack OEM# 39301

CF2111-50 50/pack OEM# 39303

Item#

CF2111HA-3

Hypoallergenic Filter

3/pack

1 1/4" x 1 3/16"

HCPCS# A7038

CF2111HA-1 1/pack

CF2111HA-2 2/pack

AirSense 11 Carry Bag
CAP1028-S11



AirMini



Item#

CF2108

Disposable Filter

2/pack • 1 1/4" x 1 3/16"

OEM# 38815

HCPCS# A7038

S9 and AirSense 10 Series



Item#

CF2107

Disposable Filter

3/pack • 2 1/8" x 1 3/8"

HCPCS# A7038

CF2107-1 1/pack OEM# 36850

CF2107-2 2/pack OEM# 36851

CF2107-50 50/pack OEM# 36853

Item#

CF2107HA

Hypoallergenic Filter

3/pack • 2 1/8" x 1 3/8"

HCPCS# A7038

CF2107HA-1 1/pack OEM# 36855

CF2107HA-2 2/pack OEM# 36856

CF2107HA-50 50/pack

OEM# 36858

Adapt SV



Item#

CF2106

Disposable Filter

3/pack • 1 1/4" x 1 1/2"

OEM# 26926

HCPCS# A7038

Item#

CF2106HA

Hypoallergenic Filter

3/pack • 1 1/4" x 1 1/2"

OEM# 26927

HCPCS# A7038

ResMed

S8 Series



Item#
CF2105
Disposable Filter
3/pack • 1³/₈" x 1³/₈"
OEM# 33916
HCPCS# A7038

CF2105-1 1/pack OEM# 33915
CF2105-2 2/pack OEM# 33973
CF2105-50 50/pack OEM# 33918



Item#
CF2105HA
Hypoallergenic Filter
3/pack • 1³/₈" x 1³/₈"
OEM# 33920
HCPCS# A7038

CF2105HA-1 1/pack OEM# 33919
CF2105HA-2 2/pack OEM# 33975
CF2105HA-50 50/pack OEM# 33922

S7 Series



Item#
CF2103
Disposable Filter
3/pack • 1¹/₄" x 1¹/₂"
OEM# 30907
HCPCS# A7038

CF2103-1 1/pack OEM# 30919
CF2103-2 2/pack
CF2103-50 50/pack OEM# 30908



Item#
CF2103HA
Hypoallergenic Filter
3/pack • 1¹/₄" x 1¹/₂"
HCPCS# A7038

CF2103HA-1 1/pack OEM# 30933

S6 Series



Item#
CF2102
Disposable Filter
3/pack • 5³/₄" x 5⁵/₈"
OEM# 21935
HCPCS# A7038

CF2102-50 50/pack
OEM# 21936

Fisher & Paykel

SleepStyle Series



Item#
CF4004-2
Disposable Filter
2/pack
HCPCS# A7038



ICON Series



Item#
CF4003
Disposable Filter
4/pack
2 1/16" x 5/8"
HCPCS# A7038



CF4003-1 1/pack
OEM# 900ICON503

CF4003-6 6/pack
OEM# 900ICON212

HC230 Series, HC600 Series



Item#
CF4002
Disposable Filter
4/pack • 2 3/4" x 7/8"
OEM# 900HC240
HCPCS# A7038
CF4002-2 2/pack



HC200, HC220, HC221



Item#
CF4001
Disposable Filter
4/pack • 5" x 3/4"
OEM# 900HC222
HCPCS# A7038
CF4001-2 2/pack



React Health/3B

3B Luna G3 and React G3X



Item# **CF8006-6**
Disposable Filter
6/pack
HCPCS# A7038
CF8006-1 1/pack
CF8006-2 2/pack

Item# **CF8006F-1**
Reusable Filter
1/pack
HCPCS# A7039

Luna and Luna II



Item#
CF8005
Disposable Filter
12/pack • 2" x 1 3/16"
HCPCS# A7038

Item#
CF8005F-1
Reusable Filter
1/pack • 2" x 1 3/16"
HCPCS# A7039
CF8005F-2 2/pack

RESmart



Item#
CF8004
Disposable Filter
3/pack • 3 5/16" x 1 5/8"
HCPCS# A7038



AEIOMed

Everest Aura, Everest III and Integrated



Item# **CF8001**
Foam Filter
2/pack • 2" x 1⁹/₁₆"
OEM# 100446
HCPCS# A7039

Beyond

ResPlus



Item# **CF9005F-2**
Reusable Filter
2/pack
HCPCS# A7039
CF9005F-1 1/pack

Breas/GE

iSleep 20 Series



Item# **CF7003**
Ultrafine Filter
2/pack • 2¹/₄" x 7⁷/₈"
OEM# 004153
HCPCS# A7038

Item# **CF7003F**
Foam Filter
Each • 2³/₄" x 1"
OEM# 004154
HCPCS# A7039

PV10



Item# **CF7002**
Ultrafine Filter
2/pack • 2³/₄" x 7⁷/₈"
OEM# 2B01975
HCPCS# A7038

Item# **CF7002F**
Foam Filter
Each • 2⁵/₈" x 7⁷/₈"
OEM# 2B01972
HCPCS# A7039

PV1



Item# **CF7001**
Ultrafine Filter
2/pack • 2¹/₂" x 1¹/₄"
HCPCS# A7038



Item# **CF7001F**
Foam Filter
3/pack • 2¹/₂" x 1¹/₄"
HCPCS# A7039

CareFusion/Viasys

PureSom



Item#
CF9003
Ultrafine Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7511
HCPCS# A7038

Item#
CF9003F
Foam Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7510
HCPCS# A7039

Orion and Pegasus



Item#
CF8002
Ultrafine Filter
6/pack • 2⁷/₈" x 1³/₄"
OEM# 774688
HCPCS# A7038



Item#
CF8002F
Foam Filter
Each • 2⁷/₈" x 1³/₄"
OEM# 774689
HCPCS# A7039

Covidien Nellcor Puritan Bennett

GoodKnight 314/318



Item#
CF6005F
Foam Filter
2/pack • 6¹/₄" x 1"
OEM# S-290609-00
HCPCS# A7039

GoodKnight 420/425



Item#
CF6003
Ultrafine Filter
2/pack • 1³/₁₆" x 3³/₁₆"
OEM# M-413950-04
HCPCS# A7038

Item#
CF6003F
Foam Filter
Each • 1³/₁₆" x 3³/₁₆"
OEM# M-413950-02
HCPCS# A7039

GoodKnight 418 A/G/P/S



Item#
CF6002
Ultrafine Filter
4/pack • 5³/₁₆" x 3³/₄"
OEM# Y-413560-01
HCPCS# A7038

GoodKnight 418



Item#
CF6001
Ultrafine Filter
4/pack • 2¹/₄" x 1¹/₂"
OEM# S-645030-00
HCPCS# A7038

Item#
CF6001F
Foam Filter
6/pack • 2⁷/₁₆" x 1¹/₂"
OEM# S-645029-00
HCPCS# A7039

DeVilbiss

IntelliPAP



Item#
CF5004
Ultrafine Filter
4/pack • 1³/₄" x 1¹/₄"
OEM# DV51D-603
HCPCS# A7038

Item#
CF5004F
Foam Filter
4/pack • 1⁵/₈" x 1³/₁₆"
OEM# DV51D-602
HCPCS# A7039

Horizon LT



Item#
CF5003
Ultrafine Filter
4/pack • 4¹/₄" x 1⁵/₁₆"
OEM# 8000D-603
HCPCS# A7038

Item#
CF5003F
Foam Filter
4/pack • 4¹/₄" x 1⁵/₁₆"
OEM# 8000D-602
HCPCS# A7039

Evo Medical

RemRest



Item#
CF8003F
Foam Filter
3/pack • 2³/₁₆" x 1¹/₁₆"
OEM# 21117
HCPCS# A7039

Comfort PAP



Item#
CF9003
Ultrafine Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7511
HCPCS# A7038

Item#
CF9003F
Foam Filter
3/pack • 1³/₈" x 1¹/₄"
OEM# 7510
HCPCS# A7039

Enhance your branding! FREE PRIVATE LABELING

Sunset's free service allows you to put your logo and contact information on our products.

Increase brand awareness and allow patients to easily contact you.

Work with our marketing team to customize the layout with your logo, phone number, website address, and reorder information.

Available on all CPAP tubes and filters.
Minimum orders apply.
Certain restrictions may apply.
Contact your sales rep for details.



Healthdyne

Quest



Item# **CF3003**
Ultrafine Filter
with Curves
5/pack • 4¹/₂" x 4³/₈"
OEM# 7302
HCPCS# A7038

Item# **CF3003F**
Foam Filter
with Curves
5/pack • 4¹/₂" x 4¹/₄"
OEM# 7301
HCPCS# A7039

Tranquility



Item#
CF3001
Ultrafine Filter
5/pack • 5¹/₈" x 1⁷/₈"
OEM# 7170
HCPCS# A7038

Löwenstein

prisma



Item# **VF008-12**
Disposable Filter
12/pack
OEM# WM29652

Item# **VF008F-2**
Reusable Filter
2/pack
OEM# WM29928

Resvent

iBreeze



Ultrafine Filter HCPCS# A7038
CF9901-1 1/pack
CF9901-2 2/pack
CF9901-6 6/pack
CF9901-50 50/pack

Invacare

Polaris LT



Item# **CF9002**
Ultrafine Filter
3/pack • 5" x 2"
OEM# 1082428
HCPCS# A7038

Item# **CF9002F**
Foam Filter
3/pack • 5" x 2"
OEM# 1082995
HCPCS# A7039

Polaris EX



Item# **CF9001**
Ultrafine Filter
3/pack • 2³/₁₆" x 1"
HCPCS# A7038

Item# **CF9001F**
Foam Filter
3/pack • 2¹/₄" x 1¹¹/₁₆"
OEM# 1124211
HCPCS# A7039

Polaris TR



Item# **CF8001**
Foam Filter
2/pack • 2" x 1⁹/₁₆"
OEM# 100446
HCPCS# A7039

ProBasics

Zzz-PAP



Item# **CF9003**
Ultrafine Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7511
HCPCS# A7038

Item# **CF9003F**
Foam Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7510
HCPCS# A7039

Respironics

Humidifier Chamber for PR System One



Item# **HUM1136502**
Humidifier Chamber for PR System
One 60 Series
Each HCPCS# A7046

Humidifier Chamber for DreamStation



Item# **HUM1122520**
Humidifier Chamber
for DreamStation
Each HCPCS# A7046

Fisher & Paykel

Humidifier Chambers for HC100 Series



Item# **HUMHC300**
Extended Life Humidifier Chamber
for the Fisher & Paykel HC100 Series
Each OEM# HC 300



Item# **HUMHC325S**
Replacement Chamber for the
Fisher & Paykel HC100 Series
Each OEM# HC 325S
HCPCS# A7046

Compatible with: • HC100
• HC105
• HC150



Humidifier Chambers for HC200 Series



Item# **HUMHC355**
Extended Life Humidifier Chamber for
the Fisher & Paykel HC200 Series
Each OEM# HC 355



Item# **HUMHC385**
Replacement Chamber for the Fisher
& Paykel SleepStyle HC200 Series
Each OEM# HC 385S
HCPCS# A7046



Compatible with: • HC200 • HC231 • HC236
• HC201 • HC232 • HC238
• HC220 • HC233 • HC242
• HC221 • HC234 • HC244
• HC254

Humidifier Chamber for HC600 Series



Item# **HUMHC365**
Replacement Chamber for the
Fisher & Paykel SleepStyle 600 Series
Each OEM# HC 365S • HCPCS# A7046

Compatible with: • HC604
• HC608



Humidifier Chamber for ICON Series



Item# **HUMHC900ICON200**
Replacement Chamber for the
Fisher & Paykel ICON Series
Each OEM# 900ICON200 • HCPCS# A7046

Compatible with: • ICON Premo
• ICON Novo
• ICON Auto



Sunset Chinstraps

Universal Neoprene Chinstrap



Item# **CS001**

Universal Chinstrap

HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows

Adjustable Neoprene Chinstrap



Size	Item#
6in	CS003S
7in	CS003R
8in	CS003M
9in	CS003L
10in	CS003XL

Compare to Puresom Adjustable Ruby

OEM# TMS-09 ADJ
HCPCS# A7036

- Latex-free
- Thin adjustable strap attaches at the top of the head. Thicker bottom strap attaches to itself at the back of the head
- Great for hard-to-fit patients

For fitting, measure along one side of the head, as indicated by the dashed line above.

Comfort Chinstrap



Item# **CS005**

Compare to ResMed Sullivan

OEM# 16015
HCPCS# A7036

- Latex-free
- One size fits most
- Easy to use —attaches to itself at the top of the head
- Soft cloth chin cup adds comfort

Basic Chinstrap



Item# **CS002**

Compare to Respireonics Basic

OEM# 302175
HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to the top strap of headgear
- Cups the chin for comfort
- Works best with any nasal mask

Deluxe Chinstrap



Size	Item#
26in	CS004S
28in	CS004L

Compare to Respireonics Deluxe

OEM# 302425
HCPCS# A7036

- Latex-free
- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort

Premium Chinstrap



Item# **CS006**


Compare to Respireonics Premium

OEM# 1012911
HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted
- Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.

Sunset Chinstraps and Headgear

Ruby Style Chinstrap



For fitting, measure completely around the head, as indicated by the dashed line above.

Size	Item#
21.5-22.5in	CS007S
22-23in	CS007M
23-24in	CS007L

Compare to Puresom Ruby
OEM# TMS-09 • HCPCS# A7036

- Latex-free
- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

Halo Style Chinstrap



Item# **CS025**
Compare to BreatheWear Halo Chinstrap
OEM# H3407
HCPCS# A7036

Ruby Style Adjustable Chinstrap



Size	Item#
Regular	CS007-ADJ
Extra Large	CS007-ADJXL

HCPCS# A7036

- Latex-free
- Fully adjustable
- Super soft material for extra comfort



Universal Mesh Headgear



Item# **HG001**
Compare to Respironics SimpleStrap Headgear
OEM# 302328
HCPCS# A7035

- Latex-free
- One size fits most
- Can be used with most 4-point nasal or full face masks

Universal Headgear



Size	Length	Item#
Regular	22.5in	HG003
Extra Large	24.5in	HG003XL

HCPCS# A7035

- Latex-free
- Can be used with most 4-point nasal or full face masks

Cleaners

Sunset Cleaning Wipes — Flowpack

Item#

CAP1003S

Flowpack of 64 wipes
12/case

- Ingredients are safe and gentle for daily cleaning of all types of CPAP and oxygen masks and equipment



Private Label CPAP Wipes

Order the flowpack wipes with your logo and contact info on a lid sticker, with a minimum quantity commitment. Talk to your sales rep for more info.



Citrus II Cleaning Wipes

Item# **CAP1003**

Canister of 62 wipes
12/case



Cleaning Wipes — Single Packs

Item# **CAP1005**

Citrus II Box of 12 single-pack wipes
12/case

- Convenient individually wrapped wipes are great for patient samples or travel use



CleanSmart CPAP Disinfectant Spray

Item#

CAP1018-L

16oz • **30/case**

CAP1018-S

20Z • **180/case**

- Kills 99.9% of viruses & bacteria
- Included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2
- Active ingredient Hypochlorous Acid does not leave a chemical residue
- EPA Registration Number 89896-2
- Formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment



Citrus II Mask Cleaning Spray

Item# **CAP1001**

1.5oz • **24/case**

Item# **CAP1002**

8oz • **12/case**

- All-natural cleanser for daily mask cleaning
- Available in 1.5oz travel size or 8oz non-aerosol spray cans




CPAP Cash Items

CPAP Tube Cleaning Brush

Item#
CAP1013
Standard 6ft
Each

- Brushes all the way through a 6ft CPAP tube without leaving debris in the middle
- Stainless steel wire and soft nylon bristles



SleepRes V-Com™

Item#
CAP1026
Each

Helping patients adjust to CPAP on their own

Adding V-Com to a patient's therapy is like Training Wheels for CPAP™. This simple device can mean the difference between your patient giving up on PAP or continuing their lifesaving treatment.

Within minutes, a respiratory therapist can show a patient what the PAP pressure feels like with and without the V-Com device inline. This feedback mechanism allows them to better understand how their airflow should feel.

As you know, comfort equals adherence and adherence equals better treatment outcomes.



90-Degree CPAP Tubing Elbow

Item#
CAP1027
Each

The Sunset 90-Degree CPAP Tubing Elbow Adapter connects standard 22mm or slim 15mm CPAP tubes to most CPAP machines at a 90°, relieving stress on the tubing connection and creating a more convenient pathway to the bed.

OEM# 37394



Handheld CPAP Pressure Meter

Item#
CAP003
Each

- Measures flow 0-18 cm/H2O
- Connects directly to the CPAP tubing without an adapter
- Lightweight device used to verify the pressure of a CPAP unit



The V-Com Difference

Improve early PAP tolerance with V-Com by providing your patients with initial comfort that doesn't compromise their treatment.

Traditionally, when a patient struggles with the constant flow of air pressure, pressure settings are reduced.

Reducing pressure settings to increase a patient's initial tolerance improves comfort, but at a cost.

Decreasing pressure settings also means decreasing expiratory positive airway pressure (EPAP) which is principally responsible for therapy.

The V-Com effect is flow-dependent with maximum comfort during the high inspiratory flows. As a result, V-Com has minimal effect on expiratory positive airway pressure.

Nasal Pad

Item#
CAP1029
5/pack



NEW

- Compatible with most CPAP masks
- Soft and gentle, sensitive skin-approved, latex- and silicone-free
- Protects against skin breakdown, redness, irritation and pressure marks
- Conveniently reusable and simple to apply, remove and clean
- Gel can improve mask stability, leak-prevention and seal for comfortable, restful sleep

CPAP Battery and Power Inverter/Charger

Portable Outlet UPS CPAP Battery

Item#

CAP1600X9

Rechargeable Battery with AC/USB Outlets & Autoswitch

One battery. Every CPAP.**One simple message for your customers.**

The Portable Outlet X9 combines universal compatibility with true Uninterruptible Power Supply (UPS) technology, giving dealers a simpler product to stock and an easier solution to recommend. Instead of managing multiple batteries and proprietary adapters, the X9 works with every CPAP, BiPAP and home ventilator using the manufacturer's original power supply.

- ✓ **Universal Compatibility for every CPAP and BiPAP**
- ✓ **True UPS automatically keeps therapy running during outages**
- ✓ **One SKU replaces multiple proprietary battery solutions**
- ✓ **Full-Night Battery Performance**
- ✓ **Engineered for Sensitive Medical Equipment**

Works with:

- Resmed
- Philips Respironics
- React Health
- Fisher & Paykel
- DeVilbiss
- 3B Medical
- Most other CPAP/BiPAP machines

**NEW**

What's Included

- Portable Outlet X9
- AC Power Cord
- Protective Travel Bag
- Owner's Manual

FSA **HSA**
ELIGIBLE
IEC **MEDICALLY**
CERTIFIED
 60601-1-2

\$10,000 **SAFE CPAP**
 POWER GUARANTEE*
See website for details

When power goes out, your CPAP stays ON!®

Unlike machine-specific batteries, the Portable Outlet X9 uses a standard AC outlet, allowing patients to use their existing manufacturer power supply instead of proprietary adapters or cables. Because it functions as a true Uninterruptible Power Supply (UPS), customers simply leave their CPAP plugged into the battery. If utility power fails, therapy continues automatically.



Specs

Operating Temperature	40°F–105°F (4°C–40°C)
Dimensions	6" W × 9.25" D × 1.5" H
Weight	3 lbs.
Warranty Policy	1-Year Limited Replacement Warranty
Contains Lithium	Yes – Rechargeable Lithium-Ion Battery
FAA Compliant	Yes – Meets TSA & FAA requirements for carry-on air travel*
Battery Capacity	159 Wh / 14.4 Ah
Output Connection Type	Standard AC Outlet + USB-A + USB-C
Number of Outlets	1 AC Outlet + 2 USB Charging Outputs
Output Voltage Capabilities	AC Output: 110V / 60Hz (matches input voltage internationally when charging); USB Output: 5V / 3.0A
Car Charging	Yes – Optional Fast Charging Cable
Solar Charging	Yes – Compatible with Portable Outlet 60W Solar Panel Charger
DC Output	USB-A & USB-C (5V / 3.0A) + DC input charging port (12–19V, 2–3.2A)
Recharge Time	Approximately 5–6 hours (wall charging)

*Always check with your airline for any additional battery requirements before flying.

DC to AC Power Inverter, Car Charger, Lighter Plug-in

Item# **CAPDC2AC150**

Recharge your Portable Outlet CPAP Battery using a car's 12V lighter outlet

- Can also power most CPAPs and other small electronics
- 150 Watts
- 3-prong outlet
- 2 USB ports



CPAP Cash Items

Transcend Micro™ 510 Auto CPAP

Each

Item#

CAPSOM507000 White

CAPSOM507001 Pink

CAPSOM507002 Gray

CAPSOM507003 Blue

- **Small:** About the size of a baseball
- **Lightweight:** Weighs less than a mobile phone
- **Works with any CPAP mask:** Use your existing mask
- **Quiet operation:** The Transcend WhisperSoft™ muffler is included

- Smallest and lightest travel CPAP device
- Compatible with any CPAP mask
- Unplug and sleep anywhere with more power options
- Adjustable GentleRise and AirRelief technology
- Stay dry, all the time
- AirMist HME – waterless humidification

Includes:

- WhisperSoft™ Muffler
- Hard Shell Case
- Power Supply
- TransFlex 6ft Hose
- AirMist HME (1 Adapter +1 Cartridge)
- Woven Filter (Inside the unit)
- Extra Hypoallergenic Filter
- 3 Year manufacturer's warranty

★★★
MADE IN
USA



SCAN FOR MORE: VIDEOS AND SPECS



Transcend Micro™ Accessories

Item#

CAPSOM505001 PowerAway® Battery

CAPSOM505049 Sun Watt Solar Panel Charger

CAPSOM505035 Plug Adapters (AU/EU/UK)

CAPSOM505010 Power Supply

CAPSOM505030 AirMist Starter Pack (2 filters and 1 holder)

CAPSOM505009 AirMist HME Pack (4 filters and 1 holder)

CAPSOM505014 AirMist HME Cartridges (4 Pack)

CAPSOM505015 AirMist HME Cartridges (8 Pack)

CAPSOM505016 AirMist HME Cartridges (12 Pack)

CAPSOM505008 SleepPak Travel Bag

CAPSOM505038 WhisperSoft Muffler Kit: 6' Hose and Muffler

CAPSOM505051 Standard Woven Filters - 2/PK

CAPSOM505052 Hypoallergenic Filters - 2/PK

Item# **CAPSOM505001**

PowerAway® Battery



Item# **CAPSOM505054**

Resupply Pack

Includes:

- Micro WhisperSoft Muffler
- 6ft AirFlex Hose
- AirMist HME Cartridges Pack (4 filters and 1 holder)
- Standard Woven filters (2 pack)



Extended Run Time

- **Long-Lasting Power:** Offers 2-3 nights or 17.7 hours of run time on a single charge at a pressure setting of 9cm, ensuring restful sleep without frequent recharging.
- **Quick Charge Time:** Recharges in less than 4 hours, so it's ready to go when you are.

Reliability and Safety

- **FAA-Approved:** Cleared for in-flight use, making it an essential for travelers
- **Backup Power:** Provides a permanent backup power source during outages, ensuring uninterrupted sleep therapy.
- **Long Life Span:** Capable of 1,000+ full charge/discharge cycles for long-term reliability.

Item# **CAPSOM505053**

Essential Pack

Includes:

- AirMist HME Cartridges Pack (4 filters + 1 holder)
- Micro WhisperSoft Muffler
- (1) 6ft AirFlex Hose
- (1) 4ft AirFlex Hose
- (2) Hose Clips
- Travel Power Supply Plug Adapters
- Standard Woven filters (2 pack)
- Hypoallergenic filters (2 pack)
- Transcend Micro Travel Pouch



Item# **CAPSOM505055**

Essential Bundle

Includes:

- Essential Pack
- White Micro 510
- SleepPak Travel Bag



Item# **CAPSOM505056**

Power Bundle

Includes:

- Essential Pack
- White Micro 510
- PowerAway Battery
- SleepPak Travel Bag



CPAP Cash Items

Snugz CPAP Mask Liners

Item#

CAP7001

Machine-Washable Liners
for Most Full Face Masks
2/pack

CAP7002

Machine-Washable Liners
for Most Nasal Masks
2/pack

- 2 per pack, 90-day supply
- Machine washable
- Fits MOST nasal or full face CPAP masks (Not recommended for ResMed Mirage Quattro or Respironics Amara View and DreamWear)
- Blend of spandex, moisture-wicking polyester and lycra fibers
- Latex-free

Snugz is a machine washable, one-size-fits-most liner for CPAP masks that is designed to help reduce noisy leaks, prevent redness and irritation, and feel great against the skin. The liner helps to grip the skin and keep the mask in place around the nose and mouth, reducing the need for a chinstrap. The improved seal is especially beneficial for high-pressure/high-flow PAP patients and can reduce nightly sleep disturbances and AHI.




Soft Headgear Strap Covers

Item#

CAP2004

- Reduces strap marks
- No bulky and uncomfortable velcro
- Reduces strap irritation on sensitive skin
- Thick padding for extra comfort
- No need for clips – just slip wraps around headgear straps
- Easy and helpful for patients with arthritis in hands or fingers
- Can be machine washed and dried
- Wraps securely grip headgear and slide to adjust
- One size fits all
- Made from 230g Coral Velvet



RemZzzs CPAP Mask Liners

Increases compliance by adding a comfortable and absorbent pad to CPAP masks. Reduces air leaks and protects skin from red marks.

Item#

CAP3001(S,M,L)

Liners for Full Face Resmed,
Respironics, DeVilbiss, and
Fisher & Paykel

OEM# RZ0-FFK1-SM32 (SMALL)
RZ0-FFK2-MD32 (MEDIUM)
RZ0-FFK3-LG32 (LARGE)

CAP3002(S,M,L)

Liners for Full Face Sunset, Fisher &
Paykel FlexiFit, Forma, and Hans Rudolph

OEM# RZ0-FFK4-SM32 (SMALL)
RZ0-FFK5-MD32 (MEDIUM)
RZ0-FFK6-LG-32 (LARGE)

CAP3003(S,M,L)

Liners for Nasal Respironics Comfort
and DeVilbiss Easy Fit

OEM# RZ0-NAK7-SM32 (SMALL)
RZ0-NAK8-MD32 (MEDIUM)
RZ0-NAK9-LG32 (LARGE)

CAP3004(XS, S,M,L)

Liners for Nasal ResMed, Respironics,
Fisher & Paykel

OEM# RZ0-NAK10A-XS32 (EXTRA SMALL)
RZ0-NAK10-SM32 (SMALL)
RZ0-NAK11-MD32 (MEDIUM)
RZ0-NAK12-LG32 (LARGE)

CAP3005(L, XL)

Liners for Respironics
FitLife Whole Face Masks

OEM# RZ0-TF6A-SM32 (SMALL AND LARGE)
RZ0-TF6B-LG32 (EXTRA LARGE)

CAP3006

Liners for Nasal Pillows, Most
Respironics, ResMed, and TapPap

OEM# RZ0-NP7A-UN32

CAP3007

Liners for Hybrid ResMed AirFit/
AirTouch, Respironics Amara View

OEM# RZ0-MF6C-UN32



Sold in
attractive
packaging

CPAP Cash Items

Pur-Sleep CPAP Aromatherapy

Item# CAP5001

VaporClear Aromatherapy

Starter Kit

OEM# VAPOR

30ml Aromatic Refills

Scent/OEM#	Item#
Black Cherry	CAP5004-BCH
Calm	CAP5004-CLM
Clear	CAP5004-CLR
Dark Chocolate	CAP5004-DCH
Deep	CAP5004-DEE
French Vanilla	CAP5004-FRV
Fresh	CAP5004-FSH
Peace	CAP5004-PCE
Spice	CAP5004-SPC
Strawberry	CAP5004-STR


Kit Includes:

- 1 aluminum diffusion stand
- 4 reusable diffusion pads
- 1.5ml bottle of Clear (peppermint, lime)
- 1.5ml bottle of Calm (eucalyptus, lavender, rosemary)

Item# CAP5002

Reusable Diffusion Pads


10/pack • OEM# DP10

Item# CAP5003

Replacement Aluminum
Diffusion Tray


Each • OEM# TRAY

- Aromatherapy products specifically for use with CPAP
- Set up the diffuser right below the CPAP filter intake, add a drop of essential oil to the pad and the fragrance is drawn into the CPAP machine



Quick Tube Release Connector

Item# CAP1016
Each



Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item# BF2200
22mm
5/pack

Attach between the machine and tube for additional protection against bacteria and allergens



Circuit Connectors

Item# RES013
22mm OD x
15mm ID/
22mm OD
10/pack



Item# RES014
22mm OD x
22mm OD
With Lip
10/pack



Oxygen Bleed-In Adapter

Item# RES019
22mm with 6mm
Rotating Stem
10/pack



RES019-CAP
With Stem Cap
10/pack

Swivel Exhalation Port

Item# VEN002
For CPAP and
Ventilator Circuits



Humidifier Tubing Adapter

Item# RES016
22mm ID x 22mm OD • with Lip
10/pack



CPAP Cash Items

CPAPMax Pillow

Item# **CAP4003**

Each



OEM# 15-501R-1-920

Traditional pillows interfere with mask performance and can cause a stiff neck, shifting, leaks, and facial discomfort.



The CPAPMax Pillow offers superb mask accommodation, excellent support, and unparalleled comfort.



Cool 3D Mesh Air-Flow
Reduces heat and perspiration



Fiberfill
Soft microfiber pillow top with hypoallergenic fiberfill

CPAPMax Replacement Cover

CAP4003RC

Each

OEM# 1-630-101-1-933

Replacement Pillow Case

CAP4003C

Each

OEM# 1-626-900-1-02

2 Comfortable Sides:



Cozy and Soft Super Plush Side
Pillow top with hypoallergenic fiberfill



NEW! Cooling Gel Side
Reduces heat and perspiration

Adjust Height to Your Comfort



10-in-1 Flip Pillow

Item# **CAP4006**

Each



The most versatile wedge and bed pillow available!

10 support positions with 1 pillow! Comfort and align your neck, shoulder and lumbar region when used as back support cushion. Use it as a knee or leg pillow wedge or even in an office chair for added support. Filled with a soft hypo-allergenic fiber-fill.

10-in-1 Flip Pillowcase Cover

CAP4006CN Navy
CAP4006CW White
Each



Leg Pillow

Item# **CAP29450R**

Each



New & Improved Patented Shape

Stays put - even when you bend your knees!



Breathable, Washable cover

Plush, soft cover easily removes for washing.



- Restores Your Natural Alignment for a Healthier Night's Sleep
- Provides perfect alignment of back and hips, relieves discomfort in hips and joints

Doctor Recommended for Back, Hip & Knee Comfort

- Helps Align Lower Back
- Stabilizes Hips
- Cushions Knees
- 25% More Leg Support
- Ventilated Cushion & Cover
- Cool Dry Comfort



New Packaging!

Clamshell design conveys high quality of product & efficient use of shelf space

CPAP Cash Items

Contour Incline Sleep Wedge

Each




Height	Item#
7 in	CAP3092407R
10 in	CAP3092410R
12 in	CAP3092412R

The Contour Incline Sleep Wedge helps patients sleep comfortably elevated! The gradual incline supports the head and chest in an upright, angled position.

The Incline Sleep Wedge will help patients sleep comfortably elevated, help aid digestion, and ease breathing and sinus pressure. The upright, angled position will provide comfort while supporting the head and chest. Plus, they'll sleep with their spine properly aligned which helps relieve muscle tension and back pain.

Features & Benefits:

- Helps aid digestion
- Also helps with other conditions often associated with Acid Reflux such as respiratory issues, nasal and sinus congestion, and snoring
- Soothes muscle tension and supports the spine
- Includes a soft, velour zippered cover for easy removal
- Cover is machine washable



Compressed In Package Expands to Full-Size Wedge



Easily keep shelves tidy and uncluttered with compressed packaging solutions

CoolPAP Pillow

Item# **CAP4007**

Each

Works with any style and brand of mask!



Rolls up for travel or storage



New retail packaging!

CoolPAP Pillow Cover
CAP4007C
Each



Deluxe Tubing Wrap

Item# **CAP2002**
6ft Nylon Tubing Wrap
Each

Compare to ResMed Tubing Wrap



- 6ft nylon wrap with fleece lining maximizes heat retention and reduces noise
- Zipper closure for ease of use



Both Tubing Wraps come in attractive retail packaging

Comfort Tubing Wrap

Item# **CAP2001**
6ft Velour Tubing Wrap
Each

Compare to Respironics Tubing Wrap



- 6ft plush velour material adds comfort while preventing rainout
- Zipper closure for ease of use

CPAP Cash Items

CPAP Bedside Storage Bag

Item#
CAP1015

- Well-made and stylish
- Designed to hold CPAP supplies such as a mask and tube
- Weighted lid allows bag to hang off top of nightstand or tabletop



Carry Bag for AirSense 11

Item#
CAP1028-S11

NEW



CPAP HCPCS and Replacement Schedule

HCPCS	Item	Recommended Replacement Schedule
A4604	Heated Circuit Tubing	1 per 3 months
A7027	Combination Oral/Nasal Mask	1 per 3 months
A7028	Replacement Oral Cushion for Combination Oral/Nasal Mask	2 per month
A7029	Replacement Nasal Pillows for Combination Oral/Nasal Mask — Pair	2 per month
A7030	Full Face Mask	1 per 3 months
A7031	Replacement Full Face Mask Cushion or Liner	1 per month
A7032	Replacement Nasal Cushion or Liner	2 per month
A7033	Replacement Pillows — Pair	2 per month
A7034	Nasal or Pillows Mask	1 per 3 months
A7035	Headgear	1 per 6 months
A7036	Chinstrap	1 per 6 months
A7037	Tubing	1 per 3 months
A7038	Filter — Disposable	2 per month
A7039	Filter — Nondisposable	1 per 6 months
A7044	Oral Interface	
A7045	Replacement Exhalation Port	
A7046	Replacement Water Chamber	1 per 6 months
A9279	Monitoring Feature/Device	
E0470	BiLevel Device without Backup Rate	
E0471	BiLevel Device with Backup Rate	
E0561	Humidifier - Non-Heated	
E0562	Humidifier - Heated	
E0601	CPAP Device	

Due to the rapid changes in reimbursement, this information is for reference only. Sunset Healthcare Solutions cannot guarantee its comprehensiveness or accuracy. Sunset urges you to contact a reimbursement specialist for up-to-date information and advice.

Oxygen Products

Oxygen	Page
Pulse Oximeter	40
SpO2 Sensors	41
Firesafe	42
Portable Concentrator	43
Concentrators	44
Masks	45
Cannula and Cannula Accessories	46-47
Tubing	47
Microdial Flowmeter	48
Regulators and Conserver	49
Oxygen HCPCS Chart	50
Cylinder Carriers	51-52
Connectors	53
Accessories	53-55
Cleaners	55
Analyzers	56
Celon Wet Bands	56
Bacteria and HEPA Filters	57
Concentrator Cabinet Filters	58-64
Caire/AirSep/SeQual	58
CHAD	59
Covidien Nellcor Puritan Bennett	59
Dynarex	59
OxyGo	59
DeVilbiss/Drive	60
Invacare	61
Rhythm, Jumao, Attendant, Lifestyle, Medacure	62
Respironics	63
Healthdyne	64
React Health/3B	64
Nidek	64
Repair Parts	65-69

Pulse Oximeter

Finger Pulse Oximeter

Item# RES5102


- Oversized 1.3" OLED display face is the largest of any finger pulse oximeter
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- CE, ISO and FDA approved
- Includes lanyard and manual



SpO2 Sensors

Disposable SpO2 Sensor Neonatal/Adult

Item# **RES7100**
Similar to Nellcor MAX-N
20/box



Neo-M-1.5 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7200**
1.5ft cord • Non-replaceable
20/box



Neo-M-3 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7201**
3ft cord • Non-replaceable
20/box




Coban Tape Roll

Item# **RES7202CT**
Use for sensor life
extension or sensor
cord stability
Each roll 1" x 177"



Infant-M-1.5 Disposable SpO2 Sensor Infant (3-20 kg)

Item# **RES7203**
1.5ft cord
20/box



Pediatric-M-1.5 Disposable SpO2 Sensor Pediatric (10-50 kg)

Item# **RES7204**
1.5ft cord
20/box



SpO2 Adapter Cable

Item#
RES7200R97-10 10ft
RES7200R97-14 14ft
Each



NEW

SpO2 Adapter Cable

Item#
RES7200R57-4 4ft
RES7200R8-10 10ft
Each



NEW

Firesafe

Firesafe Thermal Shut-Off Valve



Item#
RES010
Each

Bidirectional



- Now smaller
- Laser etched expiration date
- Human readable UDI
- Designed to last 5 years

- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

Fatalities per 100,000 patients



FireSafe Kits

Item#
RES010KIT
2 Firesafe Valves
1 Tubing Connector



Item#
RES010CANNULAKIT
2 Firesafe Valves
1 Tubing Connector
1 Adult Soft Cannula



Item#
RES010VAKIT
2 Firesafe Cannula Valves
1 Tubing Connector
1 Cannula with Female Port
1 7ft Oxygen Supply Tube



FireSafe Tubing Connector

Item#
RES010TC
100/case

RES010TC-10
10/case



FireSafe Nozzle

Item#
RES010FND-20
20/pack



Belluscura Oxygen Concentrator

X-PLOR® Portable Oxygen Concentrator

Item#

RES52051001

X-PLOR

Lightweight and compact, the X-PLOR® is ideal for patients who want to replace their heavy oxygen tanks without missing a moment. The X-PLOR is also Bluetooth®-capable and can be connected to mobile devices through Belluscura's newly developed Nomad Biometrics™ App. With a rechargeable battery, easy-to-use interface, and high oxygen output, supplementary oxygen has never been more approachable or affordable.

What's Included

- X-PLOR® Unit
- 8-Cell Battery
- User-Replaceable Sieve Bed Cartridge
- Air Intake Filter
- AC Power Adaptor
- DC Power Adaptor
- Carry Bag With Handle & Shoulder Strap

Accessories

Item#

RES50073 User-Replaceable Sieve Bed Cartridge

RES50188 8-Cell Battery

RES53012 DC Power Adaptor

RES53011 AC Power Adaptor

RES52191 Carry Bag

RES52291 Backpack



Features

- Oxygen output 82%-92% purity
- User-Replaceable Sieve Bed Cartridge
- Airgonomically™ Designed to be lightweight and compact for easy travel
- Long-Lasting, Direct Charge Battery: Up to 4.5 hours of use when fully charged
- Pulse flow settings 1 through 5 delivering 215 to 900 ml/min
- No Breath Detected (NBD) mode is a battery-saving mode that runs at Auto-pulse15 BPM at the units current setting
- FAA Compliant

Drive DeVilbiss Concentrators

Drive DeVilbiss Oxygen Concentrators

Item#

RES525DS

5 Liter

OEM# 525DS

Item#

RES1025DS

10 Liter

OEM# 1025DS



Drive DeVilbiss iG02 Portable Oxygen Concentrator

Item#

RES125D

iG02

Item#

RES125D-XB

iG02 with

Extra Battery



CAIRE Oxygen Concentrators **NEW**

CAIRE HOMESTYLE™ Aspen Stationary Oxygen Concentrator

Item#
RESSCE1000

Your O₂ Fleet Solution

Features

- Continuous flow 0.5-5 LPM
- Up to 93%±3% oxygen purity
- Sound less than or equal to 40 dB(A)*
- Power supply: 110V~, 60Hz
- Power consumption: 320VA
- Weight: 32 lbs.
- 15 inches L x 9 inches W x 23 inches H
- Built-in cord caracol
- Elegant and compact frame



- 5-year warranty on main unit and sieve beds
- 1-year warranty on accessories

CAIRE NewLife® Intensity® 10 High Flow Stationary Oxygen Concentrator

Item#
RESAS0994

Newlife Intensity 10L SOC Dual Flow
120V/60Hz

RESAS0991

Newlife Intensity 10L SOC Single
Flow 120V/60Hz



MADE IN
USA

Time-proven design

- Industry-leading reliability—greater than 99% reliability rate out of the box and less than 2%* warranty repair rate
- Comprehensive warranty—3-year warranty on the device
- Durability—engineered with an oversized compressor and sieve beds to minimize service

CAIRE FreeStyle Comfort 4-Button Portable Oxygen Concentrator

Item#
RESAS2004

Freestyle POC w/ 8-cell Battery

RESAS2001

Freestyle POC w/ 16-cell Battery

RESAS2000

Freestyle POC w/ Two 8-cell Batteries



MADE IN
USA

- 5-year warranty on main unit
- 3-year warranty on sieve beds
- 1-year warranty on accessories

Smart Technology to Ensure Saturation

- **autoSAT®** maintains a consistent pulse dose volume for every breath as the patient's respiratory rate changes
- **autoDOSE®** automatically delivers bolus of O₂ @ 20 BPM if no breath is detected within 15 secs
- **UltraSense®** offers sensitive breath detection that ensures a pulse of oxygen is delivered with each breath and when oxygen delivery is most critical
- **myCAIRE®** Telehealth app solution

Flexible flow rates to adjust to changing needs

With 5 settings and up to 1050mL of oxygen per minute, the FreeStyle Comfort offers flexibility in treatment options covers broad range of prescriptions.

CAIRE Eclipse 5® Transportable Oxygen Concentrator

Item#
RES6900LN

Eclipse 5 TOC with 2 power cartridges

RES6900BT

Eclipse 5 TOC with 1 power cartridge

MADE IN
USA



- 3-year warranty on main unit
- 1-year warranty on sieve beds and accessories

Smart Technology to Ensure Saturation

- **autoSAT®** Automatic adjustment to patient's respiratory rate
- **Adjustable Rise Time** Maximizing inhaled oxygen
- **UltraSense®** breath detection technology
- **myCAIRE®** Telehealth app solution

Continuous or pulse dose delivery to meet patients' needs

With 6 continuous flow settings up to 3 LPM and 9 pulse flow settings up to 192 mL bolus, the Eclipse 5 serves a wide range of patients.

Oxygen Masks

Simple O2 Mask

With 7ft Safety Tubing

Item# **RES2100**

Adult Size

50/case

Item# **RES2200V**

Pediatric Size

50/case

HCPCS# A4620



- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Partial Rebreather / Non-Rebreather O2 Mask with 7ft Tubing

Oxygen High Concentration
Mask with 7ft Safety Tubing
and Reservoir Bag

Item# **RES2101V**

Adult Size

50/case

OEM# 8120-7-50

HCPCS# A4620



- 5-15 LPM
- Safety vent
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Non-Rebreather O2 Mask

With 7ft Safety Tubing
& Reservoir Bag

Item# **RES2102**

Adult Size

50/case

Item# **RES2202V**

Pediatric Size

50/case

HCPCS# A4620



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape

Oxygen Cannula

Cannula with Tubing

Adult Cannula with Tubing

Length	Material	Item#	Qty
4ft	Standard	RES1104	50/case
4ft	Soft	RES1104S	50/case
7ft	Standard	RES1107	50/case
7ft	Soft	RES1107S	50/case
15ft	Standard	RES1115	50/case
25ft	Standard	RES1125	25/case
25ft	Soft	RES1125S	25/case

HCPCS# A4615

- Clean, fresh look
- 1-6 LPM
- Curved nasal prongs for added comfort
- 6-channel safety tubing designed to prevent kinks
- Soft, lightweight material is memory free for long-term comfort



Item#

RES1207V

Pediatric Cannula
with 7ft Supply Tube
50/case



RES1307V

Infant Cannula
with 7ft Supply Tube
50/case




Cannula without Tubing

Adult Cannula without Tubing

Port	Material	Item#	Qty
Male	Standard	RES1100	50/case
Male	Soft	RES1100S	50/case
Female	Soft	RES1100FS	50/case

HCPCS# A4615



Male port

Female port

Dual-Lumen (Demand) Cannula

Item# **RES4147V**

7ft

25/case

HCPCS# A4615



- 1-6 LPM
- Designed for use with dual port oxygen conservation delivery systems
- Dual tubing set provides independent tubes for oxygen delivery and sensing of patient breathing
- Smoothbore tubing

High Flow Soft Cannula

Length	Item#	Qty
4ft	RES1104SHF	25/case
7ft	RES1107SHF	25/case
15ft	RES1115SHF	25/case
25ft	RES1125SHF	25/case

HCPCS# A4615



- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

Oxygen Cannula

Humidiflo High Flow Nasal Cannula

10/case

Size	Item#	With F&P Adapter	Comparable OEM#
Extra Small	HFN001AXS	HFN001AXSA	F&P Small OPT942
Small	HFN001AS	HFN001ASA	F&P Medium OPT944
Medium	HFN001AM	HFN001AMA	F&P Large OPT946
Large	HFN001AL	HFN001ALA	N/A



F&P Adapter **HFN001-FPA**

With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerated due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O2 concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respirationics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)


See LUISA Ventilator on page



Secure Cannula Stickers

Item# **RES080**
25 pairs/case

- Hold cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



Comfort Cannula with Ear Cushions

Item# **RES1107EC**
7ft Cannula with Ear Cushions
25/case
HCPCS# A4615



Item# **RES020**
Cushions Only
50 Pairs/case

- Add cushions to existing cannula

Oxygen Tubing

Oxygen Supply Tube



Length	Item#	Qty
7ft	RES3007	50/case
15ft	RES3015	50/case
25ft	RES3025	25/case
40ft	RES3040	20/case
50ft	RES3050	20/case

HCPCS# A4616

- Crush resistant with low memory to prevent kinking
- Great durability and comfort level

Oxygen Supply Tube — Green




Length	Item#	Qty
7ft	RES3007G	50/case
15ft	RES3015G	50/case
25ft	RES3025G	25/case
40ft	RES3040G	20/case
50ft	RES3050G	20/case

HCPCS# A4616

- Latex-free
- 1-6 LPM
- 6 channels

Microdial Flowmeter

Microflow Microdial Flowmeter



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure	Inlet Pressure Range	Flow Range*
RES8160048			A
RES8160050	400 kPa (58 psi)	350 – 550 kPa	E



DISS Inlet • DISS Outlet

Item#	Inlet Pressure	Inlet Pressure Range	Flow Range*
RES8160049	400 kPa (58 psi)	345 – 500 kPa	A
RES8160051	350 kPa (51 psi)	250 – 450 kPa	E

*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
Range A	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
Range E	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

Assured, repeatable accuracy to enable safe, effective weaning.

What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.



Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned

Oxygen Regulators and Conserver

Mini Oxygen Regulator CGA 870 Style



"T" Handle for
easy mounting:
no tools needed

Barbed hose
outlet for easy
connection to
the oxygen
tubing

For use with E size or
smaller cylinders: A, B, C,
D, JD, E, M2, M4, M6, ML6,
M7, M9, M15, M22, or M24.



Item# **RES042**
Adult • Barbed
0-8 LPM • **Each**
0, .5, 1, 1.5, 2, 2.5,
3, 4, 5, 6, 7, 8 LPM

Item# **RES043**
Adult • Barbed
0-15 LPM • **Each**
0, .5, 1, 1.5, 2, 3, 4,
5, 6, 8, 10, 15 LPM

Item# **RES046**
Pediatric • Barbed
0-4 LPM • **Each**
0, $\frac{1}{32}$, $\frac{1}{16}$, $\frac{1}{8}$,
 $\frac{1}{4}$, $\frac{1}{2}$, $\frac{3}{4}$, 1,
 $1\frac{1}{2}$, 2, 3, 4 LPM

HCPSC# E1353



Pediatric models have a green dial cap

Mini Oxygen Regulator CGA 540 Style



Threaded
DISS outlet

For use with M60 and
larger cylinders: M60, M90,
M, MM, M122, M250, H, J, T
or other cylinders with CGA
540 style valve.

Item# **RES044**
Adult • DISS
0-8 LPM • **Each**
0, .5, 1, 1.5, 2, 2.5, 3,
4, 5, 6, 7, 8 LPM

Item# **RES045**
Adult • DISS
0-15 LPM • **Each**
0, .5, 1, 1.5, 2, 3, 4, 5,
6, 8, 10, 15 LPM

Item# **RES047**
Pediatric • DISS
0-4 LPM • **Each**
0, $\frac{1}{32}$, $\frac{1}{16}$, $\frac{1}{8}$,
 $\frac{1}{4}$, $\frac{1}{2}$, $\frac{3}{4}$, 1,
 $1\frac{1}{2}$, 2, 3, 4 LPM

HCPSC# E1353



Both Regulators Feature

- Pressure gauge reads to 3000 PSI
- Easy-click flow selector knob
- 2-year warranty



**Private label
laser etching
now available!**

Ask your rep for details.

Oxygen Regulator CGA 540 Style Long Stem



Item# **RES045L**
Adult • DISS
0-15 LPM • **Each**
0, .5, 1, 1.5, 2, 3, 4, 5,
6, 8, 10, 15 LPM

Bonsai Pneumatic Conserver

Item# **RES049**
Single Lumen •
Barbed Outlet
Each



- Delivers a uniform oxygen pulse in increments: 1, 2, 3, 4, 5, and 6
- Two continuous flow (CF) settings preset at 2 LPM and 4 LPM
- Delivers up to a 6:1 savings ratio from 14-40 breaths per minute
- Senses breath/delivers oxygen in first half of inspiratory cycle
- Operating pressure range: 200-3000 PSI
- Light weight and quiet operation
- Easy to operate single selector knob
- Compact design fits easily in virtually any oxygen carry bag
- Limited 2-year warranty

Bags for Regulators and Conserver



Item#
EC1012C
10×12"
1.5mil, Clear
1000 bags/box

Sample Oxygen Cylinder Specifications*

Use the handy chart below to determine which accessories to purchase for your oxygen cylinder.

* This information is intended to be used as a guide. Dimensions may vary by manufacturer. Sunset Healthcare Solutions does not sell oxygen tanks.



Older Name		A	B			C	D	JD	E		
Newer Name	M2	M4	M6	ML6	M7	M9	M15	M22	M24	M60	M/MM/M122
Diameter (in.)	2.5	3.2	3.2	4.3	4.3	4.3	4.3	5.3	4.3	7.3	8
Height (in.)	5.3	8.5	11.5	7.6	9.1	11	16.5	16.5	25.5	23	36
Empty Weight (lb)	0.7	1.6	2.2	2.8	3.3	3.7	5.3	8	7.9	22.3	39.5
Capacity (L)	42	113	164	165	198	255	425	640	680	1738	3455
Transport Method	CSBM6, CSBM6BP	CSBM6, CSBM6BP	CSBM6, CSBM6BP, CSBM6Flex, CSBM6-Carrier	CSBML6M7, CSBML6-Hardshell	CSBML6M7	CSBM9, CSBM6BP, CSBM6Flex, RES038B	CSBD, CSBDDBP, CSBD1Roller, CSBD2Roller, RES038B, RES060		CSBE, CSBD1Roller, CSBD2Roller, RES038B, RES060	Not Portable	Not Portable
Regulator Type	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES044, RES045	RES044, RES045

Oxygen HCPCS

HCPCS	Item
A4615	Cannula, nasal
A4616	Oxygen tubing, per foot
A4620	Oxygen mask, variable concentration
E0424	Stationary compressed gaseous oxygen system, rental
E0431	Portable gaseous oxygen system, rental
E0443	Portable oxygen contents, gaseous (for use only with portable gaseous systems when no stationary gas or liquid system is used), 1 month supply = 1 unit
E0445	Oximeter device for measuring blood oxygen levels non-invasively
E0560	Humidifier, durable for supplemental humidification during IPPB treatment or oxygen delivery
E1353	Oxygen supply regulator
E1354	Oxygen supply wheeled cart
E1355	Oxygen supply stand/rack
E1390	Oxygen concentrator, single delivery port, capable of delivering 85% or greater oxygen concentration at the prescribed flow rate
E1391	Oxygen concentrator, dual delivery port, capable of delivering 85% or greater oxygen concentration at the prescribed flow rate
E1392	Portable oxygen concentrator, rental
K0738	Portable gaseous oxygen system, rental; home compressor used to fill portable oxygen cylinders; includes portable containers, regulators, flow meter, humidifier, cannula or mask and tubing

Cylinder Bags

Comfort Shoulder Bags

Cylinder	Item#
ML6/M7	CSBML6M7
C/M9	CSBM9
B/M6	CSBM6
D	CSBD

- Comfort flex shoulder strap
- Double zipper for easy access
- Machine washable



CASH ITEM



CSBML6M7



CSBM9



CSBM6



CSBD

Flex O2 Bag

Item#
CSBM6Flex
Reconfigurable Bag
for M6 and M9
Cylinders
Each

- Adjusts to 5 carrying modes



CASH ITEM



Backpack

Cylinder	Item#
B/M6/C/M9	CSBM6BP
D	CSBD8BP

- Fully adjustable and comfortable backpack straps
- Machine washable
- Converts to either backpack or shoulder bag



CASH ITEM



Wheelchair/Scooter Carrier

Item# **CSBE**
For E Cylinders
Each

- Multi-part stabilizer straps for easy chair attachment
- Machine washable



CASH ITEM



Cylinder Carriers

ML6 Hardshell Tote

Item#
CSBML6-Hardshell
For Homefill Cylinders
Each



CASH ITEM



Double Tank Rolling/Shoulder Bag

Item#
CSBD2Roller
For D/E Cylinders
Each

- May also be used with shoulder strap



CASH ITEM



Oxygen Cylinder Carrier

Item#
CSBM6-Carrier
Soft Polyester Tote
for M6 Cylinders
Each

- Tote handles
- Padded shoulder strap
- Fastening strap for mobility
- Holds up to 8 M6 tanks



Cylinder Rack

Item# **RES038B**
For D, E, or M9 Cylinders —
Lightweight, Durable
Plastic Composite Rack
Each



D/E Cart

Item#
RES060
4/case

- Accommodates one D or E tank
- Sturdy yet lightweight (4.18 LBS welded metal construction)
- Adjustable height handle with PVC grip
- Large EVA Foam wheels have just the right amount of resilience for good traction and maneuverability
- Elegant, grey satin finish
- Bottom fork rest with rubber caps



Oxygen Connectors and Accessories

Christmas Tree Connector

Item# **RES001**

Green without Swivel
10/pack



Item# **RES002**

Green with Swivel
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Oxygen Tubing Y-Connector

Item# **RES017**

Fits 3/16in and 1/4in ID Tubing
(All 3 sides are the same)
10/pack



- Connects 2 oxygen sources for higher flows

Oxygen In Use Signs

Item# **RES003**

6in x 9in
With String
50/pack

- Double-sided (English and Spanish)



Washer with Rubber Ring

Item# **RES035**

Brass
25/pack



Oxygen Accessory Kit

Item# **RES1000**

RES002 Christmas Tree Connector - Swivel
RES012 Tubing Connector - Non-Swivel
RES017 Tubing Y-Connector
RES018 Tubing Connector - Swivel
RES032 Wrench
RES035 Brass Washer

Each



CUSTOMIZABLE KIT

- Order all common oxygen accessories in one easy-to-handle kit
- Contents can be customized to fit your needs



Oxygen Tubing Connector — Non-Swivel

Item# **RES012**

Fits 3/16in and 1/4in
ID Tubing
10/pack



- Barbed ends ensure secure, leak-free connections

Oxygen Tubing Connector — Swivel

Item# **RES018**

Male/Male
10/pack

RES018-25 25/pack

RES018KIT-1
Female/Male
Each

- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing



Oxygen Accessories

02 Starter Kit

Item# RES02KIT

RES1107 Cannula with 7ft Tube x2
 RES3015 15ft Supply Tube x2
 RES002 Connector — Swivel
 RES003 Oxygen In Use Sign
 RES012 Tubing Connector — Non-Swivel
 RES017 Tubing Y-Connector
 RES018 Tubing Connector — Swivel
 RES032 Wrench
 RES035 Brass Washer



CUSTOMIZABLE KIT



Wrenches and Chains

Item# RES031

With Chain **10/pack**

Item# RES032

Without Chain **10/pack**



- For use with post valve oxygen cylinders: A, B, C, D, JD or E tanks
- Made of tough, stress-resistant nylon

Item# RES031M

Metal Wrench with Chain for Post Valve Oxygen Cylinders and CGA 540 Style Regulators **5/pack**



Item# RES031-CHAIN

12in Ball Chain with Coupling Connector **50/pack**

- Add chain to an oxygen valve wrench



Item# RES031-BW

Cylinder Wrench with Two-Ring Bungee and Two Yoke Washers **10/pack**



Disposable Bubble Humidifier 500 ml

Item#

RES050

Up to 15 LPM • 6 psi • 500 ml

50/case

HPCPS# E0560

OEM# HUM001



- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient

Humidifier Connectors

Item# HUM008

L-Connector 12in • **10/pk**

Item# HUM008L

L-Connector 21in • **10/pk**

Item# HUM008S

Straight Connector 11.625in • **10/pk**



Water Traps

Item# RES037

Push On End Caps • **Each**

Item# RES037T

Threaded End Caps • **Each**



- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

Accessories and Cleaners

Homefill Tubing Connector

Item# **RES039**

10/pack

- Compatible with Invacare, Drive, and Respironics Homefill machines
- Durable, coiled tubing



Homefill Racks

Item#

RES039-RACKB

For Invacare Perfecto2

4/case



RES039-RACKW

For Invacare Platinum

4/case



RES039-RACKD

For DeVilbiss

525/1025 Series

4/case

Easily attaches onto concentrator and is designed to conveniently hold Homefill and Transfill oxygen systems.



Oxygen Liter Meter

Item# **RES034**

0-10 Liter/Minute

Each



- Oxygen pen
- Works with tanks and concentrators

Germicidal Deodorizing Cleaner

Item# **RES6100**

22oz spray bottles

12/case

- Ready-to-use spray cleans, disinfects and deodorizes
- Formulated for use on hard, non-porous surfaces
- Bactericidal, Tuberculocidal, Fungicidal, Mildewstat
- Does not contain abrasives, bleach, phosphates or phosphorous
- Effective against SARS-CoV-2 (COVID 19 virus) on hard non-porous surfaces



Germicidal Deodorizing Cleaner — Gallon

Item# **RES6102**

1 Gallon

- No-scrub, hospital strength germicidal cleaner
- Cleans, deodorizes, and disinfects surfaces, equipment and non-critical instruments



Oxygen Accessories

Smart Check Analyzer

Item#
RES7003 **NEW**
Each

Portable, handheld oxygen analyzer designed specifically for testing and maintaining oxygen concentrators.

The ultimate oxygen + flow + pressure analyzer.

The Smart Check is an all-in-one oxygen analyzer for testing continuous and pulse (POC) concentrators. Measure oxygen, flow, pressure and pulse volume with a single device. Use the Smart Check to also test CPAP, BIPAP and nebulizers.

The Smart Check is the only device that also offers a built-in flashlight, backlight, and integrated tube management system.

With its fast temperature response, SmartCheck quickly adapts to temperature changes, ensuring accurate results in less time.

Offering all these features at a very competitive price, SmartCheck delivers exceptional performance at an incredible value.



	Smart Check	UltraMax O2	Pro2 Check Elite
Measure Oxygen & Flow	✓	✓	✓
Test Pressure	✓	✓	✓
Quick Touch Pressure Check	✓	✓	✗
Calibrated for Life	✓	✓	✗
Test Pulse Concentrators	✓	✗	✗
Test CPAP + Nebulizer	✓	✗	✗
Flashlight	✓	✗	✗
Backlight	✓	✗	✗
Tube Management	✓	✗	✗
Fast Environmental Response	✓	✗	✗
Cost	\$	\$\$\$	\$\$

- Features:**
- LED Display
 - Calibration Check Button
 - Sample Inlet
 - Flashlight
 - Sample Outlet
 - Battery Cover

Celon Wet Bands

Create a perfect seal to protect oxygen tank valves

Easy to apply, safe and effective. Biodegradable and sustainable material.

- Eliminates dust, dirt & dampness
- Hand applied
- Heat-free solution
- Creates a secure and wrinkle-free seal

Item#
WB66120
66 x 120mm
for CGA 540 Valves
Bucket of 1,500



Item#
WB30540
30.5 x 40mm
for CGA 870 Valves
Bucket of 11,500



Item#
WB30525
30.5 x 25mm
for CGA 870 Valves
Bucket of 18,000



- Professional finish
- Environmentally friendly
- Biodegradable

Oxygen Bacteria And HEPA Filters

6in HEPA Capsule with Removable Sticker

Item# **BFH6S**

Threaded

Each



EverFlo Intake Filter

Item# **BF1001**

Each

Respironics OEM# 1038831



Extended Life Bacteria Filter

Item# **BF5000**

Nonthreaded

Each



Extended Life Intake Bacteria Filter

Item# **BF1103**

Each



Final Bacteria Filter

Item# **BF516**

1/4in Barbed on Both Sides

10/pack



Respironics SimplyGo Style Inlet Filter

Item#

BF8000

Each



HEPA Filter for React Health/3B Stratus

Item# **BF6020**

HEPA Filter

Each

OEM# STR1020



Intake Bacteria Filter

Item# **BF2222**

22mm fittings

5/pack



Extended Life Bacteria Filter with Muffler

Item#

OF6011

Each



Extended Life Intake Bacteria Filter

Item#

OF1025D

Each



Intake Bacteria Filter

Item# **BF7000**

Jumao JMC5A, Attendant 8X899,
Lifestyle LM5A, Medacure OC5L,
Rhythm LM5CA, Drive DECO2 5L
Concentrators

Each




Final Bacteria Filter — 1/4in to 3/8in

Item# **BF1438**

1/4in to 3/8in Stepped Barb
on Both Sides

10/pack



Final Bacteria Filter — 3/16in

Item# **BF316**

3/16 Stepped Barb and 1/4in
Straight Barb on Both Sides

10/pack



Millennium M10 Style HEPA Filter

Item# **BF1002**

Long Life Filter

Each



Caire/AirSep/SeQual

SeQual Integra



Item# **OF11005**
Foam Cabinet Filter
10/pack
75/8" x 31/2"-43/4 x 3/4"
OEM# 2185



Item# **BF5000**
Extended Life
Bacteria Filter
Each
OEM# 2033



Item# **BF516**
Final Bacteria Filter
1/4in Barbed
10/pack

AirSep Newlife, QuietLife, Intensity



Item# **OF8002**
Foam Cabinet Filter
10/pack
51/4" x 37/8" x 7/8"
OEM# F1-50-08-05



Item# **BF1438**
Final Bacteria Filter
1/4in to 3/8in
Barbed
10/pack
OEM# F1-30-03



Item# **OF8001**
Felt Filter
10/pack
2" OD x 5/8" ID x 1"
OEM# MI161-1



KIT

Item# **OF8002SP**
Newlife Kit
2 - OF8002
1 - BF1438
1 - OF8001
Each

CAIRE Companion 5



Item# **BF1103**
Extended Life Intake
Bacteria Filter
Each
OEM# 2185

SeQual Eclipse



Item# **OF11006**
Foam Cabinet Filter
10/pack
21/2" x 21/2" x 3/4"

CHAD

Total 02



Item# **BF2222**
Intake Bacteria Filter
22mm fittings
5/pack
Healthdyne OEM# 100-6631-00
Respironics OEM# H620



Item# **OF3003**
Foam Intake Filter
for use with BF2222
25/8" OD x 3/4" ID x 23/8"
10/pack
OEM# H613

Dynarex and React Health

Resp-02 and PHOENIX 5L



Item# **OF6013-10**
Foam Cabinet Filter
5.75" x 3.625"
10/pack



Item# **OF6011**
Air Intake Filter
Each
OEM# RPOC-14AIF

Covidien Nellcor Puritan Bennett

590 and 592A



Item# **OF6001**
Foam Cabinet Filter
63/4" x 43/4" x 3/8"
10/pack
OEM# 492672



Item# **BF5000**
Extended Life
Bacteria Filter
Each



Item# **BF316**
Final Bacteria Filter
3/16 Stepped Barb and
1/4in Straight Barb
10/pack
OEM# 492141

KIT

Item# **OF6001SP**
590 and 592 Kit
2 - OF6001
1 - BF5000
1 - BF316
Each

OxyGo

OxyHome



Item# **OF8006-10**
Foam Cabinet Filter
63/4" x 41/4" x 3/8"
10/pack
OEM# 1170-5401



Item# **BFH6S**
HEPA Filter
6in Capsule
Each
OEM# 2000806

DeVilbiss/Drive

SolAiris III, 303 Series



Item# OF5004
Foam Cabinet Filter
33/4" x 33/4" x 1/2"
10/pack
OEM# 30302605



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651



Item# BF5000
Extended Life Bacteria Filter
Each

SolAiris V, 505, 515 Series



Item# OF5005SP
SolAiris V Kit
2 - OF5005
1 - OF6011
1 - BF1438
Each



Item# OF5005
Foam Cabinet Filter
5 1/2" x 35/8" x 5/8"
10/pack
OEM# 50502694



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651



Item# BF5000
Extended Life Bacteria Filter
Each



Item# OF6011
Extended Life Bacteria Filter—22mm fittings
Each
Use if SolAiris was manufactured before Nov '03

1025 Series



Item# OF1025D
Extended Life Bacteria Filter
Each



Item# OF5004
Foam Cabinet Filter
33/4" x 33/4" x 1/2"
10/pack
OEM# 303DZ-605



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651


525 DS Series



Item# OF5004
Foam Cabinet Filter
33/4" x 33/4" x 1/2"
10/pack
OEM# 303DZ-605



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651



Item# OF6011
Extended Life Bacteria Filter—22mm fittings
Each

All MC Units



Item# OF5001
Foam Cabinet Filter
10/pack
6 1/8" x 3 1/2" x 1/4"
OEM# 444-514



Item# BF5000
Extended Life Bacteria Filter • **Each**



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651

Invacare

Mobilaire 3 and 5



Item# **OF9001SP**
Mobilaire Kit
2 - OF9001
1 - BFH6S
1 - BF516
Each



Item# **OF9001**
Foam Cabinet Filter
6 1/2" x 3 3/8" x 1/2"
10/pack
OEM# 2000489



Item# **BFH6S**
HEPA Filter
6in Capsule
Each
OEM# 2000806



Item# **BF516**
Final Bacteria Filter
1/4in Barbed
10/pack
OEM# P0750X030

Perfecto2



Item# **OF9004SP**
Perfecto2 Kit
2 - OF9004
1 - BFH6S
1 - BF516
Each



Item# **OF9004**
Foam Cabinet Filter
7" x 3 3/4" x 1/2"
10/pack
OEM# 1143492



Item# **BFH6S**
HEPA Filter 6in Capsule
Each
OEM# 2000806



Item# **BF516**
Final Bacteria Filter
1/4in Barbed
10/pack
OEM# P0750X030

Platinum 5 and 10



Item# **OF9002SP**
Platinum Kit
2 - OF9002
1 - BFH6S
1 - BF516
Each



Item# **OF9002**
Foam Cabinet Filter 9" x 2 7/8" x 1/2"
10/pack
OEM# 1107412



Item# **BFH6S**
HEPA Filter 6in Capsule with Sticker
Each
OEM# 2000806



Item# **BF516**
Final Bacteria Filter 1/4" Barbed
10/pack
OEM# P0750X030



Item# **HUM008S**
Straight Humidifier Connector
w/ 9in Tube
10/pack

HomeFill



Item# **OF9001**
Foam Cabinet Filter
6 1/2" x 3 3/8" x 1/2"
10/pack
OEM# 2000489



Rhythm HC, Jumao, Attendant, Lifestyle, Medacure

Rhythm LM5CA, Jumao JMC5A, Attendant 8X899, Lifestyle LM5A, Medacure OC5L

Compatible with

- **Rhythm LM5CA**



Item# **OF7001-10**
Foam Cabinet Filter
10/pack
OEM# LM209069

Compatible with

- **Jumao JMC5A**
- **Attendant 8X899**
- **Lifestyle LM5A**
- **Medacure OC5L**
- **Drive DEC02**



Item# **OF7000**
Foam Cabinet Filter
10/pack

Compatible with

- **Rhythm LM5CA**
- **Jumao JMC5A**
- **Attendant 8X899**
- **Lifestyle LM5A**
- **Medacure OC5L**



Item# **BF7000**
Intake Filter
Each
OEM# LM209066




Item# **BF516**
Final Bacteria Filter,
1/4in Barbed
10/pack
OEM# LM209068

Rhythm

Rhythm HC LM5BA



Item# **OF8000**
Foam Cabinet Filter
10/pack
OEM# LM5BA-ITF



Item# **OF6011**
Extended Life
Bacteria Filter —
22mm fittings
Each
OEM# LM5BA-ILF



Item# **BF516**
Final Bacteria Filter,
1/4in Barbed
10/pack
OEM# LM5BA-BF

Rhythm LM10A



Item# **BF5000**
Extended Life
Bacteria Filter
Each
OEM# LM10A-ILF

Rhythm P2, P2-E6, and P2-E7



Item# **OF7002-5**
Cotton Intake Filter
5/pack
OEM# P2FC-1

Respironics

EverGo



Item#
OF1005SP
EverGo Kit
2 - OF1005
1 - BF316
Each



Item# **OF1005**
Foam Cabinet Filter
3" x 2" x 3/8"
10/pack
OEM# 900-150



Item# **BF316**
Final Bacteria Filter
3/16in Stepped Barb
and 1/4in Straight Barb
10/pack
OEM# H621

SimplyGo and SimplyGo Mini



Item# **BF8000**
Inlet Filter for SimplyGo
Each
OEM# 1095248



Item# **RES1008**
Compressor Filter for
SimplyGo Mini
Each

Millennium



Item#
OF1001SP
Millennium Kit
2 - OF1001
1 - BF2222
1 - OF3003
1 - BF316
Each



Item# **OF1001**
Foam Cabinet Filter
8 3/16" x 4" x 3/8"
(1/4" tabs)
10/pack
OEM# H619



Item# **BF2222**
Intake Bacteria Filter
22mm fittings
5/pack
OEM# H620



Item# **OF3003**
Foam Intake Filter
for use with BF2222
2 5/8" OD x 3/4" ID x 2 3/8"
10/pack
OEM# H618



Item# **BF316**
Final Bacteria Filter
3/16 Stepped Barb,
1/4in Straight Barb
10/pack
OEM# H621



Item# **OF1002**
Felt Filter
2 1/8" OD x 1/2"
6/pack



Item# **BF1002**
Millennium M10 Style
Long Life Filter
Each

EverFlo



Item#
OF1003SP
EverFlo Kit
1 - BF1001
1 - BF316
Each



Item# **OF1006SP**
Cabinet Filter Kit
(No Cover)
2 Filters, 1 Casing
OEM# 1055625



Item# **OF1006**
Cabinet Filter
5 11/16" x 4 1/2" x 7/16"
10/pk



Item# **OF1006-CASE**
Foam Cabinet Filter
Casing
10/pack



Item# **BF1001**
EverFlo Intake Filter
Each
OEM# 1038831

Item#
HUM008
Humidifier
L-Connector
With 12in
Tube
10/pk



Item# **BF316**
Final Bacteria Filter
3/16in Stepped Barb
and 1/4in Straight Barb
10/pack
OEM# H621

Healthdyne

Alliance



Item# **OF3002**
Foam Cabinet Filter
8¹/₄" x 3³/₄" x 1¹/₂"
10/pack
OEM# 100-0078-10



Item# **BF2222**
Intake Bacteria Filter
22mm fittings
5/pack
OEM# 100-6631-00



Item# **OF3003**
Foam Intake Filter
for use with BF2222
2⁵/₈" OD x 3³/₄" ID x 2³/₈"
10/pack
OEM# H618



Item# **BF316**
Final Bacteria Filter
3/16in Stepped Barb
and 1/4in Straight Barb
10/pack
OEM# H621

React Health/3B

Stratus 5



Item# **OF6030**
Foam Cabinet Filter
10/pack
OEM# STR1030



Item# **BF6020**
HEPA Filter
Each
OEM# STR1020



Item# **BF516**
Final Bacteria Filter,
1/4in Barbed
10/pack
OEM# P0750X030

Nidek

Nuvo



Item# **OF11003**
Foam Cabinet Filter
9³/₈" x 3¹⁵/₁₆" x 5¹/₁₆"
10/pack
OEM# 9250-1025



Item# **BF1438**
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# 7631-1053

Item#
OF11003SP
Nuvo Kit
2 - OF11003
1 - BF1438
1 - BF5000
Each



Item# **BF5000**
Extended Life
Bacteria Filter
Each
OEM# 9250-1028

Mark IV, V, V+



Item# **BF1438**
Final Bacteria Filter
1/4in to 3/8in
Barbed
10/pack
OEM# 7631-1053

Repair Parts NEW



Invacare Platinum 5 / Platinum 10

Item# RP001
Pilot Valve and
1 O-Ring
Each
OEM# 2003052



Item# RP002
O Ring for Pilot Valve
25/pack
OEM# 1085112



Item# RP003
PE Valve
Each
OEM# 1143578



Item# RP004
On/Off Switch
Each
OEM# 1139950



Item# RP005
Hour Meter
Each
OEM# 2002939



Item# RP006
O2 Conc Hose Barb
Fitting 1/8" Conc Port
Each
OEM# 1104151



Item# RP009
Power Cord
Each
OEM# 1108764



Item# RP010
4-Way Valve
Rebuild Gasket Kit
Each



Item# RP014
Concentrator Capacitor
15MFD
Each
OEM# 2003019



Item# RP016
Concentrator Capacitor
30MFD *For Platinum 10
Each
OEM# 1128740



Item# RP019
Pressure Regulator
Each



Item# RP031
Cabinet Clips
2/pack



Item# RP044
Filter Access Door
Each
OEM# 110766



Item# RP045
Sieve Bed Filter Grid
Each
OEM# 1109520



Item# RP081
Piston Cup Seals for
Thomas Compressor
2650/2660
50/pack



Item# RP083
Brass O2 Outlet
10/pack



Invacare HomeFill

Item# RP004
Power Switch
Each
OEM# 2001134



Item# RP007
O2 Conc Fitting Receptacle Port
Each
OEM# 1104152



Item# RP008
2000psi Pressure Switch
Each
OEM# 1089467



Item# RP013
Concentrator Capacitor
10MFD
Each



Item# RP026
Cylinder Dust Cap
10/pack



Item# RP027
Compressor Dust Cap
10/pack
OEM# 1156644



Item# RP028
O-Ring for
Homefill Tank
10/pack



Item# RP031
Cabinet Clips
2/pack



Item# RP042
Male, QDV,
High Pressure, INTL
with Dust Cover
Each



Repair Parts NEW



Invacare Perfecto2 / Perfecto2V

Item# RP001
Pilot Valve and
1 O-Ring
Each
OEM# 2003052



Item# RP002
O Ring for Pilot Valve
25/pack
OEM# 1085112



Item# RP003
PE Valve
Each
OEM# 1143578



Item# RP004
On/Off Switch
Each
OEM# 1139950



Item# RP005
Hour Meter
Each
OEM# 2002939




Item# RP006
O2 Conc Hose Barb
Fitting 1/8"
Conc Port
Each
OEM# 1104151



Item# RP009
Power Cord
Each
OEM# 1108764



Item# RP010
4-Way Valve Rebuild
Gasket Kit
Each



Item# RP013
Concentrator Capacitor
10MFD
**for Perfecto2*
Each



Item# RP015
Concentrator Capacitor
22MFD
**for Perfecto2V*
Each



Item# RP018
Compressor
Rebuild Kit
**for Perfecto2V*
Each
OEM# 1191916



Item# RP019
Pressure Regulator
Each



Item# RP031
Cabinet Clips
2/pack



Item# RP036
Complete
Cabinet Set
3/case
OEM# 1153341



Item# RP043
Cabinet Access Door
**for Perfecto2*
Each
OEM# 1161831



Item# RP045
Sieve Bed Filter Grid
**for Perfecto2*
Each
OEM# 1109520



Item# RP046
PVC F-Tube
Each
OEM# 1143446



Item# RP047
Control Panel
Each
OEM# 1143469



Item# RP081
Piston Cup Seals for
Thomas Compressor
2650/2660
50/pack



Item# RP083
Brass O2 Outlet
10/pack



Item# RP101
Sticker Labels
for Perfecto2 Shroud
10/pack



Item# RP102
Sticker Labels
for Perfecto2V Shroud
10/pack



Repair Parts

NEW


Drive/DeVilbiss 525 / 1025

Item# **RP004**
On/Off Switch
Each

OEM# 1139950



Item# **RP018**
Compressor
Rebuild Kit
**Not compatible
with 1025*
Each

OEM# 1191916



Item# **RP020**
Flow Meter
Each

OEM# 505DZ-607



Item# **RP021**
Cooling Fan
Each

OEM# 525DD-634



Item# **RP022**
Exhaust Muffler
Each

OEM# 515A-705



Item# **RP023**
Circuit Breaker
Each

OEM# 525D-613



Item# **RP024**
Oxygen Outlet Port
5/case

OEM# 525DD-606



Item# **RP032**
Pressure Regulator
8.5 PSI
Each



Item# **RP048**
Power Cord
Each



Item# **RP049**
DS SSP Capacitor
**for 1025DS*
Each



Item# **RP093**
Rebuilt Drive
Rotary Valve
12/case



Item# **RP094**
Rebuild Kit for
Rotary Valve
Each



Item# **RP095**
Spacers for Rebuilt
Rotary Valves
10/pack



Respironics EverFlo (Post-2013)

Item# **RP006**
O2 Conc Hose
Barb Fitting 1/8"
Conc Port
Each

OEM# 1104151



Item# **RP012**
Flow Meter
Each

OEM# 0209184



Item# **RP035**
Pressure Regulator
Adapter Kit
Each

OEM# 1038824



Item# **RP082**
SMC Pilot and
Body Valve
Each



Respironics Trilogy

Item# **RP025**
O2 Quick Connect
L-Shape Adapter
2/pack



Item# **RP025S**
O2 Quick Connect
Straight Adapter
2/pack



RP025S is Also Compatible with:

- Löwenstein **LUISA**
- ResMed **Astral**
- React Health **V*Home**

Repair Parts NEW



Rhythm LM5BA & LM5CA, Lifestyle LM5A, Jumao JMC5A, Direct Supply Attendant 8X899, Medacure OC5L

Item# **RP004**
On/Off Switch
Each

OEM# 1139950



Item# **RP005**
Hour Meter
Each

OEM# 2002939



Item# **RP019**
Pressure Regulator
*For Rhythm LM5BA
Each



Item# **RP023**
Circuit Breaker
*for Rhythm LM5BA
Each



Item# **RP037**
Muffler
Each



Item# **RP038**
PE Valve
Each



Item# **RP039**
Control Valve
Each



Item# **RP040**
Power Cord
Each



Item# **RP041**
Casters
4/pack



Item# **RP075**
Cup Seal for
Compressor
Rebuild Kit
Each



Item# **RP079**
Pressure Regulator,
7 PSI
Each



Item# **RP080**
O-Rings for
Pressure Regulator
50/pack



Item# **RP084**
Power Cord
*for Rhythm LM5BA
Each



Item# **RP085**
Power Switch
*for Rhythm LM5BA
Each



Item# **RP086**
Cooling Fan
*for Rhythm LM5BA
Each



Item# **RP087**
Solenoid Valve
*for Rhythm LM5BA
Each



Item# **RP088**
Oxygen Sensor
*for Rhythm LM5BA
Each



Item# **RP089**
Circuit Board
*for Rhythm LM5BA
Each



Item# **RP090**
Display Screen
*for Rhythm LM5BA
Each



Item# **RP096**
Cabinet with Door for
Direct Supply Attendant
2/box



React Stratus 5

Item# **RP071**
Compressor
Each



Item# **RP073**
Regulator Valve
Each



Item# **RP074**
Cooling Fan
Each



Item# **RP076**
Power Switch
25/pack



Item# **RP072**
Compressor Seal -
Piston Ring
Each



Item# **RP077**
Pneumatic Valve
Each



Liquid Oxygen Repair Parts NEW



Caire Companion C1000 / C1000T

Item# **RP050**
Cover Side B
Each



Item# **RP051**
Side Cover Decal for Companion T
10/box
OEM# B-775499-00



Item# **RP053**
Moisture Control Pad
*for C1000
Each
OEM# B-775092-00



Item# **RP052**
Carrying Strap
Each
OEM# B-775097-00



Item# **RP054**
Top Bezel
*for C1000
Each
OEM# B-775003-00



Item# **RP055**
Cover Side A
*for C1000
Each
OEM# B-775001-00



Item# **RP058**
Tube Fill
Each



Item# **RP060**
Quick Connect Adapter
Each
OEM# B-775263



Item# **RP062**
Manifold O-Ring
Each
OEM# 775037



Item# **RP063**
Vent / Fill Assembly
*for C1000T
Each
OEM# 775037



Item# **RP065**
Lens Cap
Each



Item# **RP067**
Portable Rubber Foot
10 Sheets of 56



Caire HELIOS

Item# **RP056M**
Lip Seal - Metal
Each
OEM# B-775260-00



Item# **RP056R**
Lip Seal - Rubber
Each
OEM# B-775260-00



Puritan Bennett 41L Base Liquid Oxygen Vessel

Item# **RP057**
Rebuilt Contents Gauge
Each
OEM# B-775847-SV



Item# **RP066**
Contents Gauge Tubing
2/pack
OEM# B-775856-SV, PB21, PB41, H46



We're always adding new repair parts!

[Click here to see our full lineup](#)

Suction and Tracheostomy

Suction	Page
Suction Machine	70
Suction Canister and Kits	70-71
Open Suction Catheters	72
Open Suction Catheter Kits	73
Suction Filters	74
Suction Adapters and Connectors	74
Yankauers	74
Tracheostomy Supplies	75
Heat and Moisture Exchangers	76
Closed Suction Catheters	77-78

Suction

Suction Machine

Item#
SU100DC
With Internal
Rechargeable
Battery

SU100AC
Without Battery

SU100DPA
Includes DC Mobile
Power Adapter
and Internal
Rechargeable Battery



- Effective removal of respiratory fluids
- Durable yet lightweight
- Variable vacuum control
- Quiet and easy to use
- Two-year warranty

Comes with
Machine
Carry Bag
800cc Suction Canister
AC Adapter & Power Cord
6ft Suction Tubing
3.5in Suction Tube Connector
4.5in Suction Tube Connector
BF1438S Bacteria Filter

Accessories and replacement parts

SU100DC-PA	DC Mobile Power Adapter (for car power outlet)
SU100DC-Adapt	Replacement Power Adapter (without cord)
SU100DC-Cord	Replacement Power Cord (without adapter)
SU100DC-BAG	Carrying/Storage Bag
SU100DC-BATTERY	Internal rechargeable battery replacement

Suction Canister

Item#
RES023AW
800cc • Single Use
48/case
HCPCS# A7000
• Float valve shutoff



Covers for Suction Machines



20×24"
2mil
500 bags/roll

Item#	Color
EC2024C	Clear
EC2024T	Tan Tint

Suction

Suction Kit for SU100DC

Item# RES026S-SS

RES023AW Suction Canister
 RES024 Suction Tubing Connector 3.5in
 RES024M Suction Tubing Connector 4.5in
 BF1438S Bacteria Filter
 RES025 Suction Tubing 1/4in x 6ft

HCPCS#

CANISTER - A7000

TUBING - A7002



Suction Filter Component Kit for SU100DC

Item# RES024-MKIT

RES024 Suction Tubing Connector 3.5in
 RES024M Suction Tubing Connector 4.5in
 BF1438S Bacteria Filter



Suction Kit

Item# RES026S-MULTI

RES023AW Suction Canister 800cc
 RES024XL Suction Tube Connector 16in
 RES025 Suction Tubing 6ft x 1/4in
 BF516TS Bacteria Filter
 RES024S Suction Tube Connector 1.5in



Suction Component Kit

Item# RES024-MULTI

RES024S Connector Tube 1.5in
 RES024XL Connector Tube 16in
 BF516TS Bacteria Filter



RES026S-MULTI and RES024-MULTI are compatible with:

Precision Medical

Easy Vac	PM60
Power Vac	PM61
Easy Go	PM66 series, PM65 PM65E, PM65HG

Compass/Roscoe

Roscoe Portable Suction	50004 (AC model) 50006 (DC model)
Cliq Single Cylinder Aspirator	ASP-ROS

Medline

Heavy Duty	HCS7000
------------	---------

Invacare

Aspirator	IRC1135
Durable Suction Unit	IRC1136

Drive/DeVilbiss

Heavy-Duty Suction Machine	18600
VacuMax™ Go Suction Machine	18615

Allied Healthcare

Schuco Portable Aspirator with Sturdy Legs	S430A
Schuco Portable Suction Unit with Metal Base	S130A
Schuco Portable Aspirator with Molded Plastic Base	S330A
Gomco OptiVac® AC/DC Portable Suction Unit & 305 Aspirator	G180 G180CE G180CE-JAR G180-JAR 01-22-3005 01-22-0305
Gomco® Tabletop Aspirator Model 300	01-22-0300
Gomco® General Use Aspirators Models 3020 & 3040	01-62-3020 01-62-3040
Gomco® 405 Aspirator	01-12-0405 01-12-4005 01-23-0405
Gomco® Stand-Mounted Models 4032 & 4042	01-52-4032 01-52-4042
Gomco Table Top Gastric Aspirator 270	01-22-0270

RES026S-SS and RES024-MKIT are compatible with:

Drive/DeVilbiss

Vacu-Aide® QSU Quiet Suction Unit	7314P-D 7314P-D-EXF 7314D-D
Vacu-Aide® Suction Unit	7305D-D 7305P-D-EXF 7305D-D-EXF 7305P-D

Sunset

Suction Machine (Portable Unit)	SU100DC
Suction Machine (Stationary Unit)	SU100AC

Open Suction Catheters

50/case

Size	Length	Valve	Tip	Intended Use	Packaged	Item#
6fr	44 cm	Chimney	DeLee	Pediatric	Coiled	TRASC6-CC
6fr	44 cm	Chimney	DeLee	Pediatric	Straight	TRASC6-CS
8fr	44 cm	Chimney	DeLee	Pediatric	Coiled	TRASC8-CC
8fr	44 cm	Chimney	DeLee	Pediatric	Straight	TRASC8-CS
8fr	44 cm	Directional	DeLee	Pediatric	Straight	TRASC8-DS
10fr	58 cm	Chimney	Whistle	Adult	Coiled	TRASC10-CC
10fr	58 cm	Chimney	Whistle	Adult	Straight	TRASC10-CS
10fr	58 cm	Directional	Whistle	Adult	Straight	TRASC10-DS
12fr	58 cm	Chimney	Whistle	Adult	Coiled	TRASC12-CC
12fr	58 cm	Chimney	Whistle	Adult	Straight	TRASC12-CS
12fr	58 cm	Directional	Whistle	Adult	Straight	TRASC12-DS
14fr	58 cm	Chimney	Whistle	Adult	Coiled	TRASC14-CC
14fr	58 cm	Chimney	Whistle	Adult	Straight	TRASC14-CS
16fr	58 cm	Chimney	Whistle	Adult	Straight	TRASC16-CS
18fr	58 cm	Chimney	Whistle	Adult	Straight	TRASC18-CS

- Sterile
- Single patient use
- Latex-free

NEW




Coiled Packed

Straight Packed




Chimney
Valve



Directional
Valve

Whistle
Tip



DeLee Tip

Open Suction Catheter Kits with Cup

50/case

Size	Length	Valve Type	Tip	Intended Use	Item#
5fr	44 cm	Chimney	DeLee	Pediatric	TRAKIT5-C
6fr	44 cm	Chimney	DeLee	Pediatric	TRAKIT6-C
8fr	44 cm	Directional	DeLee	Pediatric	TRAKIT8-D
8fr	44 cm	Chimney	DeLee	Pediatric	TRAKIT8-C
10fr	58 cm	Directional	Whistle	Adult	TRAKIT10-D
10fr	58 cm	Chimney	Whistle	Adult	TRAKIT10-C
12fr	58 cm	Directional	Whistle	Adult	TRAKIT12-D
12fr	58 cm	Chimney	Whistle	Adult	TRAKIT12-C
14fr	58 cm	Directional	Whistle	Adult	TRAKIT14-D
14fr	58 cm	Chimney	Whistle	Adult	TRAKIT14-C
16fr	58 cm	Chimney	Whistle	Adult	TRAKIT16-C

- Sterile
- Includes Catheter and Gloves



Sleeved Coiled Open Suction Catheter Kits with Cup

50/case

Size	Length	Valve Type	Tip	Intended Use	Item#
8fr	44 cm	Directional	DeLee	Pediatric	TRAKIT8-SL-NG
10fr	58 cm	Directional	Whistle	Adult	TRAKIT10-SL-NG
12fr	58 cm	Directional	Whistle	Adult	TRAKIT12-SL-NG
14fr	58 cm	Directional	Whistle	Adult	TRAKIT14-SL-NG




Suction

Suction Machine Inline Filters

Item# **BF1412**

Final Bacteria Filter
1/4in to 1/2in
Stepped Barb on
Both Sides
10/pack



Item# **BF1438S**

1/4in to 3/8in
10/pack



Item# **BF3812**

3/8in to 1/2in
10/pack



Item# **BF1312**

1/3in to 1/2in
1/pack



Item# **BF516TS**

5/16in Barbed with
Threaded NPT
10/pack



Suction Tubing

1/4in

Length	Item#	Qty
18in	RES025S	50/case
6ft	RES025	10/case



3/16in

Length	Item#	Qty
18in	RES025-31618	100/case
20in	RES025-31620	100/case
6ft	RES025-31672	50/case

Suction Tube Connectors

10/pack

Size	Item#
1.5in	RES024S
3.5in	RES024
4.5in	RES024M
6in	RES024L
10in	RES024L10
16in	RES024XL



Rubber Suction Elbow Connector Drive DeVilbiss Style

Item#

RES028
6/pack



Suction Elbow Connector Precision Medical Style

Item#

RES033
10/pack



Bulb Tip Yankauers

Item# **TRAA001S**

Vented • **50/case**

Item# **TRAA002S**

Non-Vented • **50/case**

HCPCS# A4628



- Sterile
- Single patient use
- Not made with natural latex rubber

Tracheostomy

Tracheostomy Mask

Item#

RES2130

Adult Size

50/case

HCPCS# A7525

**RES2230**

Pediatric Size

50/case

Flexible Connector

50/case

Item#

TRAA004 Adult

22/15mm x 15mm

Length 4.9-7.5cm

OEM# 3222

**TRAA004P** Pediatric


15mm F x 15mm M

Length 4.9-7.5cm

OEM# 3215



Adapter for T-HME

Item# **RES027A****Each**

Flexible Trach Adapter

Item#

TRAA003

22mm Cuff

Each

OEM# 1073902, 1147898

TRAA003P Pediatric

15mm Cuff

OEM# 1073863



Tracheostomy Ties with Fastener Tabs

20/case

One-piece • Easy on


Item#

RES029P-11

1 in x 11.5 in

RES029-17

1 in x 17 in



Softer, more
breathable
100% cotton
patient-facing
material in 2026



Two-piece • Easy off

Item#

RES030P-10

0.78 in x 10 in

RES030P-14

0.78 in x 14 in

RES030P-18

0.78 in x 18 in

RES030-20

1 in x 20 in



Sunset's trach ties deliver value and comfort, with soft straps and light Velcro closures for customization. Secure, but flexible for patient movement, the ties help prevent sliding, skin irritation and accidental extubation. 60% Polyester, 20% Cotton, 20% Nylon.

Heat and Moisture Exchangers

Tracheostomy Heat Moisture Exchangers

30/case*

*RES027UOP 50/case



	T-HME RES027	T-HME with O2 Port RES027OP	HME with O2 Port RES027UOP
TV Range (mL)	150-1250	150-1250	50-1000
Moisture Output (mg H₂O/L)	36.4 at VT = 500 ml	32 mg/L	13.2 mg/L
Weight (g)	6.6	7.4	5
Resistance (cm H₂O)	0.2 @ 30 LPM	0.2 @ 30 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM
Dead Space (mL)	15	7	17
Connection (mm)	15F	15F 5.5mm O2 port	15F

Inline Heat Moisture Exchangers

30/case

*RES027FLEX 40/case




NEW



	Flex HME — Foam RES027FLEX	Inline Small Volume HME with Gas Sampling Port RES027IL	Inline Large Volume HME with Gas Sampling Port RES027ILLV	Inline HME — Neonatal RES027NE
TV Range (mL)	250-1500	70-600	150-1500	5-50
Moisture Output (mg H₂O/L)	30 mg/L	34.5 at VT = 250 ml 32.5 at VT = 500 ml	37.7 at VT = 250 ml 35.5 at VT = 500 ml 33.5 at VT = 1000 ml	32.8 at VT = 250 ml
Weight (g)	15.5	10	15	1.9
Resistance (cm H₂O)	0.5 @ 30 LPM 1.4 @ 60 LPM	0.8 @ 30 LPM 2.1 @ 60 LPM	1.0 @ 30 LPM 2.6 @ 60 LPM	0.8 @ 15 LPM 1.7 @ 30 LPM
Dead Space (mL)	90	21	34	2.2
Connection (mm)	22F, 15F	22M/15F-15M	22M/15F-15M	15M, 15F

Standard Closed Ventilation Suction System



Standard Closed Ventilation Suction System

24 Hour Adult

Endotracheal 540mm

Item#	Description
VEN46001000A	10Fr, Flex Tube
VEN46001000AM	10Fr, MDI Port
VEN46001200	12Fr, Flex Tube
VEN46001200M	12Fr, MDI Port
VEN46001400	14Fr, Flex Tube
VEN46001400M	14Fr, MDI Port
VEN46001600	16Fr
VEN46001600M	16Fr, MDI Port

Tracheostomy 305mm

Item#	Description
VEN46011000A	10Fr
VEN46011000AM	10Fr, MDI Port
VEN46011200	12Fr, Flex Tube
VEN46011200M	12Fr, MDI Port
VEN46011400	14Fr, Flex Tube
VEN46011400M	14Fr, MDI Port, Flex Tube
VEN46011600	16Fr
VEN46011600M	16Fr, MDI Port

72 Hour Adult

Endotracheal 540mm

Item#	Description
VEN56051000A	10Fr
VEN56051000AM	10Fr, MDI Port
VEN56051200	12Fr
VEN56051200M	12Fr, MDI Port
VEN56051400	14Fr
VEN56051400M	14Fr, MDI Port
VEN56051600	16Fr
VEN56051600M	16Fr, MDI Port

Tracheostomy 305mm

Item#	Description
VEN56061000A	10Fr, Flex Tube
VEN56061000AM	10Fr, MDI Port
VEN56061200	12Fr, Flex Tube
VEN56061200M	12Fr, MDI Port
VEN56061400	14Fr, Flex Tube
VEN56061400M	14Fr, MDI Port
VEN56061600	16Fr
VEN56061600M	16Fr, MDI Port



Neonatal/Pediatric

With DSE Connection

Item#	Description
VENCSS6-305	6Fr, 305mm
VENCSS8-203	8Fr, 203mm
VENCSS8-355	8Fr, 355mm
VENCSS10-203	10Fr, 203mm
VENCSS10-405	10Fr, 405mm



Ventilator

Ventilator

Page

LUISA Ventilator	79
Ventilator Filters	80-81
Adapters and Connectors	81
Ventilator Mask	81
Breathing Circuits	82-85
Heated Single-Limb	82
Non-Heated Single-Limb	83-84
Dual-Limb	85

Ventilator

Movair / Löwenstein LUISA®

Item# **VNT31390**
Each

NEW

LUISA with High Flow Therapy provides long term care with easy serviceability. Routine patient maintenance consisting of only filter changes and a device wipe down provides uncomplicated upkeep.

With simple, routine maintenance and no adapters required between circuits, the LUISA provides clinicians and patients' peace of mind with an all-encompassing system in a compact, home ventilator with the added benefit of High Flow Therapy.

Includes:

- Exhalation Module
- Leakage Circuit
- Power Supply/Power Cord
- O2 Nozzles
- Fine and Course Filters
- Protection Bag
- Accessories Bag
- USB Drive



- High Flow Therapy: all circuit systems
- Connectivity: customized patient download reports
- Uncomplicated: no adapters, intuitive operation
- Battery run time: up to 18 hours
- Monitoring: FiO2 and SpO2 measurement
- Safety: alarm displayed in two languages simultaneously
- Display: 10" angled touchscreen
- Functionality: 100 - 3,000 mL Vt without additional flow sensor
- Responsive: adjustable target volume speed
- Lightweight: ~ 8 lbs.
- Positioning: vertical or horizontal; Circuit orientation options (left or right side)

Accessories and Replacement Parts

VNT31404	Exhalation Module - Disposable
VNT31569	Power Supply
VNT27027	Power Cord
VNT30669	Oxygen Connection Nozzle
VF29652	Set of 12 Pollen/Fine Filters - 12/PK
VF29928	Set of 2 Air/Coarse Dust Filters - 2/PK
VNT31417	Protective Bag
VNT31440	Accessories Bag
VNT31414	USB Flash Drive

ResMed

Filter for Astral Ventilator

Item# **VF003-4**
4/pack

OEM# 27939



NEW

4-Filter Kit with
Filter Housing

Item# **VF003-KIT**



Item# **VF003-H**
Filter Housing
Each



02 Quick Connect Straight Adapter 1/4" Barb

Item# **RP025S**

1/4" Barb, for use with Astral
2/pack



NEW

Respironics

Air Inlet Vent Filter for V60

Item# **VF002**

For use with Philips Respironics V60
5/pack



Filters for Trilogy/Trilogy EV300 Ventilators

Item# **VF005-10**
Particulate Filter for
Trilogy Ventilator
10/pack



NEW

Item# **VF007-1**
Reusable Filter for Trilogy EV300
1/pack

VF007-10 10/pack
OEM# 1134591



Item# **CF1009F**
Reusable Filter for Trilogy
2/pack

OEM# 1029330



02 Quick Connect Straight Adapter 1/4" Barb

Item# **RP025S**

1/4" Barb, for use with Trilogy
2/pack



NEW

Breas

Vivo 45/45LS Style Air Inlet Filters

Item# **VF010-5**

Disposable filter for Vivo 45LS
5/pack

OEM# 007103



NEW

Item# **VF010F-5**
Reusable Filter for Vivo 45LS
5/pack

OEM# 007104



NEW

Item# **VF011F-5**
Reusable filter for Vivo 45
5/pack

OEM# 007105



NEW

Movair / Löwenstein

Filters for LUISA

Item# **VF008-12**
Disposable filter
12/pack

OEM# WM29652



3-Filter Kit

1 VF008 Disposable Filter
1 VF008F Reusable Filter
1 VF014F Gross Particle Filter
Item# **VF008-KIT**

Item# **VF008F-2**
Reusable Filter
2/pack

OEM# WM29928



Item# **VF014F**
Gross Particle Filter
Each



NEW

02 Quick Connect Straight Adapter 1/4" Barb

Item# **RP025S**

1/4" Barb, for use with LUISA
2/pack



NEW

Covidien/Nellcor

Ventilator Felt Filter

Item# **VF006**

For use with NPB LP6 & LP10
10/pack

OEM# L-002917-000



Universal filter

Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item#
BF2200
22mm
5/pack



Attach between the machine and tube for additional protection against bacteria and allergens

HME Filter

Bacterial/Viral HME Filter

Item#
BF104
Each



• 99.99+% BFE/VFE

Ventec VOSCN

NEW

Filters for Ventec VOSCN

Item# **VF015-10**
Internal Bacteria Filter
10/pack
OEM# PRT-00976-000



Item# **VF016-10**
External Bacteria Filter
10/pack
OEM# PRT-00850-000



Item# **VF012F-25**
Fan Filter
25/pack
OEM# PRT-00978-000



Item# **VF013F-25**
Patient Air Intake Filter
25/pack
OEM# PRT-00977-000



Ventilator Adapters, Connectors, and Mask

Valved T-Adapter — Adult

Item# **VEN001**
Inline adapter for nebulizer treatment
22F-18F-22M/15F
Each



Circuit Connectors

Item#
RES013
15mm ID x
22mm OD
10/pack



Item#
RES014
22mm ID x
22mm OD
10/pack



Exhalation Ports For Ventilator Circuits & CPAP

Item#
VEN002
Swivel
Each



Item#
VEN003
Fixed
22 mm OD • 15 mm ID
10/pack
OEM# 312149



Straight Connectors

Item#
VEN005
22mm ID,
15mm OD
Each



NEW

Item#
VEN006
22mm ID,
22mm OD
Each



NEW

Active Exhalation Valve Adapter

Item# **VEN004**
Each
OEM# 1065659

VEN004-R
Compatible with ResMed



NEW

NIV Mask with Headgear — Non-vented

Each

Item#	Size
CM038M	Medium
CM038L	Large



Heated Single Limb Ventilator Circuits without Humidifier Chamber

Active Heated Breathing Circuits

30/case

For use with MR850 Humidifier

Item#

VC1002 Active Adult

With Exhalation Valve Adapter

NEW

Includes:

- 22mm circuit
- Active exhalation valve adapter
- Optional connector

OEM# 1073222, 1073231



VC1002-R

Compatible with ResMed

Item#

VC1002P Active Pediatric

With Exhalation Valve Adapter



NEW

Includes:

- 15mm circuit
- Active exhalation valve adapter
- 15mm/22mm connector

OEM# 1073221, 1073230



VC1002P-R

Compatible with ResMed

Passive Heated Breathing Circuits

30/case

For use with MR850 Humidifier

Item#

VC1003 Adult

With Swivel Exhalation Port

Includes:

- 22mm circuit
- Swivel exhalation port
- Optional connector

OEM# 1076046, 1073227, 1073235



Item#

VC1003P Pediatric

With Swivel Exhalation Port




NEW

Includes:

- 15mm circuit
- Swivel exhalation port
- 15mm/22mm connector
- 22mm/22mm connector

OEM# 1073225, 1073234



Adult Single Limb Heated Circuit

10/case

For use with MR850 Humidifier

Item#

VC1114 Adult

OEM# RT114



Passive Ventilator Circuits available from Sunset

Item#	Heating	Exhalation Port	Filter	Water Trap	Size	Tube Diameter	OEM#
VC1000	Non-Heated	Swivel	None	No	Adult	22mm	1073228 & 1073236
VC1003	Heated	Swivel	None	No	Adult	22mm	1076046, 1073227, 1073235
VC1003P	Heated	Swivel	None	No	Pediatric	15mm	1073225, 1073234
VC1004	Non-Heated	Fixed	None	Yes	Adult	22mm	1069443 & 1069444
VC1004P	Non-Heated	Swivel	None	Yes	Pediatric	15mm	1073224 & 1073233
VC1005	Non-Heated	Swivel	Bacterial/Viral	No	Adult	19mm	1147378, 1132340, 7615
VC1005P	Non-Heated	Swivel	Bacterial/Viral	No	Pediatric	15mm	1132339
VC1006	Non-Heated	None	Bacterial/Viral	No	Adult	22mm	
VC1007	Non-Heated	Fixed	Bacterial/Viral	No	Adult	22mm	
VC1114	Heated	None	None	No	Adult	22mm	RT114
VC1100	Non-Heated	Fixed	None	No	Adult	22mm	1069467 & 1069468
VC1160	Non-Heated	Fixed	None	No	Adult	22mm	1069210

Non-Heated Single Limb Ventilator Circuits without Humidifier Chamber

Passive Non-Heated Breathing Circuit

Item#

VC1160

**Compatible with V60,
No Filters**

OEM# 1069210



Active Non-Heated Breathing Circuits

30/case

22mm cuffs

Item#

VC1001 Adult


OEM# 1073223, 1073232

VC1001-R**Compatible with ResMed**

OEM# 21321

VC1001P-R Pediatric**Compatible with ResMed**


OEM# 21312




Item#

VC1001P Pediatric

OEM# 1073220, 1073229



To use active circuits with Respiration Trilogy, simply cut ends off pressure and exhalation lines.




Non-Heated Single Limb Ventilator Circuits without Humidifier Chamber

Passive Non-Heated Breathing Circuits with Water Trap


Item#
VC1004 Adult

OEM# 1069443, 1069444



Item#
VC1004P Pediatric

OEM# 1073224, 1073233



Passive Non-Heated Breathing Circuit

For use with Respiration Trilogics
Each

Item#
VC1000
with Swivel Exhalation Port
OEM# 1073228 & 1073236

VC1000P
Pediatric • 15mm Tube Diameter
with Swivel Exhalation Port
OEM# 1052127

VC1100
with Fixed Exhalation Port
OEM# 1069467 & 1069468



Passive Circuit with Exhalation Port and Bacterial/Viral Filter

With 22mm cuffs
Each

Item#
VC1005
Adult • 19mm Tube Diameter
OEM# 1147378, 1132340, 7615

VC1005P
Pediatric • 15mm Tube Diameter
OEM# 1132339



Adult Passive Circuit with Bacteria Filter

Item# **VC1006**
Comparable to
Breas 007936
Each



Adult Passive Circuit with Bacteria Filter and Fixed Exhalation Port

Item# **VC1007**
Comparable to
Breas 007615
Each



Dual Limb Ventilator Circuits

NEW


Adult Dual Limb Active Circuit 6ft

Item#

VC1009

For use with Vyair LTV

OEM# 29657-001



Pediatric Dual Limb Active Circuit 6ft

Item#

VC1010

For use with Vyair LTV

OEM# 29658-001




Pediatric Dual Limb Active Circuit Water Trap

Item#

VC1011

For use with Vyair LTV

OEM# 29662-001




Adult Dual Limb Active Circuit with Water Trap

Item#

VC1012

For use with Vyair LTV

OEM# 29663-001



Other

Other

Page

Disinfectants	86
Cough Assist Kits	87
Breathing Masks	87
Aerosol	88
Tabletop Nebulizers	89
Handheld Nebulizers	90
Nebulizer Kits	91
Nebulizer and Anesthesia Filters	92
StrongArm Mobility and Comfort	93-94
Bath Safety	95
Maternity	96-97
Equipment Covers	98-99
Blood Testing	99

Other HCPCS

A6402	Gauze - sterile, pad size 16 sq. in. or less, without adhesive border, each dressing
A4483	Heat and Moisture Exchanger — Disposable
A4626	Tracheostomy Cleaning Brush
A4628	Oropharyngeal Suction Catheter (Yankauer)
A4629	Tracheostomy Care Kit - for established tracheostomy
A7000	Suction Canister — Disposable
A7002	Tubing for Suction Pump
A7003	Nebulizer Set - Small Volume—Disposable
A7005	Nebulizer Set - Small Volume—Nondisposable

A7013	Nebulizer Compressor Filter— Disposable
A7015	Aerosol Mask used with Nebulizer
A7525	Tracheostomy Mask
E0570	Nebulizer with Compressor
E0600	Respiratory Suction Pump - home model, portable or stationary, electronic
E2601	General Use Wheelchair Seat Cushion - width less than 22 inches, any depth
K1005	Disposable collection and storage bag for breast milk, any size, any type, each

Disinfectants

Simpli Germicidal Disposable Wipes

Item#
RES6108W
Canister of 160 wipes

- Bactericidal, Tuberculocidal, Fungicidal, Virucidal*
- Premoistened wipes containing a quaternary/alcohol-based solution.
- Made in the USA
- Kills *SARS-CoV-2, Covid 19, Mycobacterium bovis – BCG (Tuberculosis—TB) in 1 Minute
- Kills Candida Auris in 2 Minutes



CleanSmart Disinfectant

Item#
CAP1019-C
23oz Spray Bottle
6/case

- CleanSmart is EPA listed as a broad-spectrum and hospital disinfectant
- EPA Registration Number 89896-2
- CleanSmart is formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment
- CleanSmart is included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2, the virus that causes COVID-19
- Active ingredient Hypochlorous Acid does not leave a chemical residue



Cough Assist Kits

Cough Assist Kits Compatible with ABM Respiratory Care BiWaze Cough System and Philips Respironics Cough Assist T70

All kits come with 1 BF2200 Bacterial/Viral Filter and 1 RES013 Circuit Connector

Item#

CAKIT001

With 1 TUB06 6ft Tube
Each

OEM# 1090835



Item#

CAKIT002

With 1 TUB06 6ft Tube and
1 RES019 Oxygen Enrichment Adapter
Each



CAKIT003

With 1 TUB06 6ft Tube
and 1 TRAA003 Flexible
Trach Adapter
Each



ABM Respiratory Care BiWaze Cough System



Item# **OF8003-3**

Air Inlet Foam Filter
3/pack

Kits With Tube and SBC01 Series Breathing Mask

With 1 TUB06 6ft Tube

Mask Size	Item#	Comparable to OEM#
Infant	CAKIT004I	1090830
Toddler	CAKIT004T	1090831
Adult Small	CAKIT004S	1090832
Adult Medium	CAKIT004M	1090833
Adult Large	CAKIT004L	1090834



With 1 TUB04 4ft Tube

Mask Size	Item#	Comparable to OEM#
Infant	CAKIT005I	325-9237
Toddler	CAKIT005T	325-9236
Adult Small	CAKIT005S	325-9235
Adult Medium	CAKIT005M	325-9234
Adult Large	CAKIT005L	325-9217



CAKIT006

with 1 CAMP001 Mouthpiece
and 1 TUB06 6ft Tube
Each

OEM# 1090836



Breathing Mask

Air Cushion Breathing Mask with Valve – No Hooks

20/case

Mask Size	Item#
Infant	SBC01I
Toddler	SBC01T
Adult Small	SBC01S
Adult Medium	SBC01M
Adult Large	SBC01L



- Ergonomic cushion design provides an optimal seal with minimal pressure
- Clear mask permits easy patient observation
- Available in a variety of sizes for the best possible fit
- Inflation port for optimized seal

Aerosol

Aerosol Mask

Item# **RES2110**

Adult Size

50/case

Item# **RES2210**

Pediatric Size

50/case

HCPCS# A7015

- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



Peak Flow Meter

Item# **RES052**

10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**

Pediatric Mouthpiece
for Peak Flow Meter

10/case



Drainage Bag with T-Adapter

Item# **RES070**

750cc • 50/case



Drainage Bag with Y-Adapter

Item# **RES071**

750cc • 50/case



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

SperaChamber

Anti-Static Valved Inhaler Chamber
for Metered Dose Inhaler

Chamber Only

Item#

RES195-10 10/pack

RES195-50 50/pack

Mask Only

Item#

RES195MASK-S Small

RES195MASK-M Medium

RES195MASK-L Large

Chamber with Mask

10/pack

Item#

RES195MS Small (Infant)

RES195MM Medium (Pediatric)

RES195ML Large (Adult)

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free
- Dog mask appeals to kids



Tabletop Compressor Nebulizers

Compressor Nebulizers

Item#

NEB100

- Disposable Nebulizer Kit
- Foam Filters

NEB100-R

- Disposable Nebulizer Kit
- Reusable Nebulizer Kit
- Foam Filters x 5



NEB100-BAG
Carry Bag

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

Sunny the Bear

NEB300BEAR

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Foam Filters x 5

NEB300BEAR-R

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Reusable T-Piece Kit
- Foam Filters x 5



NEB300-BAG
Bag for Nebulizer Compressor
To include bag in configuration,
just add "B" to Item Number:

NEB300FISH-B
NEB300FISH-RB
NEB300BEAR-B
NEB300BEAR-RB



Gordon the Fish

NEB300FISH

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Foam Filters x 5

NEB300FISH-R

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Reusable T-Piece Kit
- Foam Filters x 5



Covers for Nebulizers



20×24"
2mil
500 bags/roll

Item#	Color
EC2024C	Clear
EC2024T	Tan Tint

Handheld Nebulizers

Handheld Mesh Nebulizer

Item#
NEB402

With internal rechargeable battery

Included

- Medication cup
- USB cord
- Mask
- Mouthpiece
- Portable bag
- Manual

- Includes internal rechargeable battery that provides up to 1 hour of continuous use after full charge
- Great value pricing and quality
- Easy to operate and transport
- FAA compliant

Particle Size MMAD < 5µm

Nebulization Rate 0.25ml/min-0.9ml/min

Medication Cup Capacity 10ml (Max)

Device Weight Approx. 1/4 lb

Device Size Approx. 1 1/2 x 2 x 3 7/8 inches



Handheld Compressor Nebulizer

Item#
NEB200

Included

- Compressor
- Medication cup
- AC power adapter
- Air hose
- Elbow
- Mouthpiece
- Mask
- 4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (Approximately 3 microns)
- No replacement filter required



Nebulizer Kits

Disposable Nebulizer Kit with T-Piece

Item# **RES091**
Disposable
Single Patient
50/case
HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
 - T-piece with mouthpiece and corrugated reservoir tube
 - 7ft supply tube



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**
Disposable
Single Patient
50/case
HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
 - Angled mouthpiece
 - 7ft supply tube



Reusable Nebulizer Kit with T-Piece

Item# **RES093**
Reusable
Single Patient
10/case
HCPCS# A7005

- Includes:**
- 10cc jet nebulizer bottle
 - T-piece with mouthpiece and corrugated reservoir tube
 - 7ft supply tube



- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°

Disposable Nebulizer Kit with Mask

Item# **RES090**
Adult Size
50/case

Item# **RES092**
Pediatric Size
50/case



HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
 - Mask
 - 7ft supply tube

Disposable single patient use for long term delivery procedures

6cc bottle is easy to assemble and disassemble



Nebulizer Filters

Plastic Intake Filter with Head

Item# **BF016**

10/pack

HCPCS# A7013



Plastic Intake Filter with Head

Item# **BF017**

5/pack

HCPCS# A7013



Foam Filter

Item# **BF018**

10/pack

HCPCS# A7013



Intake Filter for 50 PSI Compressor

Item#

BF020

3/pack

Compatible with CHAD Drive
50 PSI Compressor
OEM# 18450-FILTER



Plastic Intake Filter for Respironics InnoSpire

Item# **BF026S**

10/pack

HCPCS# A7013

- Mini
- MicroElite
- MiniElite



Foam Inlet Filter for Invacare 50 PSI

Item#

BF027-10

10/pack



Nebulizer Filter Compatibility

Caire

Medi-Mist	BF018
-----------	-------

CHAD

Drive 50psi	BF020
-------------	-------

Drive/DeVilbiss

PulmoNeb	BF017
PulmoNeb LT	BF017
Sunrise Compact	BF018
Pulmo-Aide	BF017
Pulmo-Aide Compact	BF017
PulmoMate	BF017
Power Neb Ultra 18080	BF017
CHAD 50psi	BF020

Evo

AeroMax	BF018
SportsNeb	BF018
Super SportsNeb	BF018

Invacare

Passport	BF018
Envoy	BF018
Envoy Jr.	BF018
Select	BF018

John Bunn

Nebulite EV	BF018
LT2	BF018
LX	BF018

Medline

SportsMist	BF018
------------	-------

Respironics

Inspiration 626 & 929	BF016
MisterNeb	BF016
Inspiration Elite	BF016
MiniElite	BF018
InnoSpire	BF026S

Salter

AIRE Plus	BF018
-----------	-------

Schuco

2000	BF018
3000	BF018

Sunset

NEB100 Series	BF018
NEB300 Series	BF018

StrongArm®

StrongArm® Comfort Cane

12/case

Item#	Color
BM001	Black
BM001-BL	Blue
BM001-PR	Purple

**Color pack
(one of each)**

BM001-3MC



- Ergonomic handgrip supports forearm and relieves pressure from wrist
- Offset bend centers grip over tip for maximum stability when cane is at the user's side
- Lightweight and strong
- 12 height adjustments
- Comfort Fit Handle
- Supports up to 500 lbs
- For left and right arm
- Self-standing base: Larger 'footprint' greatly improves traction and overall stability
- Safely supports user walking on grass, gravel, sand, carpet, hard surfaces, and stairs



**Now comes in
retail packaging**

StrongArm® Grabber Reacher Tool

Item#

BM002

12/case

The ultimate companion for effortless grabbing and reaching! Engineered with precision and designed for practicality, this innovative tool is perfect for those hard-to-reach items and provides extra support for heavier items. Grab with Ease and Confidence in Each Squeeze.



- Transfers weight for improved support
- Rotating jaw for vertical or horizontal use
- Integrated hook and magnetic retriever
- Durable, lightweight aluminum construction
- Folds for easy storage and portability
- 32" of added reach

StrongArm® Clear View Universal Walker Tray

12/case

Item#	Color
BM003-BL	Blue
BM003-BK	Grey
BM003-PL	Purple



Unobstructed view for safer mobility

The Clear View Walker Tray is designed to enhance the user's ability to see obstacles ahead and reduce the risk of falls. Its transparency helps maintain the user's line of sight, providing added safety and confidence while moving. Compatible with the leading brand walkers.

Transparent

See where you are walking. Avoid tripping hazards.

Universal Fit

Tray works with the leading brand walkers.

Convenience

Built-in cup holder and storage bag.

StrongArm®

StrongArm® Universal Clip-On Light

Item#
BM004
4/box

Featuring a super-bright LED light, this clip-on accessory significantly enhances visibility, providing added safety during evening strolls or walking throughout the house. Effortlessly clip onto canes, walkers, rollators, crutches, grabber reachers, knee scooters, wheelchairs, and more, ensuring powerful and reliable illumination on any mobility device.



StrongArm® LED Magnifying Glass

Item#
BM005
10/box

- Handheld 4x LED Magnifier
- 6 LED Lights
- Compact, Lightweight, and Portable
- Lens Diameter: 2.5"
- Weighs Only 1.7 oz
- Uses 2 AAA Batteries (Included)



StrongArm® Retail Starter Kit

Item#
BM006
Each

NEW

Save when you order the StrongArm Retail Starter Kit!

Includes:

- 2 **BM001** Black StrongArm Comfort Canes
- 2 **BM001-BL** Blue StrongArm Comfort Canes
- 2 **BM001-PR** Purple StrongArm Comfort Canes
- 2 **BM002** StrongArm Grabber Reacher Tools
- 2 **BM004** StrongArm Universal Clip-On Lights
- 2 **BM005** StrongArm LED Magnifying Glasses



Bath Safety

PreserveTech™ 360° Swivel Bath Chair

Item#

BSE001

Compact and customizable, the PreserveTech™ Swivel Bath Chair provides a safe transfer in and out of the tub for those with limited mobility. The seat rotates a full 360° in either direction and automatically locks at 90° intervals to help prevent skin irritation and friction wounds during transfers. The seat height can be adjusted from 15"-22", and the padded armrests can be raised out of the way.



Splash Defense™ Transfer Bench

Item#

BSE002

Splash Defense™ Transfer Bench is the first of its kind to allow the shower curtain liner to remain inside the bath tub while a user is bathing. This unique product enhancement is essential to help prevent water from spilling out onto the bathroom floor, creating a danger zone for users and caretakers. The Splash Defense™ Transfer Bench provides security for the caretaker, while giving users back their independence and allowing them to bathe safely and privately again.



PreserveTech™ 360° Universal Raised Toilet Seat

Item#

BSE003

Introducing the PreserveTech™ Universal Raised Toilet Seat. With this truly universal raised toilet seat users can set three different seat heights (2, 4 or 6") by maneuvering just two components. This seat can fit on most standard and elongated toilets and features two hygienic cutaways for easy cleaning. The pivoting armrests can be raised or removed, depending on user needs and bathroom space. This item also includes a lid with a recessed grip for effortless opening.



Maternity

Motif Roam Breast Pump with Breast Milk Collection Cups

Item#
MATAAA0026

- **Hospital-Strength Suction:** Up to 280 mmHg suction strength, giving moms the efficiency they need for every pumping session
- **Portable:** Equipped with Breast Milk Collection Cups to give moms the freedom to pump hands-free
- **Discreet:** Collection Cups fit in nursing bra and hold up to 5 oz. of milk
- **More Options:** Resupply Kits available in both flange and collection cup styles
- **Variety of Flange Sizes:** Silicone inserts included in sizes 16 mm, 19 mm, 21 mm and 24mm so mom can find the perfect fit
- **Cluster Feeding and General Massage Mode:** 2 NEW Custom Modes designed by Lactation Professionals
- **Suction levels:** Works to express milk from the breast as efficiently as possible with 9 suction levels
- **Closed System:** Backflow prevention keeps the motor clean, preventing contamination

NEW



What's Included

- | | |
|------------------------------|--|
| 1 - Motif Roam pump motor | 4 - Pairs of flange inserts (16 mm, 19 mm, 21 mm, 24 mm) |
| 2 - Collection cups | 1 - Lanyard |
| 2 - 28 mm breast shield | 1 - USB-C cable |
| 2 - Tubing w/Roam connectors | 1 - Power adapter |
| 2 - Tubing connector | 1 - Nipple Ruler |
| 4 - Duckbill valve | |
| 2 - Diaphragm | |

Motif Milk Storage Bags

Item#
MATAAC000806
100/pack

Other pack counts are available.
Ask your rep for details.

Designed to make breastmilk storage as simple as possible, these bags feature an easy-pour spout for a no-mess experience. Easily transfer milk from the bottle to the bag, seal with the leak-proof double zipper lock, label, and store in your fridge or freezer. No spills. No hassle. Made with your busy day in mind. Comes in a box of 40, 90, 100 or 300 count.

HCPCS: K1005



- Easy-pour spout design
- BPA-free
- Leak-proof with double zipper lock
- Stores up to 8 oz. of milk
- Write-on label
- Freezer safe
- Lays flat for compact freezer storage
- Self stand for refrigerator storage

Maternity

Motif Luna Double Electric Breast Pump

Item#
MATAAA0013

Sophisticated, strong, and soothing, the Motif Luna is Motif's most powerful breast pump. With a battery life of 2+ hours between charges, the Luna was built for modern motherhood. The Luna is available in both battery and non-battery versions. Lightweight and portable, the Luna weighs less than 2 pounds. With a suction strength of 280 mmHg, the Luna is proven to provide More Milk in Less Time®.

- Very quiet
- Customizable pumping sessions
- Compact size and light weight
- Less guilt, pain, & anxiety
- Convenient nightlight & digital screen



BiliSoft™ 2.0 Phototherapy System

Item#
MAT9765 BiliSoft 2.0 System – Large (10×12)
MAT11559 BiliSoft 2.0 System – Small (8×10)

The BiliSoft™ 2.0 Phototherapy System is an easy-to-use therapy device designed to treat infants with jaundice.

The right touch—and blue light—can make all the difference for infants with jaundice. BiliSoft provides research-backed phototherapy that won't interrupt precious bonding moments. Its soft cover and large pliable pad wrap babies in comfort and allow caretakers to hold, feed or swaddle the infant during treatment.

BiliSoft meets or exceeds the American Academy of Pediatrics (AAP) guidelines for intensive phototherapy to effectively treat hyperbilirubinemia and support even the tiniest patients with more comfort, better bonding and higher quality developmental care.



Accessories

Item#	
MAT9830	Large Pad Cover – 20/PK
MAT9832	Large Pad Cover – 50/PK
MAT9831	Small Pad Cover – 20/PK
MAT9833	Small Pad Cover – 50/PK
MAT13448	BiliSoft Carrying Case
MAT13447	Giraffe Phototherapy Light Meter

- **Intensive Therapy Coverage:** Large treatment area promotes developmental and family-centered care
- **Built to Last:** Durable design for extended use over time
- **Compact & Lightweight:** Easy to handle with a minimal footprint
- **Energy Efficient:** Uses a single high-powered LED instead of multiple bulbs
- **Extended LED Life:** Over 50,000 hours of effective light therapy
- **Ideal for Home Use:** Quiet operation with a long, flexible fiberoptic design

Equipment Covers

Bags for Regulators, Conservers, and Small Items



Item#
EC1012C
10×12"
1.5mil, Clear
1000 bags/box

Covers for Nebulizers and Suction Machines



20×24"
2mil
500 bags/roll

Item#	Color
EC2024C	Clear
EC2024T	Tan Tint

Respiratory Set Up Bags with Draw Tape



Item#
EC1216RB
12×16"
1mil
500 bags/box

Covers for Concentrators



25×15×30"
1.5mil
250 bags/roll

Item#	Color
EC251530C	Clear
EC251530T	Tan Tint

Covers for CPAP/Vents and Walkers



Item#
EC2430C
24×30"
1.5mil, Clear
500 bags/roll

Equipment Covers

Covers for Wheelchairs, Walkers, and Commodes



Item#	Thickness	Size	Color	Quantity
EC302035C	1.5mil	30x20x35"	Clear	200 bags/roll
EC302035T	1.5mil	30x20x35"	Tan Tint	200 bags/roll
EC5045	1mil	50x45"	Clear	250 bags/roll

- Protect equipment and materials from dust and moisture
- Packed on convenient rolls
- Low Density Polyethylene
- Packaged in dispenser carton

Mattress Covers



Item#	Vented	Thickness	Size	Color	Quantity
EC39990C	Yes	1.5mil	39x9x90"	Clear	100 bags/box
EC7252C	No	1mil	72x52"	Clear	200 bags/roll
EC7252T	No	1mil	72x52"	Tan Tint	100 bags/roll


Lancets

Pressure Activated Safety Lancet

Item# **DIALAN001**
21G, 2.0mm Depth
10/pack



NEW



Connect with us



Your sales representatives



Sean O'Donnell

Direct 312-637-6558
Email sodonnell@sunsethcs.com



Madison Gwilliam

Direct 312-533-2459
Email mgwilliam@sunsethcs.com



Evan Berger

Direct 262-788-8948
Email eberger@sunsethcs.com



John Bisceglie

Direct 312-533-2462
Email jbisceglie@sunsethcs.com



Jon Kmiec

Direct 312-533-2482
Email jkmiec@sunsethcs.com



Phil Cosens

Direct 312-533-2470
Email pcosens@sunsethcs.com



Brian Lipson

Direct 312-533-2463
Email blipson@sunsethcs.com



Jonny Deutsch

Direct 872-259-3691
Email jdeutsch@sunsethcs.com



Michelle Chojnacki

Direct 312-533-2486
Email mchojnacki@sunsethcs.com



Chris Crawford

Direct 872-258-1650
Email ccrawford@sunsethcs.com



Pavle Pantovic

Direct 872-282-0751
Email ppantovic@sunsethcs.com



Danielle Dubowski

Direct 312-533-2479
Email ddubowski@sunsethcs.com



Phone 877-578-6738
Email sales@sunsethcs.com
Visit www.sunsethcs.com

Corporate Headquarters
141 W Jackson Blvd Ste 1950
Chicago IL 60604

Central Distribution Center
2725 Alft Ln
Elgin IL 60124

Western Distribution Center
6750 W Van Buren St Ste 100
Phoenix AZ 85043

National Accounts



Gabrielle Murphy

Direct 312-638-8680
Email gmurphy@sunsethcs.com



Shannon St. Clair

Direct 312-533-2449
Email ssclair@sunsethcs.com

International Accounts



Bryan Vigna

Direct 312-533-2455
Email bvigna@sunsethcs.com