



Contracted Supplier



vizient.
Awarded Supplier





Sales Representatives



Egle Sventnickaite
872-258-1654
esventnickaite@sunsethcs.com



Mason Moss
262-788-1862
mmoss@sunsethcs.com



Owen Ruth
312-533-2451
owenr@sunsethcs.com



Cindy Geighes
312-533-2478
cgeighes@sunsethcs.com




Hunter Baudino
262-788-1855
hbaudino@sunsethcs.com



Wyatt Kanemoto
262-788-8941
wkanemoto@sunsethcs.com

Director, Strategic Accounts
Lisa Iannazzone RN, BSN
636-222-3134
lisa@sunsethcs.com

Manager, Product & Business Development
Tom Munar
312-533-2476
tmunar@sunsethcs.com



Manager, Strategic Accounts
Maggie Sinnig
562-556-7137
msinnig@sunsethcs.com

RRT Clinical Specialist
Krystle Bodine BS, RRT
kbodine@sunsethcs.com

Corporate Office

Phone 872-275-3210
141 W Jackson Blvd Ste 1950
Chicago IL 60604

Central Returns/Distribution Center

2725 Alft Ln
Elgin IL 60124

West Returns/Distribution Center

6750 W Van Buren St
Phoenix AZ 85043

Content

NIV

	Page
NIV Masks	4
CPAP/NIV Masks	5
Mask Liners	6
NIV Tubing and Connectors	7

Oxygen Regulators

	Page
Microdial Flowmeter	8
Regulators	8

Respiratory Disposables

	Page
Spirometers	9
Oxygen Cannula and Accessories	9
Oxygen Masks	10
Ventilator Filters	10
Air Entrainment Masks	11
Oxygen Supply Tubing	11
Breathing Circuits	12
Trach Ties	13
HMEs	13
Trach Masks	13
Adapters and Connectors	13

Medication Delivery

	Page
Nebulizers	14-15
Nebulizer Kits	15
Aerosol Masks	15
Inhaler Chamber	15

Closed Suction Catheters

	Page
Deluxe Switch	16
Standard Switch	17

Open Suction

	Page
Suction Catheters	18
Trach Care Kits	18
Suction Catheter Kits and Trays	18-19
Mucus/Specimen Traps	19

Laryngeal Masks

	Page
Silicone Disposable	20
Silicone Reusable	20
Flexible Silicone Disposable	20
Flexible Silicone Reusable	20
PVC Disposable	21
Flexible PVC Disposable	21
Silicone Disposable Dual Seal	21
PVC Disposable Dual Seal	21

Resuscitation Devices

	Page
Resuscitators	22
Breathing Masks	22

Anesthesia

	Page
PeDIA	23



PREMIER

Contracted Supplier

CONTRACT

KI-NS-045	Kiindo Non-Invasive Ventilation
PP-FA-2197	Medical Gas Equipment
PP-OR-2674	Disposable Anesthesia



vizient
Awarded Supplier

CONTRACT

MS1087	Active Humidification + Respiratory Disposables
MS1097	Medication Delivery
MS1049	Non-Invasive Ventilation
MS1049AB	Non-Invasive Ventilation
MS7709	PeDIA Anesthesia Induction
MS1486	Open Suction

PEDIATRIC
ALLIANCE

Ray Pediatric Vented Masks V P Comparable to Nonny masks

Full Face Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	CM061S
Medium	CM061M
Large	CM061L
Extra Large	CM061XL

Replacement Cushions

Size	Item#
Small	CUCM061S
Medium	CUCM061M
Large	CUCM061L
Extra Large	CUCM061XL

Headgear Only

Size	Item#
Small	HG061S
Large	HG061L



SCAN CODE to see
Ray Full Face Mask
on YouTube

Nasal Mask with Headgear, Tube, and Replaceable Cushions

Size	Item#
Kit w/ S&M Cushions	CM060-KITS
Kit w/ L&XL Cushions	CM060-KITL
Small	CM060S
Medium	CM060M
Large	CM060L
Extra Large	CM060XL

Replacement Cushions

Size	Item#
Small	CUCM060S
Medium	CUCM060M
Large	CUCM060L
Extra Large	CUCM060XL

Headgear Only

Size	Item#
Small	HG060S
Large	HG060L



SCAN CODE to see
Ray Nasal Mask
on YouTube

Ray Pediatric Non-Vented Masks V P


Full Face Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Large	CM061L-NV
Extra Large	CM061XL-NV

Nasal Mask with Headgear, Tube, and Replaceable Cushion(s)

Size	Item#
Large	CM060L-NV
Extra Large	CM060XL-NV
Large & XL	CM060-NVKITL

- New non-vented configuration supports use in closed-circuit ventilation systems, including ICU settings
- Designed specifically for pediatric patients, reducing reliance on ill-fitting small adult masks
- Compatible with dual-limb circuits and adaptable to single-limb systems with an exhalation valve
- May support improved therapy delivery through better fit and system flexibility
- SKU reduction by reducing the need for multiple mask types per patient



Available Size Guide

Large and Extra Large sizes help fill a critical pediatric sizing gap, with more potential to avoid intubation.

FinTek allows more adjustability to conform the mask to the face

One-piece molded silicone mask is small, light and comfortable

Headgear with 6-point velcro provides maximum comfort and versatility

Quick release clips

Dual 360 degree swivel elbow

SleepNet MiniMe 2 Pediatric Nasal Mask V P

10/pack

Vented Mask with Headgear

Size	Item#
Extra Small	CM60216
Small	CM60217
Medium	CM60218
Large	CM60221

Non-Vented Mask with Headgear

Size	Item#
Extra Small	CM60256
Small	CM60257
Medium	CM60258
Large	CM60261

novaplus+

Vented Mask with Headgear

Size	Item#
Extra Small	VCM055XS
Small	VCM055S
Medium	VCM055M
Large	VCM055L

Non-Vented Mask with Headgear

Size	Item#
Extra Small	VCM056XS
Small	VCM056S
Medium	VCM056M
Large	VCM056L


HCPCS#
MASK - A7034
HEADGEAR - A7035

Replacement Headgear

Size	Item#
Extra Small/Small	HG027S
Medium/Large	HG027L




MiniMe Pediatric Nasal Mask V



HCPCS# A7034

Headgear Size	Item#	OEM#
Small	CM004S	60213
Medium	CM004M	60214
Large	CM004L	60215

Measure around head above eyebrows for head circumference



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Pediatric SleepWeaver® CPAP Mask V

Item#
CM101325



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.

- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

KidZzzs® Pediatric Mask Liners V P

The first pediatric mask liner in the acute care space designed to prevent skin breakdown

- Made in the USA
- Latex Free



Kit: Ray Full Face Mask with KidZzzs® Mask Liner

	Size	Regular	Antimicrobial
Vented CPAP/VPAP Mask	Small	CM061SKRS	CM061SKR101
	Medium	CM061MKRM	CM061MKR102
	Large	CM061LKRL	CM061LKR103
	Extra Large	CM061XLKRXL	CM061XLKR104
Non-Vented NIV Mask	Large	CM060L-NVKRN	CM060L-NVKR503
	Extra Large	CM060XL-NVKRN	CM060XL-NVKR504

KidZzzs® for Ray Full Face Pediatric CPAP/NIV Mask

Size	25/Box		50/Box		100/Box	
	Regular	Antimicrobial	Regular	Antimicrobial	Regular	Antimicrobial
Small	CAP31KRS-25	CAP32KR101-25	CAP31KRS-50	CAP32KR101-50	CAP31KRS-100	CAP32KR101-100
Medium	CAP31KRM-25	CAP32KR102-25	CAP31KRM-50	CAP32KR102-50	CAP31KRM-100	CAP32KR102-100
Large	CAP31KRL-25	CAP32KR103-25	CAP31KRL-50	CAP32KR103-50	CAP31KRL-100	CAP32KR103-100
Extra Large	CAP31KRXL-25	CAP32KR104-25	CAP31KRXL-50	CAP32KR104-50	CAP31KRXL-100	CAP32KR104-100

KidZzzs® for Ray Nasal Pediatric CPAP/NIV Mask

Size	25/Box		50/Box		100/Box	
	Regular	Antimicrobial	Regular	Antimicrobial	Regular	Antimicrobial
Small	CAP31KRNS-25	CAP32KR501-25	CAP31KRNS-50	CAP32KR501-50	CAP31KRNS-100	CAP32KR501-100
Medium	CAP31KRNM-25	CAP32KR502-25	CAP31KRNM-50	CAP32KR502-50	CAP31KRNM-100	CAP32KR502-100
Large	CAP31KRNL-25	CAP32KR503-25	CAP31KRNL-50	CAP32KR503-50	CAP31KRNL-100	CAP32KR503-100
Extra Large	CAP31KRNL-25	CAP32KR504-25	CAP31KRNL-50	CAP32KR504-50	CAP31KRNL-100	CAP32KR504-100

KidZzzs® for Giraffe Mask

Size	25/Box		50/Box		100/Box	
	Regular	Antimicrobial	Regular	Antimicrobial	Regular	Antimicrobial
Small	CAP31KGS-25	CAP32KG101-25	CAP31KGS-50	CAP32KG101-50	CAP31KGS-100	CAP32KG101-100
Medium	CAP31KGM-25	CAP32KG102-25	CAP31KGM-50	CAP32KG102-50	CAP31KGM-100	CAP32KG102-100
Large	CAP31KGL-25	CAP32KG103-25	CAP31KGL-50	CAP32KG103-50	CAP31KGL-100	CAP32KG103-100



Increase patient compliance by reducing noisy air leaks and skin breakdown



Save time and cost by using a pre-cut liner specifically for your pediatric masks



Increase comfort by adding an absorbent liner to your patient's mask



NIV Tubing **V**

19mm with 22mm cuffs
10/case*

Item#

TUB06 6ft *Individually Packed

SUH0610 6ft

TUB010 10ft



- Sunset's most popular item
- ISO standard 22mm cuffs
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

Standard Elbow, **V** **P** Non-Vented

Item# **CAP009**
10/case



Single Limb Elbow, **V** **P** Non-Vented, AAV

Item# **CAP008**
10/case



Swivel Connector Silencer

Item#

VEN002 Each

VEN002-10 10/pack **P** **V** Contract MS1087



Single Use Exhalation Port **V** **P**

Item# **VEN003**
10/pack



Microflow Microdial Flowmeter P

Assured, repeatable accuracy to enable safe, effective weaning.

2024 Innovation Celebration Winner

What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.

Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure Range	Flow Range*
RES8160048	345 – 500 kPa	A
RES8160050	350 – 550 kPa	E



DISS Inlet • DISS Outlet

Item#	Inlet Pressure Range	Flow Range*
RES8160049	345 – 500 kPa	A
RES8160051	250 – 450 kPa	E

*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
Range A	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
Range E	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

Mini Oxygen Regulator CGA 870 Style

Item# **RES046**

Pediatric • Barbed

0-4 LPM • **Each**

0, $1/32$, $1/16$, $1/8$, $1/4$, $1/2$, $3/4$, 1, 1 $1/2$, 2, 3, 4 LPM

HCPSC# E1353

For use with E size or smaller cylinders:

A, B, C, D, JD, E, M2, M4, M6, ML6, M7, M9, M15, M22, or M24.



Barbed hose outlet for easy connection to the oxygen tubing

"T" Handle for easy mounting — no tools needed

Both Regulators Feature

- Pressure gauge reads to 3000 PSI
- Easy-click flow selector knob
- 2-year warranty

Mini Oxygen Regulator CGA 540 Style

Item# **RES047**

Pediatric • DISS

0-4 LPM • **Each**

0, $1/32$, $1/16$, $1/8$, $1/4$, $1/2$, $3/4$, 1, 1 $1/2$, 2, 3, 4 LPM

HCPSC# E1353

For use with M60 and larger cylinders:

M60, M90, M, MM, M122, M250, H, J, T or other cylinders with CGA 540 style valve.



Threaded DISS outlet

Volumetric Incentive Spirometer

12/case

Item#
RES055P

Pediatric

- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise
- Dual-sided, clear and legible identification of volumes
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration
- Compact and ergonomic design



Pediatric, Infant, and Neonate Cannula V

50/case

Item #	Size	Tubing Length	Connector	Touch	Ear Cushions
RES1100P	Pediatric	None	Standard	Standard	No
RES1107I	Infant	7 ft	Standard	Standard	No
RES1107P	Pediatric	7 ft	Standard	Standard	No
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes
RES1107SI	Infant	7 ft	Standard	Soft	No
RES1107SIU	Infant	7 ft	Universal	Soft	No
RES1107NE	Neonate	7 ft	Standard	Standard	No
RES1107SNE	Neonate	7 ft	Standard	Soft	No
RES1107SP	Pediatric	7 ft	Standard	Soft	No
RES1107SPU	Pediatric	7 ft	Universal	Soft	No
RES1114SP	Pediatric	14 ft	Standard	Soft	No
RES1114SPU	Pediatric	14 ft	Universal	Soft	No



Ear Cushions for Cannula V

Item# **RES020**
Cushions Only
50 Pairs/case



- Add cushions to existing cannula



Secure Cannula Stickers V

Item# **RES080**
25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



High Concentration Oxygen Mask V

50/case

Item #	Size	Tubing Length
RES2202I	Infant	None
RES2202I7	Infant	7 ft
RES2202P	Pediatric	None
RES2202P7	Pediatric	7 ft



Medium Concentration Oxygen Mask V

50/case

Item #	Size	Tubing Length	Connector
RES2203P	Pediatric	None	Standard
RES2100P	Pediatric	7 ft	Standard
RES2203I	Infant	None	Standard
RES2203I7	Infant	7 ft	Standard
RES2203P7	Pediatric	7 ft	Standard



- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Non-Rebreather Oxygen Mask V

50/case

Item #	Size	Tubing Length	Connector	Safety Vent
RES2102I	Infant	7 ft	Standard	No
RES2102IM	Infant	None	Standard	No
RES2102P	Pediatric	7 ft	Standard	No
RES2102PS	Pediatric	7 ft	Standard	Yes
RES2102PSM	Pediatric	None	Standard	Yes
RES2102PSU	Pediatric	7ft	Universal	Yes
RES2102PU	Pediatric	7 ft	Universal	No



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape

Filter for ResMed Astral Ventilator V

Item#
VF003-4
4/pack



Item# **VF003-H**
Filter Housing
Each



Item# **VF003-KIT**
4-Filter Kit
w/Filter Housing
OEM# 27939



Air Entrainment Mask V

With Flex Tube

50/case

Item#	Kit	Tubing Length	Connector
RES2400P	Mask Only	7 ft	Standard
RES2400PK	Yes	7 ft	Standard
RES2400PKU	Yes	7 ft	Universal
RES2400PU	Mask Only	7 ft	Universal



Oxygen Supply Tubing V

Item #	Length	Connector	Bore	Color	Qty
RES3007	7ft	Standard	Channelled	None	50
RES3007SB	7 ft	Standard	Smooth	None	50
RES3007U	7 ft	Universal	Channelled	None	50
RES3015	15 ft	Standard	Channelled	None	50
RES3015G	15 ft	Standard	Channelled	Green	50
RES3025	25 ft	Standard	Channelled	None	25
RES3025G	25 ft	Standard	Channelled	Green	25
RES3025U	25 ft	Universal	Channelled	None	25
RES3030	30 ft	Standard	Channelled	None	20
RES3030G	30 ft	Standard	Channelled	Green	20
RES3040G	40 ft	Standard	Channelled	Green	20
RES3050	50 ft	Standard	Channelled	None	20
RES3050G	50 ft	Standard	Channelled	Green	20

- 6-channel
- PHT/DEHP-free




Universal Connector

Non-Heated Breathing Circuits


For use with Respironics Trilogy
22mm cuffs
30/case

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
VC1001P	1073220 & 1073229	Active	Yes	Pediatric	15mm
VC1004P	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#
VC1001P
Active Pediatric
OEM#
1073220
1073229



To use active circuits with Respironics Trilogy, simply cut ends off pressure and exhalation lines.



Item#
VC1004P
Passive Pediatric
OEM#
1073224
1073233



Passive Heated Breathing Circuit V

Item#
VC1114P
For use with
MR850 Humidifier
10/case



Passive Heated Breathing Circuit With Swivel Exhalation Port V

Item# **VC1003P**
For use with
Philips Respironics
Trilogy
Each

OEM# 1076046,
1073227, 1073235



Includes:

- 1 VC1114P Passive Heated Breathing Circuit
- 1 VEN002 Swivel Exhalation Port

Non-Heated Dual Limb Breathing Circuits V P


30/case

REF	Humidifier Limb	Water Trap
VC2001P	No	No
VC2002P	Yes	No



Dual Limb Active Circuit with Water Trap

Item#
VC1011-10
10/pack
For use with Vyair LTV
OEM# 29662-001




Tracheostomy Ties with Fastener Tabs

20/case

One-piece

- Easy on

Item#
RES029P Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton

Two-piece

- Easy off

Item#
RES030P Pediatric



Heat Moisture Exchanger V

30/case

Item#
RES027NE-10
Inline HME — Neonatal



TV Range (mL)	15-50
Moisture Output (mg H2O/L)	30
Weight (g)	4.5
Resistance (cm H2O)	0.5 @ 10 LPM
Dead Space (mL)	2.5
Connection (mm)	15M,15F

Tracheostomy Mask with O2 Connector V

50/case

Item#
RES2230
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O2 supply tubing

Oxygen Tubing Connector — Swivel V

Item# **RES018-1**
Male/Male
1/pack



Flexible Connector V

50/case

Item# **TRAA004P**
15mm F x 15mm M
Length 4.9-7.5cm
OEM# 3215



Oxygen Bleed-In Adapter V

Item# **RES019**
22mm with 6mm
Rotating Stem
10/pack
RES019-1 1/pack



Connector with 6mm Oxygen Stem V

Item# **RES058**
22mm with 6mm Stem
10/pack



Compressor Nebulizers

Item#

NEB100 
Each

Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

NEB100-R

- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

Sunny the Bear

NEB300BEAR

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

NEB300BEAR-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



Gordon the Fish

NEB300FISH

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

NEB300FISH-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



Handheld Compressor Nebulizer

Item#

NEB200
Each

Included

- Compressor
- Medication cup
- AC power adapter
- Air hose
- Elbow
- Mouthpiece
- Mask
- 4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (approximately 3 microns)
- No replacement filter required



Handheld Mesh Nebulizer

Item#

NEB401
Each

With internal rechargeable battery

Included

- Main unit
- Medicine cup
- Disposable mouthpiece
- Disposable adult mask
- Disposable pediatric mask
- USB cord
- Carrying bag
- Instruction manual

- Small particle size: MMAD 2.2µm*
- Efficient: Low Residual Liquid
- Convenient: Small and Portable
- Quiet: ≤50dB (at one meter)
- Modern: Rechargeable Lithium Battery



Handheld Mesh Nebulizer

Item#
NEB402

With internal rechargeable battery

Particle Size
MMAD < 5µm

Nebulization Rate
0.25ml/min-0.9ml/min

Medication Cup Capacity
10ml (Max)

Device Weight
Approx. 1/4 lb

Device Size
Approx. 1 1/2 x 2 x 3 7/8 inches



Included
Medication cup
USB cord
Mask
Mouthpiece
Portable bag
Manual

Disposable Nebulizer Kit with T-Piece **V**

Item#
RES091U
Disposable
Single Patient,
Universal Connector
50/case



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

Includes:

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube

Aerosol Mask **V**

50/case

Item#
RES2210
Pediatric Size



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

Disposable Nebulizer Kit **V**

Item#
RES092
Pediatric Size
50/case

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube



SperaChamber **V**

Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler

Item#
RES195-10
10/case

RES195-50
50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



SperaChamber with Mask
10/case

Item#	Size
RES195MS	Small (Infant)
RES195MM	Medium (Pediatric)

Mask Only

Item#	Size	Qty
RES195MASK-S	Small (Infant)	Each
RES195MASK-S-10	Small (Infant)	10/box
RES195MASK-M	Medium (Pediatric)	Each
RES195MASK-M-10	Medium (Pediatric)	10/box



Closed Ventilation Suction System

Deluxe Switch (One-handed operation)

Sterile and latex-free

Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port

Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.

Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.



Closed Ventilation Suction System

Standard Switch

Sterile and latex-free

72 Hour Neonatal/Pediatric

With Swivel Connector

Item#	Description
VEN77230505	305mm, 5Fr
VEN77230506	305mm, 6Fr
VEN77230507	305mm, 7Fr
VEN77230508	305mm, 8Fr
VEN77230505M	305mm, 5Fr, MDI Port
VEN77230506M	305mm, 6Fr, MDI Port
VEN77230507M	305mm, 7Fr, MDI Port
VEN77230508M	305mm, 8Fr, MDI Port



24 Hour Neonatal/Pediatric

With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600	305mm, 6Fr
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700	305mm, 7Fr
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800	305mm, 8Fr
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000	405mm, 10Fr
VEN66001000M	405mm, 10Fr, MDI Port

With Swivel Connector

Item#	Description
VEN72430505	305mm, 5Fr
VEN72430506	305mm, 6Fr
VEN72430507	305mm, 7Fr
VEN72430508	305mm, 8Fr

Tracheostomy

With DSE

Item#	Description
VENCSS10-203	203mm, 10Fr
VENCSS10-405	405mm, 10Fr
VENCSS6-305	305mm, 6Fr
VENCSS8-203	203mm, 8Fr
VENCSS8-355	355mm, 8Fr

Swivel Straight Connector

Item#
VEN005-20


Disposable
22mm ID, 15mm OD
20/pack



Suction Catheters - Calibrated V

50/case

Item#	Size	Tip	Valve	Type	Sterile
TRASC106C	6fr	DeLee	Chimney	Coiled	No
TRASC108C	8fr	DeLee	Chimney	Coiled	No
TRASC108S	8fr	DeLee	Chimney	Straight	No
TRASC110C	10fr	Whistle	Y	Coiled	No
TRASC110S	10fr	Whistle	Y	Straight	Yes



Trach Care Kits V

50/case

Sunset Trach Care Kit

Item# **TRATK003**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- Gauze
- 2 Cotton Tips
- Brush,
- Twill Tape
- Trach Dressing

Sunset Mini Trach Kit

Item# **TRATK004**

Includes:

- 3 Pipe Cleaners
- 4 Gauze
- 2 Cotton Tips
- Brush
- Twill Tape
- Dressing

Sunset Trach Care Kit with Forceps

Item# **TRATK005**

Includes:

- 2 Gloves
- Drape
- 3 Pipe Cleaners
- 3 Gauze
- 3 Cotton Tips
- Brush
- Trach Dressing

Suction Catheter Kits with Cup and Gloves V

50/case

Item#	Catheter Size	Tip	Valve	Sterile
TRATK206	6fr	DeLee	Chimney	Yes
TRATK208	8fr	Whistle	Y	No
TRATK210	10fr	Whistle	Y	No



Pediatric Suction Catheter Kits **V**

50/case

Item#

TRATK306

Kit with 6fr Catheter

- Sterile
- Y-Valve DeLee Catheter
- Includes Gloves

TRATK308

Kit with 8fr Catheter

TRATK408

Kit with 8fr Catheter and Cup

Mini Catheter Tray with Gloves **V**

50/case

Item#

TRATK310

Kit with 10fr Catheter

- Coiled Catheter
- Calibrated


Sleeved Suction Coiled Catheter Kits with Cup **V**

50/case

Item#

TRATK410

Kit with 10fr Catheter



Mucus/Specimen Traps **V**

50/case

Item#

RES098C08

Trap with 8fr Catheter

20cc Capacity



Silicone Disposable Laryngeal Masks

Item #	Size
ALM1001	#1
ALM1015	#1.5
ALM1002	#2
ALM1025	#2.5
ALM1003	#3

- Medical grade silicone for super compatibility and patient comfort
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters
- Available with seven sizes for all patient use



Silicone Reusable Laryngeal Masks

Item #	Size
ALM1001R	#1
ALM1015R	#1.5
ALM1002R	#2
ALM1025R	#2.5
ALM1003R	#3

- Reusable Silicone Laryngeal Masks are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Serial no. is printed on each connector for easy tracking
- Additional user card for record of device usage times and patient information is supplied in each pouch



Flexible Silicone Disposable Laryngeal Masks

Item #	Size
ALM3001	#1
ALM3015	#1.5
ALM3002	#2
ALM3025	#2.5
ALM3003	#3



Flexible Silicone Reusable Laryngeal Masks

Item #	Size
ALM3001R	#1
ALM3015R	#1.5
ALM3002R	#2
ALM3025R	#2.5
ALM3003R	#3



- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness

Size	Patient Weight
#1	Up to 5 kg
#1.5	5-10 kg
#2	10-20 kg
#2.5	20-30 kg
#3	30-50 kg

PVC Disposable Laryngeal Masks

Item #	Size
ALM2001	#1
ALM2015	#1.5
ALM2002	#2
ALM2025	#2.5
ALM2003	#3

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



Flexible PVC Disposable Laryngeal Masks

Item #	Size
ALM4001	#1
ALM4015	#1.5
ALM4002	#2
ALM4025	#2.5
ALM4003	#3



Silicone Disposable Dual Seal Laryngeal Masks

Item #	Size
ALM5015	#1.5
ALM5002	#2
ALM5025	#2.5
ALM5003	#3



PVC Disposable Dual Seal Laryngeal Masks

Item #	Size
ALM6003	#3



Size	Patient Weight
#1	Up to 5 kg
#1.5	5-10 kg
#2	10-20 kg
#2.5	20-30 kg
#3	30-50 kg

Child Resuscitator

12/case

Item #	Pop-Off	Mask
RES8110H	Yes	None
RES8111	Yes	#2
RES8112	Yes	#3



Infant Resuscitator

12/case

Item #	Pop-Off	Mask
RES8120H	Yes	None
RES8121	Yes	#0
RES8122	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve

Breathing Mask

20/case

Item #	Size	Hooks
SBC01N	Neonatal (#0)	No
SBC01I	Infant (#1)	No
SBC01P	Pediatric (#2)	No
SBC02N	Neonatal (#0)	Yes
SBC02I	Infant (#1)	Yes
SBC02P	Pediatric (#2)	Yes



PeDIA - Pediatric Device for Induction of Anesthesia V P

Item#

ANEP021Y

With Mask Adapter Elbow


12/pack



2025 Innovation Celebration Winner

Use with
mask

Or use
without
mask



Turn panic into play

For children, undergoing surgery is already unsettling, and the process of induction is often the most traumatizing part. Often kids will kick and scream, believing they're being trapped by the gas or suffocated by the mask. The experiences create negative emotional scarring that can last weeks, months or even years.

Specially designed for children undergoing surgery, the Pediatric Device for Induction of Anesthesia device, or "PeDIA," for short, includes a colorful whistling balloon, providing kids with something fun to focus on and turning a stressful moment into a lighter one.

Whether they're squeezing the balloon and focusing on its whistle or directly putting themselves to sleep by breathing the gas through balloon delivery, the PeDIA balloon captivates children's attention, allowing their active participation in the induction process.

It's a simple way to add a brighter touch to the medical experience, ease anxiety and bring smiles to patients, parents and staff.

This unique approach—presenting an interactive, familiar toy-like device—not only diverts children's focus from the impending procedure but also significantly reduces the anxiety and resistance commonly associated with mask induction. The aim is to decrease the anxiety felt by over 50% of pediatrics as well as parents and staff.

The PeDIA can be used alone as a mask-free device or, for younger children, with an anesthesia mask.

SCAN FOR VIDEOS



