

Repair Parts

OXYGEN CONCENTRATORS



Invacare Platinum 5 / Platinum 10

Item# RP001

Pilot Valve and 1 O-Ring
Each

OEM# 2003052



Item# RP009

Power Cord
Each

OEM# 1108764



Item# RP033

4-Way Rebuilt Valve
Each



Item# RP002

O Ring for Pilot Valve
25/pack

OEM# 1085112



Item# RP010

4-Way Valve Rebuild
Gasket Kit
Each



Item# RP044

Filter Access Door
Each

OEM# 110766



Item# RP003


PE Valve
Each

OEM# 1143578



Item# RP011

4-Way Valve End Caps
2/pack



Item# RP045

Sieve Bed Filter Grid
Each

OEM# 1109520



Item# RP004

On/Off Switch
Each

OEM# 1139950



Item# RP014

Concentrator Capacitor 15MFD
Each

OEM# 2003019



Item# RP081

Piston Cup Seals for
Thomas Compressor
2650/2660
50/pack



Item# RP005

Hour Meter
Each

OEM# 2002939



Item# RP016

Concentrator Capacitor 30MFD
*For Platinum 10
Each

OEM# 1128740



Item# RP083

Brass O2 Outlet
10/pack



Item# RP006

O2 Conc Hose Barb
Fitting 1/8" Conc Port
Each

OEM# 1104151



Item# RP019

Pressure Regulator
Each



Item# RP098

Flow Meter
*For Platinum 10
Each



Item# RP031

Cabinet Clips
2/pack



Repair Parts **OXYGEN CONCENTRATORS**



Invacare Perfecto2 / Perfecto2V

Item# RP001
Pilot Valve and
1 O-Ring
Each
OEM# 2003052



Item# RP002
O Ring for Pilot Valve
25/pack
OEM# 1085112



Item# RP003
PE Valve
Each
OEM# 1143578



Item# RP004
On/Off Switch
Each
OEM# 1139950



Item# RP005
Hour Meter
Each
OEM# 2002939



Item# RP006
O2 Conc Hose Barb
Fitting 1/8"
Conc Port
Each
OEM# 1104151




Item# RP009
Power Cord
Each
OEM# 1108764



Item# RP010
4-Way Valve Rebuild
Gasket Kit
Each



Item# RP011
4-Way Valve End Caps
2/pack



Item# RP013
Concentrator Capacitor
10MFD
*for Perfecto2
Each



Item# RP015
Concentrator Capacitor
22MFD
*for Perfecto2V
Each



Item# RP018
Compressor
Rebuild Kit
*for Perfecto2V
Each
OEM# 1191916



Item# RP019
Pressure Regulator
Each



Item# RP031
Cabinet Clips
2/pack



Item# RP033
4-Way Rebuilt Valve
Each



Item# RP036
Complete
Cabinet Set
3/case
OEM# 1153341



Item# RP043
Cabinet Access Door
*for Perfecto2
Each
OEM# 1161831



Item# RP045
Sieve Bed Filter Grid
*for Perfecto2
Each
OEM# 1109520



Item# RP046
PVC F-Tube
Each
OEM# 1143446



Item# RP047
Control Panel
Each
OEM# 1143469



Item# RP081
Piston Cup Seals for
Thomas Compressor
2650/2660
50/pack



Item# RP083
Brass O2 Outlet
10/pack



Item# RP097
Heat Exchanger
Assembly
Each
OEM# 1186169



Item# RP101
Sticker Labels
for Perfecto2 Shroud
10/pack



Item# RP102
Sticker Labels
for Perfecto2V Shroud
10/pack



Repair Parts **OXYGEN CONCENTRATORS**



Invacare HomeFill

Item# **RP004**
Power Switch
Each
OEM# 2001134



Item# **RP007**
O2 Conc Fitting Receptacle Port
Each
OEM# 1104152



Item# **RP008**
2000psi Pressure Switch
Each
OEM# 1089467



Item# **RP013**
Concentrator Capacitor
10MFD
Each



Item# **RP026**
Cylinder Dust Cap
10/pack



Item# **RP027**
Compressor Dust Cap
10/pack
OEM# 1156644



Item# **RP028**
O-Ring for
Homefill Tank
10/pack



Item# **RP031**
Cabinet Clips
2/pack



Item# **RP042**
Male, QDV,
High Pressure, INTL
with Dust Cover
Each



Item# **RP104**
Exhaust Fan
Each
OEM# 255506,
2000817



React Stratus 5

Item# **RP071**
Compressor
Each



Item# **RP072**
Compressor Seal
—Piston Ring
Each



Item# **RP073**
Regulator Valve
Each



Item# **RP074**
Cooling Fan
Each



Item# **RP076**
Power Switch
25/pack



Item# **RP077**
Pneumatic Valve
Each



Item# **RP078**
Power Cord
Each



Item# **RP109**
Capacitor
Each

Item# **RP110**
Pneumatic Valve Normal Rebuild Kit
Each
OEM# PRTSTR0340

Item# **RP111**
Air Intake Window
Each
OEM# PRTSTR019

Item# **RP112**
Case—Front
Each
OEM# PRTSTR012

Item# **RP113**
Case—Back
Each
OEM# PRTSTR013

Repair Parts **OXYGEN CONCENTRATORS**



Drive/DeVilbiss 525 / 1025

Item# **RP004**
On/Off Switch
Each
OEM# 1139950



Item# **RP023**
Circuit Breaker
Each
OEM# 525D-613



Item# **RP092**
Power Cord Strap
10/pack
OEM# MC29D-657



Item# **RP018**
Compressor
Rebuild Kit
*For 525 Only
Each
OEM# 1191916



Item# **RP024**
Oxygen Outlet Port
5/case
OEM# 525DD-606



Item# **RP093**
Rebuilt Drive
Rotary Valve
Each



Item# **RP020**
Flow Meter
Each
OEM# 505DZ-607



Item# **RP032**
Pressure Regulator
8.5 PSI
Each



Item# **RP094**
Rebuild Kit for
Rotary Valve
Each



Item# **RP021**
Cooling Fan
Each
OEM# 525DD-634



Item# **RP048**
Power Cord
Each



Item# **RP095**
Spacers for Rebuilt
Rotary Valves
10/pack



Item# **RP022**
Exhaust Muffler
Each
OEM# 515A-705



Item# **RP049**
DS SSP Capacitor
*for 1025DS
Each



Item# **RP114**
Wiring Harness
*For 525 Only
10/pack
OEM# 525D-621



Item# **RP091**
Capacitor SSP
*For 525
Each



Inogen G3 / G5

Item# **RP103**
Filter Cage Replacement
*For G3
2/pack



Item# **RP105**
Filter Cage Replacement
*For G5
2/pack



Repair Parts **OXYGEN CONCENTRATORS**



Rhythm LM5BA & LM5CA, Lifestyle LM5A, Jumao JMC5A, Direct Supply Attendant 8X899 & A5LC-1, Medacure OC5L

Item# RP004
On/Off Switch
Each
OEM# 1139950



Item# RP005
Hour Meter
Each
OEM# 2002939



Item# RP019
Pressure Regulator
*For Rhythm LM5BA
Each



Item# RP023
Circuit Breaker
*for Rhythm LM5BA
Each



Item# RP037
Muffler
Each



Item# RP038
PE Valve
Each



Item# RP039
Control Valve
Each



Item# RP040
Power Cord
Each



Item# RP041
Casters
4/pack



Item# RP075
Cup Seal for Compressor
Rebuild Kit
Each



Item# RP079
Pressure Regulator,
7 PSI
Each



Item# RP080
O-Rings for
Pressure Regulator
50/pack



Item# RP084
Power Cord
*for Rhythm LM5BA
Each



Item# RP085
Power Switch
*for Rhythm LM5BA
Each



Item# RP086
Cooling Fan
*for Rhythm LM5BA
Each



Item# RP087
Solenoid Valve
*for Rhythm LM5BA
Each



Item# RP088
Oxygen Sensor
*for Rhythm LM5BA
Each



Item# RP089
Circuit Board
*for Rhythm LM5BA
Each



Item# RP090
Display Screen
*for Rhythm LM5BA
Each



Item# RP096
Cabinet with Door
*for Direct Supply Attendant
2/box



Item# RP106
Filter Door
*for Direct Supply Attendant
Each
OEM# O-FYM13



Item# RP107
Air Guide Cover
*for Direct Supply Attendant A5LC-1 (Prior to 2023) and Rhythm Lifestyle LM5A (2021-2022)
Each



Repair Parts **OXYGEN CONCENTRATORS**



Rhythm P2

Item# **RP117**
DC Power Cord
Each



Item# **RP118**
AC Power Cord
Each



Item# **RP100**
Compressor Rebuild Kit with
Top End and Bottom Flapper
Assembly + Cup and Sleeve
Each



Respironics EverFlo (Post-2013)

Item# **RP006**
O2 Conc Hose Barb
Fitting 1/8" Conc
Port
Each
OEM# 1104151



Item# **RP012**
Flow Meter
Each
OEM# 0209184



Item# **RP035**
Pressure Regulator
Adapter Kit
Each
OEM# 1038824



Item# **RP034**
Pressure Regulator
Adapter With Hose
Each
OEM# 1038838



Item# **RP082**
SMC Pilot and
Body Valve
Each



Item# **RP099**
Sieve Muffler
Insulation
Each



Respironics Trilogy

Item# **RP025**
O2 Quick Connect
L-Shape Adapter
2/pack



Item# **RP025S**
O2 Quick Connect
Straight Adapter
2/pack



**RP025S is Also
Compatible with:**

- Löwenstein **LUISA**
- ResMed **Astral**
- React Health **V*Home**



Respironics Millennium M10

Item# **RP108**
Flow Meter
Each
OEM# 1017847



Respironics SimplyGo Mini

Item# **RP115**
DC Power Cord
Each



Item# **RP116**
AC Power Cord
Each



Repair Parts **OXYGEN CONCENTRATORS**



CAIRE NewLife

Item# RP400

Felt F10 Filter
1.90"OD x 0.70"ID
Each
OEM# MI161-2



Item# RP402

Tie Wrap, 14.5"
Each
OEM# TW004-1



Item# RP403

Sieve Bed Assembly
Each
OEM# BE186-2S



Item# RP405

Clear Tubing 0.3125" OD
Each
OEM# TU198-064

Item# RP407

Packaging Carton
Each
OEM# PA202-1

Item# RP408

Hour Meter
Each
OEM# HM009-2



Item# RP410

1/4" Regulator
Each
OEM# RG088-1



Item# RP411

Solenoid Valve, 1/4"
Each
OEM# VA001-1



Item# RP412

Brass Fitting Adapter
Each
OEM# F0005-2



Item# RP413

Flowmeter
Each
OEM# FM056-1



Item# RP414

Solenoid Valve, 1/8"
Each
OEM# VA003-1



Item# RP415

Fan
Each
OEM# FN050-1




Item# RP416

Capacitor 20UF/370VAC
Each
OEM# CC012-8



Item# RP417

Cabinet Recepticle 120/60
Each
OEM# CA007-11



Item# RP418

Roller Cabinet Base
Each
OEM# CA007-11



Item# RP419

Neoprene Washer 1 5/8"
Each
OEM# WA025-2



NewLife 5L and 10L Parts Compatibility

Item#	Description	5L Compatible	10L Compatible
RP400	Felt F10 Filter 1.90"OD x 0.70"ID	NewLife (3 LED)	NewLife Intensity (1 LED); (3 LED)
RP402	Tie Wrap, 14.5"	Yes	Yes
RP403	Sieve Bed Assembly	No	NewLife Intensity (1 LED); (3 LED)
RP407	Packaging Carton	All LED models	All LED models
RP408	Hour Meter	Yes	Yes
RP410	1/4" Regulator	No	NewLife Intensity (1 LED); (3 LED)
RP411	Solenoid Valve, 1/4"	Yes	Yes
RP412	Brass Fitting Adapter	No	NewLife Intensity (1 LED); (3 LED)
RP413	Flowmeter	No	NewLife Intensity 10 (1 LED); (3 LED)
RP414	Solenoid Valve, 1/8"	No	(1 LED); (3 LED)
RP415	Fan	No	(1 LED); (3 LED)
RP416	Capacitor 20UF/370VAC	Yes	Yes
RP417	Cabinet Recepticle 120/60	Yes (NewLife Elite, 3 LED)	Yes (NewLife Intensity 10, 3 LED)
RP418	Roller Cabinet Base	NewLife Elite (3 LED)	NewLife Intensity (3 LED)
RP419	Neoprene Washer 15/8"	Yes	Yes

Repair Parts **OXYGEN CONCENTRATORS**




CAIRE Eclipse 5

Item# RP420
Power Cartridge
Battery
Each
OEM# 7082-SEQ



Item# RP421
Spare Loptus
Each
OEM# SP20626702




Item# RP422
DC Power Supply
Each
OEM# 4952-SEQ



Item# RP423
Universal Cart
Each
OEM# 5991-SEQ



Item# RP424
Oxygen Outlet Filter
Each
OEM# 21360366




Item# RP425
Power Supply Kit - US
Each
OEM# 208523265



Item# RP426
In-Line Filters
Each
OEM# 8069-SEQ



Item# RP427
Compressor Assembly
Each
OEM# SP5999-SEQ



Item# RP428
AC Input Cable - Type A
Each
OEM# 4997-SEQ



Item# RP429
Membrane Panel
For 3-LED Eclipse
Each
OEM# SP20618461




Item# RP430
User Manual - English
Each
OEM# 21182453




Item# RP431
Spare Membrane Kit
for 2-LED Eclipse
Each
OEM# SP21072187




Item# RP432
Power Manager and
Battery Bridge Circuit
Assembly
*for 2-LED Eclipse
Each
OEM# SP5931-2-SEQ




Item# RP433
Power Manager and
Battery Bridge Circuit
Assembly
*for 3-LED Eclipse
Each
OEM# SP5932-4-SEQ



Item# RP434
Top Case Assembly
*for 2-LED Eclipse
Each
OEM# SP21126965




Item# RP435
Control System PCB
*for 2-LED Eclipse
Each
OEM# SP21122262-31




Item# RP436
Spare Top Case Assembly
*for 3-LED Eclipse
Each
OEM# SP20623840




Item# RP437
Control System Printed
Circuit Board
Each
OEM# SP4043-9-SEQ



Item# RP438
Bottom Case Assembly
Each
OEM# SP21126771



Item# RP439
Proportional IQ Valve
Each
OEM# 5293-SEQ



Item# RP457
Preventative
Maintenance Kit
Each
OEM# 5022-SEQ



Repair Parts **OXYGEN CONCENTRATORS**



CAIRE Freestyle

Item# **RP440**
Car Charging
Power Cord
Each
OEM# CD041-1



Item# **RP441**
Carry Bag
Each
OEM# MI396-1



Item# **RP442**
Assembly Beds
Each
OEM# BE566-1S



Item# **RP443**
Gross Particle Filter —
Right
Each
OEM# FI226-1



Item# **RP444**
Gross Particle Filter —
Left
Each
OEM# FI227-1



Item# **RP445**
Compressor Intake Filter
Each
OEM# FI224-1



Item# **RP446**
O2 Outlet Filter
Each
OEM# FI225-1



Item# **RP447**
Power Supply
12V/120W
Each
OEM# PW036-1S



Item# **RP448**
16-Cell Battery
Each
OEM# BT035-1S



Item# **RP449**
O2 Outlet Filter Gasket
Each
OEM# GS049-1



Item# **RP450**
Cabinet Cover Plug —
DC Power
Each
OEM# CA361-1



Item# **RP451**
Backpack
Each
OEM# MI459-1



Item# **RP452**
Cabinet Port Cover
Each
OEM# CA359-1



Item# **RP453**
Straps
Each
OEM# MI397-1



Item# **RP454**
Top Cap Kit
Each
OEM# KI612-1S



Item# **RP455**
Replacement
Battery
Each
OEM# BT034-2S



Item# **RP456**
Accessory Bag
Each
OEM# MI320-2



Accessory Parts **OXYGEN**

Item# **RP017**
O2 Tank Toggle & Screw
10 Sets



Item# **RP029**
Filling Adaptor
CGA 870 Yoke
1/4 NPT with T/Bar
Each



Item# **RP030**
Filling Adaptor CGA
870 Yoke
1/2 NPT with T/Bar
Each

