




Caire Companion C1000




Repair Parts **LIQUID OXYGEN**




Caire Companion C1000 — Manifold

Vent Valve



Item# RP149
Vent to Fill Lever
Assembly
OEM# B-776439-00

Secondary Relief Valve




Item# RP139
Secondary Relief
Valve - 30 PSI
OEM# B-775082-00

Primary Relief Valve




Item# RP068
Primary Relief Valve
Fitting, 22 PSI
OEM# B-775078-00

Item# RP199
Primary Relief Valve
Silencer
OEM# B-776827-00



Female Fill Connector



Item# RP058
Tube Fill
OEM# B-775061-00

Item# RP198
Female QC Nut
OEM# B-775261-00

Item# RP060
Quick Connect
Adapter
OEM# B-775263


Quick Connect
Assembly



Repair Parts **LIQUID OXYGEN**



Caire Companion C1000T




Repair Parts **LIQUID OXYGEN**




Caire Companion C1000T

Warming Coil / Expansion Line




Heat Exchanger



Item# **RP063**
Vent / Fill Assembly
OEM# B-775243-00

Vent Valve






Caire Companion C1000T

Female Fill Connector

Item# **RP198**
Female QC Nut
OEM# B-775261-00



Item# **RP060**
Quick Connect
Adapter
OEM# B-775263




Item# **RP058**
Tube Fill
OEM# B-775061-00



Primary Relief Valve


Item# **RP192**
Relief Valve Fitting
OEM# B-775078-00

Item# **RP199**
Relief Valve Silencer
OEM# B-776827




Secondary Relief Valve


Item# **RP139**
Secondary Relief
Valve - 30 PSI
OEM# B-775082-00




Caire Companion C21/C41 Reservoirs




Caire Companion C21/C41 Reservoirs (continued)




Caire Companion C31 & C41—EU Reservoirs




Caire Companion C31 & C41—EU Reservoirs





Caire Hi Flow Stroller® — Side Fill




Caire Hi Flow Stroller® — Top Fill



Top Fill Manifold (Design after 1.27.2020)




Caire HELiOS Marathon H850



Caire HELiOS Plus H300

Base Bottle & Components






Caire HELiOS H36 & H46 Reservoirs



Item# RP180
Electronic Level Gauge
OEM# B-702207-LAN

Bracket Assembly


Item# RP157
Pressure Gauge Lens
OEM# B-778807-00



Item# RP167
Ogiva Clamp
OEM# B-775794-00

Item# RP167
Ogiva Clamp
OEM# B-775794-00


Item# RP190
Vent Valve Base
OEM# 20580870



Repair Parts **LIQUID OXYGEN**




Caire HELiOS U36 & U46 Reservoirs



Caire G4 Liberator® Reservoir


Configurable Parts

Design before 4/1/2020




Caire G4 Liberator® Reservoir

Economizer
Design before 4/1/2020




Economizer
Design after 4/1/2020



Caire Sprint and Stroller

Case



Repair Parts **LIQUID OXYGEN**



Caire **HELIOS**

Item# RP056M

Lip Seal - Metal
Each

OEM# B-775260-00



Item# RP056R

Lip Seal - Rubber
Each

OEM# B-775260-00



Item# RP066

Base Reservoir Contents
Gauge Tubing
2/pack

OEM# B-775856-SV, PB21,
PB41, H46



Item# RP069

Secondary Relief Valve
30 PSI
Each



Caire **C500/C550/C550i**

Item# RP148

Moisture Pad

OEM# B-777904-00

Item# RP156

Rear Side Cover

OEM# B-777607-00

Item# RP158

Front Side Cover for C-500

OEM# B-777608-00

Item# RP162

Top Bezel for C-500

OEM# B-777606-00



Puritan Bennett 21L/41L Base Liquid Oxygen Vessel

Item# RP057

41L Contents Gauge
Rebuilt
Each

OEM# B-775847-SV



Item# RP066

Contents Gauge Tubing
2/pack

OEM# B-775856-SV, PB21,
PB41, H46



Item# RP069

Secondary Relief Valve
30 PSI
Each



Accessory Parts **OXYGEN**

Item# RP017

O2 Tank Toggle & Screw
10 Sets



Item# RP029

Filling Adaptor
CGA 870 Yoke
1/4 NPT with T/Bar
Each



Item# RP030

Filling Adaptor CGA
870 Yoke
1/2 NPT with T/Bar
Each

